

A Crowning Achievement!



Premier® Implant Cement™

Designed for Implant Crown Cementation

Premier® Implant Cement™ is designed specifically for implant crown cementation. It's the ultimate solution for a versatile implant cement that provides strong, secure retention yet enables easy, non-destructive retrieval if needed at a later date. The elastomeric cement flexes under force, absorbs shock and creates a tight marginal seal. You'll appreciate the predictable and easy application with no etch, primer or bonding agent required.

Premier brings meaningful, innovative dental material and device solutions you can rely on. Every time!


premier® Inspired solutions
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • www.premusa.com

Premier® Implant Cement™ is easy to use



1 Check the fit. Complete all finishing procedures.



2 Rinse and dry abutment and interior of restoration.



3 Bleed syringe. Mount new automix tip and then bleed the syringe again. This assures proper flow of base and catalyst.




4a Using a microbrush, place a thin film of cement along the inside wall of the cervical third, making sure to completely coat the margin with cement. You have 1.5 minutes of working time at room temperature.



4b Excess cement should be avoided. Too much cement can cause future restoration retrieval to be difficult to impossible.[†]



5 Firmly seat the crown and hold in place until initial set (gel phase) occurs in 2–2.5 minutes. Safely remove excess cement before final set. It is imperative to remove all excess cement.



6 Final cement is 4 - 5 minutes, intraorally.

[†]Any excess cement (regardless of chemical make-up or brand) that remains in the sulcus around an implant could lead to peri-mucositis, peri-implantitis or implant failure if not removed immediately.

- 3001450 Implant Cement - Standard Pack
1-5ml Automix Syringe, 10 Mixing Tips and 1 Mixing Pad
- 3001455 Implant Cement - Value Pack
3-5ml Automix Syringes, 25 Mixing Tips
- 3001451 Automix Dual-barrel Syringe Tips (25)

Properties

Working Time	1.5 minutes at room temperature
Gel-phase	2-2.5 minutes at 98°F (37°C)
Set Time	4-5 minutes at 98°F (37°C) from the start of mix
Film Thickness	15 µm



Secure Retention and Retrieval
in one unique step!

"Outstanding Products"™



Made from graphite reinforced resin that will not scratch implant surfaces.



Premier® PerioWise®

- UNC 12
- 9006120 (3-pack)
- 9006121 (6-pack)
- 9006122 (12-pack)



Premier® Implant Scalers

- | | |
|--|--|
| Universal Implant Scaler
(Columbia 4L/4R) | Facial Implant Scaler
(Goldman-Fox 5) |
| 9061401 (2-pack) | 9061402 (2-pack) |
| 9061400 (5-pack) | 9061403 (5-pack) |

- | | |
|--|---|
| new! 137 Implant Scaler
(Anterior) | new! 204 Implant Scaler
(Posterior) |
| 9061421 (2-pack) | 9061423 (2-pack) |
| 9061420 (5-pack) | 9061422 (5-pack) |

premier® Inspired solutions for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • www.premusa.com • E-mail: dentalinfo@premusa.com

*An independent non-profit dental education and product testing foundation, December 2004 & 2005