



Products and Procedures for Better Dentistry



Welcome to Premier!

We are so appreciative of your choice to look to us for your important product needs. Our innovative offerings are developed and constructed with you – the dental leaders and champions of oral health in mind – and are driven by our rich family history and understanding of the global dental market place.

Your trust in Premier has allowed us to succeed with strength, and given us the opportunity to bring you technology in brands such as Enamel Pro®, Traxodent®, Big Easy® Ultralite™ and Solo.

We look forward to supplying you with our market-leading products, and for your continued partnership in the years to come.

Sincerely yours,

Julie Charlestein

President

Premier® Dental Products Company



U.S. inquiries should be directed to:

Premier® Dental Products Company

1710 Romano Drive, P.O. Box 4500, Plymouth Meeting, PA 19462, USA Customer Service: 1-888-670-6100 • Fax: 1-888-627-6160 dentalinfo@premusa.com • www.premusa.com

Overseas inquiries should be directed to:

Premier® Dental Products Company International Division

1710 Romano Drive, P.O. Box 4500, Plymouth Meeting, PA 19462, USA 610-239-6000 • Fax: 610-239-6171 • E-mail: intlinfo@premusa.com

Canadian inquiries should be directed to:

Premier® Dental Products Company

Toll Free: 1-888-773-6872 • Fax: 610-239-6171 E-mail: dentalinfo@premusa.com

Veterinary inquiries should be directed to:

Premier® Dental Products Company

Toll Free: 1-888-670-6100 • Fax: 1-888-627-6160 E-mail: dentalinfo@premusa.com

Inquiries for medical products or requests for information should be directed to:

Premier® Medical Products

1710 Romano Drive, P.O. Box 4500, Plymouth Meeting, PA 19462, USA 610-239-6000 • 1-888-PREMUSA (773-6872) • Fax: 610-239-6171 E-mail: medinfo@premusa.com




Enamelon®

Preventive Treatment Gel

Everyday Relief. Everyday Protection.

- Treats sensitivity
- · Relieves the discomfort of dry mouth tissues
- · Helps to prevent caries & gingivitis
- Non-irritating to patients with dry mouth¹
- · Contains no SLS, gluten or dyes
- 0.40% stabilized stannous fluoride (970 ppm F)
- Optimized with calcium and phosphate ions
- · Non-abrasive formula for daily use
- Non-prescription, safe for patients of all ages*



Enamelon® contains stannous fluoride along with calcium and phosphate ions that have proven effective for tubule occlusion. By blocking the exposed tubules with a combination of stannous (tin) ions along with fluoride, calcium and phosphate salts, a physical barrier is created that covers open dentinal tubules to prevent external stimuli from causing pain or discomfort.

Start your patients on Enamelon today.

9007285 Enamelon® Preventive Treatment Gel 4.0 oz (113 g), Clean Mint, 12/case

Enamelon®

Fluoride Toothpaste

Enamelon® Toothpaste contains the same clinically proven active ingredient optimized with calcium and phosphate ions as the preventive treatment gel, but in an 1150 ppm F formulation with a low abrasive for everyday use.

9007280 Enamelon® Toothpaste 4.3 oz (122 g), Mint Breeze, 12/case

9007290 Enamelon® Toothpaste 0.75 oz (21.2 g), Mint Breeze, 24/case

Not available in Canada.





Strengthen Your Smile!™







Scalers

Implant Scalers

- Thin, slightly flexible and safe for use on implants
- Autoclavable (273°F/134°C)

See full listing of implant scalers on pages 14-15.



Elevators

Cryer

• Ideal to loosen or remove dental roots, root fragments or teeth from bony sockets

See full listing of elevators on pages 71-75.

- · Paired with left and right small flag tips for excellent access
- · Large diameter handle to ensure a secure grip





Ultrasonic Inserts

Big Easy® Implant Ultrasonic Inserts

- PEEK tip is permanently attached, no special tool or replacement tips are required
- · Effectively removes biofilm and calculus
- Safe to use on implants and orthodontic brackets
- Ideal for sensitive patients
- Will not scratch metal or damage restorative surfaces
- Big Easy ergonomic comfort grips
- Compatible with most magnetostrictive hand pieces
- · Available in 25K and 30K frequencies

PEEK is considered an advanced biomaterial used in medical implants. It does not degrade with exposure to steam, and is heat, chemical and radiation-resistant. PEEK can withstand the chemical, thermal, and mechanical stresses present during the sterilization process. PEEK is non-abrasive and proven robust for long-lasting insert tips.





Big Easy® Ultrasonic Inserts

- Ergonomically designed large 0.43" (1.09 cm) diameter cushion-grip handle for maximum comfort and decreased fatigue
- Compatible with most magnetostrictive hand pieces



Resin Handle Ultrasonic Inserts

- Ergonomically designed large 0.40" (1.0 cm) diameter handle
- · Compatible with most magnetostrictive hand pieces





Hygiene/Perio



Enamel Pro®

The ONLY Prophy Paste Formulated to Deliver ACP

Strengthen their smiles!

Enamel Pro® offers innovative remineralizing technology through Amorphous Calcium Phosphate (ACP).²

- Enamel Pro prophy paste delivers 31% more fluoride uptake into the enamel than prophy paste without ACP³
- Provide a true benefit to your patients while your prophy provides a lustrous smile
- · Gluten-free



9007620	RaspberryMint Fine
9007621	RaspberryMint Medium
9007622	RaspherryMint Coarse

See full listing of Enamel Pro on page 88.



Restorative



Spot-On Etch no-slump viscosity will stay put. Only etch the areas you desire!



NEW! Spot-On™ Etch

37% Phosphoric Acid Dentin/Enamel Etch Gel

Exact placement that stays put!

Spot-On™ Etch is a 37% high viscosity phosphoric etch gel for all your routine and selective-etch procedures. It's designed with enhanced viscosity that will prevent it from dripping or slumping. Spot-On Etch easily flows through a small applicator tip, washes off readily and resists drying out.

- Ideal for routine etching and selective-etch procedures
- Blue color provides a strong contrast to enamel and dentin
- No slump formula
 Rinses clean without leaving residue
- Easy handling
 Gluten-free

3001460 Spot-On™ Etch Bulk

One 12g jumbo filled syringe with luer lock adapter, Ten 1.2ml empty syringes, Fifty 1.2ml disposable applicator tips 3001461 Spot-On™ Etch Intro

Five 1.2ml filled syringes, Twenty bent disposable applicator tips



Cleaners



NEW! Vacu Blast™

A Powerful Enzymatic Liquid Evacuation System Cleaner

Refreshing eucalyptus scent!

Its metered-dose bottle makes it easy and economical to use in each operatory to keep evacuation systems operating at peak performance.

- Compatible with amalgam separators
- · Anti-corrosive
- Non-foaming
- Neutral pH
- Use with or without an atomizer
- Environmentally friendly
- · Phosphate- and chlorine-free

9011107 12 - 32oz bottles (946ml)

Each bottle yields 64 treatments

32 oz bottle available individually through your dealer.



CE




Instruments

| Instrumente | Instrumentos | Instruments |

Diagnostic Diagnose Diagnostico Diagnostique	Page 2-9
Hygiene Hygiene Higiene Hygiène	10-21
Operative Konservierende Behandlung Operatoria Operatoire	22-47
Endodontic Wurzelinstrumente Endodoncia Instruments à Nerfs	48-53
Perio/Oral Surgery Parodontal/Mündliche Chirurgie Periodoncia/Cirugía Oral Chirurgie Parodontale et Buccale	54-84
Cassettes and Instrument Care Pflege der Kassetten & Instrumente Cassettes y Mantenimiento de Instrumente Le Soin de Cassettes et Instruments	85-86 os



Premier's state-of-the-art instrument plant is located in Philadelphia, PA USA.



Diagnostic | Diagnose | Diagnostico | Diagnostique |

Explorers Sonden Exploradores Explorateurs	Page 3-4	
Mirror Handles Mundspiegelgriffe Mangos de Espejo Manches à Mirroirs	4	
Probes Sonden Sondas Sondes	5-8	
Dressing Pliers Zahnärztliche Pinzetten Pinzas Precelles a Pansement	8-9	

| Sonden | Exploradores | Explorateurs | Explorers | available in economical 12-packs, where indicated









Explorers | Sonden | Exploradores | Explorateurs |

available in economical 12-packs, where indicated

Single-end



Mirror Handles | Mundspiegelgriffe | Mangos de Espejo | Manches à Mirroirs |



| Sonden | Sondas | Sondes | Probes




 $^{{\}it *An independent non-profit dental education and product testing foundation - January 2001.}$






Probes | Sonden | Sondas | Sondes |



| Sonden | Sondas | Sondes | Probes






CE

Probes | Sonden | Sondas | Sondes |






Dressing Pliers | Zahnärztliche Pinzetten | Pinzas | Precelles a Pansement |

Dressing



| Zahnärztliche Pinzetten | Pinzas | Precelles a Pansement | Dressing Pliers





Hygiene

| Hygiene | Higiene | Hygiène |

C	C 2	\sim	rc
.)	l a	_	רו

Page 11-13

| Zahnreinigungsinstrumente | Raspadores de Sarro | Instrument à Detartrer |

Instrument Sharpener

12

| Instrumentenschleifer | Afilador de Instrumentos | Affûteur pour Instruments |

Curettes

14-20

| Küretten | Curetas | Curettes |

Ultrasonic Inserts

21

| Ultraschall - Zahnsteinentfernung Instrumente | | Insertos Ultrasonicos | Pointes Ultrasoniques |

| Zahnreinigungsinstrumente | Raspadores de Sarro | Instrument à Detartrer | Scalers





Big Easy® Ultralite™ (BEUL) 1/2″ diam.,16g



Big Easy® (BE) 7/16" diam. 23a



Light Touch™ Round (LT) 5/16" diam. 21g




Light Touch™ Octagon (LTO)1/4" diam.
18g



^{*} Big Easy Scalers & Curettes, Townie Choice Awards, 2003.

Scalers | Zahnreinigungsinstrumente | Raspadores de Sarro | Instrument à Detartrer |



Instrument Sharpener | Instrumentenschleifer | Afilador de Instrumentos | Affûteur pour Instruments |



Now, you can quickly sharpen scalers and curettes correctly – without damaging the cutting surface.

9048200 - Autoclavable D.I.S.C. Instrument Sharpener (Metal)

9048093 - D.I.S.C. Replacement Stone

Easy-to-follow instructions for all instruments included with D.I.S.C.

Extend Instrument Life...Enhance Your Efficiency



Turn the protractor dial to the appropriate dial setting (Example: "G" for any Gracey 1-10).



Keeping the wrist stiff, use a pull stroke with moderate pressure (repeat 5 or 6 times).



Hold the instrument with the handle parallel to the Instrument Reference Lines. Position the working end towards the center.




On a curette, roll the instrument as you finish your pull stroke to maintain the rounded toe.

D.I.S.C. Sharpening stones, autoclavable and available separately, never need oil or water as is required with Arkansas stones.

ϵ

| Zahnreinigungsinstrumente | Raspadores de Sarro | Instrument à Detartrer | Scalers



Big Easy® Ultralite™ (BEUL) 1/2″ diam.,16g



Big Easy® (BE) 7/16" diam. 23g



Light Touch™ Round (LT)5/16" diam.
21g




Light Touch™ Octagon (LTO)1/4" diam.
18g



Curettes | Küretten | Curetas | Curettes |

Columbia




| Küretten | Curetas | Curettes | Curettes






Curettes | Küretten | Curetas | Curettes |

Gracey



| Küretten | Curetas | Curettes | Curettes





added visual confirmation that you

have reached the base of the pocket.

Gracey Variations

Pocket Scalers (PS) and Mini Pocket Scalers (MPS) have a slender and longer terminal shank to provide access into deeper pockets. MPS has the same slender tip design as PS, but with half the blade length. InSite PS and MPS have a depth gauge at 5mm through 10mm offering an excellent visual guide.

Big Easy Ultralite Pocket Scalers (PS) and Mini Pocket Scalers (MPS) do *NOT* feature the depth gauge found in InSite Scalers.

Rigid Gracey designs have a thicker shank than the original Graceys, but are overall slightly shorter for increased strength when removing heavy calculus.



Big Easy® Ultralite™ (BEUL) 1/2" diam.,16g




Big Easy® (BE) 7/16" diam. 23g



Light Touch™ Round (LT)5/16" diam.
21g




Light Touch™ Octagon (LTO) 1/4" diam. 18g




Curettes | Küretten | Curetas | Curettes |

Pocket Scalers



| Küretten | Curetas | Curettes | Curettes






Curettes | Küretten | Curetas | Curettes |




Big Easy® (BE) 7/16" diam. 23g



Light Touch™ Round (LT) 5/16" diam.







Ultrasonic Inserts

| Ultraschall - Zahnsteinentfernung Instrumente | Insertos Ultrasonicos | Pointes Ultrasoniques

Premier Big Easy® Ultrasonic Inserts

Big Easy cushion-grip ultrasonic inserts for improved ergonomic grip and comfort.

- Ergonomically designed large 7/16" diameter for maximum comfort and decreased fatigue
- Compatible with most magnetostrictive handpieces



Cleaning/Sterilization

Inspect, clean, and sterilize prior to each use. Do not use chemical disinfectants prior to sterilization or leave residual cleaning solutions on insert. Carefully remove gross debris, rinse and dry inserts completely prior to sterilization. Sterilize using steam autoclave or chemical vapor sterilizer. Follow sterilizer manufacturer's instructions.

- Do not use Rapid Heat Transfer (Dry Heat) sterilization method.
- Do not expose inserts to temperatures above 350°F/177°C or to chemical vapor concentrations exceeding manufacturer's recommendations.
- · Do not immerse in ultrasonic baths.
- · Do not use phenols.
- · Do not use cold disinfectant/cold sterilization solutions.
- Do not sterilize in 100% plastic sterilization pouches. (A combination paper/plastic pouch is acceptable, paper side down.)

Optimum Performance


- · Insert stacks should be straight and parallel at all times. Bent or warped inserts should be repaired or discarded.
- Bent and/or worn tips prevent effective operation of the ultrasonic insert and are subject to breakage. Inserts should be replaced when tip-wear reduces size and shape from original condition.
- · Water lines to scaler should be flushed and hand-piece should be filled with water prior to use for most efficient results.
- For optimal water spray on P style (external water flow) inserts, the water tube should be aligned with the plane tip of the insert.
- The o-ring should be carefully maintained to ensure the longest possible life. Maintenance includes water lubrication of o-ring as well as proper seating of insert within the hand-piece.
- Damaged or broken (torn) o-rings should be replaced prior to use.
- · Inserts are not indicated for use in restorative dentistry applications.

Proper Placement of Replacement O-Rings

Proper placement of o-rings will significantly extend the life of the o-rings and provide better ultrasonic scaling. During installation, the o-ring should be protected from the sharp edges of the insert stack. For 25K inserts, this is the 3/16" x 3/16" x 3.75" dark section of the insert. For 30K inserts, this is the 3/16" x 3/16" x 3.0" dark section of the insert.

- Remove and discard old o-ring.
- Place new o-ring onto a straight drinking straw (1/4" diameter), approximately 1" from the end.
- Slide the straw completely around insert stack, starting with the side of the straw that holds the o-ring.
- Slide o-ring off of the straw and place into appropriate groove. Remove and discard straw.
- · Sterilize insert before use.

Make sure you are using the supplied o-rings on your ultrasonic inserts as other o-rings may not be the appropriate size.



Operative | Konservierende Behandlung | Operatoria | Operatoire |


Cavity Prep Instruments Pa Schmelzinstrumente Instrumentos para la Preparación de Cavidades Instruments pour Preparation de Cavité	ge 23-24
Excavators Exkavatoren Excavadores Excavateurs	25-28
Amalgam Carriers/Guns Amalgam Träger/Pistolen Porta Amalgamas Portes-Amalgame/Seringues à Amalgame	29
Applicators Anwender Applicadores Applicateurs	30
Cord Packers Fadenstopfer Empacadores de Hilo Fouloirs à Fil à Écarter	30
Pluggers/Condensers Plombier/Kondensator Obturadores/Condensadores Fouloirs à Amalgame	31-33
Carvers Modellier Instrumente Talladores Sculpteurs	34-36
Amalgam Knives Amalgam-Messer Recortador de Amalgama Sculpteurs Amalgame	37
Burnishers Polier Instrumente Bruñidores Brunissoirs	37
Composite/Plastic Filling Instruments Composite Instrumente/Instrumente für Plastische Füllungen Instrumentos para Obturación Plástica y Composite Composite/Instruments pour Plombages Plastiques	38-42
Composite Instruments Composite Instrumente Instrumentos para Composite Instruments pour Plombages Plastiques et Composites	43-44
Crown Removers Kronenringabnehmer Removedores de Coronas Arrache Couronnes	44
Spatulas Spatel Espátulas Spatulas	45
Laboratory Instruments Labor Instrumente Instrumentos de Laboratorio Instruments pour Laboratoire	46-47



Cavity Preparation Instruments

| Schmelzinstrumente | Instrumentos para la Preparación de Cavidades | Instruments pour Preparation de Cavité

CE



1st # - Centigrade Scale 2nd # - Standard Scale

Black's Formula

A standard method used to describe the size, shape and angulation of cavity preparation instruments, created by Dr. Greene Vardiman Black ("the father of modern dentistry").

Hatchets, chisels and cavity prep instruments have 3 measurements: (x, y, z), based on a circular 100° (centigrade) scale, instead of the typical 360° (astronomical) scale.

- x= the width of the blade in tenths of millimeters.
- y= the length of the blade in whole millimeters.
- z = the angle the blade forms with the long axis of the handle.

A fourth number, used on select angle formers and margin trimmers, is placed between the first two numbers (x, a, y, z).

a= the degree of angle for the cutting edge to the long axis of the handle.





Cavity Preparation Instruments

| Schmelzinstrumente | Instrumentos para la Preparación de Cavidades | Instruments pour Preparation de Cavité |



| Exkavatoren | Excavadores | Excavateurs | Excavators

Black 1003225 PREMIER USA BLACK 61/62 CC (20-9-12) BC 61/62 • 1003228 1003229 PREMER USA BLACK 63/64 CC BC 63/64 • 1003229 (15-8-12) 1003230 PREMIER USA BLACK 65/66 €€ BC 65/66 • 1003230 (10-6-12) **California** 10000000 PREMIER USA CALIF 36/07 (C CC 36/37 • 1003303 (15-8-14) 1.15 1003304 PREMIER USA CALIF 38/30 (€ CC 38/39 • 1003304 (11.5-7-14) 1.10 1003432 PREMIER USA 125 (C 1.10 **125** • 1003432 1.10 1003433 PREMIER USA 165 (C 1.10 **165** • 1003433 1.10 1003434 PREMIER USA 175 CC 1.10 **175** • 1003434 1.10 1003435 PREMIER USA 185 CC 1.10 **185** • 1003435 1003436 PREMIER USA 24S CC 1.00 **245** • 1003436 **#** 1.10 1003437 PREMIER USA LS CC 1.10 **LS** • 1003437 1.10 1003411 PREMIER USA 1 (C 1.10 1(C5) • 1003411 1.10 1003413 PREMIER USA 3 CC 1.10

3(C7) • 1003413




Excavators | Exkavatoren | Excavadores | Excavateurs |



| Exkavatoren | Excavadores | Excavateurs | **Excavators**






Excavators | Exkavatoren | Excavadores | Excavateurs |





Amalgam Carriers/Guns | Amalgam Träger/Pistolen | Porta Amalgamas | Portes-Amalgame/Seringues à Amalgame |



Universal Carriers are unique placement instruments with a sleek design. They enable the operator to transport and place amalgam while minimizing finger fatigue.

REPLACEMENT PARTS FOR UNIVERSAL CARRIERS

Item	Item No.
Universal Carrier Nose Tip Straight	9062011
Universal Carrier Nose Tip Right Angle	9062013
Universal Carrier Plunger Straight	9062015
Universal Carrier Plunger Right Angle	9062017
Universal Carrier Spring (Fits all carriers)	9062021
Mini-plunger Right Angle	9062034






Applicators | Anwender | Aplicadores | Applicateurs |




Cord Packers | Fadenstopfer | Empacadores de Hilo | Fouloirs à Fil à Ècarter |



CE

Pluggers/Condensers

| Plombier/Kondensator | Obturadores/Condensadores | Fouloirs à Amalgame |






Pluggers/Condensers

| Plombier/Kondensator | Obturadores/Condensadores | Fouloirs à Amalgame |

(€

Hollenback



Ladmore



Pluggers/Condensers

| Plombier/Kondensator | Obturadores/Condensadores | Fouloirs à Amalgame |



U Detroit #1 Smooth • 1003799



Carvers | Modellier Instrumente | Talladores | Sculpteurs |

Acorn




Cleoid/Discoid



| Modellier Instrumente | Talladores | Sculpteurs | Carvers







Carvers | Modellier Instrumente | Talladores | Sculpteurs |

CE



| Amalgam-Messer | Recortador de Amalgama | Sculpteurs Amalgame | Amalgam Knives



| Polier Instrumente | Bruñidores | Brunissoirs | Burnishers





Composite/Plastic Filling Instruments | Composite Instrumente/Instrumente für Plastische Füllungen | Instrumentos para Obturación Plástica y Composite | Composite/Instruments pour Plombages Plastiques |






| Composite Instrumente/Instrumente für Plastische Füllungen | Composite/Plastic Filling Instruments | Instrumentos para Obturación Plástica y Composite | Composite/Instruments pour Plombages Plastiques





Composite/Plastic Filling Instruments | Composite Instrumente/Instrumente für Plastische Füllungen | | Instrumentos para Obturación Plástica y Composite | Composite/Instruments pour Plombages Plastiques |





| Composite Instrumente/Instrumente für Plastische Füllungen | Composite/Plastic Filling Instruments Instrumentos para Obturación Plástica y Composite | Composite/Instruments pour Plombages Plastiques

Slix™ Composite Manipulation Instruments


Nonstick - Durable - Ergonomic - Autoclavable

The tips are treated with a revolutionary protective coating that provides superior non-stick properties and a lustrous surface for long lasting esthetic appearance. No more composite sticking; dentists appreciate the easy placement and fast clean-up! The ergonomic handles are manufactured from an advance technology polymer with medical-grade silicone grips. "Surprisingly light"-only 16g - on the largest available handle (1/2" (1.27 cm) diameter) to reduce hand fatigue.





Composite/Plastic Filling Instruments | Composite Instrumente/Instrumente für Plastische Füllungen | Instrumentos para Obturación Plástica y Composite | Composite/Instruments pour Plombages Plastiques |





| Composite Instrumente | Instrumentos para Composite | Composite Instruments Instruments pour Plombages Plastiques et Composites

CCI and CMI Anodized aluminum on CMI and CCI instruments provides even greater lubricity to virtually eliminate composite sticking. The black coating does not abrade to discolor the restoration.

Composite Contact (CCI)



CCI technique




After placing composite into box, insert instrument tip into the composite, torque into position and apply gentle pressure to adjacent tooth for desired positive contact. Light-cure composite, making sure composite is not overfilled to lock in instrument. As recommended by Dr. Brian Wilk.

CMI technique



Acorn Carver (21B) is helpful in establishing anatomical contours, prior to carving.

Composite Manipulation (CMI)*



* CMI Instruments, Townie Choice Awards, 2003.


(€



Composite Instruments | Composite Instrumente | Instrumentos para Composite | Instruments pour Plombages Plastiques et Composites



Crown Removers | Kronenringabnehmer | Removedores de Coronas | Arrache Couronnes |



| Spatel | Espátulas | Spatulas | Spatulas






Laboratory Instruments

| Labor Instrumente | Instrumentos de Laboratorio | Instruments pour Laboratoire |

Wax Instruments - PK Thomas




Set of (5) PK Thomas DE Instruments • 1004388

Laboratory Instruments

| Labor Instrumente | Instrumentos de Laboratorio | Instruments pour Laboratoire |





Endodontic

| Wurzelinstrumente | Endodoncia | Instruments à Nerfs |

Excavators Exkavatoren Excavadores Excavateurs	Page 49	
Explorers Sonden Exploradores Explorateurs	49	
Specialty Spezialität Especialidad Spécialité	49-50	
Spreaders Spreizinstrumente Espaciadores Fouloirs à Canaux Radiculaires	51-52	
Pluggers Wurzelkanalstopfer Obturadores Fouloirs à Canaux Radiculaires	53	

| Exkavatoren | Excavadores | Excavateurs | Excavators





| Spezialität | Especialidad | Spécialité | Specialty



Stewart Multi-Purpose • 1003893

technique tip

As suggested by Dr. George Stewart, University of PA.

Real World Explorers have both a delicate tip to locate the canal and a rigid tip to gain access.



Tip diameter measurements in millimeters.

| Spreizinstrumente | Espaciadores | Fouloirs à Canaux Radiculaires | Spreaders







Spreaders | Spreizinstrumente | Espaciadores | Fouloirs à Canaux Radiculaires |










Perio/Oral Surgery | Parodontal/Mündliche Chirurgie | Periodoncia/Cirugía Oral | Chirurgie Parodontale et Buccale |

Scissors Scheren Tijeras Ciseaux	Page 55-59
Hemostats Arterienklemmen Pinzas Hemostaticas Pinces Hémostatiques	60
Needle Holders Nadelhalter Porta-Agujas Porte-Aiguille	61-63
Tissue Pliers Chirurgische Pinzetten Pinzas de Diseccion Pinces a Tissus	64
Spreaders/Retractors Spreizinstrumente/Wundhaken Retractores Ouvre Bouche/Ecartèurs	65
Knives Messer Cuchillos Instruments Parodontaux	66
Periotomes Periotome Periotomos Periotomes	66
Elevators - Periosteals Heber-Raspatorien für Parodontologie Periostotomos Rugines Rasp	67 patoires
Surgical Curettes Chirugische Küretten Curetas Quirúrgicas Curettes Chirurgicales	67-68
Files Feilen Limas Limes	68
Bone & Periodontal Chisels Knochenmeissel Cinceles Ciseaux	69
Biopsy Punch Hautstanze Biopsia de Piel Poincon de Biopsie	69
Rongeurs Hohlmeisselzangen Pinzas Gubias Pinces - Gouges	70
Tissue Nippers Gewebezangen Tijeras para Encías Pince a Tissus Mous	70
Aspirators Saugapparat Aspiradores Aspirateurs	71
Elevators Wurzelheber Elevadores Elevateurs	71-75
Ejectors Wurzelheber Elevadores para Raíces Elévateurs à Racines	75
Forceps Zange Fórceps Daviers	76-84
Cassettes and Instrument Care Pflege der Kassetten & Instrumente Cassettes y Mantenimiento de II Le Soin de Cassettes et Instruments	85-86 nstrumentos



| Scheren | Tijeras | Ciseaux | Scissors





Scissors | Scheren | Tijeras | Ciseaux |





| Scheren | Tijeras | Ciseaux | Scissors








Scissors | Scheren | Tijeras | Ciseaux |

Suture





| Scheren | Tijeras | Ciseaux | Scissors



Hemostats | Scheren/Zwickzangen | Pinzas Hemostaticas | Ciseaux/Instruments de Chirurgie Parodontale |





| Nadelhalter | Porta-Agujas | Porte-Aiguille | Needle Holders



Castroviejo TC • 9085159




Needle Holders | Nadelhalter | Porta-Agujas | Porte-Aiguille |





| Nadelhalter | Porta-Agujas | Porte-Aiguille | Needle Holders





Tissue Pliers | Chirurgische Pinzetten | Pinzas de Diseccion | Pinces a Tissus |









| Spreizinstrumente/Wundhaken | Retractores | Ouvre Bouche/Ecartéurs | Spreaders/Retractors



ComfortView[®] Lip & Cheek Retractor





- Retracts lips and cheeks for optimal buccal and gingival access
- Maximum visibility and patient comfort
- Suitable for use in diagnostic, preventative, therapeutic dental applications and digital impression scanning
- Autoclavable, Latex free

9061381 Universal ComfortView Lip and Cheek Retractor - 2 per box (white)

9061383 Small ComfortView Lip and Cheek Retractor - 2 per box (blue)

9061382 ComfortView Refill Cushions -10 per pack

ComfortView not available from Premier in Europe.


"Outstanding Products" ComfortView® Lip & Cheek Retractor

^{*} An independent non-profit dental education and product testing foundation, December 2006.



Knives

| Messer | Cuchillos | Instruments Parodontaux |



^{*} An independent non-profit dental education and product testing foundation, December 2005.

| Heber-Raspatorien für Parodontologie | Periostotomos | Rugines Raspatoires | **Elevators**



| Chirugische Küretten | Curetas Quirúrgicas | Curettes Chirurgicales | Surgical Curettes





Surgical Curettes | Chirugische Küretten | Curetas Quirúrgicas | Curettes Chirurgicales |

Lucas Alveolar



Files | Feilen | Limas | Limes |


Periodontal/Hirschfeld







| Knochenmeissel | Cinceles | Ciseaux | Bone & Periodontal Chisels



| Hautstanze | Biopsia de piel | Poincon de biopsie | Biopsy Punch



Uni-Punch®


- Seamless (Sharpened on both inside and outside)
- Sterile Product
- Mylar and Tyvek Sterile Peel-Pouch (Gamma radiated)



Sterile product shown in package.

Size	Item No. (Box of 25)	Item No. (Box of 50)
1.5mm •	9033509	9033519
2mm	9033502	9033512
2.5mm	9033511	9033521
3mm	9033503	9033513
3.5mm	9033510	9033520
4mm	9033504	9033514
5mm	9033505	9033515
6mm	9033506	9033516
8mm	9033508	9033518
Assorted	9033525	

Rongeurs | Hohlmeisselzangen | Pinzas Gubias | Pinces - Gouges |



Tissue Nippers | Gewebezangen | Tijeras para Encías | Pince a Tissus Mous |





| Wurzelheber | Elevadores | Elevateurs | Elevators



Elevators | Wurzelheber | Elevadores | Elevateurs |



| Wurzelheber | Elevadores | Elevateurs | **Elevators**

Coupland – Medium Handle



Cryer – Medium Handle




Cryer – Large Handle





Elevators | Wurzelheber | Elevadores | Elevateurs |

Heidbrink 13-14 • 1003396



| Wurzelheber | Elevadores | Elevateurs | **Elevators**

Potts





| Wurzelheber | Elevadores para Raíces | Elévateurs à Racines | **Ejectors**


West Root Tip





Forceps | Zange | Fórceps | Daviers |

Child Size - Maxillary









C4 • 9065054








C6 • 9065056

PREMIER MAXILLARY CHILD FORCEPS SELECTION CHART

Centrals		C150
Laterals	C4	C150
Cuspids		C150
Bicuspids		C150
First Molar	C5	C150
Second Molar	C5	C150
Third Molar		C150
Root	C6	C150



| Zange | Fórceps | Daviers | Forceps

Child Size - Mandibular



C151 • 9065058

PREMIER MANDIBULAR CHILD FORCEPS SELECTION CHART

C3 • 9065053

Centrals		C151
Laterals	C1	C151
Cuspids		C151
Bicuspids		C151
First Molar	C2	C151
Second Molar	C2	C151
Third Molar		C151
Root	C3	C151



Forceps | Zange | Fórceps | Daviers |

Maxillary





150 • 9065030 • 9065078 (serrated)



150A • 9065032



1505 • 9065034

PREMIER MAXILLARY FORCEPS SELECTION CHART

Incisors	150	150A	150S	99C	62	65				
Canines	150	150A	150S	32	99C	62				
Premolars	150	150A	150S	32	99C	62	101			
1st Molar	150S	88L	88R	53L	53R	18L	18R	24	10S	32
2nd Molar	150S	88L	88R	53L	53R	18L	18R	24	10S	32
3rd Molar	24	10S	210	32						
Roots	150	150S	65	69	URG	76S				
Periodontic	150S									

| Zange | Fórceps | Daviers | Forceps

Maxillary







88R • 9065025



53L • 9065015



53R • 9065016



18L • 9065009



18R • 9065010





Forceps | Zange | Fórceps | Daviers |

24 • 9065012

Maxillary







105 • 9065004







32 • 9065013

PREMIER MAXILLARY FORCEPS SELECTION CHART

Incisors	150	150A	150S	99C	62	65				
Canines	150	150A	150S	32	99C	62				
Premolars	150	150A	150S	32	99C	62	101			
1st Molar	150S	88L	88R	53L	53R	18L	18R	24	10S	32
2nd Molar	150S	88L	88R	53L	53R	18L	18R	24	10S	32
3rd Molar	24	10S	210	32						
Roots	150	150S	65	69	URG	76S				
Periodontic	150S									



| Zange | Fórceps | Daviers | Forceps

Maxillary










101 • 9065028

65 • 9065018





Forceps | Zange | Fórceps | Daviers |

Mandibular





151 • 9065031 • 9065077 (serrated)



1515 • 9065035



151A • 9065033



23 • 9065011 • 9065079 (serrated)



16 • 9065006

PREMIER MANDIBULAR FORCEPS SELECTION CHART

Incisors	151	151S	151A	62*	74	MD3	
Canines	151	151S	151A	62*	74	MD3	
Premolars	151	151S	151A	62*	74	MD3	
1st Molars	151S	23	16	16S	17	15	MD4
2nd Molars	1515	23	16	16S	17	15	MD4
3rd Molars	16S	222	MD4				
Roots	151	151S	151A	74	69	LRG	MD3
Pedodontic	1515	16S	62*				

^{*}Image of this item is found on page 81.



| Zange | Fórceps | Daviers | Forceps

Mandibular





17 • 9065008



222 • 9065039



Mead 4 • 9065064



74 • 9065021





Forceps | Zange | Fórceps | Daviers |

Root











Lower Root Grip • 9065044





765 • 9065045



Mead 3 • 9065063



Cassettes and Instrument Care

| Pflege der Kassetten & Instrumente | Cassettes y Mantenimiento de Instrumentos | Le Soin de Cassettes et Instruments

Instrument Cassettes

Premier stainless-steel instrument cassettes offer proven reliability and quality performance with an unparalleled satisfaction guarantee.

Unique, hingeless design offers years of trouble-free operation. Compact, double-stack feature accommodates more instruments for greater efficiency in ultrasonic cleaning, autoclaving and storing. Apertures in side walls, lid and bottom of trays permit optimum flow of cleaning solution, ultrasonic energy and greater access for rinsing.

Includes a five-year satisfaction guarantee.

9010904	Mark I (4 or 8 Instruments) 2.5" x 8.5" x 1.5"
9010908	Mark VI (8 or 16 Instruments) 4.5" x 8.0" x 1.25"
9010924	Mark IV (12 or 24 Instruments) 8.0" x 11" x 1.50"
9010901	Accessory Cassette 1.5" x 3.5" x 1.0"

"Best Products of the Year" Mark VI Cassette



Premier cassettes are easy to use!

- 1 Rinse cassette/bur block organizer under running water.
- 2 Place cassette/bur block organizer in ultrasonic unit; use Brite Shield™ to remove bio-burden and to clean & protect instruments.
- 3 Rinse cassette/bur block organizer under running water.
- 4 Allow instruments & cassette/bur block organizer to dry thoroughly.
- 5 Wrap cassettes and sterilize.





Mark I (4 or 8 Instruments) Ideal for prophy or exam setups



Mark VI (8 or 16 Instruments) Works for most operative procedures



Mark IV (12 or 24 Instruments) Accessory compartment for loose items



Accessory Cassette Holds either handpiece head or bur block

^{*} An independent non-profit dental education and product testing foundation, November - December 2007.

Cassettes and Instrument Care

| Pflege der Kassetten & Instrumente | Cassettes y Mantenimiento de Instrumentos | Le Soin de Cassettes et Instruments |

Sterilization Troubleshooting

- Don't sterilize dirty instruments. They should always be cleaned before sterilization. Otherwise, the soil and dirt will become even harder to remove because the heat, hot water, or steam bakes the contaminant onto the instrument.
- **Don't use abrasives.** Steel wool or abrasive powders/soaps should not be used to scrub stainless-steel instruments. These abrasives may wear off the protective coating of the instruments. The finish protects the instrument from oxidation and corrosion.
- Don't mix unlike metals. Sterilizing different metals together can lead to a chemical reaction which can be detrimental to the instruments.
- Don't use harsh detergents. The pH of the detergent should be neutral (7pH). Highly alkaline (high pH) detergents may produce orange-brown stains. Highly acidic (low pH) detergents will cause black pitting marks that will permanently corrode and ruin the instrument.
- Don't use chrome-plated instruments if the plating is starting to wear away. Rust and corrosion will readily occur, spreading to other instruments.
- Don't put instruments in a sterilizing unit that has a build-up of deposits in the sterilizing chamber. Always remove scale and mineral deposits before sterilization. A mild abrasive or weak solution of acetic acid (vinegar) may be used to remove these deposits.
- Don't use water with a high mineral content or water that is too soft. Always use distilled water.
- Don't exceed recommended temperatures.
 Excessive heat in the sterilizing and drying procedures can be harmful to the metallurgy of the instruments and cause "rainbow-colored" stains
- Don't use vibrating marking devices to mark instruments. Avoid the use of an "etching" device to personalize instruments. This can lead to rust/corrosion which can spread to other instruments.
- Don't open the autoclave too soon. The premature opening of the autoclave causes rapid condensation of steam on the hot instruments and will permit residue deposits.

Sterilization vs. Disinfection

Sterilization is the act or process of complete destruction of all living microorganisms, while disinfection is the act or process of destroying most disease-causing microorganisms. The two are not interchangeable! Both the Center for Disease Control (CDC) and the American Dental Association (ADA) warn that all patients should be considered potentially infected and recommend that complete sterilization procedures are used whenever possible.

Method	Conditions	Advantages
Steam Autoclave - Designed for stainless steel, textiles, liquids, heat-resistant handpieces, tubing & plastics.	20 minutes at 250°F (121°C) (15 psi) 3-5 minutes at 273°F (134°C) (30 psi)	Time efficient, good penetration,water-based liquid
Liquid Sterilants - Designed for heat-sensitive materials.	Up to 24 hours, using EPA registered solutions	Works on heat-sensitive items
Dry Heat - Designed for stainless steel, some heat-resistant hand- pieces, tubings and plastics.	2 hours at 320°F (160°C) 1 hour at 340°F (171°C)	No corrosion, can use closed containers, large capacity, immediately dry
Unsaturated Chemical Vapors - Designed for high carbon metal and stainless steel, textiles, liq- uids, heat-resistant handpieces, tubings and plastics.	30 minutes at 270°F (132°C) (20-40 psi)	Time efficient, no corrosion, items dry quickly

Steps for Proper Sterilization via Autoclaving

- 1. Pre-soaking & cleaning
- 2. Corrosion control and lubrication for hinged instruments
- 3. Packaging
- 4. Sterilization
- 5. Sterilization monitoring
- 6. Drying and cooling
- 7. Storage







Why ACP?

Amorphous calcium phosphate (ACP) technology is the first non-fluoride therapy to remineralize enamel and dentin.15

ACP - The Progressive Approach

Dental technology is ever evolving and what may have been the viewpoint a quarter century ago, the stance has now progressed to the realization that essential chemical kinetics of calcium and fluoride are beneficial in the war against dental disease.

Innovative products with ACP technology provide an arsenal of staged-potential with the use of calcium and phosphate ions.

- ACP combined with fluoride offers a multitude of valuable elements.
- ACP forms on the tooth enamel, within the dentinal tubules and when met with fluoride, provides a beneficial calcium fluoride ion reservoir on teeth, in plaque and in the saliva. Remineralization occurs when sufficient amounts of these ions are present.
- ACP products contribute to greater fluoride uptake. 1,2,10,17

Greater luster for whiter, brighter teeth^{3,8} Zygo® non-contact profilometry intensity map profiles



before prophy

polishing





polishing with leading prophy

After prophy polishing with Enamel Pro

®Zygo is not a registered trademark of Premier Dental Products Company.

Value of Fluoride Sources

Fluoride reacts with tooth mineral, forming either fluoridated apatite (firmly-bound) or calcium fluoride (loosely-bound).

Firmly-bound fluoride, incorporated onto the surface of the crystals of apatite, can reduce the solubility of tooth mineral and hence inhibit demineralization due to acids generated by plaque bacteria.17

Loosely-bound (labile F) fluoride provides a relatively slow-release form of ionic fluoride to plaque and saliva.¹⁶ A potential reservoir source of solution fluoride enhancing remineralization and retarding demineralization processes.¹⁶

Enamel Pro®

The ONLY prophy paste formulated to deliver ACP





Strengthen their smiles!

Enamel Pro® delivers 31% more fluoride uptake than conventional prophy paste²

Patient pleasing: cleans quickly, less splatter, easy to rinse, smooth consistency and gluten-free with great taste

Provide a true benefit to your patients while your prophy provides a lustrous smile

Deposits ACP

(Amorphous Calcium Phosphate)

ACP stimulates remineralization of tooth enamel1,4,5,6,7,8

ACP prevents the loss of enamel1,7,9





Box of 200 - With Fluoride

9007600	Mint Fine
9007601	Mint Medium
9007602	Mint Coarse
9007603	Mint Extra Coarse
9007604	Cinnamon Fine
9007605	Cinnamon Medium
9007606	Cinnamon Coarse
9007610	Grape Fine
9007611	Grape Medium
9007612	Grape Coarse
9007607	Strawberry Fine
9007608	Strawberry Medium
9007609	Strawberry Coarse
9007614	Bubblegum Fine
9007615	Bubblegum Medium
9007616	Bubblegum Coarse
9007617	VanillaMint Fine
9007618	VanillaMint Medium
9007619	VanillaMint Coarse

Without Fluoride

9007613 Mint Coarse

Each box of Enamel Pro® contains 200 single-use cups, packaged in a foil bag and two Comfy-Grip[™] prophy cup holders.

9007650

Comfy-grip™ cup holder (3) See page 92.

- 1 Tung, M.S., Malerman, R., Huang, S. and McHale, W.A. "Reactivity of Prophylaxis Paste Containing Calcium, Phosphate and Fluoride Salts", Journal of Dental Research, Vol. 84, Special Issue A, IADR Abstracts, 2005.
- 2 Source: "Final Report: Enamel Fluoride Uptake Study #05-106, Modified FDA Method #40," Dental Products Testing, Indiana University Emerging Technologies Center, December 2005.
- 3 Winston, A.E., Charig, A., Patel, V., McHale, W.A. and Malerman, R "Effect of Prophy Pastes on Surface of Tooth Enamel", Journal of Dental Research, Vol. 84, Special Issue A, IADR Abstracts, 2005
- 4 Tung, M.S., Eichmiller, F.C. "ACP for Tooth Mineralization". Compendium Cont. Educ Dent, Vol 25, No. 9, Suppl. 1, Sept. 2004, p 9-13. Journal Code 9600713, Journal Subset MEDJSD, ISSN 1548-8578. Corp. Author: American Dental Association Foundation.
- 5 Tung, M., Eichmiller, F., Gibson, H., Ly, A., Skrtic, D., Schumacher, G., (ADAHF Paffenbarger Research Center, NIST, Gaithersburg, MD)," Dentin Desensitization by in situ Formation of Calcium Phosphate", J. Dent Res 76 (IADR Abstracts):2985, 1997

 6 Schemehorn, B.R., Orban, J.C., Wood, G.D., et al. "Remineralization
- by Fluoride Enhanced with Calcium and Phosphate Ingredients", Journal of Clinical Dentistry, 1999; 10(1 spec. no.);13-16.

angle

Patented gearless

design

2pro®

Total Access Disposable Prophy Angle

2pro® features a patented dual-action cup and tip design to guarantee total access to all tooth surfaces.

Saves time and money

Comparably priced to other angles, yet when needed, you save time and money by removing the cup to clean remaining areas with 2pro's gentle tip

Provides total access

2pro's unique 2-in-1 design guarantees total access to all tooth surfaces

Runs smoothly and quietly

Patented gearless design ensures quiet, smooth-running, vibration-free operation with less heat build-up

Easier access to lingual areas

Ergonomic 100° design offers comfort and convenience by reaching all areas of the mouth

Latex free

Bag of 144 disposable angles

5500101	Soft/Short Cup (purple cup)†
5500102	Firm/Short Cup (blue cup)†
5500103	Soft/Long Cup (white cup)†
5500104	Firm/Long Cup (green cup)†

Bulk package of 1,000 disposable angles

5510001	Soft/Short Cup (purple cup)
5510002	Firm/Short Cup (blue cup)
5510003	Soft/Long Cup (white cup)
5510004	Firm/Long Cup (green cup)

Order your preferred cup style through your authorized dealer.

†Dealer: Available to you in a case of ten







Get Total Access!

2pro cleans and polishes all surfaces with its unique cup and tip design.

No added expense for tips or brushes!



CUP ON!

Polish and clean all facial surfaces.

Facial surfaces



Gingival/Interproximal

CUP OFF!

Perfect for interproximal access, crown and bridge, veneer margins, around implants, occlusal surfaces and orthodontic brackets.



CE


Veneer margins



Orthodontic brackets







Firm/short

blue cup,

short tip



*An independent non-profit dental education and product testing foundation, November/December 2010

- 7 Tung, M.S., Eichmiller, F.C., Paffenbarger Research Center. ADAHF. NIST. Gaithersburg, MD. "Dental Applications of Amorphous Calcium Phosphates" Journal of Clinical Dentistry, Vol.10, Issue 1, 1999.
- 8 Tung,M.S., O'Farrell, T.J., Liu, D.W.,Paffenbarger Research Center. ADAHF. NIST. Gaithersburg, MD. "Remineralization by Amorphous Calcium Phosphate Compounds" Journal of Dental Research, Vol. 72, Abstract No. 1738, 1993.
- 9 Hicks M.I. Flaitz C.M. "Enamel Caries Formation and Lesion Progression with a Fluoride Dentrifice and Calcium Phosphate Containing Fluoride Dentrifice: A Polarized Light Microscopic Study", ASDC J Dent Child. 2000 Jan.-Feb.; 67(1): 21-8, 8.
- 10 Source: "Final Report: Sound Enamel Fluoride Uptake Study #08-151, Method #40SV", Dental Products Testing Thermametric Technologies Inc., Indiana University Emerging Technologies Center,
- 11 Tung, M.S., Torres, J., "Effect of Calcium Phosphate Containing Fluoride Varnish on Dentin Permeability", J Dent Res, 86 (Spec Iss A):0985,2008 (www.dentalresearch.org)
- 12 85th General Session in New Orleans, LA.: Tung, M.S., Hwang, J., Malerman, R., McHale, W.A., "Reactivity of Varnish Containing Calcium, Phosphate and Fluoride Salts". J Dent Res 86 (Spec Iss A):
- 13 "Professionally applied topical fluoride: Evidence-based clinical recommendations", American Dental Association, Council on Scientific Affairs, J Am Dent Assoc., 2006; 137: 1151-1159

Soft/long

white cup,

long tip

- Schemehorn, B.R., Orban, J.C., Wood, G.D., et al. "Remineralization by Fluoride Enhanced with Calcium and Phosphate Ingredients", Journal of Clinical Dentistry, 1999; 10(1 spec. no.),:13-16.
- 15 Garvin, J., ACP: the next big thing. ADA News, 2007.
- 16 White, DJ and Nancollas, G.H. "Physical and Chemical Considerations of the role of firmly and loosely bound fluoride in caries prevention", J Dent Res, 1990, 69:587-594.



Firm/long

green cup,

long tip



Enamel Pro® Varnish

First advanced fluoride varnish optimized to provide valuable tooth building ions of calcium and phosphate.

Designed to create substantive forms of CaF and ACP.



Amorphous Calcium Phosphate (ACP) and calcium-fluoride (CaF) mineral on dentin surface and in tubules

Nice esthetic appearance, patient pleasing and with studies supporting fluoride uptake, fluoride release and diminished hydraulic conductance. 11,17,24

5% NaF White Varnish Box of 35

9007540	Strawberries 'n Cream 0.40ml
9007541	Bubblegum 0.40ml
9007547	VanillaMint 0.40ml
9007542	Bubblegum 0.25ml

Economical Box of 200

9007543	Strawberries 'n Cream 0.40ml
9007544	Bubblegum 0.40ml
9007545	VanillaMint 0.40ml

Each box contains single-unit doses of 5% Sodium Fluoride Varnish and an applicator brush in a foil sealed tray, one autoclavable Comfy-grip™ cup holder, Post-application patient-care instructions and doctor's directions for use.

Enamel Pro Tutorials!

Go to www.premusa.com to view online tutorials!











"Outstanding Products"*

INDICATIONS FOR USE

Enamel Pro® Varnish is a fluoride containing preparation for the treatment of dentinal hypersensitivity, and for the reduction of post operative sensitivity.

Independent studies show Enamel Pro Varnish delivers more enamel fluoride uptake and greater fluoride release than 3M Varnish[†],17,24





In a 2012 study evaluating incubation conditions on fluoride release in vitro, Enamel Pro Varnish displayed greater fluoride release than both 3M Vanish and Colgate Duraphat^{†,24}

- 18 Source: "Final Report: Sound Enamel Fluoride Uptake Study #08-230, Method #405T", Dental Products Testing, Indiana University Emerging Technologies Center, September 2008
- 19 Council on Dental Materials Instruments and Equipment. Status report: effect of acidulated phosphate fluoride on porcelain and composite restorations. J Am Dent Assoc 1988; 116:115
- 20 American Academy of Pediatric Dentistry, Guideline on Fluoride Therapy, Pediatr Dent 2008; Clinical Guidelines 121-124; www.aapd.org/media/Policies_Guidelines/G_FluorideTherapy.pdf.
- 21 "Professionally applied topical fluoride: Evidence-based clinical recommendations", American Dental Association, Council on Scientific Affairs, J Am Dent Assoc., 2006; 137: 1151-1159.
- 22 Featherstone, John D.B., The Science and Practice of Caries Prevention, J Am Dent Assoc., 2000; 131: 887-899.

^{*}An independent non-profit dental education and product testing foundation, December 2006

^{† 3}M Vanish and Colgate Duraphat are not trademarks of Premier Dental Products Company.

¹⁷ Schemehorn, B.R., Wood, G.D., McHale, W., Winston, A.E., "Comparison of Fluoride Uptake into Tooth Enamel from Two Fluoride Varnishes Containing Different Calcium Phosphate Sources" J Clin Dent, 2011; 2:251-54.



Enamel Pro® Gel

Non-acidulated Fluoride Gel



1.23% Non-acidulated Fluoride Topical Gel

- Enamel Pro® Gel with ACP technology is non-acidulated, yet provides the same fluoride uptake as acidulated phosphate (APF) gels.¹⁸
- Enamel Pro Gel offers the same advantage as neutral pH gels; safe for patients with porcelain, composite restorations and sealants
- No dulling or etching of teeth which could occurs in APF gel use¹⁹
- Use of fluorides for the prevention and control of caries is documented to be both safe and highly effective²⁰
- Fluoride gel is effective in preventing caries in school-aged children in school-aged children²¹
- ACP and Fluoride are known to inhibit demineralization and enhance remineralization^{22,23}
- Great flavors without the acidic aftertaste

16 ounce (454g) bottle - gluten free

9007200	Bubblegum	
9007201	Strawberry	
9007202	Cherry	
9007203	Orange	
Not availabl	e in Canada.	





Disposable Dual-arch Fluoride Trays with Saliva Ejector Port

E-VENT® patented fluoride tray eliminates the anxiety associated with fluoride-tray treatments. Designed by a practicing dental hygienist with you and your patients in mind

- Evacuates saliva and excess gel, foam or paste
- · Increases patient comfort and compliance
- · Virtually eliminates gagging
- · Flexible and contours to the mouth
- · Contains NO latex, phthalates or bis-phenols

Disposable Dual-arch Fluoride Trays

9007221 E-VENT® Fluoride Trays Size medium 40/box†

†Dealer: Available to you in a case of ten.

Single-patient use. U.S. Patent #5,513,986



²³ Schemehorn, B.R., Orban, J.C., Wood, G.D., et al. "Remineralization by Fluoride Enhanced with Calcium and Phosphate Ingredients", Journal of Clinical Dentistry, 1999; 10(1 spec. no.);13-16.

²⁴ Fluoride Release from Fluoride Varnishes is Affected by Incubation Conditions C. Gonzalez-Cabezas, S.E.Flannagan, R. Krell, M. Uekihara, and C.C. Niquette, Dent Res 91(Spec Iss A):1094, 2012 (www.dentalresearch.org)



Glitter®

Prophy Paste



Glitter® prophy paste is smooth, pliable and splatter-free to clean teeth quickly. Best of all, Glitter is easy to rinse off teeth. Gluten Free.

Fluoride - Box of 200 cups		
9007400	Mint Fine	
9007402	Mint Medium	
9007404	Mint Coarse	
9007406	Mint Extra Coarse	
9007410	Cherry Fine	
9007412	Cherry Medium	
9007414	Cherry Coarse	
9007420	Strawberry Fine	
9007422	Strawberry Medium	
9007424	Strawberry Coarse	
9007430	Bubble Gum Fine	
9007432	Bubble Gum Medium	
	noride - Box of 200 cups	
9007440	Mint - Coarse	

Mint - Medium



9007442

Ziroxide™

Prophy Paste



Ziroxide™ works well because it contains:

- · Zirconium Silicate
- · Cohesive Glycerine Base
- · Natural Mint Flavor

Ziroxide removes stains fast, especially tobacco and other difficult stains, yet leaves enamel with a high gloss. Gluten Free.

Fluoride - Box of 200 cups
9007173 Medium/Fine grit
9007172 Coarse grit

Fluoride - 1 lb jar
9007175 Medium/Fine grit
9007174 Coarse grit

Without Fluoride - 1 lb jar
9007171 Medium/Fine grit
9007170 Coarse grit



Comfy-grip™Prophy Cup Holders



Ergonomic cup holders for use with most prophy-paste cups and fluoride varnish applications. Autoclavable.

9007650 Comfy-grip™ cup holders (3)





ComfortView®

Lip & Cheek Retractor





"Outstanding Products"*
ComfortView® Lip & Cheek Retractor

- Retracts lips and cheeks for optimal buccal and gingival access
- Maximum visibility and patient comfort
- Suitable for use in diagnostic, preventative, therapeutic dental applications and digital impression scanning
- · Autoclavable, latex-free

9061381 Universal ComfortView Lip and Cheek Retractor - 2 per box (white)
 9061383 Small ComfortView Lip and Cheek Retractor - 2 per box (blue)
 9061382 ComfortView Refill Cushions -10 per pack

ComfortView not available from Premier in Europe.

Therma-Trol™

Potassium Oxalate Desensitizer Gel



Relieves dentinal thermal hypersensitivity in two minutes. 0.2cc MicroDose® single-patient finger applicators.

4009015 40 units

Autoclavable D.I.S.C.™

Instrument Sharpening System



The fully autoclavable D.I.S.C. sharpener is an easy-to-use, efficient way to sharpen instruments for optimum results.

The D.I.S.C. sharpener is ideal for routine sharpening as well as for chairside touch-ups. Within seconds, the D.I.S.C. will resharpen instruments to make scaling and root planing more effective, with less effort by the dental professional.

Replacement D.I.S.C. sharpening stones, autoclavable and available separately, permit convenient chair-side touch-ups. They cut faster, clean easier and never need oil as required with Arkansas Stones.

9048200

Autoclavable D.I.S.C. complete with (2) sharpening stones and Comprehensive Guide to Sharpening

9048093 Sha

Sharpening Stone (1)

™ D.I.S.C. is a trademark of Steve B. Moore/U.S. Patent # 4,821,462

See page 12 for instructions for use.



^{*}An independent non-profit dental education and product testing foundation, December 2006

[®] MicroDose is not a registered trademark of Premier Dental Products Company.



Topical Anesthetics



Topicale® Gel Pumps

- Minimal waste pump just as much as is needed
- Economical 43g package contains up to 50% more than competitors
- Reduces potential cross-contamination
- 18% benzocaine formulation

3 - 43g pumps

5 15g pun	iba
9007140	Assorted (Pina Colada,
	Raspberry, Cherry)
9007141	Cherry (3)
9007142	Raspberry (3)
9007143	Strawberry (3)
9007144	Pina Colada (3)



Topicale® Xtra

- 20% benzocaine formulation for rapid 10-second onset of surface anesthesia of soft tissue
- Three fun flavors, with no bitter after-taste
- Economical 1oz jar

0007150	Van Barri Clarin
9007152	Very Berry Cherry
9007153	Mint Freeze
9007154	Jammin' Strawberry
9007155	Assorted Package of 3



Topicale® GelPatch™

Once placed in the mouth, Topicale GelPatch dissolves gradually - releasing benzocaine to anesthetize surface tissue for up to 20 minutes. Appropriate for patient home use.

- · Ideal for denture sores
- · Prior to injection
- · Oral ulcer/irritations
- Fast acting, self-dissolving and flexible

9007149 GelPatch - (25) mint



Topicale® Gel Jar

Use on areas of injection and/or prior to scaling and curettage to alleviate patient discomfort

- · Tropical fruit flavored
- 18% benzocaine formulation

9007161 Gel - 35g / 1.25oz jar



Topicale® Ointment

Unique water-soluble base releases active ingredients on contact; permits onset of anesthesia in less than 20 seconds. Provides deeper surface anesthesia.

- Orange-raspberry flavor
- · 18% benzocaine formulation

9007151 Ointment - 35g / 1.25oz jar

Cartridge Warmer and Dispenser

Dispenses Anesthetic Cartridges at Body Temperature



Warming unit heats up to 60 cartridges to approximate body temperature to avoid anesthetic trauma from thermal shock. One cartridge is dispensed by pressing top panel. Remaining cartridges stay in case, free from environmental contamination. May be wall-mounted or placed on cabinet. Made of white, high-impact polystyrene. Complete with heating bulb.

1048021	Complete Unit (110 volt)
1048022	Replacement bulb



Cleaners







Cleaners

Vac Attak™

High Performance Evacuation System Cleaner



"Pearls For Your Practice" Dental Economics - January 2002

"Outstanding Products"*

Significantly More Effective

Unique combination of chelation chemistry and concentrated organic chlorine is significantly more effective than enzymatic cleaners in removing debris and bio-burden from evacuation lines. Vac Attak is especially effective on prophy paste and fluoride gels, the two most common causes of blocked lines. Enzymatic cleaners simply do not have enough contact time to be effective during an evacuation flush.

Fast-Acting

Vac Attak's organic chlorine formula facilitates the immediate breakdown of protein chains, found in bio-burden, into water-soluble amino acids.

Safe and Environmentally Friendly Biodegradable and non-corrosive to metals. Non-foaming - safe for use with all pumps. Prevents scale build-up in pump that otherwise forms when using hard water.

Economical and Easy To Use

One scoop of concentrated Vac Attak quickly dissolves into 32oz of warm water and then is flushed through the system. It's that easy! Up to 100 cleanings.

9011101 800g jar with scoop dispenser

Dealer: Available to you in case of six

Brite Shield™

Static Soak & Ultrasonic Cleaner



"Most Highly Rated"**

Versatile

Chemo-mechanically debrides clogged dental burs, removes alginate and stone from impression trays. Also removes set glass ionomer, carboxylate and zinc phosphate cements. Ideal for removing bio-burden from hand instruments.

Polishes and protects

A protective shield coats the instrument through a unique passivation process to restrict rust and corrosion while imparting a lustrous surface to stainless-steel instruments.

> Versatile and economical enzymatic cleaner



After

9011100 800g jar with scoop dispenser

Dealer: Available to you in case of six

Removalon® - I

Removes Set Cements



Removalon®-I effectively removes set glass ionomer, carboxylate and zinc phosphate cements from instruments, ortho bands, crowns, inlays, facings, etc. When used with an ultrasonic cleaner, Removalon-I dissolves cement away twice as fast.

Removalon- I dry pak is convenient to use and easy to store. Just extend the container and fill with warm tap water. The concentrated powder is dissolved quickly, forming a safe, effective solution.†

9011005 Dry Pak (each fills to 1/2 gallon - 1.9L)

Dealer: Available to you in case of six

†Good for two months after mix

^{*}An independent non-profit dental education and product testing foundation, December 2002

^{**}An independent non-profit dental education and product testing foundation, January 2001

^{***}An independent non-profit dental education and product testing foundation, December 2005





Rapid-Klene®

Automatic Film Processor Cleaner



Rapid-Klene® effectively and safely cleans automatic film processors in as little as 20 minutes - without noxious fumes or scrubbing.

After draining all sections, one Rapid-Klene pouch is emptied into each section to be cleaned and hot tap water is added. The transports are inserted and the machine is turned on. It's that easy.

Rapid-Klene contains no harmful chemicals, requires no rubber gloves nor protective clothing.

9018000	Case -12 cartons of 12 pouches (144 - 63g pouches)
9018101	Carton - (12 - 63g pouches)

Pro Clean™

Professional Strength Tartar and Stain Remover for Dentures, Partials, Retainers and Night Guards



Extremely effective

Easy to mix and dispense without waste Safe and environmentally friendly

Dispense for patient home use - ideal for retainers and night guards (one box lasts six months).

9011103 Box of 25 packets

9011102 Case of 10 Boxes
(25 packets/box)


"Pearls For Your Practice"
Dental Economics - November 2002

"Outstanding Products"***









Solo Diamond Sterile, Single-patient Use Diamond Instruments	Page 99-101	
Two Striper® Multi-use Diamond Instruments	102-114	
TS2000 [™] Speed-cut Spiral Diamonds	103-111	
ShortCut [™] Short-shank Diamonds	104	
Chamfer Prep Round-End Taper, Flame and Cylinder	105	
Modified Chamfer Curettage K-Style, Curettage GCP and Modified Beveled Cylinder	106	
Shoulder Prep Flat-End Taper and Cylinder	107	
KR Modified Shoulder Prep KS and KR Taper	108	
Bevel Shoulder Prep Bevel-End Cylinder	108	
Flame/Feather Edge Prep/Break Contact Flame, Cone, Interproximal Trimmer and Break Contact	109	
Occlusal & Lingual Reduction Barrel, Egg, Football, Round, Wheel and Occlusal Contour	110	
Veneer Prep Depth Cut and Depth Marker	110	
Inlay/Onlay Flat-End Cylinder, KS Cylinder, KR Cylinder and KR Taper	111	
Crown & Bridge Specialty ☐ TGE™ (Tissue Guard End-Cutting), Safe-End, Gross Reduction and SI	111 kooter™	
Endo Slow speed	111	
Operative Pear, Amalgam Removers, Inverted Cone and Diamond Versions of	112 Popular Carbides	
MicroPrep™ Micro-preparation Diamonds	112	
TSZtech™ Lithium Disilicate and Zirconia Diamonds	112	
Preparation Kits Composite Finishing Kit, Universal All-Ceramic Prep Kit	113	
Clean-A-Diamond® Mini-Squares™ Dressing Stone for Two Striper PBS® Diamonds	114	
Surgical Length Carbide Burs & Steel Burs	115	
Lab Thin-Flex® Diamonds Discs, Silverstreak™, Laboratory Diamonds, Tri-Fissure HP Carbide, Mandrels and Bur Blocks	116-117	



Get a FRESH START with every SINGLE PATIENT

- Stop guessing about worn diamonds and feel the difference cutting with a sharp new Solo diamond.
- Eliminate unnecessary bur blocks and simplify bur management with sterile, easy-to-read sealed diamond pouches in convenient boxes of 25.
- CE-marked and manufactured to ISO standards.
- Order Solo Diamonds from your favorite dental supplier today!





Ball



Solo Diamond

Single-patient use, pre-sterilized diamonds

Coarse

801016C

Gamma-sterilized and individually packaged to improve infection control

Coarse

801018C


· Convenient, fast tooth preparation with a fresh, always sharp cutting instrument

Coarse


801023C

· CE-marked, manufactured and sterilized in the U.S.A. to ISO standards











Inverted Cone

Coarse

805016C









Solo Diamond

Single-patient use, pre-sterilized diamonds

Available in 25 single-patient use diamonds per box.


ISO 9001:2000 Quality Management/ISO 13485: 2003 Medical Devices



Not all disposables are created equal!

Premier Solo single-use, pre-sterile diamonds are manufactured and sterilized in the U.S.A. Sterilization is guaranteed. Quarterly dose audits on file: validation summary reports gamma radiation (Cobalt 60)

- Bioburden Test
- Irradiation Levels
- Sterility Test
- 5-Year accelerated aging test supports the claim of "sterile until opened".







The Two Striper® Difference


Faster, Cooler Cutting and Longer Lasting

Consistent Performance

- 100% Ring Gauge Inspection for Shank Concentricity
- 100% Visual Inspection
- 100% Weight Control for Diamond Coverage
- Recommended for All Handpieces

Only the Two Striper® diamond is manufactured with the P.B.S.® diamond-particle bonding process.

Most other diamond brands are electroplated imports.



December 2000. 3) An Independent Non-Profit Dental Education and Product Testing Foundation Newsletter, Sept. 2001. 4) An Independent Non-Profit Dental Education and Product Testing Foundation Newsletter, May 2002. 5) Data on file.

P.B.S.® bonding maximizes the exposure of diamond-cutting surfaces - especially at the tips and upper circumference of the diamond instrument where most cutting occurs. Crystals, placed evenly and precisely in a uniform matrix, are fused permanently to a surgical-grade stainless-steel shank; they are guaranteed not to fall out.

Diamond particles on electroplated tools are trapped mechanically in inconsistent plating layers that can expose either too much or too little of the particle. Some particles are attached loosely and are prone to premature dislodgement. Other particles are buried completely in the plating layer and are unavailable to cut tooth structure, creating unnecessary heat and friction. Electroplated diamonds cut less consistently than Two Striper® diamonds¹ and have a shorter product life.

The combination of permanent crystal attachment, precise placement at the top and maximum crystal exposure makes Two Striper® faster, cooler cutting and longer lasting.

Natural, virgin diamond crystals on Two Striper® rotary instruments contain more corners and angles than the synthetic diamonds used on electroplated burs.



Two Striper® optimal coarse grit is preferred over supercoarse grits.

Super-coarse particles are bigger and require more plating to hold them. More plating means less diamond coverage overall. A recent study found no difference in cutting rates between electroplated coarse and super-coarse diamond grits². While matching performance, Two Striper® offer cooler cutting and are longer lasting! Two Striper® optimal coarse-grit diamond crystals are attached permanently and are guaranteed never to strip or peel.






Rotary



Also available in short-shanks! See page 104.

7\$2000th



TS2000™ provide rapid gross reduction and a finished margin without changing burs.

"Fastest Brand"

Less chair-time for you and your patient. TS2000™ saves time and money by simultaneously performing rapid gross reduction and creating a finished margin without changing instruments. The TS2000™ is truly two diamonds in one.

Consistent clinical performance.

On average, TS2000™ is 50% faster than other spiral designs – even after eight preparations⁵. This unique coarse diamond is designed to cut cooler and faster than electroplated super-coarse diamonds. The spiral action circulates water-spray more efficiently to reduce heat energy and accelerate cutting.




A finished margin is created as shown. You

is created as shown. You can easily create a chamfer, bevel, or flat shoulder by selecting one of the more than forty shapes available.



Actual Preparation



Close-up occlusal view of finished margin


"Combined lab and clinical data showed highest rated and best performing diamonds were the TS2000™ and 770.8C
Two Striper® multi-use diamonds"...










ShortCut™ Short Shank Diamonds


















































Chamfer Prep

- Round-end tapered or cylindrical diamonds create a chamfered finish line
- Recommended prep design for full coverage crowns and veneers
- Recommended prep design for PFM crowns and some ceramic crowns


			H	"L" 1.5i long ⊟ sha	ger					"L" 1.5r long sha	nm ger nk					
Round-End Taper	747.6	767.5	767.7	L767.7	767.8	767.9	770.5	770.7	770.8	L770.8	770.9	770.10	776.4	777.8	780.4	780.7
Coarse	2015058	2015452	2015453	2015652	2015454	2015455	2015456	2015458	2015457	2015650	2015460	2015459	2015461	2015972	2015462	2015463
Medium														2015930		
Fine			2015483	2015653	2015484	2015485			2015487	2015651		2015489	2015491			
Minor/Major Diameter	0.8 / 1.2	1.1 / 1.6	1.1 / 1.8	1.1 / 1.8	1.0 / 1.8	0.9 / 1.8	1.3 / 1.8	1.2 / 1.8	1.3 / 1.8	1.3 / 1.8	1.1 / 1.8	1.1 / 1.8	2.1/2.6	0.6 / 1.1	1.7/2.1	1.4/2.1


Downd Fod Town	700.0	——————————————————————————————————————	7016	70	79110		702.0	700	7026	703 0	793.10	——————————————————————————————————————			<u> </u>	
Round-End Taper	780.8	780.9	781.6	781.8	781.10	782.6	782.8	782.10	783.6	783.8	783.10	784.6	784.8	784.10	785.4	785.7
Coarse	2015464	2015467	2015465	2015478	2015466	2015423	2015424	2015425	2015426	2015427	2015468	2015429	2015430	2015431	2015470	2015469
Fine	2015494	2015497		2015498	2015496	2015433	2015434	2015435			2015438		2015440	2015441		2015499
Very Fine																2015529
Minor/Major Diameter	1.4 / 2.1	1.4/2.1	0.8 / 1.1	0.8 / 1.2	0.8 / 1.3	1.1 / 1.4	1.1 / 1.5	1.1 / 1.6	1.3 / 1.6	1.3 / 1.7	1.3 / 1.8	1.6 / 1.9	1.6 / 2.0	1.6 / 2.1	2.2 / 2.6	1.5 / 2.4



Flame	290.2
Coarse	
Medium	2015148
Fine	
Very Fine	
Minor/Major Diameter	1.6

KS Cylinder	KSO	KS1		KS2	KS3	
(Round Edge) Coarse	574.7 2015953	575.7 2015302	585.5 2015305	585.8 2015303	586.8 2015307	587.8 2015304
Fine	2015954	2015322		2015306		
Minor/Major Diameter	1.0	1.2	1.3	1.4	1.6	1.8












All-Ceramic Guide Pin

- Guide pin 0.5mm chamfer protects biological width
- Controlled tooth preparation with uniform cutting depth
- Eliminate lip or "j" margin defects




Guide Pin	SE740.8	SE259.8	SE259.2
Medium	2015966	2015964	2015993
Fine	2015967	2015965	2015994
Minor/Major Diameter	1.5/2.1	1.6/2.1	1.0/2.0

Assorted Crown Pack

(Chamfer/Shoulder)

Assorted 5-pack with our most popular Crown and Bridge Diamonds **2015500**



Curettage Modified Chamfer / Sub-gingival Preparation

- Recommended preparation design for Procera® crowns
- Provides clearly defined margins for CAD/CAM preps
- Popular for gingival curettage during crown preparation when margins are sub-gingival




	- Average of the second	-	- Francouran	
Curettage "K"	763.8	763.10	764.8	764.10
Coarse	2015977	2015979	2015978	2015980
Minor/Major Diameter	1.2	1.2	1.4	1.4

Modified Beveled Cylinder	510.6	510.8	511.8	511.10
Coarse	2015973	2015974	2015975	2015976
Minor/Major Diameter	1.0	1.0	1.2	1.2



Gold-Plated	**************************************				-	-					
Curettage (GCP™)	252.SA	252.SB	252.SC	253.SA	253.SB	253.SC	254.SA	254.SB	254.SC	257.SB	258.SB
Medium	2015603	2015605	2015607	2015609	2015611	2015613	2015615	2015617	2015619	2015623	2015627
Minor/Major Diameter	0.4 / 1.4	0.4 / 1.4	0.4 / 1.8	0.4 / 1.6	0.4 / 1.6	0.4 / 1.8	0.4 / 1.8	0.4 / 1.8	0.4/2.2	0.4/2.0	0.4 / 2.2



Premier is a registered trademark of Premier Dental Products Company. / Two Striper, P.B.S., MFS, Clean-A-Diamond and Thin-Flex are registered trademarks and TS2000, ShortCut, MicroPrep, Crown Cut, TGE, GCP, Silverstreak, Mini-Squares and Skooter are trademarks of Abrasive Technology, Inc. / U.S. Patent 5, 511, 718





Shoulder Prep/Butt-Joint Margin

 Choose flat-end tapered or flat-end cylindrical diamonds to create shoulder preparations with a 90° margin.

	-	The state of the s						"L" 1.5n long	jer 🗒	"L" 1.5n long shai	jer 🖺					
Flat-End Taper	700.5	700.6	700.8	700.9	700.11	701.5	701.7	L701.7	701.9	L701.9	702.8	703.8	703.9	703.10	704.9	708.3
Coarse	2015332	2015333	2015334	2015335	2015337	2015338	2015339	2015656	2015341	2015657	2015343	2015345	2015346	2015347	2015348	2015350
Medium																2015445
Fine				2015365	2015367		2015369		2015371		2015373	2015375		2015377		2015380
Very Fine				2015395	2015397						2015403					
Minor/Major Diameter	0.6 / 1.1	0.7 / 1.2	0.6 / 1.8	0.6 / 1.4	0.7 / 1.8	0.8 / 1.2	1.1 / 1.8	1.1 / 1.8	0.9/1.8	0.9 / 1.8	0.9 / 1.4	1.2 / 1.8	1.1 / 1.8	1.1 / 1.8	1.3 / 1.8	1.3/2.1

															-	
Flat-End Taper	708.4	721.6	721.8	721.10	722.6	722.8	722.10	723.4	723.6	723.8	723.10	724.8	724.10	725.10	726.10	727.10
Coarse	2015351	2015354	2015355	2015356	2015357	2015358	2015359	2015254	2015255	2015256	2015257	2015259	2015260	2015446	2015448	2015450
Fine				2015386			2015389			2015262				2015447		
Minor/Major Diameter	1.3 / 2.0	0.8 / 1.1	0.8 / 1.2	0.8 / 1.3	1.0 / 1.4	1.0 / 1.5	1.0 / 1.6	1.5 / 2.0	1.3 / 1.6	1.3 / 1.7	1.3 / 1.8	1.5 / 2.0	1.5/2.1	0.9 / 1.8	1.1/2.0	1.4/2.3

Flat-End Taper	798.10	Flat-End Cylinder	513.4	513.5	514.3	514.4	514.5	514.7	515.5	515.7	515.8	516.7	520.4
riat-Liiu rapei	7 30.10	riat-Liiu Cylliuei	313.4	313.3	J 17.J	717.7	317.3	317.7	313.3	313.7	313.0	310.7	320.4
Coarse	2015442	Coarse	2015271	2015272	2015273	2015274	2015276	2015278	2015275	2015277	2015279	2015985	2015281
Fine	2015443	Medium	2015288										
Minor/Major Diameter	0.5 / 1.5	Fine			2015293	2015294			2015295	2015297			2015301
	848-014	Minor/Major Diameter	1.0	1.0	1.1	1.1	1.1	1.1	1.2	1.2	1.3	1.6	1.4

TS2000												
Flat-End Taper	2000.7	2000.8	2000.9	2000.10	2001.7	2001.8	2001.9	2001.10	2002.7	2002.8	2002.9	2002.10
Coarse	2000007	2000008	2000009	2000010	2000017	2000018	2000019	2000020	2000027	2000028	2000029	2000030
Minor/Major Diameter	0.9/1.4	0.9 / 1.4	0.9/1.4	0.9 / 1.4	1.1 / 1.6	1.1 / 1.6	1.1 / 1.6	1.1 / 1.6	1.3 / 1.8	1.3 / 1.8	1.3 / 1.8	1.3 / 1.8
TS2000 Flat-End Cylinder	2013.7	2013.8	☐ ☐ 2013.10									
Coarse	2000136	2000138	2000140									
Minor/Major Diameter	1.5	1.5	1.5									










KR Modified Shoulder Prep

- Creates a margin with rounded internal angles
- Flat-end tapered or parallel diamond instruments featuring a rounded (beveled) edge









Bevel Shoulder Prep

- Recommended for porcelain-fused-to-metal or cast gold crowns
- Establish a 45° 60° angle
- Place bevel finish line



Bevel-End							
Cylinder	244.10	248.8	250.8	250.9	250.11	251.8	255.8
Coarse	2015064	2015065	2015067	2015066	2015068	2015069	2015073
Fine		2015085	2015087	2015086		2015089	2015093
Minor/Major Diameter	0.5 / 1.3	0.4 / 1.1	0.4 / 1.2	0.4 / 1.2	0.4 / 1.2	0.4 / 1.2	0.4 / 1.8

	200	
TS2000 Bevel-End	Ĥ	
Cylinder	2016.8	2016.10
Coarse	2000168	2000170
Minor/Major Diameter	1.4	1.8








Flame/Feather Edge Prep/Break Contact

- Popular design for cast-gold crowns
- Margin design for tight anterior crowns and tipped or perio-involved teeth
- Thin-shaped diamond for breaking interproximal contact


		-								"L" 1.5m long shan	jer 📗		
Flame	242.6	243.6	253.8	253.10	256.8	256.9	260.3	260.6.5	260.8	L260.8	260.10	261.8	262.6.5
Coarse	2015059	2015061	2015070	2015071	2015062	2015075	2015121	2015122	2015123	2015133	2015125	2015124	2015983
Medium									2015956				
Fine	2015060		2015090	2015091		2015095		2015152	2015153	2015163	2015155	2015154	
Very Fine								2015182	2015183		2015185		
Minor/Major Diameter	0.9/1.4	1.2 / 1.8	1.0 / 1.6	1.0 / 1.6	1.2 / 1.8	1.2 / 1.8	0.9/1.0	0.4 / 1.3	0.4 / 1.4	0.4 / 1.4	0.5 / 1.4	0.4 / 1.1	1.0 / 1.3

					A "L" 1.5m long — shan	ım A er A k —	Â		
Flame	262.8	263.8	265.6.5	265.8	L265.8	270.5	270.6.5	270.9	275.9
Coarse	2015126	2015127	2015129	2015131	2015654	2015135	2015137	2015139	2015955
Fine		2015157	2015159	2015161	2015655	2015165	2015167	2015169	
Minor/Major Diameter	1.1 / 1.7	0.4 / 1.6	0.5 / 1.4	0.4 / 1.4	0.4 / 1.4	0.5 / 1.6	0.4 / 1.9	0.5 / 1.9	0.4 / 2.2

Interproximal		
Trimmer	200.3	201.3
Fine	2015042	2015043
Very Fine	2015052	2015053
Minor/Major Diameter	0.3 / 0.8	0.3/0.8



Break Contact/Interproximal Reduction













Occlusal & Lingual Reduction (Bulk)

- An anatomically prepared occlusal surface provides adequate clearance and allows for uniform thickness of the restoration.
- Football and wheel-shape diamonds reproduce lingual concavity in crown preparation.

Football	283.4	285.5	290.4	292.3
Coarse	2015142	2015143	2015144	2015145
Fine	2015172	2015173	2015174	2015175
Very Fine	2015202	2015203	2015204	2015205
Minor/Major Diameter	2.5	0.7/2.3	0.7 / 1.9	0.4 / 1.6

Depth - mm	1.0	1.5	1.5 ", "	1.7	2.0	2.0 ",	, 2.4	2.8
	Ä		@ 1.5	mm 👸 ger 📙		(3) 1.5	mm () iger ink	
Round	115	120	L120	125	130	L130	135	138
Coarse	2015003	2015005	2015055	2015006	2015008	2015078	2015007	2015009
Fine	2015013	2015015		2015016	2015018		2015017	
Very Fine		2015025						
Minor/Major Diameter	1.2	1.6	1.6	1.8	2.3	2.3	2.6	3.0





Occlusal Contour	204.3.5
Fine	2015952
Minor/Major Diameter	3.10

Wheel	860	862	863	866
Coarse	2015543	2015544	2015545	2015546
Fine	2015553		2015555	
Minor/Major Diameter	3.5	4.1	4.7	5.7




Veneer Prep




	D	Depth of cut (in mm):							
	0.5	0.75	1.0	1.5					
D.C. Depth-Cut	DC.5	DC.75	DC1.0	DC1.5					
Coarse	2015776	2015785	2015802	2015803					
Minor/Major Diameter	1.2	1.2	1.2	1.2					












Inlay/Onlay

- Flat-end diamonds with a rounded corner create rounded internal angles for ceramic/composite inlay and onlay preparation.
- Tapered instruments establish proper divergence of axial walls.

	Exercise 1	"L" 1.5 lon sha	mm ger						"L" 1.5n long Shar	jer 📗	"L" 1.5 lon sha	·		
Flat-End Cylinder	513.4	L513.4	513.5	514.3	514.4	514.5	514.7	515.5	L515.5	515.7	L515.7	515.8	516.7	520.4
Coarse	2015271		2015272	2015273	2015274	2015276	2015278	2015275	2015797	2015277	2015798	2015279	2015985	2015281
Medium	2015288	2015289												
Fine				2015293	2015294			2015295		2015297				2015301
Minor/Major Diameter	1.0	1.0	1.0	1.1	1.1	1.1	1.1	1.2	1.2	1.2	1.2	1.3	1.6	1.4

KS Cylinder	KS0	KS1		KS2	KS3	
(Round Edge) Coarse	574.7 2015953	575.7 2015302	585.5 2015305	585.8 2015303	586.8 2015307	587.8 2015304
Fine	2015954	2015322		2015306		
Minor/Major Diameter	1.0	1.2	1.3	1.4	1.6	1.8

Modified Shoulder				
KR Cylinder	515.6	521.4	521.6	521.8
Coarse	2015986	2015987	2015989	2015990
Fine			2015988	
Minor/Major Diameter	1.2	1.4	1.4	1.4



TS2000				TS2000		Modified houlder
Flat-End Cylinder	2013.7	2013.8	2013.10	KR Cylinder	2013.8KR	
Coarse	2000136	2000138	2000140	Coarse	2000139	
Minor/Major Diameter	1.5	1.5	1.5	Minor/Major Diameter	1.0/1.6	

Specialty TGE™ **TGE**TM (Tissue Guard • Refine **End-Cutting)** 1.0 1.2 1.4 1.6 shoulder 2015836 2015838 2015840 2015842 Fine • Safe-sided 2015839 2015841 **Very Fine**

1.0

1.2

1.4

1.6

Minor/Major Diameter

Skooter™	0	1	2
Medium	TO - 2015188	T1 - 2015190	T2 - 2015192
Very Fine	FO - 2015189	F1 - 2015191	F2 - 2015193
Minor/Major Dia	ameter 2.3	2.8	3.0

Gross Reduction						HARMAN
Crown Cut™	SC5	SC8	SC10	ST6	ST8	ST11
Coarse	2015573	2015575	2015577	2015581	2015583	2015585
Minor/Major Diameter	1.7	1.7	1.8	1.4/1.8	1.4/1.8	1.4/1.8













Operative


- Featuring diamond versions of popular carbides
- Cut cooler and faster with less patient trauma
- Inverted cones for removing defective restorations

Inverted Cone	310.1	315.1.75	317.4	318.5	320.2	324.1	390.3
Coarse	2015213	2015790	2015984	2015991	2015217	2015219	2015931
Fine	2015233				2015237		
Minor/Major Diameter	0.9/1.3	1.1 / 1.6	1.6	1.8	1.3 / 1.7	1.1 / 2.1	1.2

							e de la companya de l							
Operative	1/4	1	2	6	35	169L	170L	245	330	556	557	558	701	1701
Medium	2013203	2013205	2013201	2013211	2013221	2013231	2013241	2013251	2013261	2013271	2013281	2013291	2013301	2013311
Minor/Major Diameter	0.7	0.8	1.1	1.9	1.0	0.6 / 1.1	0.7 / 1.2	0.9/1.1	0.8/0.9	0.8/0.9	0.9/1.0	1.2 / 1.3	1.0 / 1.4	1.0 / 1.2
	H1-005	H1-008	H1-010	H1-018	H2-010	H23L-009	H23L-010	H245-008	H24-008	H31-009	H31-010	H31-012	H33-012	H33R-012
Head Length mm	0.7	0.6	0.8	1.3	0.9	5.2	5.8	4.3	1.6	3.3	3.9	4.1	4.1	4.5

						"L" 1.5n long shar	jer 🖺		1.5mm longer
Pear	351.3.75	352.4	360.2	361.2.5	362.3	L362.3	364.5	365.4	L365.4
Coarse	2015229	2015230	2015794	2015796	2015792	2015799	2015227	2015225	2015228
Fine					2015793			2015245	
Minor/Major Diameter	1.1	1.3	0.8 / 1.0	1.0 / 1.2	1.1 / 1.3	1.1 / 1.3	1.0 / 1.4	1.2 / 1.5	1.2 / 1.5

	X			
Amalgam Removers	1210	1220	1240	1250
Medium	2015591	2015593	2015598	2015599
Minor/Major Diameter	0.9 / 1.3	0.9/1.6	0.9 / 1.1	0.9 / 1.3







A New Diamond Technology for Today's Advanced Materials - Lithium Disilicate and Zirconia

2015895 **TSZtech™ Diamond Assorted 5-pack** (260.8Z, 285.5Z, 703.8KRZ, 125Z, 770.8Z)

Our newest Two Striper® has sporty orange stripes indicating their proprietary diamond grit designed to cut safely through Lithium Disilicate and Zirconia crowns and gain endodontic access. TSZtech™ diamonds provide faster and cooler cutting with less trauma and no sparking. Some shapes are available in an "electric friendly" longer shank for electric handpieces with larger heads.

2	60.8Z	285.5Z 7	703.8KRZ	125Z	770.8Z
5-Pack Reorder 20	15893 2	2015891	2015894	2015892	2015890

"L" - 1.5mm longe	er shank:		
	L260.8Z	L125Z	L770.8Z
5-Pack Reorder	2015883	2015882	2015880





Composite Finishing Kit

Kit # 2013557

Contains eight diamonds and two carbides shown here. Anodyzed aluminum bur block has eight open spaces to add your own favorites.


5-pack refill numbers listed below individual shape illustrations.



Use yellow-striped (very-fine grit) and violet-striped (MFS diamonds) in a light wiping motion with copious water spray.

*No item number. Available in kit only.

Refine facial anatomy



Follow with finer grit for more subtle textures

Ħ
260.8
2015183

Refine subgingival margin

Interproximal Trimmer 200.3 Fine 2015042

Follow with finer grit for smoother margins

Interproximal Trimmer 200.3 **Very Fine** 2015052

Secondary occlusal anatomy

Football	292.3
Fine	2015175

Follow with finer grit for more refined secondary anatomy

Football	 292.3				
Very Fine	2015205				

Refine lingual anatomy

Round 130 **Fine** 2015018

Follow with finer grit for smoother lingual surfaces

MFS® Finishing Diamonds 135 MF1 (45µ) 2015850

12-fluted carbide football for primary occlusal anatomy

Carbide Football



12-fluted carbide flame for facial or occlusal anatomy

Carbide Flame H246-009*

Universal All-Ceramic Prep Kit



Kit # 2013555


Contains ten diamonds shown here. Anodyzed aluminum bur block has eight open spaces to add your own favorites.

5-pack refill numbers listed below individual shape illustrations.

All-Ceramic Prep / Shoulder

Depth Gauge / Occlusal / Lingual Reduction TS2000 Football





Axial Reduction / **Round Internal Angles**



All-Ceramic Prep / Chamfer



Taper	770.8				
Coarse	2015457				
Fine	2015487				

Seat and Finish

Margin Finishing / **Cement Removal**

Flame	263.8
Fine	2015157



Coarse

200.3
2015042

Occlusal Refinement / Adjust Bite

Football 285.5
Fine 2015173

Fine Occlusal Adjustment

Football	292.3
Fine	2015175



Rotary


Diamond Wear Factors

The unique bonding and consistent cutting action can make it difficult to determine when a Two Striper® should be discarded. Other electroplated diamonds wear by losing diamond particles and voids are easy to detect (and diminished cutting easy to feel). The Two Striper wears flat, but the diamond grains remain intact.

All diamond instruments wear preferentially at the tip. The tip also is susceptible to greater vibration and often exposed to excessive force. On very thin-tip shapes (flames, pointed cones, etc.), there is reduced space for diamond concentration so the cutting load is magnified. A clogged diamond often prompts the operator to apply more pressure, but the combination of heat and force is very detrimental.

Two Striper diamond wear can be determined by following the cleaning technique below. If rapid cutting action does not return after this protocol, the diamond crystals are worn. Wear also can be determined visually by examining the diamond tip and circumference. If the diamond's crystals remain extended above the bond, the instrument still has life. If worn flat, the instrument should not be reused.

Used Two Striper® Tip



Used Electroplated Tip



After use, you can see that the diamond crystals are still present on the Two Striper®; however, many of the crystals have pulled out from the electroplated diamond.

Two Striper® Diamond Cleaning Technique

Prolong the cutting life of your diamonds by following the simple steps below. Recommended for coarse and medium grit preparation burs, the Clean-A-Diamond® Mini-Square and Brite Shield™ enzymatic cleaner together will help keep your edge. Always use a light touch with lots of water to maximize cutting efficiency.


2014015

CLEAN - ↑ - DIAMOND®


Clean-A-Diamond® Mini-Squares™ (Package of 12)

9011100


Brite Shield™ (800g jar with scoop dispenser)



Hold handpiece with used/clogged Two Striper® diamond just above Clean-A-Diamond® Mini-Square™.



Make a plunge cut, down sideways with handpiece and the waterspray on (never tip first). Only the diamond head section should come in contact with the stone. Avoid contact between the stone and instrument neck.



Lift rotating diamond from the surface of the stone. $Clean\hbox{-} A\hbox{-} Diamond @Mini\hbox{-} Square \hbox{$^{\text{TM}}$}$ can be autoclaved after use.



Further clean and protect your Two Striper® diamonds from rust and corrosion by adding Brite Shield™ enzymatic cleaner to your ultrasonic.



For Hard Tissue Removal: Add one scoop of Brite Shield $^{\text{TM}}$ to 1.9L (64oz) of warm water.



Stir until dissolved.



Immerse instruments in solution and clean in ultrasonic unit for 1-5 minutes.



Rinse instruments thoroughly with water and dry before sterilizing.

Premier is a registered trademark of Premier Dental Products Company. / Two Striper, P.B.S., Compo-Disc, Compo-Strip, Thin-Flex, Clean-A-Diamond, Luminescence and MFS are registered trademarks and ™ Luminescence Plus, Skooter, CrownCut, TGE, MicroPrep, Mini-Squares, TS2000, GCP, SilverStreak, ShortCut and TSZtech are trademarks of the manufacturer, Abrasive Technology, Inc. • U.S. Patent 5, 511, 718 / ® Procera is a registered trademark of Nobel Biocare.



Rotary

Surgical Length Carbide Burs (FG Surgical 25mm) *Made in Germany*

1 per package

interp	#2 for :	als.În er	urface cavit ndodontics,		Cross Cut Fissure Straight To produce parallel walls and flat floors in preparations					
1 008 SFG	2 010 SFG		4 018	8 023 SFG		556 008 SFG	557 010 SFG	558 012 SFG	559 014 SFG	560 016 SFG
T	id unde	rcuts a	aper nd create rity walls	Breakthru Burs To produce undercuts		Inverted Cone For large occlusal restorations to insure rounded internal line angles				
slightly	aiverg	ern cav	nty wans	unaercus	•				iiiciii	ai
700 010 SFG	701 012 SFG	702 016 SFG	703 021 SFG	1558 014 \$F6	402 016 SF6		angles		THE THE	ui
7000010	701 012	702 016	703 021	Straight 1558	402 016	34 008	angles B	35 010		701
700 010 SFG	701 012 SFG	702 016	703 021 SFG	1558 014 \$FG	402 016 SFG	344 008 SFG	angles B S	35 010 SFG	26	
700 010 SF6	701 012 SFG	702 016	703 021 SFG	1558 014 SFG	402 016 SFG	344 000 SFG	angles B G	35 010 SFG	26 27	701
700 010 8F6 9050601 9050602	701 012 SFG	702 016	703 021 \$FG 9050607 9050608	1558 014 SF6 34 35	402 016 SFG 9050619 9050620	344 008 SFG 557 558	angles B	35 010 SFG 90506	26 27	701 702

Regular Steel Burs (HP 44.5mm, RA 22mm) *Made in Germany*


10 per package

			CONTW.			
	Round	Inv. Cone	CC Sq. Fis.	CC Tap. Fis.	LH Round	Wheel
НР	1/2, 2, 4, 6,8		558, 562	700	10, 12	
905	52053 2	9052059	8 9	9052061 10		
905	2055 4	9052083	562	9052063 12		
905	52057 6	9052086	700			
RA	1/2, 2, 3, 4, 5, 6,7,8	35, 37				14
905	52001 1/2	9052005	4	9052009 8		
905	2003 2	9052006	5	9052020 35		
905	3	9052007	6	9052022 37		

Evian Perio Burs

Surgical Length FG 25mm for Periodontal Debridement

- Blunt, tissue-kind heads
- Removes granulation tissue efficiently
- Removes deposits from root surfaces
- Easier technique than the tedious task of hand scaling



Package of	six	
9052041	2	
9052042	4	
9052043	6	
9052040	Assorted (2 of each size)	



Rotary/Laboratory





Double-sided cross-section



thin-flex®

Flexible Diamond Discs with True Edge Cutting

Designed for carving ceramic and composite materials. Thin-Flex® discs are excellent for contouring and shaping all surfaces including embrasures. Thin-Flex discs have diamond crystals wrapped around the edge of the disc to avoid "black marks" during carving procedures.

Recommended Speed - 20,000 rpm Outside Diameter - 22.2mm (7/8") Center Hole Diameter - 1.6mm (1/16")


Single-Sided - 1/pkg. 2012633 926-7 (0.10mm, 45 μ) 2012637 928-7 (0.24mm, 60 μ)

Single-Sided - 2/pkg. 2012634 926-7 (0.10mm, 45 μ) 2012638 928-7 (0.24mm, 60 μ)

Double-Sided - 1/pkg. 2012653 927-7 (0.15mm, 45 μ) 2012657 929-7 (0.28mm, 50 μ)

Double-Sided - 2/pkg. 2012654 927-7 (0.15mm, 45 μ) 2012658 929-7 (0.28mm, 50 μ)

SILVERSTREAK™ Diamond Model Trimming Wheel



The SilverStreak™ trims models faster, runs quieter and is substantially longer lasting than traditional wheels. It has an extra coarse diamond cutting surface which is patterned for aggressive, yet cool cutting. It is lightweight, durable and perfectly balanced for concentric operation.

Thickness - 6.4mm (1/4") Center Hole Diameter - 25.4mm (1")

2013401 12" Extra Coarse 2013405 10" Extra Coarse

Laboratory HP Diamonds P.B.S.® Bonded - ISO Type 2 Shank 2.35mm

	X115	X118	X120	X135	X150	X261.5	X265.10	X267.10	X285.8.5	X324.1	X350.4	X535.7	X575.7
Coarse	2012699	2012700	2012701	2012703	2012705	2012707	2012708	2012738	2012709	2012711	2012713	2012719	2012728
Fine	2012747	2012746	2012731	2012733	2012789	2012737	2012787	2012786	2012739	2012741	2012743	2012749	2012758
Reference #	801-010	801-018	801-016	801-026	801-036	860-014	863-014	863-016	862-026	805-022	807-030	836-027	830-026
Minor/Major Diameter						0.4/1.4	0.5/1.4	0.4/1.6	0.3/2.6				

	X580.6	X590.12	X701.7	X701.9	X702.10	X703.12	X767.7	X769.10	X772.10	X775.10	X776.9	X884	X889
Coarse	2012720	2012682	2012721	2012723	2012759	2012724	2012725	2012762	2012765	2012766	2012726	2012770	2012727
Medium		2012681											
Fine	2012750	2012680	2012751	2012753		2012754	2012755	2012730	2012736	2012798	2012756		2012757
Reference #	836-052	842R-012	846-018	847-016	848-020	848L-025	855-018	850-018	850-025	850L-018	856-040	825-046	825-067
Minor/Major Diameter			1.0/1.8	0.8/1.8	1.2/2.0	1.4/2.5	1.0/1.8	0.7/1.8	1.4/2.5	1.0/1.6	2.9/4.0		



Rotary/Laboratory

Lab Tri-Fissure HP Carbide

H2803-012HP

Ideal for sculpting porcelain fissure anatomy before applying stain

9050655 Package of 5



Lab Mandrels

Aluminum Bronze or Stainless Steel

These gold-color heads are made of specially formulated aluminum bronze alloy which results in minimal wear on the threads, resists rusting and lasts up to six times longer than conventional steel mandrels. The shafts are made of tough stainless steel and do not bend.

.25mm thread length



Package of 12 Aluminum Bronze 9048051 RA Shank 9048053 HP Shank

Heavy Duty Stainless Steel 2014011 HP Shank

Bur Blocks

Autoclavable, Anodized Aluminum

- Patented design features center stabilizer plate
- · Secure snap-close lid



Friction Grip (18 Hole)

9051021 Bur Block - Green **9051022** Bur Block - Blue

Right Angle/Friction Grip (17 Hole)

9051024 Bur Block 8 RA / 9 FG - Green **9051025** Bur Block 8 RA / 9 FG - Blue



14-Hole HP Bur Block 9051030

Clinical Favorites







Finishing/Polishing





Finishing/Polishing

Poli-Pro Disks™

Composite Polishing System

- Safe-Center protects composite surface
- Faster complete abrasive coverage of disk surface
- Thin and flexible conforms to tooth anatomy
- Systematic finishing and polishing single-use color-coded disks
- Two sizes:
 1/2" for greater efficiency,
 5/16" for control in gingival areas

15% more surface area than metal-centered disks!



1/2" Coarse (12mm) Refill #2019071	5/16" Coarse (8mm) Refill #2019072	1/2" Medium (12mm) Refill #2019073	1/2" Medium (12mm) Refill #2019074	5/16" Medium (8mm) Refill #2019075
1/2" Fine (12mm) Refill #2019076	5/16" Fine (8mm) Refill #2019077	1/2" Extra-Fine (12mm) Refill #2019078	5/16" Extra-Fine (8mm) Refill #2019079	RA Mandrels (4/pack) Refill #2019090
Coarse - 80 microns (grit up))			FG Mandrels (4/pack) Refill #2019091

Medium - 40 microns (grit up)

Medium - 40 microns (grit down)

Fine - 30 microns (grit both sides)

Extra-Fine - 9 microns (grit both sides)



2019070 Intro Kit (180 disks, 4 mandrels)

2019099 Starter Kit (35 disks, 1 mandrel)



Faster! Super-Charged Diamond Polishing

Single-paste System for Ceramics and Composites

Diamond Twist SCL™

Extra-Oral Polishing Kit

Glaze-like, multi-surface polishing - featuring advance open-weave fibra pre-polishers and 100% wool paste applicators - for a super-charged mirror polish

- · Pre-polish with Fibra Points (barrel, wheel or taper); open-weave design generates significantly less heat.
- Final polish with Diamond Twist Super-charged Polishing Paste applied with felt barrel, wheel or point for a glaze-like polish.







2019000

Diamond Twist SCL™ Kit

Products below shown 75% of actual size.

2019001 6g Diamond Twist

SCL Polishing Paste (unflavored)

2019033

Point.

Barrel

0470

Fibra Wheel,







Diamond Twist SCL

Super-Charged Polishing Posts















Finishing/Polishing

Diamond Twist SCO™

Intra-Oral Polishing Kit

After permanent seating and adjustments, polish crowns, inlays, onlays and nano-hybrid composites intra-orally in one easy step using Diamond Twist Super-charged Polishing Paste.

- · Use the flat or taper polishing brushes (autoclavable) for posterior occlusal surfaces.
- Use the Micro-cloth disks for polishing and final high-gloss buffing.
- · Twist syringe avoids cross-contamination and dispenses just the right amount.







2019002





2019005

2019091

2019060

Dark Blue (25/pack)

Tapered Polishing Brushes (25/pack)





1/2" Micro-Cloth Disks

2019004 Flat Polishing Brushes (25/pack)















Two Striper® MFS® **Finishing Diamonds**








MFS®

MF1 (45μ)

MF2 (20µ)

MF3 (10µ)



Skooter™	0	1	2
Medium	TO - 2015188	T1 - 2015190	T2 - 2015192
Very Fine	FO - 2015189	F1 - 2015191	F2 - 2015193
Minor/Major Dia	ameter 2.3	2.8	3.0
Reference #	905-023	905-028	905-031

Compo-Disc®

Hand-Held Autoclavable Diamond Finishing Instruments

For extremely tight contacts, use a Compo-Disc® to gain access. Safe-sided, edge and face-cutting.

Thickness: 0.1mm

2013106

Compo-Discs (Pkg. of 2 discs)



Compo-Strip®

Hand-Held Autoclavable Stainless Steel **Diamond Finishing Instruments**

Compo-Strip® diamond discs are thin, flexible, one-sided diamond abrasive strips available in two widths and three grits. They have an uncoated safe center for easy access through the contact point and are ideal for finishing Class III, IV and veneer restorations.

2013127 Package of Six Assorted Strips

Compo-Strip® **200T2** 60µ For Orthodontic Reproximation



Simultaneously strip away enamel from adjacent proximal surfaces.

Flexible, Double-sided with Safe Center

2013180 Compo-Strip 200T2 (Ortho) (Pkg. of 3 strips)



Individual Packages of Six Strips 100T **2013107** 100F 2013111 100UF 2013115

150T **2013109**

150F **2013113**

150UF **2013117**

Width
100 = 2.5mm
150 - 3.75mm

Grit	Thickness
$T = 60\mu$	0.15mm
$F = 45\mu$	0.127mm
UF = 20µ	0.1mm



CAD/CAM Laboratory Finishing Kit

Ceramic Adjustment

Packaged in Anodized Aluminum Lab Bur Block 2013552 Two Striper CAD/CAM Finishing Kit



	Ä	
Flame	X769.10	
Fine (1/pkg.)	2012730	

Contour




Round-End Taper	X590.12
Medium (1/pkg.)	2012681
Medium (1/pkg.)	2012681

Remove Sprue



Thin-flex® Disc	X929-7
Double-sided (1/pkg.)	2012657

Carve **Occlusal Anatomy**



Miniature Muslin **Buff**







Available in kit only.





LUMINESCENCE®

Single-Gel Diamond Polishing

Polish all surfaces in half the time!

Luminescence is a patented, single-gel diamond polishing system. In less than two minutes, you can polish the surface of all restorative materials: composite, glass ionomer, compomer, amalgam, precious metal and enamel. Even porcelain achieves a lustrous surface in only minutes due to the optimal concentration of micron-sized diamond particles.

2014200 Luminescence Intro Kit (3g syringe of Luminescence gel, 2 RA Mandrels, 50 applicators)

2014202 Luminescence Gel (3g)

2014204 RA Mandrels (2)


2014206 Applicators (50)






"Optimum Polishing results on the CAD/CAM Vita Mark II ceramic." 1

"Adjustment finishing with a Two Striper MF2 or MF3 alone, or as a final finishing stage after using a more coarse finishing diamond, followed by a 60-second polishing sequence with Luminescence diamond gel, was an effective procedure for optimum polishing results on the CAD-CAM Vita Mark II ceramic."



Vita Mark II Ceramic Pre-finish surface reference



Vita Mark II
Final surface appearance
after MF2 Diamond finish and
Luminescence application

LUMINESCENCE PLUS!

Diamond Polishing Paste for Patient Comfort + All-Surface Esthetics





- All-Surface Polishing composite + porcelain + enamel + glass ionomer + precious metals + amalgam
- Single-Grit/Single-Step Paste saves chair-time
- Increase Patient Comfort contains potassium nitrate topical desensitizer
- Convenient screw-tube dispenser + improved viscosity eliminates waste

2014208 Luminescence Plus[™]! Gel (3g) (Not available in Canada)

Traxodent









Dual-arch technique for 1-2 unit impressions

1

Try the appropriate Triple Tray to ensure it fits passively when the patient closes.



2

Apply adhesive to tray walls and overload both sides of the Triple Tray with a firm-setting impression material.



3

Remove cord, if used, and syringe impression material onto critical areas.



4

Place loaded Triple
Tray into the mouth
and have patient close
properly. Check contralateral side to verify
proper occlusion.



5

Inspect the final impression for detail and to ensure the patient has not closed into the sides or distal bar of the tray.



Triple Tray®

Proven Accurate In Over 90 Million Impressions

All Premier® Triple Tray® designs have very thin, loose webbing which will not distort when the patient closes.

Look for Premier on the handle, your assurance of a genuine Triple Tray.

The superior performance of Triple Tray, combined with today's very accurate materials, provides you with an ideal impression in an easy technique.

Triple Tray has been consistently rated an outstanding product by a leading independent research organization due to "least water sorption and thinnest interocclusal wafer..."

1006283	Wide Body Posterior - Box of 48
1006203	Posterior - Box of 48†
1006213	Anterior - Box of 35†
1006223	Quadrant - Box of 40†
1006243	Sideless Posterior - Box of 48 [†]
1006253	Extended Quadrant - Box of 35 [†]
1006273	Three-Quarter Arch - Box of 35 [†]
1006273	Three-Quarter Arch - Box of 35† Assorted - Box of 35
	Assorted - Box of 35
	Assorted - Box of 35 Contains (5) each: Posterior,
	Assorted - Box of 35 Contains (5) each: Posterior, Anterior, Quadrant, Sideless
	Assorted - Box of 35 Contains (5) each: Posterior, Anterior, Quadrant, Sideless Posterior, Extended Quadrant,

†Dealer: Available to you in case of ten.





Repeat Winner - 2003, 2004, 2006, 2007, 2008, 2009, 2010, 2011.





Premier's **Wide Body Posterior Triple Tray**, is 19% longer and 33% wider than the original posterior and captures detail of up to six teeth.















^{**}An independent non-profit dental education and product testing foundation, January 2001.



2007, 2008, 2009, 2010, 2011.

T-LOC™ Triple Tray®

Adhesive-free Dual-arch Impression Trays

- Patented adhesive-free design Retention tracks mechanically lock in material, eliminating the need for adhesives.
- Easier handling Large, ergonomic handles simplify loading of the tray, enable easy transfer between assistant and dentist.
- Thin, loose webbing, identical to Triple Tray, maximizes intercuspation of the teeth to ensure an accurate bite.
- Greater rigidity in posterior crossbar due to the use of a new, novel polymer.
- Ensures accuracy and saves time Low sides, thin posterior bar and wide impression areas avoid impingement when the patient closes.

Posterior - Box of 35
Anterior - Box of 28
Three-Quarter Arch - Box of 22
Posterior - Box of 35 [†]
Anterior - Box of 28 ⁺
Three-Quarter Arch - Box of 22 [†]


U.S. Patent #6,749,428 †Dealer: Available to you in case of ten.


NeoTray®

Rigid Impression Trays

- Gold Mesh[™] technology Tear-resistant material
- Exclusive KR[™] plastic Provides a rigid tray without the cost of metal disposable trays or cleaning and care of single-arch trays
- Cuts chair-time in half Simultaneously records master impressions, counter impression and bite registration
- Saves money Less chair time, impression material and corrections
- Increases accuracy Avoids costly errors and chairside adjustments
- Assortment of seven sizes available in fun colors that patients love.

2206NT	Introductory Kit - Box of 45 - Contains: 10 each of Wide Posterior, Sideless Posterior, 5 each of Anterior, Posterior, Full Arch, Quadrant and Extended Quadrant
2204NT	Anterior - Box of 30
2205NT	Full Arch - Box of 30
2202NT	Sideless Posterior - Box of 48
2208NT	Wide Posterior - Box of 48
2201NT	Posterior - Box of 48
2207NT	Extended Quadrant - Box of 35
2203NT	Quadrant - Box of 40







Wide Posterior



Traxodent®

Hemodent® Paste Retraction System

Provides predictable hemostasis and soft-tissue management in minutes!

Easy, effective hemostasis and retraction The sleek syringe with bendable tip permits easy application of Traxodent directly into the sulcus. After two minutes, it is rinsed, leaving an open, retracted sulcus.

Traxodent is gentle, absorbent and fast. The soft paste produces gentle pressure on the sulcus while it absorbs excess crevicular fluid and blood. The aluminum chloride creates an astringent effect without irritating or discoloring surrounding tissue. Traxodent provides predictable hemostasis and retraction in less time and with greater comfort.

Quick effective retraction

About 2 minutes - no fumbling with hemostatic solutions - save chair time

Superior hemostasis

Stops bleeding and crevicular seepage while retraction occurs - easily absorbs excess crevicular fluid

Increased patient comfort

Avoid tearing gingiva - less pressure on sensitive tissue

Convenient delivery system

No cumbersome and expensive guns or cartridges - ergonomic disposable syringe and bendable syringe tip provides excellent reach

Maximum freshness

Each syringe is packaged in resealable foil pouches

Excellent value

Ready to use, no hidden cost - kits include syringes and disposable tips

9007091 Traxodent Value Pack (25 syringes and 50 tips)

9007092 Applicator Tips (50)

9007093 Traxodent Starter Pack (7 syringes and 15 tips)

(E 0470



Traxodent applied



Photos courtesy of Dr. Dean Elledge, DDS, MS

Versatile!
Use alone, with cord or with retraction caps.

Dentists* have their say!

98% will recommend to a colleague

90% found Traxodent was easy to rinse

97% found Traxodent provided sufficient isolation and adequate hemostasis



Traxodent is highly effective for all your tissue management challenges.

Hemostasis and moisture control prior to composite restorations

Tissue management prior to veneer or crown cementation

Hemostasis after removal of a temporary crown

Post-op hemostasis after laser troughing

Before a CAD-CAM scan













Bend tip for optimal access.



traxodent

Position tip parallel to axial plane of the tooth.



traxodent.

Dispense sufficient material to displace tissue.



Place Premier® Retraction Cap on tooth structure.



Have the patient bite and maintain pressure on the cap. Wait 2 minutes.



Rinse away, leaving sulcus open and accessible.

Optional Cord Technique Traxodent eliminates the need for a second cord.



Place a single strand of retraction cord.



Dispense sufficient material to displace tissue. Wait about 2 minutes.



Rinse away, leaving sulcus open and accessible.





Premier® Retraction Caps

Anatomic formed cotton caps work synergistically with Traxodent to quicken and assist retraction and hemostasis.

- Easily guides Traxodent into the sulcus
- · Compression assists and quickens hemostasis
- · Cotton absorbs fluids
- Anatomic form allows easy placement

9048255	Size 5 (60)
9048253	Size 3 (120)
9048251	Size 1 (120)
9048260	Assorted - (30) #5, (15) #1, (15) #3





Hemodent®

Buffered Aluminum Chloride Epinephrine-free Hemostatic Liquid



For over 50 years, the world's most sought after hemostatic product.

Hemodent® effectively stops minor gingival bleeding. It contains no epinephrine to help avoid cardiac reactions. Hemodent is stable and offers a long shelf life.

Recommended to arrest bleeding during crown preparations, before impression taking, placing fillings, cement bases, or cementing inlays or crowns; also after minor surgery, curettage, gingivectomies, gingivoplastys and deep scaling.

Packag	e of	10
--------	------	----

900/0/1	Tucc bottle*
9007072	20cc bottle*
9007073	40cc bottle*

^{*} Individual bottles available through your dealer.



Knit-Pak™

Knitted Retraction Cord Non-impregnated



Knit-Pak™ gingival retraction cord is a non-impregnated knitted cord made from 100% cotton.

Knit-Pak is ideal to use in combination with Hemodent® hemostatic agent or Traxodent for tissue displacement and management prior to taking an impression. The knitted cord contains numerous interlocking chains that carry greater quantities of hemostatic solution compared to twisted or braided cord. Knit-Pak easily packs into sulcus by compacting under pressure. It will not fray and resists popping out of the sulcus. If touched by a bur, Knit-Pak will cut clean and not get entangled.

Knit-Pak's 6 different sizes are color coded for quick and easy identification.

Each Knit-Pak contains 100 inches (254 cm) of cord.

	Size	Color
9007551	000	Green
9007552	00	Brown
9007553	0	Purple
9007554	1	Blue
9007555	2	Orange
9007556	3	Yellow





A stainless steel cutting blade inside the cap automatically cuts the cord upon closing.

FS Hemostatic™

Ferric Sulfate Epinephrine-free Hemostatic Liquid



A 15.5% ferric sulfate liquid solution offers instantaneous coagulation for difficult hemostatic control. It contains no epinephrine so it is safe for use on all patients.

9007057	1oz dropper bottle (29.5ml)
9007058	2oz refill bottle (59ml)



W.O.W.™ & W.O.W.™2

Accurate Disposable Articulator Quick and Simple

Self-articulated and precise models in two fast pours!

The W.O.W. and W.O.W.2 are innovative single-use articulators. Now you can articulate the case exactly as the dentist took the impression. Use for all types of cases, dual-impression trays, pre-poured models and conventional impressions.

Intelligent design provides accurate dies and models.

- Clear plastic view allows for accurate pin placement.
- Pin placement design eliminates rotation.
- Shallow retention grooves minimize dimensional change from stone expansion.
- Fixed hinge allows for positive occlusal relationship between upper and lower models and appropriate lateral movement.
- Remount bite, if necessary, with heat in less than a minute.

W.O.W.2 features Flash Air Dam™ technology

- Keeps unwanted debris from lodging in pin holes
- Keeps cases neat and delivers professional presentation

0069	W.O.W. Intro Kit - Contains: 30 Posterior, 20 Anterior/Full Arch and 300 dowel pins
0065	W.O.W. Anterior/Full Arch Box of 50 and 300 dowel pins
0065AD	W.O.W.2 Anterior/Full Arch Box of 50 and 300 dowel pins
0072	W.O.W. Posterior Box of 50 and 300 dowel pins
0072AD	W.O.W.2 Posterior Box of 50 and 300 dowel pins
0067	Dowel pins Bag of 300
0020	Die Maker Bag of 50







Die Maker



Posterior Bite Tray



Especially useful in close-bite registration cases as the patient will not bite through the material

Box of 50 9048076 Posterior Trays

Prosthetic



Alfa™ Triple Tray®

For Full-arch Closed-mouth Alginate Impressions




- · Reduces gagging
- Cuts chair time in half simultaneously captures upper and lower impressions
- Saves money requires less alginate and no clean-up
- Increases efficiency closed-mouth technique frees operator from holding tray
- Ideal for alginate impressions


U.S. Patent #6,749,428

"Pearls For Your Practice" Dental Economics - June 2003

Box of 24		
1006313	Small (purple)	
1006323	Medium (white)	
	. , ,	
1006333	Large (green)	
	<u> </u>	
1006343	1006343 Assorted Contains:	
10005 15	Assorted Contains.	
	(5) Small, (14) Medium, (5) Large	
	(3) Sindii, (17) Medidiii, (3) Laige	



Two Impressions:Mandibular and Maxillary





Perfecta® REV!® Hydrogen Peroxide Tooth-whitening Gel	Page 134	
Perfecta® Bravó® Hydrogen Peroxide Tooth-whitening Gel	135	
Perfecta® Carbamide Peroxide Tooth-whitening Gel	136-137	
Perfection White™ In-Office Hydrogen Peroxide Tooth Whitening	138	
Alfa™ Triple Tray® Full-arch, Closed-mouth Alginate Impressions	139	
Tray Magic™ Soft Tray Trimmer	140	
PREMA® Enamel Microabrasion System	141	

Hydrogen Peroxide Take-home Tooth Whitening

Hydrogen Peroxide is a powerful oxidizing agent that penetrates the enamel and dentin. The length of treatment depends on the level of discoloration and the whitening product used. The higher the concentration of hydrogen peroxide, the shorter the wear-time.

Stable and Strong Gel

Take-home tooth whitening is faster, easier and more effective with Premier's hydrogen peroxide products. Offered in strengths of 9% or 14% hydrogen peroxide, the patented water-based gel was developed with the patient in mind, to minimize sensitivity and assure patient compliance. Premier's products meet the demands of today's patient; fast results, short wear-times and reduced sensitivity.

Assure Patient Compliance

The leading cause of treatment failure with take-home whitening is patient non-compliance. Even the busiest schedule can accommodate the Perfecta® Bravó® treatment time of 30 minutes or Perfecta® REV!® with only 15 minutes of treatment time.

Perfecta® REV!®

14% Hydrogen Peroxide

Hydrogen Peroxide
Take-home Tooth Whitening

Only 15 minutes, once a day to your patients' whitest smiles!

Short Wear Time - 1/4 the time of the leading strips.^{1,2}

Less Sensitivity - 66% less sensitivity than the leading strips ^{1,2}

Superior Whitening - Patented technology with proven results¹

REV! Finishing Rinse - Specially buffered formula immediately accelerates the breakdown of the gel into its active components, providing an added boost to the final whitening effect.

Rinsing effect safely removes residual gel from teeth and any soft tissue.

Contains an advanced form of highly stable Vitamin C which is known to protect soft tissue from any available free radicals and is known to promote healthy teeth and gingiva.

Contains Pentasodium Triphosphate, a sequestering ingredient that helps to "lift and remove" enamel stain.

Alcohol-free, cool mint flavor which leaves the mouth clean and refreshed after every treatment.

Suggest a REV! Refresher Pak for all your previous whitening patients. Get them back to your office and keep them smiling white.

 Clinical Evaluation of Bleaching Efficacy of Two Different Systems, Sheila Koh, D.D.S., F.A.G.D. University of Texas Health Science Center at Houston - Dental Branch, December 2004.

2 Crest® Whitestrips® Supreme

REV! Patient Pak Contains:

2 - 3cc syringes with 4oz (118.3ml) REV! finishing rinse. You add trays and tray cases.

Each Patient Pak is a complete 2-week treatment.

4000142 REV! Patient Pak (6)

Available individually through your dealer.

REV! Refresher Pak Contains:

1 - 3cc syringe with 4oz (118.3ml) REV! finishing rinse. You add trays and tray cases.


Ideal for touch-ups.

4000141 REV! Refresher Pak (6)

Available individually through your dealer.









Perfecta® Promotional Material Available

Contact Customer Service at 888-670-6100 or 610-239-6053





Perfecta® Bravó®

9% Hydrogen Peroxide



Hydrogen Peroxide Take-home Tooth Whitening

30 Minutes, Once A Day!

Bravó® is ideal for your patients.

30 minutes, once a day - improves patient compliance

Short wear time - only 7 hours is required to complete the 2-week take-home treatment.

Strong - stable $9\% H_2O_2$ gel

Avoid sensitivity - patented water-based formulation hydrates teeth as they whiten to avoid sensitivity.

Lasting results - 2 to 3 years, depending on patient habits

Bravó Patient Pak Contains:

4 - 3cc syringes with dispensing tips, patient instructions with shade guide, doctor's instructions, cosmetic patient pouch with mirror, tray material and tray case. Sufficient material for 2-week patient treatment.

4000092 Bravó Patient Pak

Bravó Value Pak Contains:

24 - 3cc syringes - enough whitening material for six patients, dispensing tips, patient instructions with shade guide, doctor's instructions and cosmetic patient pouches. You add trays and tray cases.

4000093 Bravó Value Pak

Bravó 50-Pak Contains:

50 - 3cc syringes with dispensing tips, patient instructions with shade guide and doctor's instructions. You add trays and tray cases.

4000094 Bravó 50-Pak

Bravó Touch-up Pak Contains:

1 - 3cc syringe with dispensing tip and patient instructions with shade guide.

4000091 Bravó Touch-up Pak



Carbamide Peroxide Take-Home Tooth Whitening

Avoid Dehydration

Perfecta® carbamide peroxide gels are formulated with a predominantly glycol base for minimal dehydration and sensitivity. Other gels use glycerin as a base. Glycerin tends to dehydrate the tooth, thought to be a leading cause of thermal hypersensitivity.

Avoid Thermal Hypersensitivity

A small percentage of patients may experience thermal sensitivity associated with whitening. By following a few simple guidelines, you can minimize or eliminate this side effect.

Start Slowly

Patients should wear the tray for only one hour the first day, then one hour twice daily for the first week. Overnight wear is not recommended due to a greater incidence of thermal sensitivity and tissue irritation.

Select The Right Gel Strength

Perfecta 16%: Can remove a full range of stains including those resulting from tobacco, dark beverages and even tetracycline. Wear time can be adjusted for light stains and for those patients with a history of thermal sensitivity.

Perfecta 21%: Tetracycline stains, fluorosis. Dispense only to compliant patients with no history of thermal sensitivity. Give good instructions and frequently monitor treatment progress.





Perfecta® 16%

Carbamide Peroxide



16% Patient Pak Contains:

6 - 3cc syringes, 6 dispensing tips, tray case, tray material, patient instructions with shade guide, doctor's instructions and cosmetic patient pouch with mirror.

4000161 Unflavored 16% Patient Pak 4000162 Mint 16% Patient Pak

16% Value Pak Contains:

36 - 3cc syringes, 36 dispensing tips, patient instructions with shade guide, doctor's instructions and cosmetic patient pouches. You add trays and tray cases.

Enough whitening material to create six 2-week patient treatments or use for touch-up treatments.

 4000163
 Unflavored 16% Value Pak

 4000164
 Mint 16% Value Pak

16% 50-Pak Contains:

50 - 3cc syringes, 50 dispensing tips, patient instructions with shade guide, 25 prescription bags and doctor's instructions.

 4004160
 Unflavored 16% (50-Pak)

 4004165
 Mint 16% (50-Pak)

Perfecta® 21%

Carbamide Peroxide



21% Patient Pak Contains:

6 - 3cc syringes, 6 dispensing tips, tray case, tray material, patient instructions with shade guide, doctor's instructions and cosmetic patient pouch with mirror.

4000211 Unflavored 21% Patient Pak 4000212 Mint 21% Patient Pak

21% Value Pak Contains:

36 - 3cc syringes, 36 dispensing tips, patient instructions with shade guide, doctor's instructions and cosmetic patient pouches. You add trays and tray cases.

Enough whitening material to create six 2-week patient treatments or use for touch-up treatments.

 4000213
 Unflavored 21% Value Pak

 4000214
 Mint 21% Value Pak

21% 50-Pak Contains:

50 - 3cc syringes, 50 dispensing tips, patient instructions with shade guide, 25 prescription bags and doctor's instructions.

4004210 Unflavored 21% (50-Pak) 4004215 Mint 21% (50-Pak)

With Value Paks and 50-Paks, order Tray Material (4008305) and Tray Cases (4008311). See page 138.



Standard Tubes

2-arch Supply Per Tube



A tremendous value with 2 ounces (59cc) of whitening gel in a squeeze tube. Choose 11%, 16%, or 21% carbamide peroxide; fresh mint flavor. Glycol-based for minimal dehydration-induced sensitivity. Ideal package for long-term cases (tetracycline), or for several members in the same family to whiten simultaneously.

Mint only - case of 24

4007115 11% Standard Tube

4007165 16% Standard Tube

4007215 21% Standard Tube

Available individually through your dealer.

Recommend ordering with Tray Material (4008305) and Tray Cases (4008311)

Standard Tube Kits

2-arch Kit



Standard Kits offer the economy of tube delivery with the efficacy of Perfecta® gels. Each kit contains enough gel for a two-arch treatment as well as a touch-up treatment.

Contains:

2 ounces (59cc) of whitening gel in a squeeze tube, tray case, tray material, patient instructions with shade guide and doctor's instructions.

Mint only - case of 12

4006115 11% Standard Kit

4006165 16% Standard Kit

Available individually through your dealer.

In-office Perfection White™

Four Patient - 8-arch Chairside 35% Hydrogen Peroxide Kit



Perfection White™ is a unit-dosed chairside whitening system sufficient to treat 8 arches. This kit should be used only with a rubber dam (not included). A gingival protectant gel is included for use with a rubber dam. Light-activation is not required. Only one hour of chair-time needed per treatment.

Components:

- (8) 4cc ampules of 35% hydrogen peroxide
- (4) 2cc ampules of gingival protectant gel
- (4) activator powder cups (mix with peroxide prior to application to form gel)
- (4) mixing spatulas, four applicator brushes
- · Proprietary powder formula potentiates the chemical whitening action of the hydrogen peroxide

Perfection White - Four Patient 4008104

Accessories



Perfecta® Cosmetic Patient Pouches

Great for enhancing the presentation of any whitening treatment and helps with patient compliance.

4008315 Cosmetic Patient Pouches (6)



Perfecta® Block-Out

Perfecta® Block-Out is a light-cure resin used to block out undercuts on dies or to create whitening-gel reservoirs (optional with Perfecta products). 20-second cure.

4008300 Perfecta Block-Out

2-3cc syringes, 20 tips

Trav Material

Clear tray material for whitening trays. Optimum comfort and adaptability.

4008305 Tray Material (25 sheets)

.040", 5"x 5"

Tray Cases

Perfecta tray cases can store one or two trays. Ultra-trim to slip easily into pocket or purse.

4008311 Tray Cases (24)



Cosmetic/Whitening

Alfa™ Triple Tray®

For Full-arch Closed-mouth Alginate Impressions




- · Reduces gagging
- Cuts chair-time in half simultaneously captures upper and lower impressions
- Saves money requires less alginate and no clean-up
- Increases efficiency closed-mouth technique frees operator from holding tray
- · Ideal for alginate impressions

U.S. Patent #6,749,428

"Pearls For Your Practice"
Dental Economics - June 2003

Box of 24	
1006313	Small (purple)
1006323	Medium (white)
1006333	Large (green)
1006343	Assorted
	Contains : (5) Small,
	(14) Medium, (5) Large





Two Impressions: Mandibular and Maxillary

Tooth-whitening Trays - One Step. One Hour. One Visit.

Combine take-home tooth whitening with your routine hygiene or restorative appointments.



1. Initiate discussion on tooth whitening by recording patient's tooth shade.



5. Using fast-set stone, pour one side, allow to set.



2. Load both sides of Alfa Triple Tray; insert as you would a maxillary single-arch tray.



6. Pour other side and allow to set.



3. Patient closes in normal comfortable position.



7. Fabricate trays with vacuum former & easily trim with Tray Magic™ Trimmer.



 Simultaneously capture maxillary/mandibular impression in as little as four minutes.



8. Patient leaves with whitening trays in hand.



Cosmetic/Whitening

Tray Magic™

Soft Tray Trimmer

Tray Magic™ is designed especially for soft tray material including soft denture relines, athletic mouthguards and snore guards. Improved handle design provides greater blade control. The curved, double-edge, heated tip is perfect for cutting trays and is specially designed for interproximal cuts. Smoothly cuts edges for a comfortable fit with an adjustable heat-setting and durable metal housing. One-year limited warranty.

4008200	Tray Magic Trimmer (110V)
4008207	Tray Magic Trimmer (220V)
4008201	Tray Magic replacement
	cord, wand and blade



Tray Magic technique



Model Should have no palate, flat base, and no over-extended base. Apply lubricant spray to model. Block-out is not required with Perfecta® whitening gels.



Vacuum Forming Cast should rest flat on the vacuum former. Allow tray material to sag only 1/2" to 1". Run vacuum motor for 30-60 seconds. Check for tight adaptation along gingival margin.



Trimming Use heated Tray Magic blade to rough-cut and remove excess material.



Start at most posterior molar and guide the heated tip along the free gingival margin. Apply only gentle pressure. Heat should cut tray material, leaving a smooth edge.



Keep the Tray Magic blade perpendicular to the facial surface. Trim 1/2 to 1mm short of tissue. To finish, turn down the trimmer heat and gently glide warm blade over any rough edges to smooth.

PREMA® Enamel Microabrasion System



PRĒMA® enamel microabrasion compound is a patented acid-abrasive formulation used to correct superficial white, brown or multicolored "dysmineralization" enamel defects. PRĒMA creates a plaque-resistant surface as a result of its enamel glazing ability.³

PREMA compound is applied with a 10:1 gear reduction contra angle, on a low-speed handpiece, to safely remove a controlled, insignificant amount of superficial enamel.

Enamel microabrasion does not recognizably alter enamel surface contours; however, it does permanently eliminate many developmental enamel coloration defects which do not penetrate deeply into the enamel surface. The compound is an ideal blend of specially formulated abrasive powder and a mild concentration of hydrochloric acid.

Enamel microabrasion is different from bleaching techniques and is much more conservative than preparing the enamel and restoring the tooth with composite or veneers.

9008010	PRĒMA Mini Kit
9008011	(5) Mandrels
9008015	(25) Mandrel tips and (2) Mandrels
	Developed by Dr. Ted Croll

0470

Cosmetic/Whitening

PREMA case study

1

17-year-old with chief complaint of "white spots on my front teeth."



2

PREMA® slurry placed upon tooth surface, after rubber dam application.



3

PREMA® slurry compressed with side flutes of polishing cup at very slow speed.



4

Tip of PREMA® polishing cup completes slow speed microabrasion of both tooth surfaces.



5

PREMA® slurry rinsed away with water spray.



0

Enamel Pro® fluoride topical gel used to fluoridate microabrased surface. Enamel Pro® fluoride topical gel saturation of tooth surfaces using cotton swab




1

Immediately after enamel microabrasion using PREMA® microabrasion slurry. Post treatment smile view.




Photography courtesy of Theodore Croll, D.D.S.













CompCore™ AF Dual-Cure

Fluoride-Releasing, Automix Composite Core Build-up Material

Advanced Formula, Dual-Cure CompCore™
Provides exceptional performance and the
enhanced light-curing properties save
considerable time. If light-curing is not an
option – allow CompCore™ AF to
self-cure in 4.5 minutes intraorally!

Improved natural appearance

The greater translucency of the enhanced CompCore AF presents a more natural tooth appearance which minimizes the chance of shadows under ceramic crowns- yet provides good contrast and excellent radiopacity.

Outstanding handling for ease of use CompCore AF Stack will stay put – even in difficult maxillary core build-ups, yet flows under pressure for excellent adaptation. Preps like dentin without ditching.

Superior durability for peace of mind CompCore AF composite core build-up material utilizes exclusive Hyperbranched Technology™ and nano-fillers for superior mechanical properties and long term clinical success - available in automix SyringeMix™ and cartridge dispensers.















"Pearls For Your Practice"
Dental Economics - April 2005

SyringeMix STACK: Dual-Cure/Automix

3001512	Shade White	
3001513	Shade A3	
3001514	Shade Blue	
Contains two	Og automiy syringos	15 miy tin

Contains two 9g automix syringes, 15 mix-tips, 15 standard & 10 extra-fine intraoral tips

Cartridge STACK: Dual-Cure/Automix

3001430	Intro Kit / Shade White
3001431	Intro Kit / Shade A3
Contains o	ne 50g cartridge, type-25 gun,
10 ACCOR®	matrices, 4ml each IntegraBond
& Auto-cur	e Activator, 25 mix-tips,
25 intraora	l tips

Cartridge Refills

3001432	Shade White
3001433	Shade A3
Contains or	ne 50g cartridge, 10 ACCOR°
matrices, 25	mix-tips, 25 intraoral tips

CE 0470

Twist STACK: Auto-Cure/Syringe Hand-Mix

3001490	Shade White	
3001491	Shade A3	
	one 6g catalyst & base,	
mixing pa	ds, 20 spatulas	

Accessories

3001418	SyringeMix Mix-Tips & Intraoral Tips (35)
3001468	SyringeMix Extra-Fine
	Intraoral Tips (50)
3001416	Applicator Gun, Type-25
3001417	Cartridge Mix-Tips &
	Intraoral Tips (35)
9048081	ACCOR® Matrices, 50 multi-sized
	core matrix forms



See pages 151-153 for full listings of matrix bands and wedges, including ACCOR core matrix forms.



CompCore™ AF

Advanced Formula, Dual-Cure Composite Core Build-up Material









Bonding Agent

(use only bonding agents compatible with dual-cure composites, such as IntegraBond $^{\mathtt{m}}$)

- 1 Clean and dry prepared tooth surface.
- 2 Apply 37% phosphoric acid to the cavo surface for 10 to 15 seconds.
- 3 Rinse and dry cavo surface, do not desiccate the tooth.
- 4 Apply first coat of IntegraBond™ on prepared cavo surface. Allow to sit for 15 seconds, and then air dry. Apply 2 to 3 more coats; without waiting, dry again after last application.
- 5 Light-cure for 20 seconds









Core Build-Up

- 1 Bleed syringe. Assemble intraoral tip to mix-tip, attach to syringe and bleed tip.
- 2 Place intraoral tip into the base of the prepared tooth and dispense. Let material rise around the tip while withdrawing to eliminate air bubbles.
- 3 Light-cure for 30 seconds. If the core is larger than 8mm, light-cure for at least an additional 10 seconds.*

Check for complete hardness. The preparation can be finished immediately. Alternately, the core can be chemical-cured, final set-time is 4.5 minutes.

Discript Duranting		
Physical Properties		
Compressive Strength:	250 MPa	

Light-cure Setting Times


Working Time: 2.5 minutes at 72° (22°C)

* Final Set Time	Shades		
Depth of cure	А3	White	Blue
6mm	30 sec	30 sec	30 sec
8mm	30 sec	40 sec	60 sec

Chemical-cure Setting Times

 Working Time:
 2.5 minutes at 72° (22°C)

 Final set:
 4.5 minutes at 98° (37°C)



ACCOR® is a registered trademark of ACCOR Inc. / 1 Evaluation of Hyperbranched Polymers Effect on Dental Composite Properties, Lizenboim,K., Dodiuk-Kenig, H., Eppelbaum, I, Shenkar College for Engineering and Design, Ramat-Gan, Israel. / 2 Independent laboratory study from Creighton University, data on file. Dent. Research (1153,2002)

CompCore™ AF (Advanced Formula) is CE marked. CompCore Tubs are not CE marked.

Restorative

CompCore[™]

Smooth Flowing, Syringeable, Composite Core Build-up Material with Fluoride - Hand Mix



Ideal Handling - Saves Time
Premier® CompCore™, available in white and natural shades, is easily syringed, permitting precise placement without slumping and

self-cures in only 3 minutes to save chair-time.

CompCore is radiopaque and contains fluoride. Its unique composite formulation provides outstanding compressive and tensile strength, yet easily trims like dentin. Create a strong, unitary core, bonding CompCore to tooth structure with Bond Boost SE™ and IntegraBond™.

Outstanding Value - Saves Money Satisfaction Guaranteed!

3001405 Shade Natural (28g) 3001401 Shade White (50g)	3001402	Shade White (28g)	
	3001405	Shade Natural (28g)	
3001404 Shade Natural (50g)	3001401	Shade White (50g)	
300 1404 Shade Natural (309)	3001404	Shade Natural (50g)	



IntegraCem™

Dual-cure Resin Cement



Easy Technique

IntegraCem's™ low viscosity and thin film thickness allows excess cement to vent evenly for precise marginal fit. After a 10-sec. light-cure (halogen), you easily will remove cement flashing from seated crowns, bridges, inlays and onlays.

Versatile

Illustrated cards guide you through all step-by-step techniques, including PFM crown repair, cementation of Maryland bridges, posts, ortho brackets and bands.

Ideal for "Bonded" Amalgam IntegraCem has superior bond-strength to all surfaces, especially metal. Formulated to chemically bond to fresh amalgam, IntegraCem will eliminate sensitivity^{3,4} and permanently seal all restorative margins.

Dual-Cure with Anaerobic Accelerator⁵ Posts can be stabilized immediately by light-curing; IntegraCem's anaerobic accelerator guarantees complete curing even at the bottom of the post preparation.

3001410 Base & Catalyst- 4g each



IntegraBond™

Universal Adhesive System





Over 10 years Proven Performance & Reliability⁶⁻¹⁷

Premier® IntegraBond™ is universal. Unlike other brands, IntegraBond - fifth generation bonding agent - works with all composite materials and resin cements: light-cured, dual-cured and chemical-cured.

Use IntegraBond after acid-etching.

IntegraBond is easy to use. By priming and bonding in a single application, you save time while sealing dentin tubules to eliminate sensitivity and protect vital tooth surface.

IntegraBond is economical. Up to 50% less cost than other brands.

IntegraBond is strong. In every application, IntegraBond offers proven, reliable bond strength due to its nano-filler technology and unique formulation.

3001421 **IntegraBond Adhesive Kit**Contains: 4ml IntegraBond, 4ml Auto-Cure
Activator, (2) 1.2ml Premier Etch Syringes,
50 Microbrushes, 20 Disposable Etch Tips,
Mix Wells

3001411	IntegraBond - 7ml	
3001412	Auto-Cure Activator - 7ml	
for use wh	en light-curing is not possible	
3001413	IntegraBond - 4ml	









Bonding Agent

- 1 Clean and dry prepared tooth surface.
- 2 Apply 37% phosphoric acid to the cavo surface for 10 to 15 seconds.
- 3 Rinse and dry cavo surface, do not desiccate the tooth.
- 4 Apply first coat of IntegraBond™ on prepared cavo surface. Allow to sit for 15 seconds, and then air dry. Apply 2 to 3 more coats; without waiting, dry again after last application.
- 5 Light-cure for 20 seconds

When Light-Curing is not an option: Endodontic Post Cementation

In step 5, add equal drops of Auto-cure Activator with IntegraBond in a clean mixing well/dappan dish and combine for 2 seconds. Continue with step 5 using the IntegraBond/Activator mixture. Use a paper point to wet the deepest portions of the canal. Avoid pooling of the adhesive. Proceed to post cementation.

Note: The use of IntegraBond with Auto-cure Activator can accelerate the working time of the self- and dual-cure cements.

3. Hadi Z, Rosenstiel SF, Rashid RG. A comparison of post-treatment sensitivity using adhesive cavity liner vs conventional cavity varnish (abstract 1335). J Dent Res 1998;77(special issue A):272 / 4. Davis RD, Overton JD. Efficacy of bonded and nonbonded amalgam in the treatment of teeth with incomplete fractures. JADA 2000;131:469–78. / 5. I. Eppelbaum, H Dodik, S Kenig, B. Zalsman, A. Valdman, and R. Pilo. 'The Role of an Anaerobic Accelerator in Dental Adhesives' J. Adhesion Sci., Technl., Vol 10, No 10, pp 10075-1087 (1996) / 6-17. Full bibliography available.



NexTemp™

The Next Generation in Temporary Cement

"Pearls For Your Practice" Dental Economics - May 2006

"Best Products of the Year"*



Greater Retention, **Less Cement!** Just dispense a small bead of cement around the inside margin.

Firm Retention/Easy Removal 4X greater retention than most zinc-oxide cements

Eugenol-Free - Non-irritating to soft tissue. Will not inhibit permanent cementation Fluoride Release

Potassium Nitrate - Known to reduce sensitivity and act as a sedative to the pulp Chlorhexidine - A proven anti-bacterial agent **Resin-Matrix Formulation** Good marginal seal with no wash-out

Two-Stage Curing - In approximately 2 minutes, initial gel-set stabilizes temporary, allowing easy removal of excess. Final set occurs one minute later.

Automix Syringe - Saves application time; guarantees consistent mix

Neutral Shade - Blends esthetically with temporary restoration / no show-through

"Outstanding Products"**

3001470	Bulk - 4 Automix Syringes (5ml each)
3001471	Standard - 1 Automix Syringe (5ml) & 10 Mix Tips
3001451	Automix Dual-barrel Syringe Tips (25)





Premier® Implant Cement™

Secure Retention & Retrievability for Implant Retained Crowns



NexTemp Temporary Cement-5ml

Non-Eugenol Temporary Resin Cement with Fluoride, m Nitrate and Chlorhexidine

> **Secure retention** and retrievability in one unique cement!

High retention- Long-term cementation Retrievability - Off when you need it Radiopaque - Improves detection of excess cement

Low solubility - Good marginal seal/ no washout

Virtually no taste or odor - Increased patient comfort

Easy handling - Two-stage cure features initial gel-phase in 2-2.5 minutes for easy removal of excess cement. Rigid final set ensures marginal seal and firm retention.

Saves time - Automatic, consistent mix and controlled dispensing

Easy clean-up

3001450 **Implant Cement Standard Pack** 1-5ml Automix Syringe, 10 Mixing Tips and

1 Mixing Pad

3001455 **Implant Cement Value Pack** 3-5ml Automix Syringes, 25 Mixing Tips

Automix Syringe Tips (25) 3001451

Properties

Working Time 1.5 minutes at room temperature Gel-phase 2-2.5 minutes at 98°F (37°C) Set Time 4-5 minutes at 98°F (37°C)

Film Thickness 15 µm







IntegraPost™ is easy to carry and place using the IntegraPost carrier.



CompCore™ is syringed through the patented **Flow-Thru-Head®** providing a void-less stress-free mechanical interlocking of post and core.



IntegraPost's rounded head eliminates stress points and the Flow-Thru-Head eliminates voids, ensuring complete post and core integration.

IntegraPost®

Ultra-passive Integrating Core Retainer



CE

Safe and Secure

The IntegraPost features an innovative round Flow-Thru-Head® and a parallel passive design to optimize the adhesive buffer zone, conserve tooth structure and minimize stress. 18,19,20

Balances Strength with Elasticity
Made from bio-compatible titanium alloy,
the IntegraPost is stronger than resin fiber
yet more elastic than stainless steel.²¹

Easy to Use

Six post sizes and corresponding drills make post selection a simple process. IntegraPost heads are optimally sized in relation to shank length. The unique post-carrier instrument facilitates handling. Passive design minimizes stress points on root, reducing potential for root fracture. IntegraDrill is ideally sized in relation to post diameter, permitting 2.5X more cement into post space. This creates a "buffer zone" to absorb potentially stressful forces and the greater cement volume increases resistance to vertical dislodgement.²²

Innovative Design

Midflange and shoulder are designed to dissipate stress away from canal. The safe, non-threaded parallel shank minimizes potential root fracture. Anti-rotational knurling locks into cement.

In a 2006 ADA Professional
Product Review on Endodontic
Posts, which included metal, fiber and a
ceramic post, IntegraPost received
the highest laboratory
rating of 92% with a
buyer's satisfaction of excellent!²³





18-22. Full bibliography available.

23. American Dental Association Professional Product Review™
Publication by the ADA Council on Scientific Affairs, Vol. 1, Spring
2006. / IntegraPost designed by Endodontist, Dr. Phil Klein.



IntegraPost®

Mini Intro Kit



3001001 Mini Intro Kit

Contains: 6 posts

(2) each #3, #4, #4.5, (3) drills - (1) each #3, #4, #4.5,

IP Carrier



IntegraPost®

Refills



Refill	s - (10	pack)

3001105	Size 3	
3001110	Size 4	
3001115	Size 4.5	
3001125	Size 5	
3001130	Size 6	
IntegraDrill	•	
3001205	Size 3 (vellow)	

3001205	Size 3 (yellow)	
3001210	Size 4 (red)	
3001215	Size 4.5 (purple)	
3001225	Size 5 (blue)	
3001230	Size 6 (green)	

IntegraDrills (5 pack)

3001235 Assorted

Integra Instruments

3001301	IP Carrier
3001305	IP Probe

Cure-Thru® IntegraPost®

Light-Transmitting Zirconia Glass Fiber Post

A simple technique for post and core utilizing Zirconia glass fibers for the right balance of strength and elasticity



3001600 Intro Kit

Contains: 20 posts, (5) each #4, #4.5,#5, #6 and (4) corresponding drills

3001601 Mini Intro Kit

Contains: 10 posts, (5) each #4, #5 and (2) corresponding drills

Refills (10 pack)

3001610 Size 4 (red)

3001611 Size 4.5 (purple)

3001612 Size 5 (blue)

3001613 Size 6 (green)







Premier® Etch

37% Phosphoric Acid Dentin/Enamel Etch Gel



Premier® Etch contains phosphoric acid for etching tooth enamel and dentin when preparing surfaces for the application of composites.

3001422

37% Dentin/Enamel Etch Kit Contains: one 12g jumbo syringe, with luer lock adapter, ten empty 1.2ml syringes and (50) disposable applicator tips for the ten 1.2ml syringes



Porcelain Etchant

9.6% Hydrofluoric Acid Gel



Yields superior retention. Each 0.2cc finger applicator unit is sufficient for etching two to four porcelain facings.

4009001 Porcelain Etchant - (22 units)

MicroDose®
"Outstanding Products"*

38% Etchant - MicroDose®

Fast Acting, Silica-Free Dentin/Enamel Etchant



15-20 second etch time. Each 0.2cc unit is sufficient for etching two to four teeth. Bright blue in color.

4009051 Refill (50) 0.2cc
MicroDose pipettes
4009099 MicroDose Syringe Pack
(3) autoclavable syringes

Porcelain Silane

Advanced Silane Chemistry



Each 0.2cc finger applicator unit is sufficient for etching two to four facings, crowns or inlays. A clear liquid.

4009020 Porcelain Silane - (18 units)

Porcelain Bonding Kit



4009025 Porcelain Bonding Kit - 18 packs Contains: one unit each of

Porcelain Etchant and Porcelain Silane

^{*}An independent non-profit dental education and product testing foundation, December 2002

^{**}An independent non-profit dental education and product testing foundation, January 2001



Sycamore Wood Wedges

The Original Anatomically Carved

- Fits interproximally to support restoration and provide proper contact
- Slightly compressible, yet will not splinter, even in thin sections.
- Teeth can be separated easily to compensate for the thickness of the matrix band.
- After wedge is removed, teeth return to their original position. A perfect contact point is restored.

Eight color-coded sizes (400 count)

13T (Green)

13 (Yellow)

15 (Blue)

16 (Red)

17 (Neutral)

Mini (Light Blue)

UT - Ultra Thin (Orange)

13XT - Extra Thin (Neutral)

9061109

9061106

9061100

9061101

9061102

9061103 9061104


9061105

9061108

R	1	•
\		

Ultra-thin is half the width of regular wedges, yet high through the embrasure.

"Most Highly Rated"**



Number designates length of wedge in mm. (13=13mm)



Premier's anatomically shaped wedges fit interproximally to support your restoration and provide proper contour.



Other wedges fail to match the natural anatomy between teeth.

technique tip

Presoaking Wedges in water or Hemodent® just prior to use causes them to swell while in place, creating greater wedge efficiency.

Assorted (UT, 13, 13T & 15) Presoaking W Hemodent® jus them to swell w greater we

Sycamore Wood Wedges Have 1st Separate Wood Wedges Have 1st Separate Productive Have 1st Separate Productive Have 1st Separate Productive Have 1st Separate Productive Have 1st Separate Research Have 1st Separat

Wood Wedge Organizer

9061115 Six-compartment Wedge Organizer


Multi-purpose Wedge Positioner

PREMIER CO

9065146 Multi-purpose Wedge Positioner



Cure-Thru® Reflective Curing Wedges



- Clear, anatomically-shaped wedge with a unique, soft underside provides superior adaptability to tooth structure
- The clear, acrylic body and the white core reflect up to 90% of the curing light towards the restoration, resulting in a better contact
- Ensures superior adaptation on proximal root surfaces. Helps to prevent crevicular fluid from contaminating the restoration.

9061216 Intro. System
Contains: (40) Ultra-Thin, (30)
Small and (30) Medium
Wedges, (30) each Blue
Cure-Thru Premolar and
Molar Matrix Bands

9061221 100 Small
9061223 100 Medium
9061225 100 Ultra-Thin

Cure-Thru® Clear Cervical Matrices

Compress restorative material as you lightcure. Flexible to the exact contour of the tooth for dense restorations with excellent marginal integrity.

9061011 Assorted box of 275

Contains: (50) Anterior/Premolar and Molar/Universal, (35) all other shapes, positioning instrument, gingival retraction instrument

Available in boxes of 150

9061012	719CT
9061013	720CT
9061014	721CT
9061015	722CT
9061016	723CT
9061017	Anterior/Premolar
9061018	Molar/Universal
9061340	Positioning Instrument

"Pearls For Your Practice " Dental Economics - March 2000







Actual	Size	_		_		
719	720	721	722	723	Anterior/ Premolar	Molar/ Universal

Cure-Shield

Hand-held Visible Light Shield

Premier Cure-Shield effectively blocks out potentially harmful ultra-violet and blue light at wave-lengths shorter than 525 nanometers which, after prolonged exposure, may cause photochemical burn of the retina.

- Protects the entire staff
- Can be held by the assistant, the patient or the free hand of the operator
- Is not cumbersome or expensive like eye glasses
- · Does not restrict your visibility



9006166 Pack of three

Pat.Nos.Des.294,639; 4,640,685; 4,737,104






ACCOR®

Multi-sized Matrices for Core Build-up





Easy to use

One cut across the Premier ACCOR® Matrix produces basic gingival contour. No lubrication required, nearly impossible to overseat.

Versatile

Can be used on anterior and posterior teeth, with glass ionomers, composite, pattern acrylic or softened wax.

Disposable

Eliminates cross-contamination possible with trial-and-error fitting of individually sized forms.

9048081 50 Matrices

ACCOR® is a registered trademark of ACCOR Inc.



Blue Cure-Thru® Matrix Bands

Transparent Contoured Matrices 0.075mm



- Blue color increases the contrast between the matrix and tooth without compromising polymerization of the composite
- Thin, yet strong enough to pack against
- Simplifies posterior composite placement

9061231 Premolar Bands (100) 9061233 Molar Bands (100)

Pre-contoured Metal Bands



- Anatomical shape achieves better contact points
- Ideal for large Class II restorations with incremental curing technique
- · Finished edges ensure patient comfort
- 0.035mm thickness

9061070 30 Pre-contoured Metal Bands

Restorative

Stop Strip

Mylar Anterior Bands with Integrated Stopper

Turns a three-handed procedure into a two-handed one. Light-cure Class III and IV restorations with ease.



Unique integrated stopper secures the matrix band between adjacent teeth.



Pull the band taut; you have improved access for light-curing.

9061370	Assorted, 8mm & 10mm
	(75 each)
9061375	Straight, 10mm (100) (blue)
9061380	Straight, 8mm (100) (red)

"Pearls For Your Practice"

Dental Economics, April 2001

RC-Prep





RC-Prep® For Chemo-mechanical Preparation of Root Canals	Page 156	
RC-Prep® MicroDose® Cleans & Lubricates for Chemo-mechanical Preparation of Root Canals	157	
K-Files, Hedstrom Files, Reamers Stainless Steel Endodontic Instruments	158-159	
Root Fillers Carrier for Root-canal Sealer	160	
Gates-Glidden Drills Shape & Taper the Canal	160	
Barbed Broaches Flexible Stainless Steel Avoids Breakage	161	
Finger Spreaders For Lateral Condensation	161	
Finger Pluggers For Vertical Condensation	161	
Gutta Percha Points Machine Rolled	162	
Two Striper® Endo/Access Diamonds	162	
Surgical-length Carbide Burs	163	
Absorbent Points Sterile	163	
Peeso Drills Prepare Post Space in Endodontically Treated Teeth	163	
, , , , , , , , , , , , , , , , , , , ,		



RC-Prep® Technique



RC-Prep is easily introduced into the canal directly via the MicroDose Pipette or a reamer, file or root-canal filler. As RC-Prep cleans and lubricates, instruments cut freely without binding.



When RC-Prep reacts with sodium hypochlorite solution, pulp debris is actually lifted out of the canal. The pulp chamber is brightened, helping to locate additional canals.



With RC-Prep and proper instrumentation, the canal is thoroughly debrided. Dentin tubules and lateral canals are opened and cleaned, essential for successful treatment.

RC-Prep®

For Chemo-mechanical Preparation of Root Canals

Over 40 years of clinical use and numerous publications demonstrate the clinical efficacy of RC-Prep for chemo-mechanical preparation of root canals.

RC-Prep's unique formulation of glycol, urea peroxide and EDTA in a special water-soluble base helps remove calcifications and lubricates the canal to permit more efficient instrumentation. Whether you prefer finger or engine-driven instruments, RC-Prep allows reamers and files to move easily without binding. RC-Prep is also excellent for use with apex locators as it permits consistently reliable readings.

9007131 18g Jar Dealer: Available to you in pack of ten.	9007129 RC-Prep Syringe Kit (5) 3cc syringes and (50) tips
9007132 18g Pump	9007139 Metal Syringe Tip Refill (50)
Dealer: Available to you in pack of ten.	9007135 9g Syringe with 1 Plastic Tip (package of 2)
9007133 227g Jar	Dealer: Available to you in pack of ten.
	9007127 Plastic Syringe Tip Refill (50)







RC-Prep® MicroDose®

Cleans and Lubricates for Chemo-mechanical Preparation of Root Canals

Direct Application Into the Canal

RC-Prep® MicroDose® is convenient, easy to dispense and eliminates the risk of cross-contamination. Each MicroDose pipette contains the ideal amount of RC-Prep (0.2cc) for more efficient instrumentation during root-canal procedures. You can apply RC-Prep directly into the pulp chamber or into a dappen dish. After each use, the empty pipette is discarded and the MicroDose syringe can be autoclaved.


9007138	RC-Prep MicroDose Intro Kit (100) Pipettes and (2) Autoclavable Syringes
9007137	RC-Prep MicroDose Refill (100) Pipettes
4009099	MicroDose Syringe Pack (3) Autoclavable Syringes











Standard feature on all K-Files and Hedstrom Files, sizes 06-40

K-Files ■

Premier K-Files are used to enlarge canal after debridement with a reamer. After reamer, use a file one size smaller than the reamer so that the fragments of dentinal shavings can be carried away from the root-canal walls. Use in clockwise direction (1/4 turn only), then withdraw.

Safe-End Tip on sizes 06 - 40

Stainless Steel Endodontic Instruments

With Tactile Plastic Handles

ISO-Coded

Premier's tactile plastic handles, on all Reamers, K-Files and Hedstrom Files, have the ISO symbols and number of the instrument clearly and indelibly marked on the color-coded handle. The dentist or the assistant can easily locate the proper instrument in any set-up.

Firm Grasp

These peanut-shaped handles provide the utmost in tactile feel with a positive grip. They enable the dentist to negotiate torturous and curved canals easily and effectively. The safety eye can be easily threaded for additional protection.

Virtually Break-proof

These highly flexible instruments resist fracture. They can be adapted to curved root canals, while practically eliminating breakage from binding.

Extremely Sharp

Premier Endodontic Instruments stay sharp and cut efficiently time after time because they are made of the finest stainless steel in the world.

Precise to AAE and ISO Standards
Exact tolerances and uniformity are built
in by craftsmen who specialize in making
endodontic instruments exclusively, through
precision quality controlled machinery.

Sterilizable

These heat-proof plastic handles can be sterilized in an autoclave, dry heat, or chemicals without fear of deformation.

Premounted Silicone Stop Saves TimeAll Premier K-Files, Hedstrom and Reamers come with a premounted silicone stop at no additional charge.

Rubber Stop Color Coding for Files & Reamers Red/21mm. blue/25mm and black/31mm.





Hedstrom Files ●

Premier Hedstrom Files are used for more rapid enlargement of root canals. Note that the cutting edges are inclined backwards. The insertion is made easily since there is no cutting taking place on the insertion stroke. When the instrument is withdrawn, cutting takes place and debris is carried away from the apex. Use with an up-and-down motion.

Safe-End Tip on sizes 06 - 40



Premier Reamers are used for initial debridement of root canals. Start with smallest numbers and use each next larger size in turn to prevent jamming and straying from course of root canal. Spiral cutting edges carry dentinal tissue away from apex. Use in clockwise direction (1/4 turn only).

Safe-End Tip on sizes 06 - 40



21mm (6)	per box						
9053399	08	9053404	30	9053409	55	9053416	15-40
9053400	10	9053405	35	9053410	60	9053417	45-80
9053401	15	9053406	40	9053411	70		
9053402	20	9053407	45	9053412	80		
9053403	25	9053408	50	9053413	90		
25mm (6) per box							
9053424	08	9053429	30	9053434	55	9053441	15-40
9053425	10	9053430	35	9053435	60	9053442	45-80
9053426	15	9053431	40	9053436	70		
9053427	20	9053432	45	9053437	80		
9053428	25	9053433	50	9053438	90		
31mm (6) per box							
9053876	08	9053881	30	9053886	55	9053894	45-80
9053877	10	9053882	35	9053887	60		
9053878	15	9053883	40	9053889	80		
9053879	20	9053884	45	9053890	90		
9053880	25	9053885	50	9053893	15-40		



21mm (6) per box						
06	9053005	25	9053010	50	9053015	90
08	9053006	30	9053011	55	9053018	15-40
10	9053007	35	9053012	60		
15	9053008	40	9053013	70		
20	9053009	45	9053014	80		
per box						
06	9053035	25	9053040	50	9053045	90
08	9053036	30	9053041	55	9053046	100
10	9053037	35	9053042	60	9053048	15-40
15	9053038	40	9053043	70	9053049	45-80
20	9053039	45	9053044	80		
31mm (6) per box						
06	9053095	25	9053100	50	9053107	110
08	9053096	30	9053101	55	9053111	15-40
10	9053097	35	9053102	60	9053112	45-80
15	9053098	40	9053103	70		
20	9053099	45	9053106	100		
	06 08 10 15 20 per box 06 08 10 15 20 per box 06 08 10 15 20 per box	06 9053005 08 9053006 10 9053007 15 9053008 20 9053009 per box 06 9053035 08 9053036 10 9053037 15 9053038 20 9053039 per box 06 9053095 08 9053096 10 9053097 15 9053098	06 9053005 25 08 9053006 30 10 9053007 35 15 9053008 40 20 9053009 45 per box 06 9053035 25 08 9053036 30 10 9053037 35 15 9053038 40 20 9053039 45 per box 06 9053095 25 08 9053096 30 10 9053097 35 15 9053098 40	06 9053005 25 9053010 08 9053006 30 9053011 10 9053007 35 9053012 15 9053008 40 9053013 20 9053009 45 9053014 per box 9053035 25 9053040 08 9053036 30 9053041 10 9053037 35 9053042 15 9053038 40 9053043 20 9053039 45 9053044 per box 06 9053095 25 9053100 08 9053096 30 9053101 10 9053097 35 9053102 15 9053098 40 9053103	06 9053005 25 9053010 50 08 9053006 30 9053011 55 10 9053007 35 9053012 60 15 9053008 40 9053013 70 20 9053009 45 9053014 80 per box 06 9053035 25 9053040 50 08 9053036 30 9053041 55 10 9053037 35 9053042 60 15 9053038 40 9053043 70 20 9053039 45 9053044 80 per box 06 9053095 25 9053100 50 08 9053096 30 9053101 55 10 9053097 35 9053102 60 15 9053098 40 9053103 70	06 9053005 25 9053010 50 9053015 08 9053006 30 9053011 55 9053018 10 9053007 35 9053012 60 15 9053008 40 9053013 70 20 9053009 45 9053014 80 oper box 06 9053035 25 9053040 50 9053045 08 9053036 30 9053041 55 9053046 10 9053037 35 9053042 60 9053048 15 9053038 40 9053043 70 9053049 oper box 9053044 80 80 9053049 9053044 80 oper box 9053095 25 9053100 50 9053107 08 9053095 25 9053101 55 9053111 10 9053096 30 9053101 55 9053111 10 9053097







Premier SensiPast root fillers have a miniature clutch in the mandrel that allows them to be used on any kind of right-angle hand-piece without breaking. *Guaranteed!*

Package of 4				
21mm				
9053571	25			
9053572	30			
9053575	Assorted (25-40)			
25				
25mm				
9053576	25			

Root Fillers



Premier Root Fillers are used with a contraangle at low speed to carry root-canal sealer to all recesses and walls of root canal. Numbers of Root Fillers correspond with numbers of root-canal instruments. Always use a Root-canal Filler of at least one size smaller than last instrument used for debridement of root canal.

Package of 2	1
9053556	25
9053557	30
9053558	35
9053559	40
9053560	Assorted (25 - 40)
25mm	
9053561	25
9053562	30
9053563	35
9053564	40
9053565	Assorted (25 - 40)

Gates-Glidden Drills



Premier Gates-Glidden Drills are used to open the root canal and to further shape and taper the canal. They are designed so that, should they bind while drilling in the root canal, they will break at the top of the shaft so they can be retrieved easily. Shafts are indented for easy size identification.

Package of 6				
Package of	0			
9055001	No. 1			
9055002	No. 2			
9055003	No. 3			
9055004	No. 4			
9055005	No. 5			
9055006	No. 6			
9055009	Assorted (1-6)			





Barbed Broaches



Premier plastic-handle Barbed Broaches are made of extremely flexible stainless steel to avoid breakage. Each instrument has more than 40 barbs, spiraling around the core, which change their angulation when introduced into narrow canals - to avoid breaking.

The color-coded plastic handle easily identifies size. It can be autoclaved without affecting color.

Package of 144				
9053631	25			
9053632	30			
9053633	35			
9053634	40			
9053635	50			
9053636	60			
9053637	Assorted (25-60)			

Individual 12 packs available through your dealer.

Finger Spreaders



Premier Finger Spreaders are precisely sized to provide excellent lateral condensation of gutta percha toward canal walls. Premier Finger Spreaders are available in 25mm lengths and ISO sizes 15-40. Each finger spreader is made of surgical-grade flexible stainless steel with a highly polished working surface. The handles provide optimum tactility and are sterilizable up to 200°C/390°F.

Package of 6				
21mm				
9053814	15			
9053816	20			
9053818	25			
9053820	30			
9053822	35			
9053826	Assorted (15-40)			
25mm				
9053830	15			
9053832	20			
9053834	25			
9053836	30			
9053838	35			
9053840	40			
9053842	Assorted (15-40)			

Finger Pluggers




Premier plastic-handle Finger Pluggers are used to condense gutta percha for the final sealing of the canal and dentinal tubules. They are flexible and precision made to match ISO-sized reamers and files. Sterilizable up to 200°C/390°F.

Package of 6	5
21mm	
9053771	15
9053772	20
9053773	25
9053774	30
9053788	Assorted (15-40)
25mm	
9053791	15
9053792	20
9053793	25
9053794	30
9053795	35
9053796	40
9053808	Assorted (15-40)





Gutta Percha Points



Consistency and "Feel"

Specially formulated Premier Gutta Percha is neither too soft nor too brittle. These points offer the ideal consistency for plugging or for lateral condensation. Gives excellent radiopacity. Premier Gutta Percha is softened easily and dissolves readily in chloroform.

A.A.E. Sizes - Color-Coded

All A.A.E. Gutta Percha are machine-made (not hand-made) to ensure an identical shape and size to the corresponding Reamer and/or File. A.A.E. sizes of Gutta Percha, in boxes of 6 vials.

Box of 60	
9056013	90
9056014	100
9056017	130
9056021	Assorted 90-140
Box of 120 9056001	15
9056002	20
9056003	25
9056004	30
9056005	35
9056006	40
9056007	45
9056008	50
9056009	55
9056010	60
9056011	70
9056012	80
9056019	Assorted 25-40
9056020	Assorted 45-80

Gutta Percha Accessory Points




Regular Taper

These Gutta Percha Points are firm enough for secure placement in the canal and can be compacted easily to fill all tubules under the pressure of lateral or vertical condensation.

Box of 60 9056107	Medium/Large
Box of 120	
9056101	Extra Fine
9056102	Fine/Fine
9056103	Medium/Fine
9056104	Fine
9056105	Fine/Medium
9056106	Medium
9056116	Assorted (XF, FF, MF, F, FM, M)

Two Striper®

Endo/Access Diamonds



Gaining access through porcelain, metal, enamel and dentin is accomplished easily using a coarse, round Two Striper® Diamond.

S -79 has a safe-end for enlarging the entrance of the pulpal chamber without disrupting the chamber floor.

The small, round surgical-length diamonds help in locating calcified, hard-to-find canal orifices. The flame-shaped diamond then is used to flare out the opening, making easy access for further instrumentation.

Endo/Access Diamonds				
2013194	1C (each)			
2013197	291.4C (each)			
Popular FG	Access Diamonds			
2015003	115C (5-pack)			
2015005	120C (5-pack)			
2015006	125C (5-pack)			
2015008	130C (5-pack)			
2015835	S-79F (5-pack)			



Surgical-length Carbide Burs

Long Shank for Deep Pulp Chamber



These long-shank burs enable penetration deep into the pulp chamber without interference from cusps or crowns. Made of the finest tungsten carbide in friction-grip stainless-steel surgical shanks. 1 per package.

9050601	1	
9050602	2	
9050603	4	
9050604	6	
9050605	8	
9050607	34	
9050608	35	
9050612	402	
9050615	1558	
9050618	556	
9050619	557	
9050620	558	
9050621	559	
9050622	560	
9050625	700	
9050626	701	
9050627	702	
9050628	703	

See page 140 for IntegraPost™

Absorbent Points

Sterile



Five quality Absorbent Points in each easy-opening cell. Maintain sterile condition until each cell is open. Consistent in size; consistent in texture.

Box of 200	- 40 cells of 5
9055101	XX-Fine
9055102	X-Fine
9055103	Fine
9055104	Medium
9055105	Coarse
9055106	Extra Coarse
9055109	Assorted (XXF, XF, F, M, C)
AAE Sizes	- Box of 200
9055123	15
9055124	20
9055125	25
9055127	35



Peeso Drills

These engine-driven instruments are designed to prepare post space in endodontically treated teeth. To protect the integrity of the tooth and root-canal filling, gutta percha should be removed with hot pluggers to the desired depth before use of Peeso drills.





Package of	6	
9055021	No. 1	
9055022	No. 2	
9055023	No. 3	
9055024	No. 4	
9055025	No. 5	
9055026	No. 6	
9055029	Assorted (1-6)	



Diagnostic

Hygiene

Operative

Endodontic

Premier Instruments

2-9

10-21

22-47

48-53

Perio/Oral Surgery	54-84
Cassettes and	85-86
	03-00
Instrument Care	
204 Series Scalers	13
A	
-	20
Amalgam Carriers/Guns	
Angle Former/Cavity Prep	
Apical Picks	
Applicators	
Articulating Paper Forcep	
Aspirators	
Aspirators	
В	
Barnhart Curettes Big Easy®	6, 20-21, 66 1-16, 18-19 69 31
С	
Carriers - Amalgam	29
Carvers	
Cassettes	85
Cattoni Scalers	13
Cavity Finder™ Explorers	3-4
Cavity Preparation Instruments	23-24
Cement Applicator	30
Child-Size Forceps	76-77
Chisels - Operative	23
Chisels - Bone & Periodontal	
Clear-View™ Periodontal Probes	7-8
Columbia Curettes	14
ComfortView®	65, 93
Composite Instruments	43-44
Composite Placement	34, 38-44
Composite/Plastic Filling Instruments	38-42
Condensers/Pluggers	
Cord Packers	
Coupland Elevators	73
Crown & Collar Scissors	
Crown Removers	44

Curettes Periodontal. 14-2 Surgical 67-6	
D	
Dressing Pliers	
E	
Ejectors - West Root Tip	75 19
Files	
Bone	58 14 77 33 31
Frazier suction rubes	' '
G	
Gingival Cord Packers. 3 Glicks 5 Gracey Curettes 16-2	50
Н	
Handles - Mirror	24
Implant Scalers	
I	
Jacquettes Scalers	3
K	
Knives - Amalgam	
L	
Laboratory Carvers	17 15 20 93 58
Mallets	59
Margin Trimmers/Cavity Prep 2 McCall Curettes 1 Mini-Pocket Scalers 19-2	14

N
Needle Holders 61-63 Ni-Ti Spreaders 51
P
Packers - Cord
Periodontal Chisels 69
Periodontal Elevators
Periodontal Files (Hirschfeld)
Periodontal Probes5-8
Periosteals
Periotomes
PerioWise® Probes/Screening System5
Pin and Point Setting Forceps 50
Placement
Plastic Filling
Pliers
Dressing
Endodontic
Pluggers/Condensers31-33, 53
Pocket Scalers
Potts Elevators
Probes5-8
R
Real World Endo Explorers
Retractors
Rongeurs
Root Canal Instruments 49-53
Rubber Dam Clamp & Punch
S
Scalers
Scalers 11-13 Scissors 55-59
Scalers
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58
Scalers
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralite™ - Big Easy® 11-16, 18-19
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralite™ - Big Easy® 11-16, 18-19 Ultrasonic Inserts 21
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralite™ - Big Easy® 11-16, 18-19
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralitet™ - Big Easy® 11-16, 18-19 Ultrasonic Inserts 21 Uni-Punch® Biopsy Punch (Disposable) 69
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralitet™ - Big Easy® 11-16, 18-19 Ultrasonic Inserts 21 Uni-Punch® Biopsy Punch (Disposable) 69
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralitet™ - Big Easy® 11-16, 18-19 Ultrasonic Inserts 21 Uni-Punch® Biopsy Punch (Disposable) 69 Universal Carriers 29
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralite™ - Big Easy® 11-16, 18-19 Ultrasonic Inserts 21 Uni-Punch® Biopsy Punch (Disposable) 69 Universal Carriers 29
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralite™ - Big Easy® 11-16, 18-19 Ultrasonic Inserts 21 Uni-Punch® Biopsy Punch (Disposable) 69 Universal Carriers 29 W Wax Instruments
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralite™ - Big Easy® 11-16, 18-19 Ultrasonic Inserts 21 Uni-Punch® Biopsy Punch (Disposable) 69 Universal Carriers 29 W Wax Instruments P. K. Thomas 46
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T Tissue Nippers 70 Tissue Pliers 64 U Ultralite™ - Big Easy® 11-16, 18-19 Ultrasonic Inserts 21 Uni-Punch® Biopsy Punch (Disposable) 69 Universal Carriers 29 W Wax Instruments P. K. Thomas 46 Wax Spatulas 45 Wedge Positioner 151
Scalers 11-13 Scissors 55-59 Sharpener (D.I.S.C.™) 12, 93 Slix™ Composite Manipulation Instruments 41-42 Spatulas 45 Spreaders 51-52, 65 Stewart Probe 49 Suture Scissors 58 T T Tissue Nippers 70 Tissue Pliers 64 U Ultralite™ - Big Easy® 11-16, 18-19 Ultrasonic Inserts 21 Uni-Punch® Biopsy Punch (Disposable) 69 Universal Carriers 29 W Wax Instruments P. K. Thomas 46 Wax Spatulas 45



Premier Products

Λ	L
Absorbent Points 163 ACCOR® Matrices 153 ACP Formula Products 88, 90-91 Alfa™ Triple Tray® 131, 139 Anesthetic - Topicale® 94 Articulator 131 Auto-Cure Activator 146	Enamel Pro® Gel. 91 Enamel Pro® Prophy Paste 88 Enamel Pro® Varnish. 90 Endo Finger Instruments 158-159, 161 Etchant 150 Evacuation System Cleaner 96 E-VENT® 91 Evian Perio Burs 115
Barbed Broaches161	
Biopsy Punches (Disposable) .69 Bite Tray/Registration	Files 158 K-Files 159 Hedstrom 159 Finger Pluggers 161 Finger Spreaders 161 Finishing Diamonds 113, 122 FS Hemostatic™ 130
	Gates Glidden Drills160
Surgical Length Carbide	Glitter® 92 Gutta Percha Points 162
CAD/CAMI also and a Fig. 11 1/2/2	Gutta Percha Accessory Points162
CAD/CAM Laboratory Finishing Kit122 Carbide Burs Surgical Length115, 163	H
Cartridge Warmer & Dispenser94	Hedstrom Files
Cassettes85	Hemodent®130
Cements	HP Rotary Organizer117
Cervical Matrices	
Clean-A-Diamond®	I
Cleaners 114 Brite Shield™ 96, 114 Pro Clean™ 97 Rapid-Klene® 97 Removalon®-I 96	Implant Cement™ .147 Impression Trays/Dual-arch .131, 139 Alfa™ Triple Tray® .127 NeoTray® .127
Vac Attak™96	T-LOC TM Triple Tray®
ComfortView®	Triple Tray®
	Instrument Sharpener D.I.S.C.™
Comfy-grip™ Prophy Cup Holders92	IntegraBond™146
CompCore™145	IntegraCem™146
CompCore™ AF	IntegraPost®
Compo-Disc®	IntegraPost®, Cure-Thru®
Compo-Strip®	
Compo-Strip® 200T2122	K
Composite Finishing122	K-Files
Cure-Shield152	K-riles
Cure-Thru®	NIIIL-FdK''''130
Blue Transparent Contoured Matrix Bands153	т
Clear Cervical Matrices152	L
IntegraPost®149	Laboratory Diamonds116
Reflective Curing Wedges152	Lab Tri-Fissure HP Carbide
	Lip & Cheek Retractor
D	Luminescence TM
	Luminescence Plus! TM
Diamonds	Lummescence Plus!""123
Diamond Discs116	M
Diamond Trimming Wheel - SilverStreak \dots . 116	M
Diamond Twist SCL™ & SCO™121	Mandrels117
D.I.S.C.™ Instrument Sharpener	Matrices
Disposable Diamonds 99-101	Cervical Matrices
•	cervical iviatrices152

Hygiene/Perio	87-94
Cleaners	95-97
Rotary	98-117
Finishing/Polishing	118-123
Prosthetic	124-131
Cosmetic/Whitening	132-141
Restorative	142-153
Endo	154-163

Matrix Bands	
Blue Cure-Thru®	153
Pre-contoured Metal	153
Stop Strip	153
MicroDose®	
38% Etchant (Enamel)	150
Porcelain Bonding Kit	150
Porcelain Etchant	150
Porcelain Silane	150
RC-Prep®	157
N	
NeoTray®	127
NexTemp™	
P	
Peeso Drills	163
Perfecta® Tooth Whitening	
Accessories	130
Bravó®	
Perfecta® 16%, 21%	
REV!®	
Standard Tubes	
Perfection White™ (In Office)	
Points	
Absorbent	163
Gutta Percha	
Poli-Pro Disks™	
Polishing Paste	
Porcelain Bonding Kit	
Porcelain Etchant	
Porcelain Silane	
Posts	
Cure-Thru® IntegraPost®	140
IntegraPost®	
PREMA® Enamel Microabrasion System .	
Premier Etch®	
Premier® Implant Cement™	
Pro Clean™	



Premier Products

Hygiene/Perio	87-94
Cleaners	95-97
Rotary	98-117
Finishing/Polishing	118-123
Prosthetic	124-131
Cosmetic/Whitening	132-141
Restorative	142-153
Endo	154-163

P (continued)

Prophy Angles	
2pro®	89
Prophy Paste	
Enamel Pro®	
Glitter®	
Ziroxide™	92
R	
Rapid-Klene®	97
RC-Prep®15	6-157
Reamers	159
Retraction Caps	129
Removalon®-I	96
REV!®	134
Root Fillers	160
S	
SensiPast®	160
Short Shank (ShortCut™) Diamonds	104
$Silver Streak^{\text{TM}} - Diamond \ Trimming \ Wheel \ \dots$	116
Solo Diamonds	9-101
Stop Strip	153
Sycamore Wood Wedges	151
T	
T-LOC™ Triple Tray®	127
Temporary Cement	147
Therma-Trol™	93
Thin-Flex® Diamond Discs	116
Tooth Whitening	
Accessories	
Perfecta® Cosmetic Patient Pouches	138
Parfacta® Block-Out	130

 Tray Material
 138

 Tray Cases
 138

Tooth Whitening (continued) Perfecta®
Bravó®135
Perfecta® 16%, 21%
REV!®
Standard Tubes
Topicale®
Topicale® Gel94
•
Topicale® GelPatch™94
Topicale® Ointment94
Topicale® Xtra94
Traxodent®
Tray Magic™ Trimmer140
Triple Tray®
2pro® Disposable Prophy Angles89
Two Striper® Diamonds
TS2000™ Diamonds 103-108, 110-111
Two Striper® All-Ceramic Prep Kit113
Two Striper® Diamonds102-113, 122, 162
Two Striper® Endo/Access Diamonds111, 162
Two Striper® MFS® Diamonds113, 122
Two Striper® MicroPrep™ Diamonds 104, 112
Two Striper® ShortCut™
(Short Shank) Diamonds104
Two Striper® Skooter™ Diamonds111
Two Striper® TSZtech™ Diamonds112
U
Ultrasonic Inserts
Ultrasonic Solution96
Uni-Punch® Biopsy Punch (Disposable)69
V
•
Vac Attak™96
•
Vac Attak™96 Varnish
Vac Attak™ 96 Varnish Enamel Pro® 90
Vac Attak™ 96 Varnish Enamel Pro® 90 W
Vac Attak™ 96 Varnish 90 W W Wedge Positioner 151 Wedges 151
Vac Attak™ 96 Varnish 90 Enamel Pro® 90 W 90 Wedge Positioner 151 Wedges 152 Cure-Thru® Reflective Curing 152
Vac Attak™
Vac Attak™ .96 Varnish
Vac Attak™
Vac Attak™ .96 Varnish