Alfa™ Triple Tray®

US/UK Alginate Full-Arch Closed-Mouth Impression Tray

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
One Step - Two Impressions

Simultaneously obtain maxillary and mandibular impressions - plus bite registrations!

Inspired solutions
for daily dentistry®

Manufacturer:
Premier® Dental Products Company
1710 Romano Drive,
Plymouth Meeting, PA 19462 U.S.A.
888-670-6100 • 610-239-6000 •
premierdentalco.com


Made in USA
Fabriqué aux États-Unis
Hecho en EE. UU.
Hergestellt in den USA
Fabbricato negli USA

U.S. Pat. Nos. D464,138; D413,386; D403,771
DPU# 10067/0221040 Rev8 PMM

EC

REP

Alfa™ Triple Tray® US/UK

Intended Use: Premier® Dual-Arch Trays® are intended to carry, confine and control impression material, while simultaneously capturing master/counter impression and bite registration.

Indications: All alginate preliminary impressions - ideal for fabrication of custom trays, tooth whitening trays, temporaries, study models, mouth guards/night guards, removable prosthetics and orthodontics. Alfa Triple Tray can also be used with alginate substitute vinyl polysiloxane (VPS) material. Do not use a rigid setting or heavy body VPS material.

Device Description: Premier® Dual-Arch Trays® are designed with thin transparent webbing or mesh so it can allow the capture of occlusal details during impressioning. The mesh also offers firm support of impression material before and after curing/setting. The combine dual-arch tray and impression material captures the three-dimensional shape of the teeth and oral tissue in the maxillary, mandibular arches and the bite registration.

Contents:

| Description | |
|---|---------|
| Alfa™ Triple Tray® Small (Purple) – Box of 24 | 1006313 |
| Alfa™ Triple Tray® Medium (White) – Box of 24 | 1006323 |
| Alfa™ Triple Tray® Large (Green) – Box of 24 | 1006333 |
| Alfa™ Triple Tray® Assorted – Box of 24 | 1006343 |
| • Five(5) –Small (Purple) • Fourteen(14) –Medium (White) • Five(5) –Large (Green) | |

Contraindications: Do not use Alfa Triple trays under the following conditions:

- Dental impression trays should not be used on patients with mouth/oral cancer.
- Patient is bleeding intraorally or has crevicular fluid seepage.
- If there is a missing tooth directly above or below the tooth being prepped.
- If there is a missing or absent tooth adjacent to the tooth being prepped on both (mesial and distal) sides.
- The tray has become damaged or distorted.
- Patient has lost or compromised dentition or restorations.

Precautions

- Tray should fit comfortably prior to use to ensure there is no impingement on the tray.
- Proper hygiene and tissue management shall be performed prior to taking an impression.
- Premature removal from the mouth.
- Not allowing enough working time.
- Only use with algine or alginate substitute impression materials.

Warnings

- The reuse of a single-use device can expose patients and staff to cross-contamination.
- Do not use with crown and bridge impression materials.

Adverse Events

- Incorrect removal technique may tear gingival tissue upon removal from mouth.
- If a serious incident has occurred in relation to the device, please report the incident to the Manufacturer and/or competent authority of the Member State where the device was used.

Materials:

| Alfa Triple Trays - Materials | |
|--|----------------------------|
| Trays – HDPE (High Density Polyethylene) FDA 21 CFR Part 177.1520 (Olefin Polymers) | Mesh – Cerex Nylon type 23 |

Directions for Use

- Select appropriate size tray. Medium fits most adult patients. Premier logo on top of handle ensures proper maxillary/mandibular orientation.
- Try in selected tray by inserting into the patient's mouth. It should fit passively, with the crossbar at the back of the tray distal to the most posterior teeth in arch. As with any impression making technique, ensure patient can breathe comfortably through nose with mouth closed.
- Prepare 2-3 scoops of alginate, depending on tray size (two for the small, three for medium or large). Load alginate into tray along outer walls to ensure it flows into retentive undercutts. Remove excess; alginate should extend 1-2mm just beyond tray wall. It is recommended that the maxillary anterior section be slightly overfilled to record the buccal vestibule.
- Place loaded tray into the mouth. Retract the cheek with a dental mirror to provide easier access. The tray should be centrally placed and not impinge on tissue. Have the patient place their tongue into the posterior notch of the tray.
- If necessary, assist the patient by pulling their lips beyond and over tray wall as they close into a natural bite.

6. Remove the set impression alginate impression by placing your finger on either purchase point and pressing down and out as you ask the patient to open. Once the suction is broken, the tray will be easily removed. Disinfect and inspect the impression. Pour the stone as per the manufacturer's directions.

Storage: Device storage and handling temperature should not exceed (37.7° C) 100° F.

Shelf Life: Non applicable

Disposal: Device should be disposed as biohazardous medical waste in accordance to facility or institutional guidelines/procedures.

For Further Information: To obtain SDS or IFU, visit premierdentalco.com or call Premier at 610-239-6000.

Lab Techniques UK/US

- Mix the stone and pour up one arch, making sure that no bubbles are incorporated into the stone.
- If needed, apply tray to vibrator to eliminate bubbles from stone.
- Wrap the tray in a damp towel to preserve the integrity of the alginate while mixing stone for the second arch.
- Once the first side is set, pour the second arch without removing the first arch from the tray. If needed, apply tray to vibrator to eliminate bubbles from stone. Allow stone to set.
- Remove tooth models from Alfa Triple Tray. Trim as needed.

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SUGGESTIONS: Consider using cold water when mixing alginate to increase working time. Apply a separating medium to the impression to facilitate easier removal of models from tray. Use model formers to create a model base. To save time, use a fast-set stone when pouring up the models.

If the stone from the second pour wraps around the tray and joins the first pour, separate the two arches using a lab bur/diamond to remove them from the tray.

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