# IntegraPost® Ultra-passive Endodontic Post


- FR Tenon Radiculaire Ultra Passive
- ES Poste Endodontic Ultra-Passivo
- DE Extra Passiv Endodontische Wurzelstift
- IT Perno Endodontico ultra-passivo
- TR Ultra Pasif Endodontik Post

Safe and Secure  
Easy to Use  
Innovative Design

## Instructions For Use

The primary function of the **IntegraPost®** is to provide retention for the overlying core build-up and final restoration, when insufficient coronal tooth structure remains.

**1 Preparing Post Space**  
Establish post space length and width by superimposing the enclosed acetate overlay size guide over a radiograph of the treated tooth. Select the smallest diameter post with sufficient coronal head length. After choosing the post size, select the corresponding drill size. Integra drills are sized slightly larger than the corresponding posts to maximize the adhesive and elastic properties of resin cements in the prepared canal space. **Sterilize before use.**

**2** Begin with the smallest Integra drill (yellow size 3), incrementally enlarge the canal, ending with the drill size that matches the post size selected in step 1. The cutting length (fluting) for all drills is 9mm. Maintain 4-5mm of Gutta Percha in the canal as an apical seal.

**4** Try-In the IntegraPost with the IP Carrier instrument. The shoulder (lower flange) should sit level with the top of the post space. If necessary, adjust the post to ensure proper seating. The length from the tip of the shank to the shoulder of the post should equal the radicular post space. For example, if the shoulder is sitting 2mm coronally to the top of the post space, remove 2mm from the tip of the shank. To adjust, remove from carrier and grasp post firmly and trim shank apically with a rotary burr or cutting pliers.

**3** Clean and rinse post space of debris. Residual Gutta Percha can be removed with a heated root canal plugger.

**5 Post Cementation (Etch-Bond-Lute)**  
Apply a dentin/enamel etchant such as Premier Etch, to the post space and coronal tooth structure for 10-15 seconds. Rinse with water for 5-10 seconds. Dry dentin in the canal until no visible moisture is apparent on the paper point.

**8** Remove solvents with a vigorous blast of air for at least 15 seconds or until no visible movement of the resin is observed. Apply multiple coats to the cavo surface and blow off solvents after each coat.

Light-cure cavo surface for 20 seconds.

**9** Cementation with a resin luting cement is recommended for maximum retention and compressive strength of the post and core restoration. If using IntegraCem, dispense equal amounts of base and catalyst onto a mixing pad and mix for 20 seconds.

**10** For cementation with resin cement, hold the IntegraPost in the carrier, and coat shank. Next apply cement to the post space with a paper point, plugger or Lentulo Spiral. Seat the coated post in the coated canal.

**Tip: The use of an activator with bonding agent can accelerate curing times of resin cements.**

**11** IntegraCem is dual-cured. Stabilize the post by light-curing the coronal aspect for 20 sec. IntegraCem's anaerobic accelerator ensures set at the bottom of the post space. Oral set time is approximately 3.5 minutes.

## Premier IntegraPost® System

The **IntegraPost** is a passive parallel design that is safe and easy to use.

Sterilize posts and drills before use.

The IntegraPost comes in six sizes and is manufactured from bio-compatible, non allergenic Titanium Alloy.

The micro-knurled shaft is pre-sandblasted, and designed to vent cement.

The unique round, flow-through head mechanically integrates core material for maximum retention and force distribution.

When handled in the IP carrier, the post is easy to position and pass.

Integra drills correspond numerically to the post sizes, but their diameters are sized slightly larger to maximize the adhesive and elastic properties of resin luting cements.

IntegraPost sizes 4 and 4.5 are most common and compare to Flexi-Post sizes 0 and 1. IntegraPost diameters generally correspond to ParaPost diameters. IntegraDrill numbers (sizes) correspond to ParaPost Drill numbers (sizes).

## Integra Probe Instrument

Included in the Intro Kit with Hot Box, the IP Probe (with blue stops) provides a quick and accurate method of determining post length to ensure that the shoulder is properly positioned at the top of the canal with the shank fully seated. Place the long end of the probe into the post space. The tip should go to the bottom. Move the rubber stop to a convenient reference point on the tooth, place the other end of the probe into the post space. The tip will rest on the top of the post space. Move the stop to the same reference point. Measuring both ends of the probe with a finger ruler, the difference will be the "post space length".

## Integra Drill

Technique Tip: Always start with the smallest drill and incrementally increase drill size to corresponding post size.

Size 3 Yellow	Size 4 Red
Size 4.5 Purple	Size 5 Blue
Size 6 Green	

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