

## 1. PURPOSE



These instructions are recommended to ensure cleaning, disinfection, and sterilization of select Premier® Endodontic devices (see section 2) before first use for non-sterile devices and before each reuse for reusable devices. This document aims to help healthcare professionals handle Premier® Endodontic devices in a safe manner as well as (re)process and maintain them in an appropriate manner, in accordance with the requirements of ISO 17664.

## 2. SCOPE OF APPLICATION

\*These instructions apply to the Premier® Endodontic Devices/Instruments listed here:

Instrument Description	Applicable SKUS
Premier® Barbed Broach	9053631, 9053633, 9053634, 9053636
Premier® Endo Drills	9055001, 9055002, 9055003, 9055004, 9055005, 9055006, 9055009, 9055021, 9055022, 9055023, 9055024, 9055025, 9055029
Premier® Finger Pluggers	9053791, 9053792, 9053793, 9053794
Premier® Finger Spreaders	9053814, 9053816, 9053818, 9053820, 9053826, 9053830, 9053832, 9053834, 9053836, 9053842
Premier® Hedstrom Files	9053399, 9053400, 9053401, 9053402, 9053403, 9053404, 9053405, 9053406, 9053407, 9053416, 9053417, 9053424, 9053425, 9053426, 9053427, 9053428, 9053429, 9053430, 9053431, 9053432, 9053433, 9053434, 9053435, 9053441, 9053442, 9053877, 9053878, 9053879, 9053893
Premier® K-Files	9053200, 9053201, 9053202, 9053203, 9053204, 9053205, 9053206, 9053207, 9053208, 9053209, 9053210, 9053211, 9053212, 9053218, 9053219, 9053230, 9053231, 9053232, 9053233, 9053234, 9053235, 9053236, 9053237, 9053238, 9053239, 9053240, 9053241, 9053242, 9053243, 9053244, 9053248, 9053249, 9053290, 9053291, 9053292, 9053293, 9053294, 9053295, 9053296, 9053297, 9053298, 9053299, 9053300, 9053301, 9053302, 9053304, 9053311, 9053312
Premier® Reamers	90053002, 90053003, 90053004, 9053030, 9053031, 9053032, 9053033, 9053034, 9053035, 9053036, 9053037, 9053039
Premier® Root Fillers	9053556, 9053560, 9053561, 9053562, 9053565, 9053571

Refer to the information on the label to determine the processing applicable to the device(s):

Sterility	Single-use device	Processing required before first use	Processing required after each use
		Yes	No
	No		Yes

## 3. WARNINGS AND PRECAUTIONS

### Warnings and precautions for the user:

- The devices covered by these instructions are intended for use in medical or hospital environments by qualified healthcare professionals.
- Use a dental dam when using the device(s) to avoid, for example, aspiration or ingestion by the patient.
- For your own safety, use personal protective equipment required during processing of the devices.
- For your own safety, wear surgical masks, gloves, and safety goggles.
- Carefully read the label or marking on the packaging to ensure you are using the correct device.

#### Warnings and precautions for the processing of devices:

- Use approved cleaning and disinfecting agents (e.g., approved by the VAH/DGHM or FDA, or bearing the CE marking) and use them according to the recommendations in their respective instruction manual.
- It is the user's responsibility to check the devices before each use in order to identify any possible defects.
- Cracks, deformations, signs of corrosion, loss in color or marking are signs that the device is no longer able to achieve the required performance level and should be discarded.
- Do not exceed a sterilization temperature of 135 °C.

#### In case of incident:

All serious events occurring in relation to the product must be reported to the manufacturer and the competent authority according to local regulations.

## 4. LIMITATIONS ON REPROCESSING

Generally speaking, any device showing visible signs of wear and tear or damage should be discarded (see sections 8 and 12).



#### Single use devices:

Devices labeled for single use only must not be reprocessed for reuse, as they are not designed to perform as intended after first use. Changes in mechanical, physical or chemical characteristics occurring with repeated use and/or (re)processing may compromise the integrity of the design and/or material, thus reducing the safety, performance and/or compliance of the device(s). When single use devices are supplied non-sterile and require sterilization before use, the relevant sections of these instructions are applicable.

#### Reusable devices:

These devices can be reused and reprocessed up to 8 times depending on the complexity of the canal to be treated. Where necessary, refer to the instructions in the device's instructions for use.

## 5. INITIAL PROCESSING AT THE POINT OF USE FOR REUSABLE DEVICES

After use, follow the steps below:

1. **Disassembly:** Remove the endo stop(s) from the instrument(s).
2. **Pre-cleaning:** Within a maximum of 30 minutes after use, remove excess soiling from the device(s) with disposable, lint-free wipes or a soft brush. Immerse the device(s) in a solution of water and neutral detergent.
3. **Rinsing:** Thoroughly rinse the device(s) with plenty of running water for at least 1 minute.

## 6. PREPARATION BEFORE CLEANING

#### Precautions:

- The device(s) should be reprocessed as soon as possible after use.
- The user should observe the concentrations and soaking times indicated in these instructions. An excessive concentration may cause corrosion or other defects on the devices.
- The disinfectant solution should not contain aldehyde so as to avoid fixation of blood residue.
- Do not use a disinfectant solution containing phenol, aldehyde or substances not compatible with the devices.
- The washer/disinfector must comply with ISO 15883 and undergo regular maintenance and calibration.

## 7. CLEANING/DISINFECTION

Follow one of the two methods described below (manual or automated) for cleaning and disinfection:

- Manual and mechanized devices before first use (if applicable, see section 2) and before each reuse (if applicable, see section 2).
- Endo stops before first use and before each reuse.

#### **Manual cleaning/disinfection:**

Equipment: Cleaning/disinfectant solution, brush, ultrasonic bath, purified running water, absorbent cloth.

1. Place the device(s) in a container, limiting any contact between the parts as much as possible.
2. Immerse the device(s) in the recommended cleaning/disinfectant solution. If necessary, use a soft nylon brush to gently scrub the device(s) until all visible soiling has been removed. If needed, use ultrasonic equipment as well.
3. Remove the device(s) from the solution and container and thoroughly rinse them under purified running water for at least 1 minute.
4. Dry the device(s) with single use absorbent cloth.

#### **Automated cleaning/disinfection:**

Equipment: Washer/disinfector, purified water, cleaning/disinfectant solution:

- Washing: Neodisher® Mediclean Forte (0.5 % concentration)
  - Thermal disinfection: Neodisher® Mediklar Special (0.03 % concentration)
1. Place the device(s) in a washer/disinfector basket, limiting any contact between the parts as much as possible.
  2. Process using a standard washer/disinfector cleaning cycle for at least 10 minutes at 93 °C or A<sub>0</sub> value > 3000 and complete with a hot air drying cycle for at least 15 minutes at 110 °C.

## **8. INSPECTION AND MAINTENANCE**

1. Before sterilization, discard any device(s) that has/have the following defects:
  - Plastic deformation
  - Bent device
  - Untwisted device
  - Damaged or blunt cutting edges
  - No marking
  - Corrosion
  - Discoloration
  - Other visible defects
2. Reassemble the endo stop(s) on the appropriate device(s).
3. Thoroughly inspect each device to check that all visible contamination has been eliminated. In case of contamination being observed, repeat the cleaning/disinfection process described above.

## **9. PACKAGING**

#### **Precautions:**

- Check the use-by date of the sterilization pouch stated by the manufacturer.
- Use packaging that can withstand temperatures up to 141 °C and complies with ISO 11607 and EN 868.
  1. The device(s) should be packed in a medical grade sterilization pouch (compliant with ISO 11607-1).  
Limit any contact between the devices and seal the pouches in accordance with the manufacturer's recommendations.

## 10. STERILIZATION

### Precautions:

- Autoclave (moist heat) sterilization using a pre-vacuum (forced air removal) cycle is recommended.
- Place the pouches in the sterilizer in accordance with the recommendations of the sterilizer's manufacturer.
- The autoclaves should comply with the requirements of the applicable standards (EN 13060 and EN 285) and should be approved, maintained and checked in accordance with these standards and the manufacturer's recommendations.
- Before any sterilization cycle, make sure that the maximum load indicated by the sterilizer's manufacturer is not exceeded.

Device class	Class B
Exposure time	Min. 3 minutes. The exposure time can be extended to 18 minutes to comply with the recommendations of the World Health Organization (WHO), the Robert Koch Institute (RKI), etc. These medical devices are able to withstand such sterilization cycles.
Temperature	134 °C
Drying time	Recommended: 20 minutes (minimum, in chamber)
Visual inspection	Check the device(s) in accordance with section 8 and verify proper performance of the sterilization cycle (packaging integrity, no humidity, color change of sterilization indicators, physical and chemical integrators, and digital records of various cycle parameters).





## 11. STORAGE














### Precautions:

- If the packaging has been opened, damaged or become wet, the sterile state of the devices inside the packaging is not guaranteed. Perform a new complete (re)processing cycle or discard the device(s).
  1. Store the device(s) in sterile packaging in a well-ventilated area, protected from dust, moisture, insects and temperature/humidity extremes, and at the temperature specified by the paper-plastic pouch by the manufacturer of the steam sterilizer.
  2. The packaging of the sterile devices should be carefully examined before opening (packaging integrity, no humidity, and expiry date) to ensure that the packaging's integrity has not been compromised during storage.

## 12. DISPOSAL

When a device reaches the end of its life, make sure that it is discarded in accordance with the applicable laws and regulations.

	Symbol	Designation	Description	Reference / ISO registration number
General symbols	 <a href="http://www.premierdentalco.com">www.premierdentalco.com</a>	Operating instructions	Indicates the need for the user to consult the applicable operating instructions and/or the processing instructions.	ISO 15223-1 (ISO no. 7000-1641)
	<b>Rx Only</b>	Prescription device	Caution : Federal law restricts this device to sale by or on the order of a "dentist/physician" licensed by the law of the State in which he/she practices to use or order the use of the device	21 CFR 801.109
Product identification & other information specific to each batch		Catalog number	Indicates the manufacturer's catalog number so that the medical device can be identified.	ISO 15223-1 (ISO no. 7000-2493)
		Batch code	Indicates the manufacturer's batch code so that the batch can be identified.	ISO 15223-1 (ISO no. 7000-2492)
		Date of manufacture	Indicates the date on which the medical device is manufactured. Date format: YYYY-MM-DD.	ISO 15223-1 (ISO no. 7000-2497)

Symbol	Designation	Description	Reference / ISO registration number	
Other symbols relating to product identification and specifications		Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1 (ISO no. 7000-1051)
		Clockwise rotation	Indicates that an instrument is intended to be used in a clockwise direction. This symbol is accompanied by a rotational speed expressed in revolutions per minute [rpm].	ISO 21531 (ISO no. 7000-0258)
		Material symbol	Indicates the material from which the device is made:  Nickel-titanium alloy  Stainless steel  Silicone  Plastic  Aluminum	ISO 9687 (ISO no. 7000-2793)
		Tip diameter	Indicates the tip diameter of an instrument, expressed in hundredths of a millimeter.	Not applicable
		Taper	Indicates the taper of an instrument, expressed in millimeters per millimeter of length (e.g., .02 corresponds to 2 %)	Not applicable
		Length	Indicates the usable length of the device.	Not applicable
		Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.	ISO 15223-1 (ISO no. 7000-2609)
		Sterilizable in a steam sterilizer (autoclave) at temperature specified	Indicates that the medical device is sterilizable in a steam sterilizer (autoclave) at the temperature specified.	ISO 9687 (ISO no. 7000-2868)

