Cure-Thru® IntegraPost®

Light-Transmitting Zirconia
Glass Fiber Post

US/UK

Premier® Cure-Thru® IntegraPost®
Light-Transmitting Zirconia Glass Fiber Post

Product Description:

Cure-Thru IntegraPosts are made of zircon-rich, AR-glass fibers for radiopacity and strength. These glass fibers are embedded in a special translucent epoxy resin matrix making the post compatible with light and/or dual curing composites (Premier CompCore AF). Cure-Thru IntegraPosts are factory silanated and chemically adhere to composite core build-up materials and resin cements resulting in a homogeneous all-composite restoration. The tip is tapered for an anatomical fit.

All post sizes are 19mm in length and color-coded by diameter size.
Size 4, Red, diameter 0.96mm / taper 4mm

Size 4.5, Purple, diameter 1.08mm / taper 4mm

Size 5, Blue, diameter 1.21mm / taper 6mm

Size 6, Green, diameter 1.46mm / taper 6mm

Cure-Thru Fiber Post posts can be cleaned with the usual decontamination products typically used in a dental office and are autoclavable up to 275°F / 135°C.

NOTE: To achieve maximum radiopacity use the lowest X-ray intensity and time exposure. Too much exposure can compromise visibility on X-rays.

INDICATIONS FOR USE:

The Cure-Thru Integra Post is indicated where less than half of the coronal tooth structure is remaining subsequent to the completion of endodontic therapy. The purpose of the post is to retain the core build-up material and reinforce the remaining coronal tooth structure.

DIRECTIONS FOR USE:

I Prepare canal for post:

1. Select post size and corresponding drill size by using a radiograph to plan canal depth. Cure-Thru IntegraPosts and drills are color-coded for easy size selection. All posts are 19mm in length with a tapered tip.

2. Remove obturation material by preferred method (Peeso drill, post drill, or hot plugger). Maintain a minimum of 4-5mm of gutta percha in canal as an apical seal.

3. Complete post space preparation and try in post for fit. Cut to length outside of mouth with a diamond disk. Decontaminate post with ethyl or isopropyl alcohol after try-in.

4. Clean or clean and disinfect the root-canal preparation, being sure to remove any residual eugenol sealer or cement. Follow manufacturer's instructions for etching all canal and coronal tooth surfaces. Thoroughly rinse and dry post-hole preparation prior to final post cementation.

II Cementing the post:

Carefully read and follow the specific instructions for use of the chosen adhesive and composite cement.

1. Apply selected bonding agent to the post and prepared root canal. Prepare the resin composite cement. Apply cement to paper point, lentulo spiral or plugger to deliver into canal preparation. Coat the portion of the post to be inserted into the root canal. Be sure to work quickly to avoid premature setting of the cement prior to inserting post.

2. Insert the post into the canal. If using more than one post, coat and insert first post before coating and placing subsequent post(s).

3. Light-cure cement through and around post.

III Core build-up: Premier CompCore AF or other composite core material recommended.

1. The final restoration must be circumferentially supported by sound tooth structure, not the core build-up (**Ferrule Effect**).

2. If using a matrix core-form like the Premier ACCOR® matrix, cut and fit.

3. Coat both the exposed post and the coronal tooth structure with selected bonding agent, being sure to follow manufacturers's recommendation for adhesive.

4. Syringe composite core material into remaining cavo surface and around exposed post. Light-cure core material according to manufacturers's instructions.

The Cure-Thru IntegraPost and core restoration is now complete. Preparation for final crown placement can be done upon complete curing of core build-up material.

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