Traxodent® Unit Dose

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Traxodent® Unit Dose

Hemostatic Retraction Paste System

Intended use / Indications for use:

Traxodent® has been specifically formulated to provide gingival retraction and hemostasis. It is intended for use prior to taking an impression, for cavity preparation or when hemostasis and retraction is required.

Product Description: Traxodent® is a product designed for dental hemostasis and/or retraction. The two main components consists of clay and 15% aluminum chloride. When the clay is in the sulcus, its bulk provides retraction. The aluminum chloride is an effective and proven dental astringent-hemostatic agent.

Contents:	
Item No.	Product Configurations
9007094	Traxodent Unit Dose Pack (24)
9007090	Traxodent Unit Dose Pack (50)

Accessories:	
Item No.	Product Configurations
9007089	Traxodent Unit Dose Syringe Dispenser

Contraindications: Do not use Traxodent in cases of:

- Do not use on persons with known allergy to aluminum chloride
- Patients with significant periodontal disease or furcation involvement.
- An intraoral site that has active infection or inflammation

Precautions:

- Unit dose capsules are designed for single patient use only. Do not reuse unit dose capsules
- Do not use if the unit dose dispenser or capsule is broken/damaged.
- Do not insert capsule dispensing tip into the sulcus.
- Avoid contamination of product. Do not return dispensed material (paste) to their original container. Dispose of leftover material properly.
- Avoid sunlight or nearby heat sources. Shelf life will be reduced if not stored properly.
- Prescription device, Rx only.

WARNINGS:

- DO NOT swallow retraction paste.
- Ensure all paste is removed by rinsing prior to dental impression, to avoid inhibition of polymerization of impression material
- Do not use on persons with known allergies to Aluminum Chloride
- Aluminum Chloride is irritating to skin/eyes.
- Do not use after the shelf life period indicated on product labeling/packaging.
- Do not freeze

Adverse Reactions:

- Adverse reactions may include but are not limited to:
 - Recession of the gums from prolonged use or excessive force in sulcus
 - Additional procedures may be necessary if impression material was inhibited due to residual paste.
- Allergic type of reactions such as swelling, inflammation, irritation, and pain could occur in persons with known allergy to Aluminum chloride.
- Cross-contamination and infection within the gums is possible if single use capsules are reused with multiple patients.

If a serious incident has occurred in relation to the device, please report the incident to the Manufacturer and/or competent authority of the Member State where the device was used.

Materials: 15% Aluminum Chloride Hexahydrate Epinephrine-free hemostatic agent in a clay matrix.

Directions for Use (Please follow carefully) FOR UNIT-DOSE

- Remove a Unit Dose Traxodent retraction capsule from the individually wrapped foil pack.
- Firmly seal the Unit Dose retraction capsule in the Premier Unit Dose syringe dispenser
- Remove the sealing cap from the Unit Dose capsule tip (immediately before application only) and discard.
- Note: To prevent sharp edges, do not modify the capsule tip.
- For initial use, dispense a small aliquot of the Unit Dose Traxodent paste onto a mixing pad.
- Rinse and air dry prepared tooth. Extrude Traxodent paste slowly into the sulcus while maintaining the dispensing tip just above the sulcus and align approximately parallel to the axial plane of the tooth preparation. Be careful not to jam dispensing tip into sulcus.
- Ensure the extrusion of sufficient volume of material into the sulcus and around the prepared tooth to achieve adequate tissue retraction.
- Allow Traxodent paste to remain in the sulcus for two to three minutes. Often, the marginal gingiva will grow pale when retraction is achieved.
- After hemostasis has been achieved, Traxodent Paste should be removed thoroughly by rinsing with a mixture of air and saliva ejector or suction tip.
- Examine treatment site to verify complete removal of the Traxodent prior to taking the impression.
- Dispose of unit dose capsule as per local laws

Reprocessing: (ONLY for Syringe Dispenser)

Cleaning: Inspect, clean and sterilize before each use. Do not use phenols or iodophors containing cleaning agents and do not heat above 177°C (350°F).

- Place each syringe dispenser in a separate paper/plastic steam sterilization pouch. If using a sterilizing cassette, ensure that the sterilizer's maximum load is not exceeded.
- Place bagged syringe dispenser into a steam autoclave per the autoclave manufacturer's instructions.
 - Pre-vacuum Steam Sterilization - Full Cycle: 132°C (270°F) for 4 minutes
 - To dry, use the drying cycle of the steam autoclave. Set cycle for 20 to 30 minutes.
 - Do not exceed 137°C (278°F).

Storage:

- Do not store at temperature over 28°C (82°F).
- Do not freeze.
- Store in a dry place

Shelf Life:

- Traxodent® has a shelf life of three (3) years (36 months) from the date of manufacturing.
- Do not use after product has expired. See expiry date on the product label.

For Further Information: To obtain an SDS or IFU visit www.premierdentalco.com or call Premier at 610-239-6000.

Symbols: A full explanation of symbols used on Premier Dental packaging and literature is located at premierdentalco.com/symbol-glossary.

Symbols	Symbol Meaning
	Expiry Date
	Date and Country of Manufacture
	Manufacturer
	Batch number
	Catalog Number
	Caution: See instructions for use
	Indicates the upper limit of temperature to which the device can be safely exposed.
	Irritant (Skin and Eye) Skin Sensitizer.
	Keep away from sunlight
	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician
	Single Use
	Medical Device
	European Union authorized representative
	Swiss Union authorized representative

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