

Priva Rotary Instruments™

Diamond Burs

Devices: Priva Rotary Instruments™ Diamond Burs

WARNINGS

General

1. Thoroughly dry the devices after all steps, if they will sit for more than 15 minutes prior to the next step.
2. All solutions used must be compatible with stainless steel dental burs (alkaline, acidic, or aldehyde solutions may cause damage).
3. Follow all applicable solution or equipment manufacturer's instructions.
4. All applicable solutions and equipment should meet national regulatory requirements.
5. **These devices are NOT supplied sterile.**
Clean and sterilize prior to use, as follows.

Processing Limitations

None, other than diamond wear, to be determined by the user.

INSTRUCTIONS

Initial Treatment at the Point of Use

Reprocess the devices as soon as possible after use. Do not allow embedded material to dry. Keep the devices immersed in solution from the time they are used until they are cleaned – use enzymatic cleaner, cleaning/disinfecting, cleaning, or water (in order of preference).

Preparation Before Cleaning

Not applicable

Cleaning: Automated

1. Rinse the devices under cold, heavy running tap water.
2. Ultrasonically clean the devices in enzymatic cleaner (Premier Brite Shield™ Enzymatic Static Soak and Ultrasonic Cleaner, or equivalent, dissolved in tap water) for at least five minutes. Avoid the devices in contact with each other during cleaning (do not pile).
3. Rinse thoroughly with tap water.

Cleaning: Manual

Not applicable

Disinfection

Not applicable

Drying

Dry by wiping with absorbent, disposable wipes, dry heat, or filtered, pressurized air. Devices that are not completely dry may corrode.

Maintenance, Inspection, and Testing

Visually examine the devices for contamination. Magnification is recommended. **Repeat the above steps for any devices that are still contaminated.**

Packaging

Sterilization pouch

Sterilization

Steam sterilize: pre-vacuum, average pressure of 30 psi (207 kPa), minimum of 4 minutes at 132°C, and a minimum dry time of 22 minutes.

Storage

Store the devices in the sterilized wrap.

Additional Information

Some surface discoloration or dulling may occur after moist heat sterilization. This does not affect performance.

Manufacturer Contact: +1 740-548-4100

The instructions provided above have been validated by the manufacturer as capable of preparing a medical device for use/reuse. It remains the responsibility of the processor to ensure that processing, as actually performed using equipment, materials, and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

Date Issued: 1/2022

Designed & Manufactured by:



Distributed by:

Premier® Dental Products Company
1710 Romano Drive Plymouth Meeting, PA 19462 U.S.A.
(888) 670-6100 • (610) 239-6000 • premierdentalco.com



premierdentalco.com/promise

IFU# 30074 • 1025054 Rev0 PDF