## T-LOC™ Triple Tray®

US/UK Instructions for Use – Adhesive-free Dual-Arch Impression Tray

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
UK Responsible Person: Parkway House, Palatine Rd  
Northenden, Wythenshawe  
Manchester M22 4DB  
United Kingdom

**CH REP** MDSS CH GmbH · Laurenzenvorstadt 61  
5000 Aarau, Switzerland

**EU REP** MDSS GmbH · Schiffgraben 41  
30175 Hannover, Germany

U.S. Pat. No. 6,749,428  
DFU# 10063 0426025 Rev8 PDF

### T-LOC™ Triple Tray® – US/UK

**Intended Use:** Premier® Dual-Arch Trays® are intended to carry, confine and control impression material, while simultaneously capturing master/counter impression and bite registration.

**Indications:** T-LOC Triple Trays are indicated for crown and bridge procedures - for one or two unit cases per quadrant. Also used in the creation of veneers, inlays, onlays, provisional restorations, single tooth implant crowns and bridges.

**Device Description:** Premier® Dual-Arch Trays® are designed with thin transparent webbing or mesh so it can allow the capture of occlusal details during impressioning. The mesh also offers firm support of impression material before and after curing/setting. The combine dual-arch tray and impression material captures the three-dimensional shape of the teeth and oral tissue in the maxillary, mandibular arches and the bite registration.

**Contents:**

Description	
T-LOC™ Triple Tray® Posterior – Box of 35	<b>1006204</b>
T-LOC™ Triple Tray® Posterior – Case of 10	<b>1006204CS</b>
T-LOC™ Triple Tray® Anterior – Box of 28	<b>1006214</b>
T-LOC™ Triple Tray® Anterior – Case of 10	<b>1006214CS</b>
T-LOC™ Triple Tray® Three-Quarter Arch – Box of 22	<b>1006274</b>
T-LOC™ Triple Tray® Three-Quarter Arch – Case of 10	<b>1006274CS</b>

**Contraindications:** Do not use T-LOC Triple trays under the following conditions:

- Dental impression trays should not be used on patients with mouth/oral cancer.
- Patient is bleeding intraorally or has crevicular fluid seepage.
- If there is a missing tooth directly above or below the tooth being prepped.
- If there is a missing or absent tooth adjacent to the tooth being prepped on both (mesial and distal) sides.
- The tray has become damaged or distorted.
- Patient has lost or compromised dentition or restorations

**Precautions**

- Tray should fit comfortably prior to use to ensure there is no impingement on the tray.
- Triple Tray with a firm setting material may be used with hard-setting Heavy Body/Light Body one-step techniques, not two-step techniques.
- Proper hygiene and tissue management shall be performed prior to taking an impression.
- Premature removal from the mouth.
- Not allowing enough working time.

**Warnings**

- The reuse of a single-use device can expose patients and staff to cross-contamination.
- Ensure patient has no allergies to adhesives or impression materials prior to use.

**Adverse Reactions**

- Incorrect removal technique may tear gingival tissue upon removal from mouth.

If a serious incident has occurred in relation to the device, please report the incident to the Manufacturer and/or competent authority of the Member State where the device was used.

**Materials:**

Materials	
Trays – HDPE (High Density Polyethylene)/ABS (Acrylonitrile Butadiene Styrene) FDA 21 CFR Part 177.1520 (Olefin Polymers)	Mesh – Cerex Nylon type 23

**Directions for Use**

- Try-in unloaded tray. The crossbar of the posterior tray must be distal to the last tooth. Center the tray buccally/lingually.
- Have patient “close” on the unloaded tray and check opposite side of patient’s mouth to make sure that the patient is closed into maximum intercuspation.
- Overload both sides of tray with firm-setting impression material. It may be used with hard-setting Heavy Body/Light Body one-step techniques. Two-step techniques are contraindicated.
- Remove retraction cord, if used. Syringe impression material around important cervical, interproximal and edentulous areas.
- With syringed impression material in place and still unset, place loaded tray into proper position. (Make sure crossbar is correctly placed.)
- Have patient close into maximum intercuspation.
- When set, gently remove from mouth. Inspect final impression. No plastic pressure areas on the tray walls should be visible.
- Disinfect the impression tray using your preferred impression disinfectant before sending it to the dental lab.

**Storage:** Store in a cool, dry place away from direct sunlight and excessive heat sources.

**Shelf Life:** Non applicable .

**Disposal:** Device should be disposed as biohazardous medical waste in accordance to facility or institutional guidelines/procedures.

**For Further Information:** To obtain SDS or IFU, visit premierdentalco.com or call Premier at 610-239-6000.

Symbols	Symbol Meaning
	Manufacturer / Hersteller / Fabricante / Fabricant / Fabbricante / Tillverkare / Producent

	Recyclable / Recyclebar / Reciclabile / Recyclable / Riciclabile / Återvinningsbar / Genanvendelig
--	--

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**Lab Techniques** – **UK/US**

Evaluate impression for any visible plastic pressure areas which may affect the accuracy of the impression.

- Trim excess impression material from the borders of the impression.

**Helpful hint:** Before pouring, mark the handle of the upper arch. This will eliminate articulating the model upside down.
- Box impression containing die with wax and position dowel pin. Seal off any holes that stone may leak through. First, the working side is poured, adding stone or epoxy in small increments. Next, the counter side is poured. Allow to set. Pour base and allow to set. Remove boxing wax.
- Articulate model and let plaster set.
- Separate casts by pulling apart carefully.
- Continue with normal crown and bridge procedures.

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