AeroPro™

Cordless Prophy System

REF 5500520 Autoclavable Outer Sheath

EN

Instructions for use: Please refer to the AeroPro™ Directions for Use manual for proper operating instruction for the Outer Sheath. The Outer Sheath is designed for use only with the Premier® AeroPro™ Cordless Prophy System and should not be used with any other brands of cordless prophy handpiece.

The enclosed Outer Sheath is not sterile. Upon receipt, it must be steam autoclave sterilized prior to use in accordance with the following instructions:

Instructions for steam autoclave sterilizing: The following instructions are for use ONLY on the Outer Sheath. Please refer to Section 4.8 Infection Prevention Procedures of the AeroPro Directions for Use manual for proper cleaning, disinfecting or sterilizing methods for all other AeroPro components

Disinfection: Do not use chemical disinfectants prior to sterilization or rapid deterioration of the materials may occur.

Packaging: Place each Outer Sheath in a separate paper/plastic steam-sterilizing pouch. If using a sterilization cassette ensure that the sterilizer's maximum load is not exceeded.

Steam Autoclave Sterilization: Place bagged Outer Sheath into a steam autoclave, per the autoclave manufacturer's instructions.

- **Gravity Steam Sterilization**
 - o Full Cycle: 135°C (275°F) for 10 minutes
- **Pre-vacuum Steam Sterilization**
 - o Full Cycle: 132°C (270°F) for 4 minutes

Alternate method: Place non-bagged Outer Sheath into the steam autoclave and run at the above listed cycles.

NOTE: Outer Sheath should remained bagged until ready for use.

Drying: To dry, use the drying cycle of the steam autoclave. Set cycle for 20 to 30 minutes. Do not exceed 137°C (278°F).

Cold liquid disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.

Manufacturer contact:

Call Premier Dental Customer Experience at 888-670-6100 or 610-239-6000.

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1710 Romano Drive, Plymouth Meeting, PA 19462 U.S.A. • (610) 239-6000