IntegraCem™

Dual-Cure Resin Cement

Premier® Dental Products Company
1710 Romano Drive, Plymouth Meeting, PA 19462 U.S.A.
(888) 670-6100 • (610) 239-6000 • www.premusa.com

Made in Israel
Technology in partnership
BJM Laboratories Ltd. /

Indications for Use: Premier® IntegraCem™ is a resin, dual-cure universal dental adhesive cement. IntegraCem can be used for the cementation of endodontic posts, restorations to enamel, dentin or metal, amalgam bonding, porcelain repair, and orthodontic attachment bonding.

Physical Properties:

- Working time is 3 minutes. While on the pad, ambient air inhibits setting. When placed orally, heat and anaerobic conditions accelerate set. Initial set is 1.5 minutes after mix. Final set is 3.5 minutes.
- Dual-cure adhesives can shorten working time of dual-cure resin cements.

Contraindications: Do not use on patients with known sensitivity or allergy to resins. Wash exposed area completely to remove material. If an allergic reaction is observed, consult a physician immediately and discontinue use with this patient.

To obtain a SDS, visit Premier Dental at www.premusa.com or call 610-239-6000.

Storage: Do not store IntegraCem over 75°F (24°C). Refrigerate at 50°F (10°C) for maximum shelf-life. Allow material to warm to room temperature before use.

Shelf life: 24 months from date of manufacture. Refrigerate for maximum shelf-life.

Rx ONLY

DFU# 33140 • 1217060 Rev3 BJM

1 Cementation of Endodontic Post

- Remove all debris from the post space.
- Apply Premier® Etch to the post space and cavo surface for 10 to 15 seconds. Rinse to remove etchant with a spray of oil-free water for 5 to 10 seconds.
- Dentin in the canal should be immediately and thoroughly dried until no visible moisture is apparent on the paper point.
- Mix equal number of drops of IntegraBond™ and Premier Auto-Cure Activator in a clean mixing well for 2 seconds. If using another manufacturer's bonding agent, follow directions for dark-curing the canal space.
- Apply generous amounts of mixed IntegraBond and Activator with a clean brush to the moist cavo surface. Use a paper point to coat the deepest part of the canal. Leave undisturbed for 15 seconds.
- Remove solvents with a vigorous blast of air for at least 15 seconds or until no visible movement of the resin is noted.
- Dispense equal amounts of IntegraCem™ Base and Catalyst onto a mixing pad and mix for 20 seconds until a smooth, uniform paste results. Working time at room temperature is at least 3 minutes.
- Apply mixed cement on to post and into post space using Premier Sensi-Past® Root Fillers or alternately apply cement with paper points. Seat the post.
- Stabilize post by light-curing coronal portion of post for 20 seconds. Oral set time is approximately 3.5 minutes. Continue with core build-up after cement sets.

2 Metal Restoration to Dentin or Enamel

Adhesive Cementation of Crown & Bridge, Inlays, and Maryland Bridges

- Apply Premier® Etch to tooth for 10 to 15 seconds. Caution should be used when etching thin dentin.
- Rinse with a spray of oil-free water for 5-10 seconds to remove etchant. Dentin should be blotted until there is no pooling of water and surface is left moist and glistening.
- Apply a generous coat of IntegraBond™ to the prepared tooth surface; leave undisturbed for 15 seconds.**
- Remove solvents with a vigorous blast of oil-free air for at least 15 seconds. Apply 2 to 3 layers of IntegraBond; dry with a full blast of oil-free air after the last application. Look for light-reflection (shiny surface) as an indication of sufficient coating. If necessary, apply additional coats of IntegraBond until light-reflection is achieved. Light-cure for 20 seconds.
- Dispense equal amounts of IntegraCem™ Base and Catalyst onto a mixing pad and mix for 20 seconds until a smooth, uniform paste results. Working time at room temperature is at least 3 minutes.
- Apply mixed cement to the restoration.
- Seat restoration; apply pressure until IntegraCem reaches initial set, 1.5 minutes after mix. For easy removal of excess, light-cure all sides 2-3 seconds, no more than 10 seconds total. Check for firmness and remove excess.
- Following removal of flash, light-cure margins for 20 seconds from each direction. Oral set time is approximately 3.5 minutes.

** For the thinnest possible film thickness (i.e. inlays), and / or when light curing IntegraBond is not convenient, IntegraBond can be self-cured by utilizing the Premier Auto-Cure Activator. Mix equal drops of IntegraBond™ and Premier Auto-Cure Activator in a clean mixing well for 2 seconds and proceed with step 3. Avoid light curing bonding agent.

3 Ceramic Restoration to Dentin/Enamel

Adhesive Cementation of All Ceramic Crown & Bridge, Inlays, and Onlays

- Apply Premier® Etch to tooth for 10 to 15 seconds. Caution should be used when etching thin dentin.
- Rinse with a spray of oil-free water for 5-10 seconds to remove etchant. Dentin should be blotted until there is no pooling of water and surface is left moist and glistening.
- Etch internal surface of porcelain restoration with MicroDose™ Porcelain Etch for 60 seconds. Rinse and dry the surface well with a blast of oil-free air. The etched porcelain should have a frosted white appearance.
- Brush MicroDose™ Porcelain Silane on etched porcelain surface and let dry for 3 minutes.
- Apply a generous coat of IntegraBond™ to the moist prepared tooth surface; leave undisturbed for 15 seconds.**
- Remove solvents with a vigorous blast of oil-free air for at least 15 seconds. Apply IntegraBond 2 to 3 more times; dry with a full blast of oil-free air after the last application. Look for light-reflection (shiny surface) as an indication of sufficient coating. If necessary, apply additional coats of IntegraBond until light-reflection is achieved. Light-cure for 20 seconds.
- Dispense equal amounts of IntegraCem™ Base and Catalyst onto a mixing pad and mix for 20 seconds until a uniform paste results. Working time at room temperature is at least 3 minutes.
- Apply mixed cement to restoration. Seat restoration. Apply pressure until IntegraCem reaches initial set, 1.5 minutes. For easy removal of excess, light-cure all sides 2-3 seconds, no more than 10 seconds total. Check for firmness and remove excess.
- Following removal of flash, light-cure margins for 20 seconds from each direction. Oral set time is approximately 3.5 minutes.

** For the thinnest possible film thickness (i.e. inlays), and / or when light curing IntegraBond is not convenient, IntegraBond can be self-cured by utilizing the Premier Auto-Cure Activator. Mix equal drops of IntegraBond™ and Premier Auto-Cure Activator in a clean mixing well for 2 seconds and proceed with step 3. Avoid light curing bonding agent.

4 Amalgam Bonding

To Form a Chemical Bond between Tooth Surface, Cement, and Amalgam

- Apply Premier® Etch to tooth for 10 to 15 seconds. Caution should be used when etching thin dentin.
- Rinse with a spray of oil-free water for 5-10 seconds until all etchant is removed. Dentin should be blotted until there is no pooling of water and surface is left moist and glistening.
- Mix equal number of drops of IntegraBond™ and Premier Auto-Cure Activator in a clean mixing well for 2 seconds.
- Apply generous amounts with a clean brush onto the moist cavity walls; leave undisturbed for 15 seconds.
- Remove solvents with a vigorous blast of air for at least 15 seconds or until no visible movement of the resin is noted.
- Dispense equal amounts of IntegraCem™ Base and Catalyst onto a mixing pad and mix for 20 seconds until smooth, uniform paste results. Working time at room temperature is at least 3 minutes.
- Apply mixed cement to the cavity walls. Use microbrush to absorb excess cement that might be pooling in corners of preparation. Proceed immediately with placing and condensing of amalgam. NOTE: The matrix band should be waxed or covered with separating medium to avoid cement adhesion. Carving and finishing of the restoration may proceed immediately following the placement of the amalgam.

5 Intraoral Porcelain Repair

Veneering Metal with Composite - Use of a Rubber Dam is Mandatory

- Roughen exposed metal surface with a small round Two Striper® Diamond bur, Solo™ Disposable Diamond bur or Micro Etcher clinical sandblaster.
- Rinse the sand-blasted or roughened metal and dry.
- If metal surface is partly covered with porcelain, etch exposed edge of porcelain with MicroDose™ Porcelain Etch for 60 seconds, rinse, and dry well. Apply MicroDose Porcelain Silane, let dry for 3 minutes. Brush on a layer of IntegraBond to exposed edge of porcelain. Remove solvents with a vigorous blast of air for at least 15 seconds and light-cure for 20 seconds.
- Dispense equal amounts of IntegraCem™ Base and Catalyst onto a mixing pad and mix for 20 seconds until a smooth, uniform paste results. Working time at room temperature is at least 3 minutes.
- Apply a thin, uniform layer of cement, using a plastic instrument or injection system, onto exposed metal surface only (avoid contacting exposed porcelain edge with cement) and light-cure for 20 seconds. The IntegraCem opaques the metal surfaces and bonds to it.
- Apply composite resin of your choice and finish.

6 Metal Restoration to Metal

Adhesive Cementation of Crown and Bridge to Cast Core Implant Prosthesis

- Roughen internal surface of the crown with a small round Two Striper® Diamond bur, Solo™ Disposable Diamond bur or Micro Etcher clinical sandblaster. The same surface treatment can be used on a cast core or implant abutment, prior to placement in mouth.
- Rinse the sand-blasted or roughened metal and dry.
- Dispense equal amounts of IntegraCem™ Base and Catalyst onto a mixing pad and mix for 20 seconds until smooth, uniform paste results. Working time at room temperature is at least 3 minutes.
- Apply mixed cement to bonding surface of crown.
- Seat restoration; apply pressure until IntegraCem reaches initial set 1.5 minutes after mix. For easy removal of excess, light-cure all sides 2-3 seconds, no more than 10 seconds total. Check for firmness and remove excess.
- Following removal of flash, light-cure margins for 20 seconds from each direction. Oral set time is approximately 3.5 minutes.

7 Orthodontic Attachment Bonding

Bonding Metal, Plastic or Ceramic Orthodontic Attachments on Teeth, Acrylic and Ceramic Crowns

- Clean teeth with plain pumice using a rotating bristle brush or rubber cup.
- Rinse, isolate and air-dry to completely dry the bonded surfaces. For anterior teeth, add etching is not required.
- Dispense IntegraBond™ directly onto a microbrush.
- Apply a generous coat of IntegraBond onto the enamel or porcelain surfaces. Remove solvents with a vigorous blast of air for at least 15 seconds and light-cure for 20 seconds.
- Dispense equal amounts of IntegraCem™ Base and Catalyst onto a mixing pad and mix for 20 seconds until a smooth, uniform paste results. Working time at room temperature is at least 3 minutes.
- Apply small amount of cement on bracket base. Cement should be added to one bracket at a time.
- Press bracket to tooth with slight rotating action. A small amount of cement should exude around the entire bracket base periphery. Final bracket adjustment must be made before cement begins to cure (approximately 60 seconds).
- For snap setting, IntegraCem can be light-cured for 20 seconds from each direction.
- Keep teeth isolated for 3 minutes after last bond.
- Arch wires may be placed 10 minutes after last bond.

Ordering Information

3001410	IntegraCem™ (4g each of base and catalyst)	CompCore™ AF SyringeMix™ (Two 9g automix syringes, 15 mix-tips, 15 standard and 10 extra-fine intraoral tips)
3001421	IntegraBond™ Adhesive Kit (4ml IntegraBond, 4ml Auto-Cure Activator, two 1.2ml Premier Etch Syringes, 50 Microbrushes, 20 Disposable Etch Tips, Mix Wells)	3001512 Shade White 3001513 Shade A3 3001514 Shade Blue
3001411	IntegraBond™ (7ml)	CompCore™ AF Cartridge - Intro Kits (One 50g cartridge, type-25 gun, 4ml each IntegraBond and Auto-Cure Activator, 25 mix-tips, 25 intraoral tips)
3001413	IntegraBond™ (4ml)	3001430 Shade White 3001431 Shade A3
3001412	Auto Cure Activator (7ml)	CompCore™ AF Cartridge - Refills (One 50g cartridge, 25 mix-tips, 25 intraoral tips)
3001422	Premier® Etch 37% (One 12g jumbo syringe, with luer lock adapter, ten empty 1.2ml syringes and 50 disposable applicator tips for the ten empty 1.2ml syringes)	3001432 Shade White 3001433 Shade A3
3001461	Spot-On™ Etch Intro 5.12ml filled syringes 20 bent disposable applicator tips	4009001 MicroDose™ Porcelain Etch (22 units; 0.2cc each)
3001460	Spot-On™ Etch Bulk 1.12g jumbo filled syringe with luer lock adapter 10 1.2ml empty syringes 50 1.2ml disposable applicator tips	4009020 MicroDose™ Porcelain Silane (18 units; 0.2cc each)

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