**premier**

## IntegraPost® Ultra-passive Endodontic Post

FR Tenon Radiculaire Ultra Passive  
ES Poste Endodontic Ultra-Passivo  
DE Extra Passiv Endodontische Wurzelstift  
IT Perno endodontico ultra-passivo  
TR Ultra Pasif Endodontik Post

**Safe and  
Secure**


**Easy to Use**

**Innovative  
Design**

### Instructions For Use

The primary function of the **IntegraPost®** is to provide retention for the overlaying core build-up and final restoration, when insufficient coronal tooth structure remains.

#### Premier IntegraPost X-ray Overlay

- 1 Preparing Post Space**  
Establish post space length and width by superimposing the enclosed acetate overlay size guide over a radiograph of the treated tooth. Select the smallest diameter post with sufficient coronal head length. After choosing the post size, select the corresponding drill size. Integra drills are sized slightly larger than the corresponding posts to maximize the adhesive and elastic properties of resin cements in the prepared canal space.


**Sterilize before use.**

**5 Post Cementation (Etch-Bond-Lute)**


Apply a dentin/enamel etchant such as Premier Etch, to the post space and coronal tooth structure for 10-15 seconds. Rinse with water for 5-10 seconds. Dry dentin in the canal until no visible moisture is apparent on the paper point.


- 6** If bonding with Integra System products, mix equal drops of IntegraBond™ and Auto-Cure Activator in a clean mixing well for two seconds. If using another manufacturer's adhesive, follow their instructions for dark-curing the canal space. If IntegraBond/Activator mix is being used with dual-cure or self-cure composite build-ups like CompCore™ AF, apply multiple coats to the cavo-surface and light-cure the bond/activator mixture.

- 7** Apply generous amounts of mixed IntegraBond and Activator to all tooth surfaces, including the cavo surface for bonding to the composite core material. Multiple coats of the bond/activator mixture on the cavo surface are highly recommended. Use a paper point or tufted applicator to coat the deepest part of the canal. Leave undisturbed for 15 seconds.


- 8** Remove solvents with a vigorous blast of air for at least 15 seconds or until no visible movement of the resin is observed. Apply multiple coats to the cavo surface and blow off solvents after each coat.

- 9** Cementation with a resin luting cement is recommended for maximum retention and compressive strength of the post and core restoration. If using IntegraCem, dispense equal amounts of base and catalyst onto the cavo surface and blow off solvents after each coat.

- 10** For cementation with resin cement, hold the IntegraPost in the carrier, and coat shank. Next apply cement to the post space with a paper point, plugger or Lentulo Spiral. Seat the coated post in the coated canal.

- 11** IntegraCem is dual-cured. Stabilize the post by light-curing the coronal aspect for 20 sec. IntegraCem's anaerobic accelerator ensures set at the bottom of the post space. Oral set time is approximately 3.5 minutes.

- 12 Core Build-Up**  
Once the post is firmly cemented, begin the core build-up. If using a core form like Premier's Accor Matrix, trim and seat matrix. Syringe Premier CompCore™ AF first against the floor of the pulp chamber and then around the post flange. In the same motion syringe the core material through the top hole in the post head and continue stacking to form the finished core.

- 13 The IntegraPost acts as an integrated core retainer for the final restoration. The rounded flow-thru head mechanically locks the core material to the post. The combination of this head design with undercuts provided by the flanges, creates anti-rotational stability and helps eliminate voids and reduce stress points for the life of the post and core restoration.**

### Premier Integrapost® System

The **IntegraPost** is a passive parallel design that is safe and easy to use.

Sterilize posts and drills before use.

The Integrapost comes in six sizes and is manufactured from bio-compatible, non-allergenic Titanium Alloy.

The micro-knurled shaft is pre-sandblasted, and designed to vent cement.

The unique round, flow-through head mechanically integrates core material for maximum retention and force distribution.

When handled in the IP carrier, the post is easy to position and pass.

### Integra Probe Instrument

Included in the Intro Kit with Hot Box, the IP Probe (with blue stops) provides a quick and accurate method of determining post length to ensure that the shoulder is properly positioned at the top of the canal with the shank fully seated. Place the long end of the probe into the post space.

The tip should go to the bottom. Move the rubber stop to a convenient reference point on the tooth, place the other end of the probe into the post space. The tip will rest at the top of the post space. Move the stop to the same reference point. Measuring both ends of the probe with a finger ruler, the difference will be the "post space length".

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	REF
Posts	
Intro Kit / Shade White.....	3001430
Intro Kit / Shade A3.....	3001431
Intro Kit contains one 50g cartridge, type-25 gun, 4ml each IntegraBond & Auto-cure Activator, 25 mix-tips, 25 intraoral tips	
Cartridge Stack: Dual-Cure	
IntegraPost Mini Intro Kit.....	3001001
IntegraPost Refill (10) Size 3.....	3001105
IntegraPost Refill (10) Size 4.....	3001110
IntegraPost Refill (10) Size 4.5.....	3001115
IntegraPost Refill (10) Size 4.5L.....	3001120
IntegraPost Refill (10) Size 5.....	3001125
Cartridge Refills	
CompCore AF Refill Shade White.....	3001432
CompCore AF Refill Shade A3.....	3001433
Refill contains one 50g cartridge, 10 ACCOR® matrices, 25 mix-tips, 25 intraoral tips	
CompCore Auto-Cure/Tub Hand-Mix (not CE marked)	
Shade White (28g).....	3001402
Shade Natural (28g).....	3001405
Shade White (50g).....	3001401
Shade Natural (50g).....	3001404
Tub contains catalyst & base, mixing pads, spatulas	
Accessories	
SyringeMix Mix-Tips & Intraoral Tips (35).....	3001418
SyringeMix XXF Extra-fine Intraoral Tips (50).....	3001468
Aplicator Gun, Type-25.....	3001403
Cartridge Mix-Tips & Intraoral Tips (35).....	3001417

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