## T-LOC™ Triple Tray®

Adhesive-free

Dual-arch Impression Trays

### Directions For Use

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
 Manufacturer:  
Premier® Dental Products Company

1710 Romano Drive, Plymouth Meeting, PA 19462 U.S.A.  
(888) 670-6100 • (610) 239-6000 • premierdentalco.com

U.S. Pat. No. 6,749,428  
REF 1006204 Posterior (35 trays/box)  
REF 1006214 Anterior (28 trays/box)  
REF 1006274 Three-Quarter Arch (22 trays/box)

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### US/UK

#### T-LOC™ Triple Tray® Technique

After Completion of Tooth Preparations and Gingival Retraction

1. Try-in unloaded tray. The crossbar of the posterior tray must be distal to the last tooth. Center the tray buccally/lingually.
2. Have patient "close" on the unloaded tray and check opposite side of patient's mouth to make sure that the patient is closed into maximum intercuspsation.
3. Overload both sides of tray with firm-setting impression material. It may be used with hard-setting Heavy Body/Light Body one-step techniques. Two-step techniques are contraindicated.
4. Remove retraction cord, if used. Syringe impression material around important cervical, interproximal and edentulous areas.
5. With syringed impression material in place and still unset, place loaded tray into proper position. (Make sure crossbar is correctly placed.)
6. Have patient close into maximum intercuspsation.
7. When set, gently remove from mouth. Inspect final impression. No plastic pressure areas on the tray walls should be visible. Disinfect.

**Caution:** Using a polyether and / or firm setting impression material with the Three-Quarter Arch T-LOC™ Triple Tray® may cause excessive retention and be difficult to remove.

#### Laboratory technique for pouring Triple Trays

*Evaluate impression for any visible plastic pressure areas which may affect the accuracy of the impression.*

1. Trim excess impression material from the borders of the impression.

**Helpful hint:** Before pouring, mark the handle of the upper arch. This will eliminate articulating the model upside down.

2. Box impression containing die with wax and position dowel pin. Seal off any holes that stone may leak through. First, the working side is poured, adding stone or epoxy in small increments. Next, the counterside is poured. Allow to set. Pour base and allow to set. Remove boxing wax.

3. Articulate model and let plaster set.

4. Separate casts by pulling apart carefully.

5. Continue with normal crown and bridge procedures.

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