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premier®

Premier® Ultrasonic Inserts

For use in
magnetostrictive
ultrasonic scalers.

premier®

Manufacturer:
Premier® Dental Products Company
1710 Romano Drive,
Plymouth Meeting, PA 19462, U.S.A.
610-239-6000
premierdentalco.com

For intended use by dental
professionals only.

US/UK

For use in magnetostrictive ultrasonic units

Indications:

The ultrasonic scaler and inserts are intended for use by medical professionals during dental cleaning and periodontal therapy to remove calculus and tartar deposits and stains from teeth using the application of an ultrasonically vibrating tip.

Contraindications:

Use of any ultrasonic device is not recommended in cases where the operator or patient has been fitted with a cardiac pacemaker.

Caution:

- Inserts are intended for dental hygiene (cleaning) of dental implants and prostheses only as a preventive procedure

- Not intended to be used as tool to assist other restorative, surgical, operative, or similar dental procedures

- Rx Only

Cleaning/ Sterilization:

Inspect, clean, and sterilize prior to each use. Do not use chemical disinfectants prior to sterilization or leave residual cleaning solutions on insert. Carefully remove gross debris, rinse and dry inserts completely prior to sterilization. *Premier® Dental has independently validated their instruments to be sterile when using a pulsed pre-vacuum steam sterilization process for five minutes at 134°C (273°F).*

- Do not use Rapid Heat Transfer (dry heat) sterilization method.

- Do not expose inserts to temperatures above 350°F/177°C or to chemical vapor concentrations exceeding manufacturer's recommendations.

- Do not use phenols.

- Do not use a bleaching agent.

- Do not use cold disinfectant/cold sterilization solutions.

- Do not sterilize in 100% plastic sterilization pouches.

(A combination paper/plastic pouch is acceptable, paper side down.)

Optimum Performance:

- Insert stacks should be straight and parallel at all times. Bent or warped inserts should be repaired or discarded.

- Worn tips prevent effective operation of the ultrasonic insert and are subject to breakage. Inserts should be replaced when tip-wear reduces size and clinical access from original condition. See wear guide.

- Water lines to scaler should be flushed and handpiece should be filled with water prior to use for most efficient results.

- The o-ring should be carefully maintained to ensure the longest possible life. Maintenance includes water lubrication of o-ring as well as proper seating of insert within the hand-piece.

- Damaged or broken (torn) o-rings should be replaced prior to use.

Warranty:

Premier® ultrasonic inserts are warranted for 90 days against manufacturer's defect.

Proper Placement of Replacement O-Rings

Proper placement of o-rings will significantly extend the life of the o-ring and provide better ultrasonic scaling. During installation, the o-ring should be protected from the sharp edges of the insert stack.

- Remove and discard old o-ring.

- Place new o-ring onto a straight drinking straw (1/4" diameter), approximately 1" from the end.

- Slide the straw completely around insert stack, starting with the side of the straw that holds the o-ring.

- Slide o-ring off of the straw and place into appropriate groove. Remove and discard straw.

- Sterilize insert before use.

Make sure you are using the supplied o-rings on your ultrasonic inserts as other o-rings may not be the appropriate size.

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