ES

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**premier®** ComfortView®

US/UK Lip and Cheek Retractor • Latex Free

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Manufactured for:  
Premier® Dental Products Company  
1710 Romano Drive  
Plymouth Meeting, PA 19462 U.S.A.  
888-670-6100 • 610-239-6000  
premierdentalco.com • info@premierdentalco.com

Made in Czech Republic.  
United States Patent No. 7,077,652

**Rx ONLY**

IFU# 90613 • 1220054 Rev8 NKP

### Instructions for Cleaning and Sterilizing

#### Cleaning:

1. Device must be disassembled for effective cleaning.
2. Inspect the device for debris. Remove visible debris or organic matter from the device with a disposable wipe or a soft brush.
3. Using a mild detergent and warm tap water, clean and rinse device.
4. Dry device.

#### Sterilizing:

The retractor and cushions are autoclavable. Sterilize in a steam autoclave using distilled water at temperature of 134°C /273°F.

- Re-assemble device after sterilization.
- Cold sterilization is not recommended.

### Instructions for Use

Figure 1

Figure 2

Figure 3

*Helpful Hint: Apply non-petroleum jelly or lip balm to patient's lips prior to inserting lip and cheek retractor.*

1. With the pink cushions facing the patient, contract the vertical lip holders to ease placement (Figure 1).
2. Place one side of the lateral lip holder against the inside of the patient's cheek (Figure 2).
3. Verify that the pink cushions are resting inside the patient's anterior vestibular region. The upper cushions should be resting adjacent to the frenulum. Lift patient's lips, if needed.
4. Once securely in place, gently guide the other side of the retractor in the other cheek by slightly bending it.
5. Verify that the patient's lips are comfortably in place (Figure 3).

#### Reorder Information: ComfortView®

- |            |         |   |
|------------|---------|---|
| <b>REF</b> | 9061381 | Lip and Cheek Retractor with Mounted Cushions - UNIVERSAL (2 per box) |
| <b>REF</b> | 9061383 | Lip and Cheek Retractor with Mounted Cushions - SMALL (2 per box)     |
| <b>REF</b> | 9061382 | Refill Cushions (10 per pack)   |

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