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Brite Shield™
Instrument Cleaner and Protectant

Removes Protein (Blood, Saliva, Bone, Tissue, etc.), Light Stains, Set Cements (non-resin reinforced) and Protects Instruments from Rust and Corrosion

Enzymatic Cleaner - Environmentally friendly
Directions for Use:
1. For daily cleaning and rust protection, place one scoop of Brite Shield powder and 64 oz. (1894 ml) of warm water in ultrasonic unit.
• For heavy contaminated instruments, use two scoops of powder.
• For cement/stain removal, place two scoops of Brite Shield powder and 16 oz. (473ml) of warm water in ultrasonic unit.
2. Stir until dissolved.
3. Immerse instruments in solution and clean in ultrasonic unit for 1-5 minutes.
4. Rinse instruments thoroughly with water and dry before sterilizing.

Read SDS before using product. For SDS, visit www.premusa.com or call Premier at 888-670-6100 or 610-239-6000.

For optimal results, change solution daily.
• Should not be used on anodized aluminum instruments or bur blocks
• Also protects stainless steel instrument cassettes
• Not for internal use
• Ideal for removing protein and debris from diamond-cutting tools and burs
• Removes alginate from metal trays
• May be used without ultrasonic machine; soak 15-20 minutes

Premier® Dental Products Company
1710 Romano Drive, Plymouth Meeting, PA 19462 U.S.A.
(888) 670-6100 • (610) 239-6000 • www.premusa.com
Made in U.S.A.

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