


Ensure accuracy and save time



T-LOC™ Triple Tray®

Adhesive-Free Dual-Arch Impression Trays

- ✓ **Patented adhesive-free design** - Retention tracks mechanically lock in material, eliminating the need for volatile, unpleasant smelling and expensive adhesives. You save both time and money.
- ✓ **Easier handling** - Large, ergonomic handles simplify loading of the tray, enable easy transfer between assistant and dentist, and facilitate tray insertion and removal from the mouth.
- ✓ **Accurate bite** - Thin, loose webbing maximizes intercuspation of the teeth and exhibits "...least water sorption & thinnest occlusal wafer"¹ to ensure an accurate bite.
- ✓ **Greater rigidity in posterior** due to the use of a new, novel polymer.
- ✓ **Ensure accuracy and save time** - Low sides, thin posterior bar and wide impression areas avoid impingement when the patient closes. Premier® Triple Tray® is proven² accurate - as you simultaneously capture master impression, counter impression and bite registration.


U.S. Patent #6,749,428

T-LOC™ is the most significant advance in closed-mouth impression making since the development of the Premier® Triple Tray®

- 1006274 Three-Quarter Arch (22 trays/box)
- 1006204 Posterior (35 trays/box)
- 1006214 Anterior (28 trays/box)



dental
townie
choice
8 Time Winner



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



premier® Inspired solutions
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • www.premusa.com

Proven Accurate in Over 100 Million Impressions

Triple Tray®

Dual-Arch Impression Trays

✓ Award-winning for over a decade

Premier® Triple Tray® is the professionals' most trusted name for accurate dual-arch impressions. The thin, adaptable mesh supports your impression material, ensures complete intercuspation, resists tearing and is proven to have "...least water sorption"¹, reducing possibility of distortion.

- 1006283 Wide Body Posterior (Box of 48)
- 1006243 Sideless Posterior (Box of 48)
- 1006203 Posterior (Box of 48)
- 1006253 Extended Quadrant (Box of 35)
- 1006213 Anterior (Box of 35)
- 1006273 Three-Quarter Arch (Box of 35)
- 1006223 Quadrant (Box of 40)
- 1006233 Assorted (5 each of the styles above)



Townie Choice
14 years and
counting!



"Can't live
without!"

Alfa™ Triple Tray®

Full-Arch Alginate Impressions

✓ Revolutionizing full-arch alginate impression making


Alfa™ Triple Tray® improves patient comfort while helping to prevent gagging. Alfa cuts chair time in half and is easy to use for most full-arch alginate impressions. Use for custom trays, whitening trays, night guards, athletic guards, study models, and as a matrix for temporaries.

View online tutorial at www.premusa.com

Alfa™ Triple Tray® - Box of 24

- 1006313 Small (purple)
- 1006323 Medium (white)
- 1006333 Large (green)
- 1006343 Assorted (5 small, 14 medium, 5 large)

U.S. Patent #6,749,428



premier® Inspired solutions
for daily dentistry®



premusa.com/promise

Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA

Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • www.premusa.com • E-mail: dentalinfo@premusa.com

© Premier and Triple Tray are registered trademarks and Alfa and T-LOC are trademarks of Premier Dental Products Company. / ©2005 Premier Dental Products Company. / 1. Independent Research Organization, July 2000. / 2. Complete Triple Tray® Bibliography available upon request. / 3. 2003, 2004, 2006, 2007, 2008, 2009, 2010, 2011, 2012, 2013, 2014, 2015, 2016, 2017 Townie Choice Award, DentalTown Magazine. / 4. The Dental Advisor, November 2002. / 5. Independent Research Organization, July 1997 and July 2004.

8698741 0418067 Rev6