

Unique Knitted Microfibers with Superior Absorbency!



Aluminum Chloride (AlCl₃)
0.5 ± 0.1mg/inch

Knit-Pak™+ Impregnated Knitted Cord

Knit-Pak™+ is an impregnated retraction cord engineered from unique knitted microfibers to provide unsurpassed performance. The innovative fibers are distinguished by providing high durability and outstanding retention properties.

These attributes create a retraction cord with superior absorbency and a fray-resistant, pliable structure for efficient procedures and great patient outcomes.

- Offers exact-dosage medicinal loading (AlCl₃· 6H₂O 0.5mg± 0.1mg/inch)
- Easily packs into sulcus
- Cuts clean – will not entangle in burs
- Lint and particle-free
- 100 inches (254cm) per container



A ceramic cutting blade inside the cap automatically cuts the cord upon closing firmly.




The cord piece is held by the cap until removed, preventing it from dropping. Protects cord piece and container from cross-contamination.

Ordering Information

9007560	000.0	Dark Green	9007563	0	Purple
9007561	000	Green	9007564	1	Blue
9007562	00	Brown	9007565	2	Orange

Knit-Pak™

Knitted Retraction Cord



Knit-Pak™ is ideal to use in combination with Hemodent® hemostatic agent or Traxodent® for tissue displacement and hemostasis prior to taking an impression.

- Carries greater quantities of hemostatic solution than twisted and braided cord
- Easily packs into sulcus
- Will not fray and resists popping out of the sulcus
- Cuts clean – will not entangle in burs
- 100% cotton – color coded for easy identification
- 100 in (254cm) per container



Ordering Information

9007551	000	Green	9007554	1	Blue
9007552	00	Brown	9007555	2	Orange
9007553	0	Purple	9007556	3	Yellow

Premier® Retraction Caps

Cotton caps work with Traxodent to quicken retraction and hemostasis.

- Easily guides Traxodent into the sulcus when patient closes
- Compression assists and quickens hemostasis
- Anatomic form allows easy placement

- 9048251 Size 1 Small (120)
- 9048253 Size 3 Medium (120)
- 9048255 Size 5 Large (60)
- 9048260 Assorted - (30) #5, (15) #1, (15) #3



Hemodent®

Buffered Aluminum Chloride Hemostatic Solution

- Effectively stops minor gingival bleeding
- Epinephrine-free to avoid cardiac reactions
- Stable and offers a long shelf-life

- 9007071 10cc bottle
- 9007072 20cc bottle
- 9007073 40cc bottle



Premier® Cord Packers

Premier offers a selection of cord packing instruments –

- Balshi and Sherman packers, with or without serrations
- Round or taper ends in angle and straight configurations, with or without serrations



Please contact your dealer concerning product availability in your country.

premier® Inspired solutions for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
 Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • www.premusa.com • E-mail: dentalinfo@premusa.com

8610325 0818035 Rev5 NKP10M