



ComfortView®

Providing greater comfort while maximizing access and visibility



Small



Universal

ComfortView®

Lip & Cheek Retractor

- Retracts lips and cheeks for optimal buccal and gingival access
- Maximum visibility and patient comfort
- Ideal for both short and long procedures
- Anatomic design with cushioned lip-supports that position lips away from anterior teeth and gums
- Suitable for use in diagnostic, preventative, therapeutic dental applications, digital impression scanning, in-office tooth whitening and intra-oral photography
- Autoclavable, Latex free

9061381 Universal ComfortView Lip and Cheek Retractor - 2 per box (white)

9061383 Small ComfortView Lip and Cheek Retractor - 2 per box (blue)

9061382 ComfortView Refill Cushions -10 per pack

ComfortView not available from Premier in Europe.

- 1 Anatomically shaped lip holders for superior access to buccal and gingival areas – maximum access and visibility
- 2 Lip holders with opening for the frenulum for maximum comfort – enhanced patient comfort for short- and long-term procedures
- 3 Cushioned lip-supports position lips away from anterior teeth and gums – allow for visual and instrument access
- 4 Lateral flaps gently hold cheeks and contract for easy placement/removal – ensure patient comfort and assist in keeping patient's mouth open.
- 5 ComfortView frame material and cushions are autoclavable.