


Surprisingly practical.
Get a 2-fer with 2pro!



2pro® Total Access Disposable Prophyl Angle

Award-winning **2pro** is the versatile 2-in-1 disposable prophyl angle that helps you easily access hard-to-reach areas of the mouth quickly. One cup. One tip. No extra cost or time to switch. Which means a better cleaning experience, happier patients and less chair time! 2pro is your all-in-one polishing angle. Every time.

- ✓ Patented gearless, jam-free design - quiet, smooth polishing
- ✓ Angled neck for easier access - makes hard-to-clean areas a breeze
- ✓ Advanced 2-in-1 design - saves time, improves workflow
- ✓ Convenient access brush – efficiently cleans around ortho bracket



Top Prophyl Angle
2013, 2014 & 2015



premier® Inspired solutions
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • www.premusa.com

Get Total Access!

2pro cleans and polishes all tooth surfaces with its unique cup and tip design!

No added expense for tips or brushes!


Use the cup to polish all tooth surfaces. Use the tip for gingival and proximal areas, occlusal surfaces, veneer margins, around implants and orthodontic brackets.

Cup On



Facial surfaces


Cup Off



Veneer margins



Orthodontic brackets



Gingival/proximal



Soft/short
purple cup,
short tip

Firm/short
blue cup,
short tip

Soft/long
white cup,
long tip

Firm/long
green cup,
long tip

“ I absolutely love the Premier® 2pro® total access soft disposable prophyl angle. The cup shape reduces splatter and is easy to use. Being able to remove the cup and use the rubber point is great. This allows access around the orthodontic brackets and wires. ”

- Sandra B., RDH, Pennsylvania

	Bag of 144	Bulk Package (1,000 angles)
Soft/Short Cup (purple cup)	5500101	5510001
Firm/Short Cup (blue cup)	5500102	5510002
Soft/Long Cup (white cup)	5500103	5510003
Firm/Long Cup (green cup)	5500104	5510004

Order through your authorized dealer.

premier® Inspired solutions
for daily dentistry®



premusa.com/promise

Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • www.premusa.com • E-mail: dentalinfo@premusa.com

U.S. Patent # 7,104,796, U.S. Patent # 7,153,133.



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



8625555 0318062 Rev6 NKP5M