


# Split Personality

Cool.



Voted top prophylaxis angle...*really cool!*

**2pro**  
Total Access  
Disposable Prophylaxis Angle



Award-winning 2pro® provides total access to all tooth surfaces for **all patients** by combining a prophylaxis cup and tip in one angle! Improve your efficiency and reduce cost by using 2pro!

See how your office can reduce its cost!  
**Take the 2pro Challenge and  
Discover Your Savings!**

[www.premusa.com/try2protoday](http://www.premusa.com/try2protoday)



Patented gearless design



100° angle



Disposable Prophy Angle

## Get Total Access!

2pro cleans and polishes all tooth surfaces with its unique cup and tip design!

### No added expense for tips or brushes!

Use the cup to polish all tooth surfaces. Use the tip for gingival and proximal areas, occlusal surfaces, veneer margins, around implants and orthodontic brackets

### Saves time and money

Comparably priced to other angles, yet when needed, you save time and money by removing the cup to clean remaining areas with 2pro's gentle tip

### Provides total access

2pro's unique 2-in-1 design guarantees total access to all tooth surfaces

### Runs smoothly and quietly

Patented gearless design ensures quiet, smooth-running, vibration-free operation with less heat build-up

### Easier access to lingual areas

Ergonomic 100° design offers comfort and convenience by reaching all areas of the mouth

### Latex free

**QUICK TIP!**



Use Around Orthodontic Brackets!



**Soft/Short Purple Cup**  
Short Tip



**Firm/Short Blue Cup**  
Short Tip



**Soft/Long White Cup**  
Long Tip



**Firm/Long Green Cup**  
Long Tip

**Bag of 144**

REF 5500101

REF 5500102

REF 5500103

REF 5500104

**Bulk Package (1,000 angles)**

REF 5510001

REF 5510002

REF 5510003

REF 5510004

Order through your authorized dealer.



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, U.S.A.

Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • www.premusa.com • E-mail: dentalinfo@premusa.com