

# Innovate with Premier® Dental

*The name you know, the quality you expect*



---

Endodontic Products & Techniques

## RC-Prep®

Root Canal Preparation Cream for Chemo-Mechanical Preparation of Root Canals

#1  
Selling

Root Canal Preparation Cream in America!<sup>1</sup>



### Direct Application Into the Canal

RC-Prep's unique formulation of glycol, urea peroxide and EDTA in a special water-soluble base helps to remove calcifications and lubricates the canal to permit more efficient instrumentation. RC-Prep is also excellent for use with apex locators, as it permits consistently reliable readings.

- 9007132 18g pump (case of 10)<sup>†</sup>
- 9007133 227g jar
- 9007131 18g jar (case of 10)<sup>†</sup>
- 9007135 2 - 9g filled syringes (case of 10)<sup>†</sup>
- 9007129 Syringe kit (5) 3cc syringes and 50 tips
- 9007127 Plastic tip refill (50) (for use with 9g syringe)
- 9007139 Metal Syringe tip refill (for use with 3cc syringe)
- 9007138 RC-Prep MicroDose intro kit (100) pipettes and (2) autoclavable syringes
- 9007137 RC-Prep MicroDose refill (100) pipettes
- 4009099 MicroDose syringe pack (3) syringes

CE  
0470

<sup>†</sup> Individual units available through your dealer



### RC-Prep Technique



RC-Prep is easily introduced into the canal directly via the MicroDose Pipette or a reamer, file or root canal filler. As RC-Prep cleans and lubricates, instruments cut freely without binding.



When RC-Prep reacts with sodium hypochlorite solution, pulp debris is actually lifted out of the canal. The pulp chamber is brightened, helping to locate additional canals.



With RC-Prep and proper instrumentation, the canal is thoroughly debrided. Dentin tubules and lateral canals are opened and cleaned, essential for successful treatment.

1. SDM 2016.



www.premusa.com

# Endodontic Instruments




**Glick 1** • 1003850



**Stewart Probe** • 1003893



**Explorer 16** • 1003844



**Plugger 5/7** • 1003863



**Excavator 33L** • 1003885


*technique tip*



The **Stewart Multi-Purpose Probe** has rigid tips to break through pulp stones.

*For full line of instruments, see Premier's catalog*

## Spreaders (Stainless Steel & Nickel Titanium- Niti)




**D11** • 1003879

**Ni-Ti D11** • 1003905



**D11T** • 1003872

**Ni-Ti D11T** • 1003906



**25S** • 1003869

**Ni-Ti 25S** • 1003908

*"Outstanding Product"<sup>2</sup>*  
Premier® Ni-Ti Spreaders



**Endodontic Plier** • 9065086




**Locking Plier (Grooved)** • 9065088



## Real World Endo™ Explorers

RWE, an endodontic educational company, offers a comprehensive course offered by Dr. Ken Koch and Dr. Dennis Brave for the general practitioners focusing on ways to perform endodontics faster and more predictably. For more information, visit [realworldendo.com](http://realworldendo.com)

*"Outstanding Product"<sup>2</sup>*  
Real World Endo™ Explorers



**Posterior** • 1003892


**Anterior** • 1003891

2. An independent non-profit dental education and product testing foundation, December 2004.

## Files & Reamers

### Stainless-Steel Endodontic Instruments

- Extremely Sharp - Efficient Instrumentation
- Pre-mounted Silicone (Non-Latex) Stop
- All Sizes - Same Low Price
- Firm Grasp - Peanut Shaped Handle
- ISO Symbols on Color-Coded Handle
- Extremely Flexible - Virtually Break-Proof
- Precise to AAE and ISO Standards
- Sterilizable
- Swiss-Made



## Premier® Safe-End Files Avoid Ledging

Standard feature on all K-Files  
& Hedstrom Files, sizes 06-40

### K-Files

21mm (6) per box													
9053200	06	9053203	15	9053206	30	9053209	45	9053212	60	9053215	90	9053219	45-80
9053201	08	9053204	20	9053207	35	9053210	50	9053213	70	9053216	100		
9053202	10	9053205	25	9053208	40	9053211	55	9053214	80	9053218	15-40		
25mm (6) per box													
9053230	06	9053233	15	9053236	30	9053239	45	9053242	60	9053245	90	9053248	15-40
9053231	08	9053234	20	9053237	35	9053240	50	9053243	70	9053246	100	9053249	45-80
9053232	10	9053235	25	9053238	40	9053241	55	9053244	80	9053247	110		
31mm (6) per box													
9053290	06	9053294	20	9053298	40	9053302	60	9053306	100	9053310	140		
9053291	08	9053295	25	9053299	45	9053303	70	9053307	110	9053311	15-40		
9053292	10	9053296	30	9053300	50	9053304	80	9053308	120	9053312	45-80		
9053293	15	9053297	35	9053301	55	9053305	90	9053309	130				

### Hedstrom Files

21mm (6) per box													
9053399	08	9053402	20	9053405	35	9053408	50	9053411	70	9053416	15-40		
9053400	10	9053403	25	9053406	40	9053409	55	9053412	80	9053417	45-80		
9053401	15	9053404	30	9053407	45	9053410	60	9053413	90				
25mm (6) per box													
9053424	08	9053427	20	9053430	35	9053433	50	9053436	70	9053441	15-40		
9053425	10	9053428	25	9053431	40	9053434	55	9053437	80	9053442	45-80		
9053426	15	9053429	30	9053432	45	9053435	60	9053438	90				
31mm (6) per box													
9053876	08	9053879	20	9053882	35	9053885	50	9053889	80	9053894	45-80		
9053877	10	9053880	25	9053883	40	9053886	55	9053890	90				
9053878	15	9053881	30	9053884	45	9053887	60	9053893	15-40				

### Reamers

21mm (6) per box											
9053000	06	9053003	15	9053006	30	9053009	45				
9053001	08	9053004	20	9053007	35	99053018	15-40				
9053002	10	9053005	25	9053008	40						
25mm (6) per box											
9053030	06	9053034	20	9053038	40	9053042	60				
9053031	08	9053035	25	9053039	45	9053045	90				
9053032	10	9053036	30	9053040	50	9053048	15-40				
9053033	15	9053037	35	9053041	55	9053049	45-80				
31mm (6) per box											
9053090	06	9053094	20	9053098	40	9053111	15-40				
9053091	08	9053095	25	9053099	45	9053112	45-80				
9053092	10	9053096	30	9053101	55						
9053093	15	9053097	35	9053102	60						

#### Rubber Stop Color Coding for Files & Reamers

Red = 21mm  
Blue = 25mm  
Black = 31mm

Features a Safe-End Tip  
on sizes 06-40



## Sensipast®

### Root Fillers

Premier Sensipast root fillers have a miniature clutch in the mandrel that allows them to be used on any kind of right-angle handpiece without breaking.

Package of 4  
21mm

9053571 25  
9053572 30

25mm

9053576 25



## Root Fillers

### Carrier for Root Canal Sealer

Package of 4  
21mm

9053556 25  
9053557 30  
9053558 35  
9053559 40  
9053560 Assorted (25 - 40)

25mm

9053561 25  
9053562 30  
9053563 35  
9053564 40  
9053565 Assorted (25 - 40)



## Gates-Glidden Drills

### Shape & Taper the Canal

Package of 6

9055001 No. 1	9055005 No. 5
9055002 No. 2	9055006 No. 6
9055003 No. 3	9055009 Assorted
9055004 No. 4	


## Peeso Drills

### Prepare Post Space in Endodontically Treated Teeth

To protect the integrity of the tooth and root-canal filling, gutta percha should be removed with hot pluggers to the desired depth before use of Peeso drills.

Package of 6

9055021 No. 1  
9055022 No. 2  
9055023 No. 3  
9055024 No. 4  
9055025 No. 5  
9055026 No. 6  
9055029 Asst.



## Barbed Broaches

### Flexible Stainless Steel Avoids Breakage

Package of 144

9053631 25  
9053632 30  
9053633 35  
9053634 40  
9053635 50  
9053636 60  
9053637 Assorted (25-60)

Individual 12 packets available  
through your dealer.

## Finger Spreaders

### For Lateral Condensation

Premier Finger Spreaders are precisely sized to provide excellent lateral condensation of gutta percha toward canal walls.

Package of 6  
21mm

9053814 15  
9053816 20  
9053818 25  
9053820 30  
9053822 35  
9053824 40  
9053826 Assorted (15-40)

25mm

9053830 15  
9053832 20  
9053834 25  
9053836 30  
9053838 35  
9053840 40  
9053842 Assorted (15-40)

## Finger Pluggers

### For Vertical Condensation

Package of 6

21mm  
9053771 15  
9053772 20  
9053773 25  
9053774 30  
9053788 Assorted (15-40)

25mm

9053791 15  
9053792 20  
9053793 25  
9053794 30  
9053808 Assorted (15-40)

## Gutta Percha Points


### Machine Rolled / Color Coded

Box of 120

9056001 15  
9056002 20  
9056003 25  
9056004 30  
9056005 35  
9056006 40  
9056007 45  
9056008 50  
9056009 55  
9056010 60  
9056011 70  
9056019 25-40  
9056020 45-80

Box of 60

9056021 90-140



## Gutta Percha Accessory Points

Box of 120

9056101 Extra Fine  
9056102 Fine/Fine  
9056103 Med./Fine  
9056104 Fine  
9056105 Fine/Med.  
9056106 Medium  
9056116 Asst.

## Absorbent Points

### Sterile - 40 cells of 5



Box of 200


9055101 XX-Fine  
9055102 X-Fine  
9055103 Fine  
9055104 Med  
9055105 Coarse  
9055106 X-Coarse  
9055109 Asst.



## Cure-Thru® IntegraPost®

Light-Transmitting Zirconia  
Glass Fiber Post

A simple technique for post and core utilizing Zirconia glass fibers embedded in composite for the right balance of strength and elasticity.



Length 19mm

Safe and conservative parallel-sided passive post.

Uses IntegraDrills listed below.



3001600 Cure-Thru IntegraPost Intro Kit  
Contains: 20 posts, (5) each #4, #4.5, #5, #6 and (4) corresponding drills

3001601 Cure-Thru IntegraPost Mini Intro Kit  
Contains: 10 posts, (5) each #4, #5 and (2) corresponding drills

### Refills (10 pack)

3001610 Size 4 (red) 0.96mm

3001611 Size 4.5 (purple) 1.08mm

3001612 Size 5 (blue) 1.21mm

3001613 Size 6 (green) 1.46mm

## IntegraPost®

Ultra-passive Integrating Core Retainer

IntegraPost's rounded head eliminates stress points and the Flow-Thru-Head eliminates voids, ensuring complete post and core integration.

### Refills - Package of 10

3001105 Size 3

3001110 Size 4

3001115 Size 4.5

3001125 Size 5

3001130 Size 6



3001001 Mini Intro Kit

Contains: 6 posts (2) each #3, #4, #4.5, (3) drills - (1) each #3, #4, #4.5, IP Carrier

### IntegraDrills (3 pack)

3001205 Size 3 (yellow)

3001210 Size 4 (red)

3001215 Size 4.5 (purple)

3001225 Size 5 (blue)

3001230 Size 6 (green)

### IntegraDrills (5 pack)

3001235 Assorted



## Two Striper® TSZtech™

Diamond Technology for Today's Advanced  
Materials - Lithium Disilicate and Zirconia

Assorted 2015895  
(Includes 260.8Z, 285.5Z,  
703.8KRZ, 125Z, 770.8Z)



### Reorder (5-pack)

260.8Z 2015893

285.5Z 2015891

703.8KRZ 2015894

125Z 2015892

770.8Z 2015890



### "L" - 1.5mm longer shank (5-pack)

L260.8Z 2015883

L125Z 2015882

L770.8Z 2015880



# #1

Top Performing Diamond  
to adjust, remove or access  
all-ceramic restorations!

## Endodontic Access



125Z Round

## Endo/Access


Slow Speed Diamonds

### Reorder (1-pack)

1C 2013194 Round

2C 2013195 Round

291.4C 2013197 Pointed Football



## Perfecta® REV!®

14% Hydrogen Peroxide

Non-vital whitening for  
endodontically treated teeth

4000142 REV! Six Patient Paks

Available individually through dealer.

**Each Patient Pak is a complete 2-week treatment.**

Contains: 2 - 3cc syringes with 4oz (118.3ml) REV! finishing rinse and whitening toothpaste. You add trays and tray cases.

4000141 REV! Six Refresher Paks

Available individually through dealer. **Ideal for touch-ups.**

Contains: 1 - 3cc syringe with 4oz (118.3ml) REV! finishing rinse. You add trays and tray cases.



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA

Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • www.premusa.com • dentalinfo@premusa.com

®Two Striper is a registered trademark and TSZtech is a trademark of the manufacturer, Abrasive Technology, Inc. / U.S. Patent 5,511,718