

Make a colorful impression

Rigid Tray & Tear-Resistant Gold Mesh™ Design



Neotray™
Impression Trays

Gold Mesh™ technology

Tear-resistant material

Exclusive KR™ plastic

Provides a rigid tray without the cost of metal disposable trays or cleaning and care of single-arch trays

Cuts chair-time in half

Simultaneously record master impressions, counter impression and bite registration

Saves money

Less chair time, impression material and corrections


Increases accuracy

Avoid costly errors and chairside adjustments

**Assortment of seven sizes available
in fun colors that patients love.**

NeoTray™

Impression Trays




- #2206NT** Assortment
Contains: 10 each of Wide Posterior, Sideless Posterior, 5 each of Anterior, Posterior, Full Arch, Quadrant and Extended Quadrant
- #2204NT** Anterior - Box of 30
- #2205NT** Full Arch - Box of 30
- #2202NT** Sideless Posterior - Box of 48
- #2208NT** Wide Posterior - Box of 48
- #2201NT** Posterior - Box of 48
- #2207NT** Extended Quadrant - Box of 35
- #2203NT** Quadrant - Box of 40

Impressions are as easy as....



1.
Try-in the appropriate tray to ensure a correct fit and proper closure.



2.
Load both sides of tray with firm-setting impression material.




3.
Syringe material onto critical areas and place loaded tray in mouth; have patient close correctly.



4.
Inspect final impression for good detail and no impinging on tray walls.
Send disinfected impression to the lab.



Premier Dental Products Co. • Premier Dental (Canada)
1710 Romano Drive, Plymouth Meeting, PA 19462, USA • Phone 888-670-6100
610-239-6000 • Fax 610-239-6171 • www.premusa.com • E-mail: dentalinfo@premusa.com



MDSS
Burckhardtstr. 1
30163 Hannover, Germany

