

EXCELLENT  
VISIBILITY  
& PATIENT  
COMFORT

## UNC 12

1mm increments up to 12mm

3-6-9-12mm


Original 3-5-7-10mm

# PerioWise®

## "The Friendly Probe®"

*PerioWise® periodontal probes are extremely easy to read and gentle to the tissue.<sup>1</sup> Slightly flexible, with a rounded tip, PerioWise probes offer greater patient comfort and are safe around implants. They can be autoclaved and used repeatedly.*

The vibrant color lines against the white probe provide excellent contrast. At first glance, you can determine the depth of the pocket quickly and accurately.



*The vibrant colors on PerioWise offer excellent visibility.*



*Metal probes can be difficult to read.*

PerioWise, multi-colored periodontal probes are available with (3-5-7-10) and (3-6-9-12) markings. A green mark from 0-3mm indicates healthy tissue; readings in the probe's red markings indicate the need to perform a comprehensive periodontal examination. The PerioWise multi-colored probes are ideal for quick periodontal screening or complete examinations.

PerioWise UNC 12 provides the same comfort and ease of use as our original multi-colored PerioWise probes in the popular UNC 12 style. The intense purple markings provide excellent contrast, indicating 1mm increments up to 12mm. Recommended for complete periodontal exams.



# PerioWise® "The Friendly Probe®"

**Plastic, autoclavable periodontal probes for excellent readability, patient comfort and acceptance.**

- Slightly flexible for greater patient comfort
- Clearly visible when using intra-oral camera

**Fast and Accurate Screening System<sup>2</sup>**

The patient Screening Report helps with patient education and compliance.



**FREE**

**Screening System Components**

- Screening Records for patient chart
- Screening Reports given to patient
- Chairside Patient Education Guide

**Safe For Use Around Implants**  
Restored implants can be difficult to identify. PerioWise is safe; it will NOT scratch implants.


**Green - Go**      **Red - Stop**

Easy-to-read 3mm green band indicates periodontal health.

Reaching a red mark indicates the possible presence of periodontal disease.

For your **FREE** Screening System Components  
Call: 1-888-670-6100  
or E-mail: dentalinfo@premusa.com



**Repeatedly Autoclaved and Reused**  
After dozens of uses, the non-toxic, FDA-approved coloring on the tip of PerioWise gradually wears away. That's when it's time to reorder.



To learn more about the PerioWise Screening System, visit [www.periowise.com](http://www.periowise.com)

## Premier® Implant Scalers Engineered for reliable implant maintenance

Fiber reinforced graphite that will not scratch implants. The thin and sharp well-defined tips allow access into tight areas for effective and safe scaling. The handle has been redesigned for greater comfort and control. Reusable, autoclave at (273°F/134°C).



Premier® Dental Products Company • 888-670-6100 • 610-239-6053 • Fax 610-239-6171 • Email dentalinfo@premusa.com • 1710 Romano Drive, Plymouth Meeting, PA 19462 • www.premusa.com

1. "Improved Diagnostics: Clinical Evaluation of a Color-Coded, Polymeric Periodontal Probe." CLINICAL PREV. DENT.; Vol. 14, No. 4; July/August 1992. / 2. "A Patient-Centered Approach to Periodontal Disease Detection." JOURNAL OF PRACTICAL HYGIENE; Vol. 8, Number 3; May/June 1999. / © Premier is a registered trademark of Premier Dental Products Company/® PerioWise and The Friendly Probe are registered trademarks of PerioWise Inc.

8690608 111343 Rev3 PG10M