

Spot-On™ 37% Phosphoric Etch Gel

Spot-On™ Etch is a 37% high-viscosity phosphoric etch gel for all your routine and selective-etch procedures. It's designed with enhanced viscosity that will prevent it from dripping or slumping. The deep blue color provides excellent contrast to enamel and dentin. Spot-On Etch easily flows through a small applicator tip, washes off readily and resists drying out.

Key benefits

- · Ideal for routine etching
- Selective-etch procedures
- Blue color provides a strong contrast to enamel and dentin
- · Easy handling
- Rinses clean without leaving residue

Spot-On Etch no-slump viscosity will stay put. Only etch the areas you desire!








Selective-Etch






Selective-etch enamel before...

- Using self-etch adhesives
- Etching dentin; reducing the chance of over etching and hypersensitivity
- Composite or porcelain repairs
- · Cementing ortho brackets



SEM untreated Incisor – 2000X



Spot-On™ Etch treated Incisor- 2000X

Distinctive change in the micromorphological appearance of enamel

Ordering Information: Spot-On™ High viscosity 37% phosphoric acid dentin/enamel etch gel

3001461 Spot-On™ Etch Intro
Five 1.2ml filled syringes
Twenty bent disposable applicator tips

3001460 Spot-On™ Etch Bulk
One 12g jumbo filled syringe
with luer lock adapter
Ten 1.2ml empty syringes
Fifty 1.2ml disposable applicator tips

