

#1

Top performing diamond to adjust, remove or access all-ceramic restorations!

Recommended for
Zirconia and
Lithium Disilicate

TSZTECH™

The P.B.S.® Brazed Advantage

Faster, Cooler Cutting & Longer Lasting!



Highest rated all-ceramic cutting for today's new generation of tougher ceramic materials.

Those materials include but are not limited to IPS e.max®, BruxZir®, BruxZir® Now, Lava™ Plus, Cercon®, Zirlux® FC, Obsidian™, IPS Empress, Procera® and Lava™ Ultimate.

TSZtech Contraindications – Will not cut or separate metal.

Use light touch and copious spray: 160,000 rpms recommended

2015895 **TSZtech™ Diamond Assorted 5-pack**
(260.8Z, 285.5Z, 703.8KRZ, 125Z, 770.8Z)



5-pack Reorder Item Numbers:



260.8Z
2015893
Minor/Major
Diameter
0.3/1.3



285.5Z
2015891
Minor/Major
Diameter
2.2



703.8KRZ
2015894
Minor/Major
Diameter
1.1/1.8



125Z
2015892
Minor/Major
Diameter
1.7



770.8Z
2015890
Minor/Major
Diameter
1.2/1.7

"L" - 1.5mm Longer Shank
Available in 5-packs:



L260.8Z
2015883
Minor/Major
Diameter
0.3/1.3



L125Z
2015882
Minor/Major
Diameter
1.7



L770.8Z
2015880
Minor/Major
Diameter
1.2/1.7

Laboratory Diamonds
Available one per pack:



X590.12Z
2012847
Minor/Major
Diameter
1.2



X769.10Z
2012845
Minor/Major
Diameter
0.7/1.8

Intra-oral Adjustments



285.5Z **Football**

Atraumatic Crown and Bridge Removal



770.8Z **Round-End Taper** Section and Remove

Endodontic Access



125Z **Round**



Learn More About Crown Removal
Scan the QR Code or search "TSZtech" on YouTube.



FASTER!

Super-Charged Diamond Polishing Single-paste System for Ceramics and Composites


Recommended for
Zirconia and
Lithium Disilicate

Diamond Twist SCL™

Extra-Oral Polishing Kit

Glaze-like, multi-surface polishing - featuring advance open-weave fibra pre-polishers and 100% wool-paste applicators - for a super-charged mirror polish.

- Pre-polish with Fibra Points (barrel, wheel or taper); open-weave design generates significantly less heat.
- Final polish with Diamond Twist Super-Charged Polishing Paste applied with felt barrel, wheel or point for a glaze-like polish.



2019000
Diamond Twist SCL™ Kit

Products below shown 75% of actual size.

2019001
6g Diamond Twist
SCL Polishing Paste
(unflavored)



2019032
Fibra Wheel,
Thin
(6/pack)



2019021
Felt Wheel,
Thin
(10/pack)



2019033
Fibra Wheel,
Thick
(4/pack)



2019020
Felt Wheel,
Thick
(5/pack)



2019031
Fibra Point,
Taper



2019022
Felt
Point



2019030
Fibra
Point,
Barrel



2019023
Felt
Barrel



2019051
SurGrip
Mandrel



2019050
Standard
Mandrel



Diamond Twist SCO™

Intra-Oral Polishing Kit

After permanent seating and adjustments, polish crowns, inlays, onlays and nano-hybrid composites intra-orally in one easy step using Diamond Twist Super-Charged Polishing Paste.

- Use the flat or taper polishing brushes (autoclavable) for posterior occlusal surfaces.
- Use the Micro-Cloth disks for polishing and final high-gloss buffing.
- Twist syringe avoids cross-contamination and dispenses just the right amount.



2019002
Diamond Twist SCO™ Kit

Products below shown 75% of actual size.

2019003
3g Diamond Twist SCO
Polishing Paste
(mint flavored)



2019005
Tapered Polishing
Brushes (25/pack)



2019004
Flat Polishing
Brushes
(25/pack)



2019091
FG Mandrels (4/pack)
not included in kit



2019090
RA Mandrels
(4/pack)



2019060
1/2" Micro-Cloth
Disks
Dark Blue (25/pack)



2019061
5/16" Micro-Cloth
Disks Light Blue
(25/pack)



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, U.S.A.

Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • www.premusa.com • E-mail: dentalinfo@premusa.com

© Premier is a registered trademark of Premier Dental Products Company. / ®Two Stripper and P.B.S. are registered trademarks and TSZtech is a trademark of the product manufacturer, Abrasive Technology, Inc.
® IPS e.max, BruxZir, Cercon, ZirLux FC, IPS Empress, Procera are not registered trademarks and ™ Lava Plus, Lava Ultimate and Obsidia are not trademarks of Premier Dental Products Company.
"Intra-Oral Adjustment" image courtesy of Shalom Mehler D.M.D. "Traumatic Crown Removal" image courtesy of Abdi Sameni D.D.S.