

**NEW!**

# Premier® Instrument Cassettes


**Premier stainless steel instrument cassettes offer proven reliability and quality performance with an unparalleled satisfaction guarantee.**

Unique, hingeless design offers years of trouble-free operation.

Compact, double-stack feature accommodates more instruments for greater efficiency in ultrasonic cleaning, autoclaving and storing.

Apertures in side walls, lid and bottom of trays permit optimum flow of cleaning solution, ultrasonic energy and greater access for rinsing.

5-year satisfaction guarantee.



**9010904**  
Mark I (4 or 8 Instruments)  
2.5" x 8.5" x 1.5"  
(6.35cm x 21.59cm x 3.81cm)



**9010908**  
Mark VI (8 or 16 Instruments)  
4.5" x 8.0" x 1.25"  
(11.43cm x 20.32cm x 3.18cm)



**9010924**  
Mark IV (12 or 24 Instruments)  
8.0" x 11" x 1.50"  
(20.32cm x 27.94cm x 3.81cm)



**9010901**  
Accessory Cassette  
1.5" x 3.5" x 1.0"  
(3.81cm x 8.89cm x 2.54cm)

## Premier Cassettes Are Easy To Use

**1)** Rinse cassette /bur block organizer under running water. **2)** Place cassette /bur block organizer in ultrasonic unit; use Brite Shield™ to remove bio-burden and to clean & protect instruments. **3)** Rinse cassette/bur block organizer under running water. **4)** Wrap cassettes and sterilize.

**NEW!**

**Premier Hi-Q™  
Chrome Bur Block**

*Metal Bur Block for diamond  
and carbide burs*

The Hi-Q bur block features an advanced material design that out-performs standard aluminum blocks. The bright chrome finish maintains a clean appearance through repeated sterilizations, and will not fade or discolor like anodized blocks.

Suitable for all methods of sterilization, including cold soak, and compatible with ultrasonic cleaners, the Hi-Q Chrome bur block organizes and presents multi-use dental rotary instruments in a professional and functional manner.

*This chrome-plated, anodized aluminum bur block measures 55mm x 25mm x 27mm, and is available in two configurations: 15 hole FG and 10 FG & 5 RA.*



**9051000**

Hi-Q Chrome Bur Block 15 Hole - FG

**9051001**

Hi-Q Chrome Bur Block 15 Hole 10FG/5RA

**Premier Brite Shield™**

*The most versatile and economical enzymatic cleaner on the market*

**Cleans up to 35% faster** - Powerful protease enzymes rapidly break down proteins and bio-burden for greater cleaning efficiency.

**Polishes and Protects** - A protective shield coats the instrument through a unique passivation process to restrict rust and corrosion while imparting a lustrous surface to stainless steel instruments.

**Costs Less** - This economically priced, super-concentration provides up to 34 gallons at a significant savings.

**User Friendly** - Brite Shield™ is non-hazardous, biodegradable and has no offensive odors.

**Versatile** - Brite Shield will remove set glass ionomer, zinc phosphate and polycarboxylate cement and is also useful for removing caked-on stone and alginate in impression trays.

*"...89% rated it excellent or good and worthy of trial by colleagues."*



**9011100**


**Brite Shield** case of six, 800gm containers with scoop dispenser.

*Available individually through your dealer.*

**Premier Big Easy™**

*Stainless Steel Scalers and Curettes with thicker, softer grips for better control*

**Sterilizable - Will Not Crack!**



*Under the medical-grade silicone grips, Big Easy has a light-weight stainless steel hollow handle for improved tactile sensitivity and outstanding durability.*

*Rated best of products tested in 1999.<sup>2</sup>*

Thicker, softer, sculpted medical-grade silicone grips cushion your fingers.

Cushion grips are resilient, providing passive feedback so you use less pressure.

Thin, sharp blades and precise angles offer optimal efficiency.

*Big Easy handles available on all popular designs.*



Premier Dental Products Co. Premier Dental (Canada) 1710 Romano Drive, Plymouth Meeting, PA 19462, USA

888-670-6100 610-239-6000 Fax 610-239-6171 www.premusa.com E-mail: dentalinfo@premusa.com

