

Introducing!



Sleeve-It™

Disposable Barrier Sleeves

Help Prevent Cross-Contamination

Premier's Sleeve-It™ disposable barrier sleeves provide proven protection against cross-contamination. Made from high-grade durable material, Sleeve-It barriers have passed extensive testing¹ to ensure patient safety.

Easy to use, Sleeve-It barriers slide over syringe dispensers and feature a pre-cut opening to pass through the dispensing tip.

Ideal for multiple-use dispensers that cannot be reprocessed using heat-sterilized or high-level disinfectants, Sleeve-It barriers can be conveniently used for each new patient per FDA, CDC and OSHA recommendations.^{2*}



Sleeve-It™

Disposable Barrier Sleeves

- ✓ Protection against cross-contamination
- ✓ Passed extensive testing¹ to ensure patient safety
- ✓ Easy to use



100 SERIES

9011501 100 series
(Fits 0.7 -3.0ml syringe)
Contains: 300 barrier sleeves



200 SERIES

9011503 200 series
(Fits 5ml dual-barrel syringe)
Contains: 300 barrier sleeves



Vac Attak™ GREEN

Safe, Effective and EPA-Compliant Evacuation System Cleaner

- ✓ Powerful enzymatic, EPA-compliant cleaner
- ✓ Compatible with amalgam separators
 - pH neutral
 - No chlorine
 - No oxidizing agents
 - Non-foaming
 - Non-corrosive
- ✓ Up to 80 cleanings per jar
- ✓ Economical and easy-to-use powder
- ✓ Fresh citrus scent

9011105 Six, 800g jars with scoop dispenser



Brite Shield™

Instrument Cleaner and Protectant

- ✓ Static soak and ultrasonic cleaner
- ✓ Removes glass ionomer, carboxylate, zinc phosphate cements and bio-burden
- ✓ Unique passivation process coats instruments restricting ruts and corrosion
- ✓ Removes alginate and stone from impression trays
- ✓ Economical and easy to use.
One jar makes up to 34 gallons

9011100 Six, 800g jars with scoop dispenser



premier® Inspired solutions
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com

¹ Viral penetration: ASTM F 1671
Synthetic blood penetration: ASTM F 1670
Resistance to puncture: ASTM F 1342
Tear resistance: ASTM D 1004
Tensile testing: ASTM D 882
Biocompatibility: ISO 10993-5, ISO 10993-10, ISO 10993-12

2. Data on file.
* IPAC Checklist for Dental Practices-Core Elements, 1st Revision: July 2019
follows FDA recommendation regarding use of a new barrier before each
patient on multiple-use dental dispenser devices that cannot be sterilized.

©2020 Premier Dental Co. | All Rights Reserved

8610427 0120064 Rev2 NKP10M