

Introducing!



## Premier® Disposable Face Shields

Protect your team!

In this 'new normal,' Premier is pleased to offer premium disposable face shields to help protect your face from aerosols, sprays and splatters.

Optically-clear with full peripheral vision, each face shield has an effective anti-fog coating on both the front and back sides for reliable visibility.

In addition to function, Premier's face shields also enable maximum comfort with a 1.5-inch foam padding that lifts the shield away from the face. The extra clearance allows room for most eyewear while an elastic band holds the face shield firmly in place to prevent slippage.

### The Premier Advantage:

- ✓ Helps protect face from aerosols, sprays & splatter
- ✓ Optically-clear with full peripheral vision
- ✓ Anti-fog coating on both sides for reliable clarity
- ✓ Comfortable fit with foam padding
- ✓ Fits over most eyewear (glasses & loupes)
- ✓ Elastic band holds shield firmly in place
- ✓ Single-use, disposable, non-sterile product



### Available now!

1070000-PK20 20 shields/box  
1070000-PK50 50 shields/box

The Premier Disposable Face Shields is compliant with the International and American Test Standard: ANSI/ISEA Z87.1-2015. American National Standard for Occupational and Educational Personal Eye and Face Protection Devices. Test results on file and available upon request.



[premierdentalco.com/promise](https://premierdentalco.com/promise)

**premier** Inspired solutions  
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • [premierdentalco.com](https://premierdentalco.com)

©2020 Premier Dental Co. | All Rights Reserved 8610705 0620057 Rev0