

**Reliability & Efficiency**

## Want reliable cements?

Premier® has you covered!

### **NexTemp®**

**Non-Eugenol Temporary Resin Cement with Fluoride**

For your provisional restorations, NexTemp provides firm retention that lasts between visits yet can easily be removed when the time comes. Now available in clear shade!

### **ZR-Cem™**

**Self-Adhesive Universal Resin Cement**

Specially formulated for zirconia restoration and an excellent choice for all C&B substrates. Contains a BPO/amine-free initiation system for long lasting color stability.

### **Premier® Implant Cement™ Plus**

**Non-Eugenol Temporary Resin Cement for Implant Retained Crowns**

Strong, secure retention enables easy, non-destructive crown retrieval if necessary. Unsurpassed radiopacity<sup>†</sup> enables clear visibility on radiographs to help you deliver better patient outcomes.

*Choose Premier cements for reliable and predictable results you can count on – every time.*

### **Premier® Universal Primer**

Zirconia and Glass Ceramics, Metal and Composite Primer

An ideal primer to increase retention when cementing in less retentive preps (short or over-tapered). It contains two coupling agents, MDP and silane methacrylate for stronger adhesion to all substrates – *now that's peace of mind!*

**premier®** Inspired solutions  
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • [premierdentalco.com](http://premierdentalco.com)



Now in Clear Shade!

## NexTemp®

Non-Eugenol Temporary Resin Cement with Fluoride

Available in original opaque and new clear shade!

- ✓ Lasts until the next scheduled appointment – yet can easily be removed when needed
- ✓ Keep patients happy – predictable retention minimizes premature dislodgement
- ✓ Eugenol-free formula – no interference with bonding agent and cementation of final restoration
- ✓ Contains fluoride, potassium nitrate and chlorhexidine

### NexTemp Standard Pack

One 5ml Automix Syringe, 10 Mix Tips

3001471 Shade Opaque

3001472 *new!* Shade Clear

### NexTemp Bulk Pack

Four 5ml Automix Syringes

3001470 Shade Opaque

3001473 *new!* Shade Clear

### Automix Syringe Tips

3001451 Package of 25



## ZR-Cem™

Self-Adhesive Universal Resin Cement

Latest Technology = Better Patient Experience

- ✓ Innovative MDP formula enhances adhesion & bond to zirconia substrates
  - ✓ Lasting color stability
  - ✓ Virtually no post-operative sensitivity ensures patient comfort
- Easy to Use**
- ✓ No pre-treatment is necessary, just dispense the cement and seat the restoration
  - ✓ Optimal viscosity, easy handling and a fast & easy clean-up
  - ✓ Room temperature storage – no need for refrigeration

### ZR-Cem Standard Pack

One 8g Syringe, 10 Mix Tips, 10 Intraoral Tips

3001560 Shade White

3001561 Shade A2

3001562 Shade Translucent

### ZR-Cem Refill Tips

3001418 35 mix tips & 35 intraoral tips



## Premier® Implant Cement *Plus*

Non-Eugenol Temporary Resin Cement for Implant Retained Crowns

- ✓ Easy to use – no etch, primers or adhesive required
- ✓ Gel stage in just 2 minutes – excess easy to remove
- ✓ Unsurpassed radiopacity – excellent visibility on radiographs
- ✓ Secure retention – for long-lasting hold
- ✓ Retrievable – without damage to the crown, abutment or implant
- ✓ Available in 2 shades – white and pink

### Premier Implant Cement *Plus* Standard Pack

One 5ml Automix Syringe, 10 Mix Tips

3001550 Shade Original (Pink)

3001551 Shade White

### Premier Implant Cement *Plus* Value Pack

Three 5ml Automix Syringes, 25 Mix Tips

3001555 Shade Original (Pink)

3001556 Shade White

### Automix Syringe Tips

3001451 Package of 25



## Premier® Universal Primer

Zirconia and Glass Ceramics, Metal and Composite Primer

3001568 5ml Bottle

New Premier® Universal Primer, is a single-component adhesive primer that increases the retention between restorative substrates and resin cements.

- ✓ Ideal retention for zirconia, alumina, metal, glass ceramics + composite materials
- ✓ Dual coupling agents provides versatility and work with all substrates
- ✓ Easy to use – just apply, wait 60 seconds, and dry



Available to order now at  
[premierdentalco.com/cements](http://premierdentalco.com/cements)

**premier®** Inspired solutions  
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA  
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • [premierdentalco.com](http://premierdentalco.com) • [info@premierdentalco.com](mailto:info@premierdentalco.com)

©2020 Premier Dental Co. | All Rights Reserved / †When compared to competitive implant cements – data on file.

8621271 1120023 Rev0 PDF