

## AeroPro® Delivers Prophy Peace of Mind!

Great Ergonomics, Less Aerosol



## AeroPro®

The new prophy standard your practice needs

AeroPro's light-weight, ergonomic cordless and pedal-less handpiece gives you the mobility you want, the choice of prophy angles you like, and the control you need. Fully CDC-compliant, AeroPro's programmed speeds effectively polish teeth yet generate less volume of aerosol, less aerosol spread, and less splatter vs. other leading brands according to a recent UCI study.\* Why risk over-spinning with a rheostat-driven handpiece?

AeroPro's innovative ergonomic design helps reduce hand fatigue by allowing you to rotate the angle 360° for improved intraoral access and comfort. Its quick-touch button enables fingertip control for easy switching between low and high speeds. Its quiet yet powerful motor maintains a consistent speed during polishing—no extra pressure required. ChargeSMART™ Battery Technology gives you long-lasting power to get through the workday.



**premier** Inspired solutions  
for daily dentistry®

# AeroPro®

## Cordless Prophy Handpiece

### Less Aerosol, Less Splatter vs. Other Leading Brands\*

- ✓ AeroPro provides constant, effective cleaning
- ✓ Generates less aerosol and less splatter vs. other leading brands.  
*Why risk overspinning with a rheostat-driven handpiece?*
- ✓ CDC-compliant with a convenient infection prevention protocol



### Ergonomics Promote Hand & Wrist Health

- ✓ Easy 360° angle rotation - improved intraoral access and comfort
- ✓ Quick-touch button - enables easy switching between speeds
- ✓ Powerful motor maintains a consistent speed  
- no extra pressure required

### Long-Lasting Power + 3 Year Warranty

- ✓ ChargeSMART™ Battery Technology gets you through the workday
- ✓ 3 Year Warranty - Includes handpiece + battery

#### THE RESULTS ARE IN: AeroPro® generates less aerosol than the competition\*



- 5500510 AeroPro® Complete Package**  
(3) Autoclavable Outer Sheaths,  
1 each: Motor Component, Charging Station,  
AC Adapter, Stand, Box of 500 Disposable Barriers,  
Accessories
- 5500500 AeroPro® Basic Kit**  
(1) Autoclavable Outer Sheath,  
1 each: Motor Component, Charging Station,  
AC Adapter, Stand, Box of 500 Disposable Barriers,  
Accessories
- 5500520 AeroPro® Autoclavable Outer Sheath**
- 5500530 AeroPro® Disposable Barriers - Box of 500**

**premier®** Inspired solutions  
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA  
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com

\*Effect of Handpiece and Polishing Speed on Volume of Aerosol Generated and Distance of Aerosol Spread Using No Suction.  
Based on a University of California Irvine study, 2021. Data on file. / Freedom and Infinity are not trademarks of Premier.

8610327 0620061 Rev5 PDF