

# AeroPro® New & Even Better!

“ AeroPro is lightweight  
and takes the strain off  
my wrist that I didn't even  
realize was going on!

- Jessica I, RDH



## AeroPro® Cordless Prophylaxis Handpiece

AeroPro® just got even better. The enhanced quick-touch button provides better ergonomics and precise control to remove stains and polish your patients' teeth. Utilizing new sealing technology, AeroPro can be reprocessed safely and quickly with increased durability. AeroPro generates less aerosol, less splatter vs. other leading handpieces according to a University of California Irvine study<sup>1</sup>.

1. Based on a University of California Irvine study, 2021.



Enhanced quick-touch button  
*Better ergonomics, control & protection*



Lightweight, 100% mobile design  
*Promotes hand & wrist health*



New sealing technology  
*Increased durability for safe reprocessing*



Less aerosol & less splatter  
*Precise finger-tip speed control*



[premierdentalco.com/promise](https://premierdentalco.com/promise)

**premier** Inspired solutions  
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • [premierdentalco.com](https://premierdentalco.com)

©2023 Premier Dental Co. | All Rights Reserved

# AeroPro®

## Cordless Prophylaxis Handpiece

### Ergonomics Promote Hand & Wrist Health

- ✓ Easy 360° angle rotation - improved intraoral access and comfort
- ✓ Quick-touch button - enables easy switching between speeds
- ✓ Powerful motor maintains a consistent speed  
- no extra pressure required

### Long-Lasting Power + 3 Year Warranty



- ✓ ChargeSMART™ Battery Technology gets you through the workday
- ✓ 3 Year Warranty - Includes handpiece + battery

### Less Aerosol, Less Splatter vs. Other Leading Brands\*

- ✓ AeroPro provides constant, effective cleaning
- ✓ Generates less aerosol and less splatter vs. other leading brands.  
*Why risk overspinning with a rheostat-driven handpiece?*
- ✓ CDC-compliant with a convenient infection prevention protocol

#### THE RESULTS ARE IN:

AeroPro® generates less aerosol than the competition\*



**Accepts most disposable angles**  
Rotates 360° for improved intraoral access

**New! Enhanced quick-touch button**  
Easily switch speeds

**New & Even Better!**

**New! Sealing technology**  
Increased durability for safe reprocessing

**Dimpled silicone grip**  
Improves tactile feel while preventing slippage

**Quiet and powerful German-made motor**  
Constant speed without extra pressure

**ChargeSMART™ Battery Technology**  
Charges quickly and provides lasting power for all day use

**5500510 AeroPro® Complete Package**  
(3) Autoclavable Outer Sheaths,  
1 each: Motor Component, Charging Station,  
AC Adapter, Stand, Box of 500 Disposable Barriers, Accessories

*Best value!*

**5500500 AeroPro® Basic Kit**  
(1) Autoclavable Outer Sheath,  
1 each: Motor Component, Charging Station,  
AC Adapter, Stand, Box of 500 Disposable Barriers, Accessories

**5500520 AeroPro® Autoclavable Outer Sheath**

**5500530 AeroPro® Disposable Barriers - Box of 500**

**premier®** Inspired solutions  
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA  
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com

\*Effect of Handpiece and Polishing Speed on Volume of Aerosol Generated and Distance of Aerosol Spread Using No Suction.  
Based on a University of California Irvine study, 2021. Data on file. / Freedom and Infinity are not trademarks of Premier.

8610327 0123056 Rev6 NKP10M