



Reliability & Efficiency

Want reliable cements?

Premier® has you covered!

NexTemp®

Non-Eugenol Temporary Resin Cement
with Fluoride

For your provisional restorations, NexTemp provides firm retention that lasts between visits yet can easily be removed when the time comes. Now available in clear shade!

ZR-Cem®

Self-Adhesive Universal Resin Cement

Specially formulated for zirconia restoration and an excellent choice for all C&B substrates. Contains a BPO/amine-free initiation system for long lasting color stability.

Premier® Implant Cement™ Plus

Non-Eugenol Temporary Resin Cement
for Implant Retained Crowns

Strong, secure retention enables easy, non-destructive crown retrieval if necessary. Unsurpassed radiopacity[†] enables clear visibility on radiographs to help you deliver better patient outcomes.

Choose Premier cements for reliable and predictable results you can count on – every time.

Premier® Universal Primer

Zirconia and Glass Ceramics, Metal and Composite Primer

An ideal primer to increase retention when cementing in less retentive preps (short or over-tapered). It contains two coupling agents, MDP and silane methacrylate for stronger adhesion to all substrates – *now that's peace of mind!*

premier® Inspired solutions
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • premierdentalco.com



Now in Clear Shade!

NexTemp®

Non-Eugenol Temporary Resin Cement with

Available in original opaque and clear shade!

- ✓ Lasts until the next scheduled appointment – yet can easily be removed when needed
- ✓ Keep patients happy – predictable retention minimizes premature dislodgement
- ✓ Eugenol-free formula – no interference with bonding agent and cementation of final restoration
- ✓ Contains fluoride, potassium nitrate and chlorhexidine

NexTemp Standard Pack

One 5ml Automix Syringe, 10 Mix Tips

3001471 Shade Opaque

3001472 Shade Clear

NexTemp Bulk Pack

Four 5ml Automix Syringes

3001470 Shade Opaque

3001473 Shade Clear

Automix Syringe Tips

3001451 Package of 25



ZR-Cem®

Self-Adhesive Universal Resin Cement

Latest Technology = Better Patient Experience

- ✓ Innovative MDP formula enhances adhesion & bond to zirconia substrates
- ✓ Lasting color stability
- ✓ Virtually no post-operative sensitivity ensures patient comfort

Easy to Use

- ✓ No pre-treatment is necessary, just dispense the cement and seat the restoration
- ✓ Optimal viscosity, easy handling and a fast & easy clean-up

ZR-Cem Standard Pack

One 8g Syringe, 10 Mix Tips, 10 Intraoral Tips

3001560 Shade White

3001561 Shade A2

3001562 Shade Translucent

ZR-Cem Refill Tips

3001418 35 mix tips & 35 intraoral tips



Premier® Implant Cement Plus

Non-Eugenol Temporary Resin Cement for Implant Retained Crowns

- ✓ Easy to use – no etch, primers or adhesive required
- ✓ Gel stage in just 2 minutes – excess easy to remove
- ✓ Unsurpassed radiopacity – excellent visibility on radiographs
- ✓ Secure retention – for long-lasting hold
- ✓ Retrievable – without damage to the crown, abutment or implant
- ✓ Available in 2 shades – white and pink

Premier Implant Cement Plus Standard Pack

One 5ml Automix Syringe, 10 Mix Tips

3001550 Shade Original (Pink)

3001551 Shade White

Premier Implant Cement Plus Value Pack

Three 5ml Automix Syringes, 25 Mix Tips

3001555 Shade Original (Pink)

3001556 Shade White

Automix Syringe Tips

3001451 Package of 25



Premier® Universal Primer

Zirconia and Glass Ceramics, Metal and Composite Primer

3001568 5ml Bottle

New Premier® Universal Primer, is a single-component adhesive primer that increases the retention between restorative substrates and resin cements.

- ✓ Ideal retention for zirconia, alumina, metal, glass ceramics + composite materials
- ✓ Dual coupling agents provides versatility and work with all substrates
- ✓ Easy to use – just apply, wait 60 seconds, and dry



Available to order now at
premierdentalco.com/cements

premier® Inspired solutions
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • info@premierdentalco.com

©2023 Premier Dental Co. | All Rights Reserved / †When compared to competitive implant cements – data on file.

8621271 0523074 Rev1 NKP5M