



Introducing!

New! Two Striper® Perforated Diamond Finishing Strips



Perforated Strip



Combo-Strip



Two Striper®
— AN  abrasive
technology BRAND —

Two Striper® Perforated Diamond Finishing Strips are designed for interproximal reduction, shaping, and contouring. The perforations assist in debris removal.

Available in two options:

- ✓ Perforated Strip - Features an uncoated safe center for easy access through the contact point
- ✓ Combo-Strip - Includes a serrated center to facilitate access in tight spaces between teeth

Available through your authorized dealer.

Available in 10-packs

2 widths
300=2.5mm
350=3.75mm

Perforated

2013132 300 Fine
2013142 300 Medium
2013152 300 Ultra-Fine
2013162 350 Fine
2013172 350 Medium
2013182 350 Ultra-Fine

2013192 Assorted

*Contains: (2) 300 Fine, (2) 300 Medium,
(1) 300 Ultra-Fine, (2) 350 Fine,
(2) 350 Medium, (1) 350 Ultra-Fine*

Combo-Strip

2013161 350 Fine
2013171 350 Medium
2013181 350 Ultra-Fine
2013191 Assorted

*Contains: (4) 350 Fine, (3) 350 Medium,
(3) 350 Ultra-Fine*

Two Striper® is a registered trademark of the product manufacturer, Abrasive Technology, LLC / Lit# 8621303 0823089 Rev0 PDF



Distributed by

premier®

Inspired solutions
for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • premierdentalco.com

©2023 Premier Dental | All Rights Reserved

Introducing!

New! Two Striper® Bulk Packs

Your favorite Two Striper® diamonds are now available in economic packs of 50!



- ✓ Cuts faster – *less chair time, happier patients*
- ✓ Sharper diamonds that stay put – *Proprietary P.B.S.® brazing process*
- ✓ Lasts longer for multiple preps – *Use fewer burs than 'bargain' diamonds*



	245	330	557	2005.8C	230C	250.8C	251.8C	260.10VF	260.8C	260.8F
5-Pack	2015143	2015457	2015467	2015459	2015123	2013281	2013251	2015345	2015455	2015173
Bulk 50-Pack	2050143	2050457	2050467	2050459	2050123	2050557	2050245	2050345	2050455	2050173

	260.8VF	261.8C	283.4C	285.5C	285.5F	285.5SC	285.5VF	515.7C	702.8C	703.8C
5-Pack	2015183	2015124	2015142	2015143	2015173	2015700	2015203	2015277	2015343	2015345
Bulk 50-Pack	2050183	2050124	2050142	2050143	2050173	2050700	2050203	2050277	2050343	2050345

	767.7C	767.9C	770.10C	770.8C	770.8F	770.8SC	780.9C	782.8C	790.8C	863C
5-Pack	2015453	2015455	2015459	2015457	2015487	2015701	2015467	2015424	2015471	2015545
Bulk 50-Pack	2050453	2050455	2050459	2050457	2050487	2050701	2050467	2050424	2050471	2050545

Two Striper® and P.B.S.® are registered trademarks of the product manufacturer, Abrasive Technology, LLC

Lit# 8621304 0823092 Rev0 PDF



Distributed by
premier®
 Inspired solutions
 for daily dentistry®

Premier® Dental Products Company • 888-670-6100 • premierdentalco.com

©2023 Premier Dental | All Rights Reserved