

Inspired Innovation with Premier®



**new! VeneerNow™**  
Injectable Matrix System  
*See page 1.*

## Premier® Restorative Products

Every procedure you perform shapes your patients' impression of you. We've got you covered with solutions engineered to deliver reliable results. Every time.

**new! Traxodent®**  
Unit Dose Value Pack  
*See page 4.*



**new! PremierAir®**  
Titanium Implant Scaler  
*See page 12.*



*new!*  
**VeneerNow™**  
Injectable Matrix System



**Simplify your anterior composite restorations!**

Elevate your dental practice and improve your workflows with the VeneerNow™ Injectable Matrix System – your solution for quickly creating chair-side attractive, natural-looking composite veneers and anterior restorations. Designed for ease of use, it enables composite injection with minimal clean-up.

VeneerNow also works with packable composite, offering versatility. This essential addition to your operatory toolkit makes anterior restorations fast with predictable results.

**VeneerNow™** Fast, Predictable Results!



**Full anatomical design**  
Replicates natural tooth shape for superior results



**Side wings**  
Enhance stability & minimize interproximal clean-up



**Closed system matrix**  
Provides a void-free smooth composite layer for easier finishing and polishing



**Natural cervical line**  
For a precise anatomical emergence profile



**Injection stalk for flowable composite**  
Provides easy injection for optimal placement



**Versatile**  
Use with flowable or packable composite

## Ordering information

### 9069237 Starter Kit

Contains:

32 Small/Medium - SM Teeth 5-12, Palmer 4 - 4 (4 ea.)

16 Medium/Large - ML Teeth 5-12, Palmer 4 - 4 (2 ea.)

### 9069244 Mini Kit

Contains:

8 Small/Medium - SM Teeth 7-10, Palmer 2 - 2 (2 ea.)

4 Medium/Large - ML Teeth 7-10, Palmer 2 - 2 (1 ea.)



### Refills:

#### Small/Medium - SM

Contains: 6 matrices

|         |          |                 |        |
|---------|----------|-----------------|--------|
| 9069268 | Tooth 5  | Palmer <u>4</u> | FDI 14 |
| 9069275 | Tooth 6  | Palmer <u>3</u> | FDI 13 |
| 9069282 | Tooth 7  | Palmer <u>2</u> | FDI 12 |
| 9069299 | Tooth 8  | Palmer <u>1</u> | FDI 11 |
| 9069305 | Tooth 9  | Palmer <u>1</u> | FDI 21 |
| 9069312 | Tooth 10 | Palmer <u>2</u> | FDI 22 |
| 9069329 | Tooth 11 | Palmer <u>3</u> | FDI 23 |
| 9069336 | Tooth 12 | Palmer <u>4</u> | FDI 24 |

### Refills:

#### Medium/Large - ML

Contains: 6 matrices

|         |          |                 |        |
|---------|----------|-----------------|--------|
| 9069343 | Tooth 5  | Palmer <u>4</u> | FDI 14 |
| 9069350 | Tooth 6  | Palmer <u>3</u> | FDI 13 |
| 9069367 | Tooth 7  | Palmer <u>2</u> | FDI 12 |
| 9069374 | Tooth 8  | Palmer <u>1</u> | FDI 11 |
| 9069381 | Tooth 9  | Palmer <u>1</u> | FDI 21 |
| 9069398 | Tooth 10 | Palmer <u>2</u> | FDI 22 |
| 9069404 | Tooth 11 | Palmer <u>3</u> | FDI 23 |
| 9069411 | Tooth 12 | Palmer <u>4</u> | FDI 24 |



Refill




View VeneerNow™  
Instructional Videos

For additional video instruction by indications, visit [premierdentalco.com/VeneerNow](http://premierdentalco.com/VeneerNow) or scan QR code.



## Traxodent®

Hemostasis Retraction Paste System



Traxodent applied



Post-op

## The #1 Choice for hemostasis & retraction<sup>1</sup>

Award-winning Traxodent® is the #1 choice<sup>1</sup> among dental professionals because it works fast, helps capture accurate impressions and rinses easily.

And since Traxodent doesn't irritate or discolor surrounding tissue, patients leave happier too.

*Count on Traxodent for fast hemostasis & retraction!*



**Effective Hemostasis  
& Retraction**



**Fast & Easy  
Handling**



**Rinses Clean  
With No Residue**



**Versatility**




View Traxodent® Videos

Scan QR or visit  
[premierdentalco.com/Traxodent](https://premierdentalco.com/Traxodent)

*"The product provided excellent  
retraction and is easy to use  
and gentle to the tissue."*

-Mark P. Cohen, DDS  
Roselle Park, NJ

Ergonomic,  
autoclavable unit dose  
syringe dispenser



Sleek syringe with bendable tips

1. SDM Northcoast, North America 2024.





Unit Dose Value Pack  
(50 Capsules)

## Routine Hemostasis



Old restoration is removed  
then Traxodent is placed.



Rinse after 2 minutes. Bleeding  
is stopped and field is isolated.

Photos courtesy of  
Richard H. Nagelberg, DDS, PC

## Ordering Information

- 9007094 Traxodent Unit Dose Pack  
Contains: 24 capsules
- 9007090 **NEW!** Traxodent Unit Dose Value Pack  
Contains: 50 capsules
- 9007089 Unit Dose Syringe Dispenser (3)
- 9007091 Traxodent Value Pack  
Contains: 25 syringes and 50 tips
- 9007093 Traxodent Starter Pack  
Contains: 7 syringes and 15 tips
- 9007092 Applicator Tips (50)

## Premier® Retraction Caps

Anatomic formed cotton caps work synergistically with Traxodent to quicken and assist retraction and hemostasis.

- ✓ Easily guides Traxodent into the sulcus
- ✓ Compression assists and quickens hemostasis
- ✓ Cotton absorbs fluids
- ✓ Anatomic form allows easy placement

- 9048251 Size 1 - Small (120)
- 9048253 Size 3 - Medium (120)
- 9048255 Size 5 - Large (60)
- 9048260 Assorted - (15) #1, (15) #3, (30) #5



## Crown and Bridge with Retraction Cap



Syringe Traxodent and place retraction cap  
on prep; have patient bite down on cap for  
2 minutes.



Tooth is prepared for  
a clean impression.

Photos courtesy of  
Dean Ellledge, DDS, MS

## Crown and Bridge with Retraction Cord



A single retraction cord is packed  
around the preps.



Place Traxodent over the cord.  
Wait 2 minutes.

Photos courtesy of  
Dean Ellledge, DDS, MS

## Loose Tissue



Cord is placed.



Place retraction cap  
on prep; have patient  
bite down on cap for  
2 minutes.



Impression with  
well-defined margins

Photos courtesy of  
William Crinzy, DDS, MS

## IntegraPost®

Ultra-passive Integrating Core Retainer




- ✓ Safe and Secure - The IntegraPost features an innovative round Flow-Thru-Head® and a parallel passive design to optimize the adhesive buffer zone, conserve tooth structure and minimize stress.<sup>2,3,4</sup>
- ✓ Balances Strength with Elasticity - Made from bio-compatible titanium alloy, IntegraPost is stronger than resin fiber, yet more elastic than stainless steel.<sup>5</sup>
- ✓ Easy to Use - Five post sizes make post selection a simple process. The unique post carrier facilitates handling.

In an **ADA Professional Product Review** on  
Endodontic Posts, which included metal, fiber, and a ceramic post,  
**IntegraPost** received the **highest laboratory**  
**rating of 92% with a buyer's satisfaction of excellent!**<sup>6</sup>

| IntegraPost Refills (10 each) |          |         |            |
|-------------------------------|----------|---------|------------|
| 3001105                       | Size 3   | 3001301 | IP Carrier |
| 3001110                       | Size 4   | 3001305 | IP Probe   |
| 3001115                       | Size 4.5 |         |            |
| 3001125                       | Size 5   |         |            |
| 3001130                       | Size 6   |         |            |

## Cure-Thru® IntegraPost®

Light-transmitting Zirconia Glass Fiber Post



- ✓ Safe and Simple to Use – Composite Zirconia Glass Fiber composition won't fracture root and provides easy non-traumatic removal.
- ✓ Aesthetic – White color of post guarantees no grey shadow or show through.
- ✓ Strong and Flexible – Elastic modulus and flexural strength more like tooth structure allowing the post to respond like the tooth to traumatic stresses. (E-modulus 52 GPa/Flexural Strength 1650 MPa)
- ✓ Factory Silanated for Secure Adhesion – Post is chemically treated to bond with both resin cement and composite core build-up material.
- ✓ Light Transmitting – Photo curing through the post assures complete polymerization of cement around and down the post length.
- ✓ Four Color Coded Post Sizes – Multiple sizes for easy fit in all teeth and canal space. Color coding make size selection easy.
- ✓ Tapered Tip Design – Parallel body with slightly tapered tip conforms to apical root anatomy for safe and effective adaptation of post.



| Refills (10 pack) |          |          |
|-------------------|----------|----------|
| 3001610           | Size 4   | (red)    |
| 3001611           | Size 4.5 | (purple) |
| 3001612           | Size 5   | (blue)   |
| 3001613           | Size 6   | (green)  |

Not available in Canada.

2-5. Full bibliography available. / 6. American Dental Association Professional Product Review™  
Publication by the ADA Council on Scientific Affairs, Vol. 1, Spring 2006. / IntegraPost designed by  
Endodontist, Dr. Phil Klein.



Non-impregnated  
Knitted Weave


## Knit-Pak™

### Knitted Retraction Cord

Use in combination with Hemodent® or Traxodent® for tissue displacement and hemostasis.

- ✓ Absorbent knitted cotton can take in an abundant amount of hemostatic solution
- ✓ Compacts under pressure, making it easier to pack
- ✓ Will not get entangled in a bur
- ✓ Measures 100 inches / 254 cm

|         | Size | Color  |
|---------|------|--------|
| 9007551 | 000  | Green  |
| 9007552 | 00   | Brown  |
| 9007553 | 0    | Purple |
| 9007554 | 1    | Blue   |
| 9007555 | 2    | Orange |
| 9007556 | 3    | Yellow |



A cutting blade inside the cap automatically cuts the cord upon closing firmly.



Impregnated  
Knitted Weave

## Knit-Pak+™

### Impregnated Knitted Cord

Unique Knitted Microfiber with Superior Absorbency

Knit-Pak+ is knitted impregnated retraction cord manufactured from unique microfibers.

- ✓ Aluminum chloride impregnated cord (AlCl<sub>3</sub>·6H<sub>2</sub>O)
- ✓ High durability and retention properties
- ✓ Resists fraying and easily packs into sulcus
- ✓ Measures 100 inches / 254 cm

|         | Size  | Color      |
|---------|-------|------------|
| 9007560 | 000.0 | Dark Green |
| 9007561 | 000   | Green      |
| 9007562 | 00    | Brown      |
| 9007563 | 0     | Purple     |
| 9007564 | 1     | Blue       |
| 9007565 | 2     | Orange     |



## Hemodent®

### Buffered Aluminum Chloride Hemostatic Solution

- ✓ Fast & effective
- ✓ Unique formulation
- ✓ Wide range of uses
- ✓ 3-Year shelf life



|         |             |
|---------|-------------|
| 9007071 | 10cc bottle |
| 9007072 | 20cc bottle |
| 9007073 | 40cc bottle |

## IntegraBond™

### Universal (5th gen) Adhesive System

Clinically proven IntegraBond<sup>7-15</sup> - provides outstanding performance for all restorative challenges!

- ✓ Offers outstanding shear bond strength to all popular light, dual- or auto-cure composites.

Use IntegraBond after acid-etching.

|         |   |
|---------|---|
| 3001411 | IntegraBond - 7ml   |
| 3001412 | Auto-Cure Activator - 7ml for use when light-curing is not possible |
| 3001413 | IntegraBond - 4ml   |



## Premier® Universal Primer

### Zirconia and Glass Ceramics, Metal and Composite Primer

- ✓ Dual coupling agents works with all substrates
- ✓ Easy to use – just apply, wait 60 seconds, and dry
- ✓ No pretreatment of restoration required
- ✓ Can be used with adhesive and self-adhesive resin cement to achieve maximum bond strength



|         |            |
|---------|------------|
| 3001568 | 5ml bottle |
|---------|------------|



Not available in Canada.

7-15. Full bibliography available.


Not available in Canada.

## CompCore™ AF

Dual-Cure Composite Core Build-up Material

**Faster & stronger**  
**A better core experience!**

Build cores with confidence! CompCore™ AF resin composite is easy to place, shape and finish. The composite cuts like dentin with excellent visibility on radiographs while its superior mechanical properties and fluoride release deliver long-term clinical success.



Proprietary Polymer Strengtheners



**Outstanding  
Handling**



**Optional  
Dual-Cure**



**Superior Mechanical  
Properties**



**Minimum Shrinkage  
of < 2.0%**



### CompCore™ AF SyringeMix STACK:

9g Automix Syringes/Dual-Cure

3001512 SyringeMix STACK / Shade White

3001513 SyringeMix STACK / Shade A3

3001514 SyringeMix STACK / Shade Blue

Contains: two 9g automix syringes, 15 mix tips,  
15 standard and 10 extra-fine intraoral tips.

### Refill Kit

3001432 Shade White

3001433 Shade A3

Contains: 50g cartridge of CompCore AF Dual-Cure,  
25 mix tips, 25 intraoral syringe tips

### Accessories

3001418 SyringeMix Mix-Tips & Intraoral Tips (35 each)

3001468 SyringeMix Extra-Fine Intraoral Tips (50)


3001417 Cartridge Mix-Tips & Intraoral Tips (35 each)





## Bonding Agent

(use only bonding agents compatible with dual-cure composites, such as IntegraBond™)



1. Clean and dry prepared tooth surface.
2. Apply 37% phosphoric acid to the cavosurface for 10 to 15 seconds.
3. Rinse and dry cavo surface, do not desiccate the tooth.
4. Apply first coat of IntegraBond™ on prepared cavosurface. Allow to sit for 15 seconds, and then air dry. Apply 2 to 3 more coats; without waiting, dry again after last application.
5. Light-cure for 20 seconds.

## Core Build-up



1. Bleed syringe. Assemble intraoral tip to mix-tip, attach to syringe and bleed tip.
2. Place intraoral tip into the base of the prepared tooth and dispense. Let material rise around the tip while withdrawing to eliminate air bubbles.
3. Light-cure for 30 seconds. If the core is larger than 8mm, light-cure for at least an additional 10 seconds.
4. Check for complete hardness. The preparation can be finished immediately with a Two Striper® bur. Alternately, the core can be chemical-cured, final set-time is 4.5 minutes.

## NexTemp®

Non-Eugenol Temporary Resin Cement

**Your go-to temporary cement**  
*Stays Put, Removes Easily*



NexTemp is the reliable temporary cement you can count on. That's because NexTemp creates an excellent marginal seal with firm retention to help prevent premature dislodgment or microleakage.

What's more, NexTemp is the only eugenol-free resin cement containing potassium nitrate, chlorhexidine and fluoride – which studies show help mitigate post-op sensitivity.<sup>16</sup>



**Temporary  
Restorations Stay Put**



**Eugenol-  
Free**



**Proprietary  
Formula**



**Affordably  
Priced**



### Opaque shade


3001470 Bulk pk.: 4-5ml automix syringes

3001471 Standard pk.: 1-5ml automix syringe, 10 mix tips

### Clear shade

3001473 Bulk pk.: 4-5ml automix syringes

3001472 Standard pk.: 1-5ml automix syringe, 10 mix tips




### **Greater Retention, Less Cement!**

*Just dispense a small bead of cement  
around the inside margin.*


*"This is such a superior product  
to all the other temp cements  
we have tried over the years."*

– Peter Schott, DMD  
Middletown, RI


### Premier® NexTemp® Temporary Cement is easy to use




Dry and isolate the  
prepared tooth.




Bleed the syringe.  
Mount a new automix  
dispensing tip and bleed  
the syringe again.




Dispense a small bead  
of cement around the  
inside margin.



Seat the temporary  
and hold in place. In  
approximately 2 minutes,  
excess cement reaches a  
distinct gel phase and can  
be removed safely.



Final set occurs  
between 4-5 minutes  
after seating of  
temporary restoration.†



Use a PerioWise®  
Periodontal Probe  
for post-treatment  
periodontal evaluation.

\*An independent non-profit dental education and product testing foundation, December 2006 & 2007 / †Final set time is determined intraorally. Ambient set times vary and are not representative of clinical conditions. / 16. Jalalian E, Meraji N, Miraei M. A comparison of the efficacy of potassium nitrate and Gluma desensitizer in the reduction of hypersensitivity in teeth with full-crown preparations. J Contemp Dent Pract. 2009 Jan 1;10(1):66-73.

## ZR-Cem®

**A better zirconia cementation solution is here!**  
New ZR-Cem® Self-adhesive resin cement

ZR-Cem® is an innovative self-adhesive universal resin cement specially formulated for zirconia restorations. Dual-cured ZR-Cem enables a strong bond to all ceramic materials, dentin and enamel - ensuring superior retention and marginal integrity.



**Innovative  
MDP Formula**



**Lasting  
Color Stability**



**Virtually No  
Post-op Sensitivity**



**No  
Pre-treatment**



### Standard Pack

1 - 8g syringe, 10 mix-tips and 10 intraoral tips

3001560 Shade Opaque White

3001561 Shade A2

3001562 Shade Translucent

### Refill Tips

3001418 35 mix-tips & 35 intra-oral tips

### Light cure

**2-3 seconds tack cure - easily remove excess cement**

*For high-translucency zirconia and other light-transmittable restorations (material thickness ≤ 2.5mm), light cure all surfaces including margins for 20 seconds to accelerate final set.*

### Chemical cure

**72°F (22°C) temperature**

Working time: 1.5 - 2.5 minutes.

*Higher temperatures accelerate set.*

**98.6°F (37°C) temperature**

Gel phase: approximately 1-1.5 minutes

Final set: 3.0 - 4.0 minutes.

*Higher temperatures accelerate set.*

### Physical Properties

Compressive strength: 150 MPA

Flexural strength: 100 MPA

Film thickness: 10 µm

Radiopacity: 250% Al

Shelf life: 2 years




3-second tack cure



Easy clean-up

Photos courtesy of Timothy M. Bizga, DDS

### ZR-Cem: superior bond strength vs. other leading self-adhesive cements to zirconia†



MaxCem Elite, RelyX Unicem 2 and Calibra Universal are not trademarks of Premier Dental.  
† Latta, M. Creighton University School of Dentistry Study 2019.

*Not available in Canada.*

## Premier® Implant Cement™ Plus


The ultimate solution for implant crown cementation.

- ✓ Secure retention – *for long-lasting hold*
- ✓ Retrievable – *without damage to the crown, abutment or implant*
- ✓ Easy to use – *no etch, primers or adhesive required*
- ✓ Unsurpassed radiopacity\* – *excellent visibility on radiographs*

Unsurpassed  
Radiopacity




### Premier® Implant Cement™ is easy to use



1 Check the fit. Complete all finishing procedures.

2 Rinse and dry abutment and interior of restoration.


3 Bleed syringe. Mount new automix tip and then bleed the syringe again. This assures proper flow of base and catalyst.



4a Using a microbrush, place a thin film of cement along the inside wall of the cervical third, making sure to completely coat the margin with cement. You have 1.5 minutes of working time at room temperature.

4b Excess cement should be avoided. Too much cement can cause future restoration retrieval to be difficult to impossible.\*


5 Firmly seat the crown and hold in place until initial set (gel phase) occurs in 2–2.5 minutes. Safely remove excess cement before final set. It is imperative to remove all excess cement.



6 Final cement is 4–5 minutes, intraorally.

\*Any excess cement (regardless of chemical make-up or brand) that remains in the sulcus around an implant could lead to peri-mucositis, peri-implantitis or implant failure if not removed immediately.

Premier Implant Cement Plus offers unsurpassed radiopacity\* with excellent visibility on radiographs.



- 1 Premier® Implant Cement™ Plus - Shade Original (Pink)
- 2 Premier® Implant Cement™ Plus - Shade White
- 3 Improv®
- 4 Implantlink® Semi Xray
- 5 Retrieve™ DC

\*As compared to other implant cements. Improv, Implantlink Semi Xray and Retrieve DC are not trademarks of Premier® Dental Products Co.

### Standard Pack

1-5ml Automix Syringe, 10 Mix Tips  
3001550 Shade Original (Pink)  
3001551 Shade White

### Value Pack

3-5ml Automix Syringes, 25 Mix Tips  
3001555 Shade Original (Pink)  
3001556 Shade White

### Accessories

3001451 Automix Dual-barrel Syringe Tips (25)

Not available in Canada. Ask your Canadian dealer about Premier Implant Cement.



## Triple Tray®

Proven Accurate in Over 185 Million Impressions

- ✓ Three techniques in one  
– *Saves time, money and materials.*
- ✓ Adaptable, ultra-thin webbing  
– *Provides accurate bite registration and low water sorption<sup>17</sup> reducing the possibility of distortion*
- ✓ Disposable – *No cleaning involved, aseptic*
- ✓ Seven configurations – *Designed to address all clinical cases*

|         |   |
|---------|---|
| 1006203 | Posterior - Box of 48   |
| 1006243 | Sideless Posterior - Box of 48                                  |
| 1006283 | Wide Body Posterior - Box of 48                                 |
| 1006213 | Anterior - Box of 35  |
| 1006223 | Quadrant - Box of 40  |
| 1006253 | Extended Quadrant - Box of 35                                   |
| 1006273 | Three-Quarter Arch - Box of 35                                  |
| 1006233 | Assorted - Box of 35<br>Contains 5 each of designs listed above |



17. Independent Research Organization, July 2000.

## T-LOC™ Triple Tray®

Adhesive-free Dual-arch Impression Trays

- ✓ Patented adhesive-free design
- ✓ Easier handling
- ✓ Accurate bite
- ✓ Greater rigidity in posterior
- ✓ Ensure accuracy and save time

|         |                                |
|---------|--------------------------------|
| 1006204 | Posterior - Box of 35          |
| 1006214 | Anterior - Box of 28           |
| 1006274 | Three-Quarter Arch - Box of 22 |



## Spot-On™ Etch

### 37% Phosphoric Etch Gel

Spot-On™ Etch is a 37% high viscosity phosphoric etch gel for all your routine and selective-etch procedures. It's designed with enhanced viscosity that will prevent it from dripping or slumping. Spot-On Etch easily flows through a small applicator tip, washes off readily and resists drying out.

- ✓ Ideal for routine etching and selective-etch procedures
- ✓ Blue color provides a strong contrast to enamel and dentin
- ✓ No slump formula
- ✓ Easy handling
- ✓ Rinses clean without leaving residue



Spot-On™ Etch no-slump viscosity will stay put. Only etch the areas you desire!

Spot-On™ Etch clinical image courtesy of Shalom Mehler DMD, Teaneck, NJ.



- 3001460 Spot-On™ Etch Bulk  
One 12g jumbo filled syringe with luer lock adapter  
Ten 1.2ml empty syringes  
Fifty 1.2ml disposable applicator tips
- 3001461 Spot-On™ Etch Intro  
Five 1.2ml filled syringes  
Twenty bent disposable applicator tips

## Premier® Etch

### 37% Phosphoric Acid Dentin/ Enamel Etch Gel

Premier® Etch contains phosphoric acid for etching tooth enamel and dentin when preparing surfaces for the application of composites.



- 3001422 37% Dentin/Enamel Etch Kit  
Contains: one 12g jumbo syringe, with luer lock adapter, ten empty 1.2ml syringes and 50 disposable applicator tips for the ten 1.2ml syringes

## PremierAir®

Titanium Implant Scalers

4 designs crafted specifically for supra and subgingival maintenance



Scan to view full line of PremierAir® instruments!




Barnhart 5/6 Ti • 9061460



NV4 Ti • 9061461



Langer 3/4 Ti • 9061462



Neb 128B/Langer 5 Mini Ti • 9061463

## Premier® Implant Scalers

Made from graphite reinforced resin that will not scratch implant surfaces.




137 Anterior • 9061421 (2-pack) • 9061420 (5-pack)



204 Posterior  
• 9061423 (2-pack) • 9061422 (5-pack)




Facial (Goldman-Fox)  
• 9061402 (2-pack) • 9061403 (5-pack)



Universal (4L/4R)  
• 9061401 (2-pack) • 9061400 (5-pack)

## Amalgam Carrier




DE Regular (2.00mm) / Large (2.75mm) • 9062201  
DE Large (2.75mm) / Jumbo (3.0mm) • 9062202


## Premier® PerioWise®

"The Friendly Probe"®


- ✓ Easy to read
- ✓ Safe around implants
- ✓ Autoclavable



Original 3-5-7-10mm • 9006107 (3-pack) • 9006102 (6-pack) • 9006105 (12-pack)



3-6-9-12mm • 9006108 (3-pack) • 9006104 (6-pack) • 9006106 (12-pack)



1mm increments  
up to 12mm

UNC 12 • 9006120 (3-pack) • 9006121 (6-pack) • 9006122 (12-pack)

## Premier X5 Sectional Matrix System®

Reliability. Convenience. Value.

*A superior sectional matrix system is here.*

Premier X5 – an innovative, complete, high-quality, 5 component system designed to achieve accurate Class 2 composite restorations.



### 9061150 X5 Intro Kit

Contains:

- (1) Matrix Pin Holder
- (1) Matrix Ring Forceps
- (10) Matrix Rings (gray)
- (25) each

#### Wedges

- Small (yellow)
- Medium (green)
- Large (blue)

#### Matrices

- 3.5mm
- 4.5mm
- 5.5mm
- 6.5mm

### 9061150 X5 Essential Kit

Contains:

- (10) Matrix Rings (gray)
- (25) each

#### Wedges

- Small (yellow)
- Medium (green)
- Large (blue)

#### Matrices

- 3.5mm
- 4.5mm
- 5.5mm
- 6.5mm


### Reorder

-  9061190 X5 Matrix Ring (35)
-  9061170 X5 Wedge - Small (yellow) (100)
-  9061171 X5 Wedge - Medium (green) (100)
-  9061172 X5 Wedge - Large (blue) (100)
-  9061160 X5 3.5mm Matrix (25)
-  9061161 X5 4.5mm Matrix (80)
-  9061162 X5 5.5mm Matrix (80)
-  9061163 X5 6.5mm Matrix (80)

- 9061180 Matrix Pin Holder
- 9061181 Matrix Ring Forceps







**Strong**



**Contouring/Contacts**



**Compatible**



**Easy**




**Priced Right**

***There's never been a better time to make the switch to X5 and save without compromising quality results.***




**1a**  
**Placement of Matrix.** Select matrix based on the approximate occlusogingival height. Place prongs of Pin-Holder through corresponding holes in the matrix with the gold side visible to easily disengage the instrument once the matrix is in place.


When placing the matrix before the wedge, place a finger at the occlusal aspect of the matrix to hold it in place & prevent it from dislodging, or use Pin-Holder to secure matrix, proceed with wedge placement.



**1b**  
Fold the matrix tab towards you so you can insert the matrix apically into the gingival crevice & adjust matrix to the proper height to obtain ideal benefit from the contour of matrix.




**2a**  
**Placement of Premier X5™ Ring.** Use the ring forceps to grasp and slightly expand the ring. Carefully place the ring as low as possible with the tines of the ring close to the gingival margin and straddling the wedge(s).




**2b**  
When you are ready to release the ring from the Ring Forceps, place your finger and press down on the matrix tab to prevent any unwanted movement.




**2c**  
Premier X5™ Ring in the ideal position from a buccal view.



**2d**  
**Dual-Wedging:** The Premier X5™ system allows for the insertion of a second wedge, to prevent any leakage, to be placed from the opposite side without removing the ring. The second wedge should be inserted underneath the first wedge.




**2e**  
**Mesial & Distal Ring Placement:** The ring may be placed mesially and/or distally for better access such as the case with MOD restoration.



**3**  
**Dual-Mesial Ring Placement:** Premier X5™ ring may be stacked over each other for better access such as the case with MOD restoration.

**Perform Restoration:** Proceed to fill entire cavity preparation with desired composite material.




**4**  
**Premier X5™ Removal:** Remove ring using Ring Forceps, remove wedge(s) using Pin-Holder, separate matrix from composite using a thin bladed instrument by pushing the instrument into the embrasure. Fold matrix wings backwards and cure composite again from buccal and lingual surfaces. Remove matrix with Pin-Holder; grip one of the matrix holes with the Pin-Holder and gently wiggle the matrix out.

If needed, you may use another Pin-Holder on the "opposite" side of the matrix and move the matrix back and forth, or upward, until loose.


## Stop Strip

Anterior Bands with  
Integrated Stopper

- ✓ Turns a three-handed procedure into a two-handed one
- ✓ Light-cure Class III and IV restorations with ease



Unique integrated stopper  
secures the matrix band  
between adjacent teeth.



|         |                   |
|---------|-------------------|
| 9061380 | 8mm (red) (100)   |
| 9061375 | 10mm (blue) (100) |

## Sycamore Wood Wedges

The Original Anatomically Carved

- ✓ Round tip easily passes over gingiva
- ✓ Concave sides conform matrix band to tooth's anatomy
- ✓ Square end is easily inserted and removed



Wedge length in mm  
(13=13mm)  
Mini = 11mm  
UT = 13mm

### Eight color-coded sizes (400 count)

|         |                                |
|---------|--------------------------------|
| 9061109 | Mini (Light Blue)              |
| 9061106 | UT - Ultra Thin (Orange)       |
| 9061100 | 13XT - Extra Thin (Neutral)    |
| 9061101 | 13T (Green)                    |
| 9061102 | 13 (Yellow)                    |
| 9061103 | 15 (Blue)                      |
| 9061104 | 16 (Red)                       |
| 9061105 | 17 (Neutral)                   |
| 9061108 | Assorted<br>(UT, 13, 13T & 15) |

## Cure-Thru® Clear Cervical Matrices

Transparent Pre-contoured Matrix Forms

- ✓ Flexes to the exact contour of the tooth for dense restorations with excellent marginal integrity
- ✓ Compresses restorative material as you light-cure
- ✓ Helps eliminate air inhibited layer
- ✓ Reduces time to final polish


Actual Size



719 720 721 722 723 Anterior/  
Premolar Universal

### Available in boxes of 150

|         |   |
|---------|---|
| 9061012 | 719CT   |
| 9061013 | 720CT   |
| 9061014 | 721CT   |
| 9061015 | 722CT   |
| 9061016 | 723CT   |
| 9061017 | Anterior/Premolar   |
| 9061018 | Molar/Universal   |
| 9061340 | Positioning Instrument - double-ended, autoclavable   |
| 9061011 | Assorted Box - Contains: 275 matrix forms, positioning instrument, gingival retraction instrument |



## Pre-Contoured Metal Bands

- ✓ Anatomical shape achieves better contact points
- ✓ Ideal for large Class II restorations with incremental curing technique
- ✓ Finished edges ensure patient comfort
- ✓ Measures 0.035mm thickness



9061070 Box of 30

## BioGREEN™ Micro-Applicators

Bendable. Economical. Eco-Friendly.

- ✓ Available in regular, fine and super-fine sizes
- ✓ Eco-Friendly
- ✓ Gluten-free & Lint-free brush tips

3005100 Regular (green)

3005101 Fine (blue)

3005102 Super-Fine (yellow)

Contents 100 micro-applicators



## ComfortView® Lip & Cheek Retractor

- ✓ Retracts lips and cheeks for optimal buccal and gingival access
- ✓ Maximum visibility and patient comfort
- ✓ Suitable for use in diagnostic, preventative, therapeutic dental applications, in-office whitening, digital impression scanning and intra-oral photography
- ✓ Autoclavable, latex-free

9061381 Universal ComfortView Lip and Cheek Retractor  
2 per box (white)

9061383 Small ComfortView Lip and Cheek Retractor  
2 per box (blue)

9061382 ComfortView Refill Cushions  
10 per pack



## NexTemp®

Non-Eugenol Temporary Resin Cement

See page 9.



## ZR-Cem®

Self-Adhesive Universal Resin Cement


See page 10.



## Premier X5 Sectional Matrix System®

Reliability. Convenience. Value.

See page 15.



**premier®** Inspired solutions  
for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA  
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • [premierdentalco.com](http://premierdentalco.com) • [info@premierdentalco.com](mailto:info@premierdentalco.com)

©2025 Premier Dental Co. | All Rights Reserved

8610307 0325031 Rev16 PDF