

Stacks when you want,

Flows when you need.



StackFlow™

Injectable Light-Cured Composite

StackFlow™ is a versatile flowable composite with highly thixotropic behavior that offers clinicians the best of both worlds: control when you want it and adaptability when you need it.

A proprietary hydrophobic silane surface-treatment applied to nano fillers improves dispersion and strengthens the resin–filler interface, supporting excellent mechanical integrity and smooth, controlled handling.

StackFlow places precisely and polishes to a natural, long-lasting gloss. The result is a streamlined workflow and durable, aesthetic restorations across a wide range of indications.

Order through your authorized dealer.



Versatile Flowable Composite
Suitable for a wide range of indications



Thixotropic Handling
Holds its shape when dispensed yet flows on demand for efficient adaptation



Hydrophobic Silane Surface Treatment
Promotes uniform filler dispersion for enhanced strength & predictable handling



Strength Equivalent to a Conventional Composite
Delivers durability without sacrificing ease of placement



Zero-Bubble Syringe Tip
Ensures smooth placement while helping minimize bubbles and oozing

StackFlow™

Injectable Light-Cured Composite

StackFlow™ Standard Pack

Contains: 2 x 2.5g Syringes, 20 Applicator Tips

3005270 Shade: Light


3005271 Shade: Neutral

3005272 Shade: Warm





3005273 Shade: Bleached

3005275 Applicator Tips

Contains: 100 Zero-Bubble Applicator Tips
(20 Gauge)



Strong enough to replace a conventional composite!



premier® Inspired solutions for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA

Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com

8640011 0426023 Rev0 NKP10