

Introducing!

## Two Striper® EURO Diamond Burs




NEW!



### Trusted Two Striper® Performance in Euro-numbered Shapes!

Two Striper® EURO is the newest addition to the renowned Two Striper® family. This convenient offering is designed to make it easy for clinicians to select their preferred diamond shape by featuring familiar Euro-shape numbering – for use with their favorite handpiece.

*Two Striper EURO follows the industry standard in shapes and nomenclature.*



Euro-shape convenience  
*Popular European-numbered shapes simplify the selection process when desired*



Grit indication  
*Grit-indicating grooves located by the neck of the shank for easy identification*



Two Striper® Performance  
*Cuts faster, lasts longer for multiple preps*



Sharper diamonds that stay put  
*Proprietary P.B.S.® brazing process*

Distributed by


**premier**® Inspired solutions  
for daily dentistry®











Premier® Dental Products Company • 888-670-6100 • premierdentalco.com

©2026 Premier Dental | All Rights Reserved











# Two Striper® EURO Diamond Bars

5-packs  
 Available in over 50 grits/shapes.  
 \*856-018ZX available in 3-pack  
 Order through your authorized dealer.













										
	132-008	133-010	134-014	135-014	247-009	368-023	379-018	379-023	392-016	801-023
Super Coarse						2015737		2015741		
Coarse						2015735	2015738	2015740		2015745
Fine	2015725	2015727	2015729	2015731	2015733	2015736	2015739	2015742	2015744	
Very Fine	2015726	2015728	2015730	2015732	2015734			2015743		
Minor/Major/Length (mm)	0.83 / 1.05 / 3.01	1.03 / 1.25 / 4.22	1.42 / 1.65 / 6.01	1.42 / 1.65 / 9.01	0.92 / 1.15 / 3.22	- / 2.30 / 5.00	- / 1.80 / 3.40	- / 2.30 / 4.20	0.33 / 1.60 / 5.00	- / 2.30 / 1.92

										
	811-033	811-037	830-012	835-010	835-012	837-012	847-016	850-014	850-016	856-012
Super Coarse	2015747	2015749								
Coarse	2015746	2015748	2015750	2015752	2015753	2015754	2015755	2015756	2015757	2015758
Medium			2015751							
Minor/Major/Length (mm)	- / 3.30 / 4.30	- / 3.70 / 7.00	- / 1.20 / 2.70	- / 1.00 / 4.00	- / 1.20 / 4.00	- / 1.20 / 8.00	1.05 / 1.60 / 8.00	0.73 / 1.40 / 10.00	0.93 / 1.60 / 10.00	0.80 / 1.20 / 8.00

										
	856-014	856-016	856-018	856-021	856L-016	856L-018	859-010	859-014	862-012	909-040
Super Coarse		2015762	2015764		2015767	2015769				
Coarse	2015759	2015760	2015763	2015765	2015766	2015768			2015772	2015774
Medium										
Fine		2015761					2015770	2015771	2015773	
ZX			2013687*							
Minor/Major/Length (mm)	0.87 / 1.40 / 8.00	1.08 / 1.60 / 8.00	1.29 / 1.80 / 8.00	1.60 / 2.10 / 8.00	1.01 / 1.60 / 9.00	1.21 / 1.80 / 9.00	0.41 / 1.10 / 10.00	0.41 / 1.40 / 10.00	0.45 / 1.20 / 8.00	- / 4.00 / 1.52

Distributed by  
**premier** Inspired solutions for daily dentistry®



Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, USA  
 Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • E-mail: info@premierdentalco.com