The Diamond Difference



Premier®

Rotary Products

From everyday procedures to specialized needs, Premier's rotary lineup delivers the performance, precision, and efficiency you rely on.



Two Striper®
Ultra-Premium Diamond Burs

See page 3.



Priva Rotary Instruments®
Diamond Burs

See page 27.



Solo Diamond®
Pre-sterilized Diamond Burs


See page 31.

premier® Inspired solutions for daily dentistry®

The Diamond Difference

Premier's rotary lineup - from the unmatched performance of **Two Striper®** to the smooth, efficient operation of **Priva Rotary Instruments®** Diamond Burs, and the convenience of pre-sterilized **Solo Diamond®** – is designed to deliver exceptional results with every procedure.







	Two Striper® Ultra-Premium Diamond Burs	Priva Rotary Instruments® Diamond Burs	Solo Diamond® Pre-sterilized Diamond Burs
	KEY POINTS OF DIFFERENTIATION		
Overview	Setting the Industry Standard for over 50 Years, Two Striper Premium Diamond Burs are the product of world-class engineering expertise in material grinding and abrasion. The P.B.S.® Brazed bonding system maximizes the exposure of diamond-cutting surfaces.	High-performance electroplated diamond burs featuring SmartEdge® technology, designed to optimize every procedure and enhance your diamond portfolio.	Pre-sterilized, individually packaged electroplated diamond burs.
Technology	The proprietary P.B.S.® braze bonding technology features permanently attached sharp natural diamond crystals.	High-performance electroplated process featuring SmartEdge® technology – intentionally positioned, uniformly shaped engineered diamond crystals.	Standard electroplating process with engineered diamond crystals.
Features & Benefits	P.B.S.® technology features permanently attached natural diamond crystals to cut tooth structure faster and stay sharp longer. This advanced technology helps create an optimal cutting edge enhancing efficiency and performance during extended use.	SmartEdge® technology optimizes rotational accuracy by creating an innovative cutting surface – maximizing diamond-to-tooth contact. This innovative technology helps provide fast cutting, reduced vibration, efficient tooth surface prep and prevent diamond dislodgement.	Pre-sterilized diamond burs conveniently packaged in single-unit pouches providing a 5-year sterile shelf-life.
Diamond Crystals	Carefully sized natural diamond crystals ideal for fastest cutting in tooth structure.	Carefully sized and uniformly shaped engineered diamond crystals.	Engineered diamond crystals.



Two Striper®
Ultra-Premium Diamond Burs

Deliver faster, cooler, more precise cutting - with proprietary P.B.S.® braze technology.

See page 3.



Priva Rotary Instruments®
Diamond Burs

Featuring SmartEdge® technology, designed to enhance your rotary portfolio.

See page 27.



Solo Diamond®
Pre-sterilized Diamond Burs

Solo diamonds are engineered to be sharp and provide reliable results.

See page 31.

Index

Two Striper® Difference 3-4

Two Striper® EURO 5-6

TSZtech™ 7-8

TS2000™ Speed-Cut Spiral 9

ShortCut™ - Short Shank Diamonds 10

Chamfer Prep 11
Round-End Taper, Flame, Cylinder and KS System

Deep Chamfer, Guide Pin & Assorted Crown Pack 12

Shoulder Prep 13
Flat-End Taper and Cylinder

KR Modified Shoulder Prep 14
KS Cylinder and KR Taper

Bevel Shoulder Prep 14
Bevel-End Cylinder

Flame/Feather Edge Prep/Break Contact 15
Flame, Cone, Interproximal Trimmer and Break Contact

Occlusal & Lingual Reduction 16
Barrel, Egg, Football, Round, Wheel and Occlusal Contour

Veneer Prep 16
Depth Cut and Depth Marker

Inlay/Onlay 17
Flat-End Cylinder, KS Cylinder, KR Cylinder and KR Taper

Crown & Bridge Specialty 17
TGE™ (Tissue Guard End-Cutting), Safe-End, Gross Reduction

Operative 18
Pear, Amalgam Removers, Inverted Cone and Diamond Versions of Popular Carbides

MicroPrep™, MFS® Diamonds, Skooter™ 18

Two Striper® Bulk Packs 19

Planmeca Planmill 20
Milling Tools

Laboratory 21
Thin-Flex® Discs, Silverstreak™ Diamond Model Trimming Wheel

Laboratory 22
HP Lab Diamonds, TSZTech™

Finishing & Polishing 23
Two Striper® Perforated Diamond Finish Strips

Finishing & Polishing 24
Two Striper® Diamond Finishing Strips, Two Striper® 200T2

Diamond Twist SCL™/SCO™ 25
Ceramic/Composite Polishing

Finishing & Polishing 26
Luminescence Plus™ and Luminescence®,

Poli-Pro Disks™ 26
Composite Polishing

Priva Rotary Instruments® 27-30

Solo Diamond® 31-34

Two Striper®
AN abrasive technology BRAND

All Diamonds are Not Created Equal


Experienced dentists know that not all diamond burs are created equal. That's why they've preferred Two Striper® for over 50 years.

Two Striper® delivers faster, cooler, more precise cutting than electroplated brands. What's more, Two Striper® stays sharp longer for multiple preps – and that means a more efficient workflow, less chair time, and happier patients.


Two Striper® is the only bur made from a select grade of natural diamonds with the proprietary P.B.S.® brazing process – the best bonding system in the industry. You can literally feel the difference.

- ✓ Cuts faster – *less chair time, happier patients*
- ✓ Sharper diamonds that stay put – *Proprietary P.B.S.® brazing process*
- ✓ Lasts longer for multiple preps – *Use fewer burs than 'bargain' diamonds*

P.B.S.® - Bonded Two Striper



Electroplated Brands




Cuts faster, stays sharp longer!

Independent study proves Two Striper superiority[2]

PRODUCT EVALUATION SNAPSHOT

CRITERIA BASED ON AVERAGE SCORE (OUT OF 5)



TWO STRIPER DIAMONDS
DISTRIBUTED BY PREMIER DENTAL

FINAL SCORE
AVERAGE OF SECTIONS A AND B

4.7


Dental PRODUCT SHOPPER BEST PRODUCT

Read Evaluation


Visit:
premierdentalco.com/ts-eval

1. Compendium, February 1997. An Evaluation of Operator Preference of Diamond Burs in Coronal Tooth Preparation. / 2. J Prosthet Dent 2014;111(1):53-63. PubMed.gov ID# 24246370. Used with permission. / 3. ZX Grit - Data on file. / 4. Features proprietary P.B.S.® brazing process. / 5. An Independent Non-Profit Dental Education and Product Testing Foundation Newsletter, Sept. 2001. / 6. An Independent Non-Profit Dental Education and Product Testing Foundation Newsletter, May 2002. / 7. American J. Dentistry, Feb. 2000. Post adjustment polishing of CAD/Cam ceramic with Luminescence Diamond Gel. W.J. Finger. / ®Premier is a registered trademark of Premier Dental Products Company. / Two Striper®, P.B.S.® and MFS® are registered trademarks and, LTech™ Skooter™, Crown Cut™, TGE™, MicroPrep™, TS2000™, GCP™, ShortCut™, TSZtech™ are trademarks of the manufacturer, Abrasive Technology, LLC. / ®Brasseler, Neo, SSW Piranha and SS White are not trademarks of Premier.

Distributed by premier®



"Because advancement rate and load are interrelated (in general, the application of higher loads should lead to more rapid advancement), combining the two parameters provides an overall indicator for the ease of advancement of the rotary cutting instrument inside the substrate. Indeed, a significant difference was found between the Two Striper and the other rotary cutting instrument groups."



Rate of advancement of 5 different groups shown for 8 cuts. Rate showed decreasing trend with increase in cut number when averaged over all groups.

Two Striper® — SS White — Neo — SSW Piranha — Brasseler

The P.B.S.® Braze Difference is a bonding system that maximizes the exposure of diamond-cutting surfaces. Crystals, placed evenly and precisely in a uniform matrix, are fused permanently to a surgical-grade stainless-steel shank.

Electroplated diamonds cut less consistently than Two Striper® diamonds[1] and have a shorter product life. That's because electroplated tools are trapped mechanically in inconsistent plating layers that can expose either too much or too little of the particle. Some particles are attached loosely and are prone to premature dislodgement. Other particles are buried completely in the plating layer and are unavailable to cut tooth structure, creating unnecessary heat and friction.

The combination of permanent crystal attachment, precise placement at the top and maximum crystal exposure makes Two Striper® cut faster, cooler and stay sharp longer.



Better Diamonds... Better Cut, Control & Consistency

The quality of dental diamond instruments starts with the quality of the diamond crystals themselves. And with Two Striper®, the difference is crystal clear. We use natural, virgin diamond crystals - which have more corners and angles than synthetics - for superior abrasion. The individual crystals are selected by exact size and shape to ensure high diamond exposure and a uniform cutting surface. We manufacturer to exacting standards and rigorously test to ensure the highest levels of quality.

Feeling is Believing

Two Striper® burs just feel right. Carefully crafted to create a more consistent surface for each bur, Two Striper® reduces unwanted vibration, skipping and bur chatter. This allows for a smoother surface finish. The result is faster, more efficient cutting, less chair time and happier patients.

Color-Coded Stripes On Shank Indicate Grit


FG Shapes:

- Super Coarse (150µ)
- Coarse (120µ)
- Medium (85µ)
- Fine (60µ)
- Very Fine (45µ)

MFS® Finishing Diamonds:

- MF1 (45µ)
- MF2 (20µ)
- MF3 (10µ)

Instrument Coding




Two Striper®
AN abrasive technology BRAND

Two Striper® EURO

Diamond Burs






Trusted Two Striper® Performance in Euro-numbered Shapes!

Two Striper® EURO is the newest addition to the renowned Two Striper® family. This convenient offering is designed to make it easy for clinicians to select their preferred diamond shape by featuring familiar Euro-shape numbering – for use with their favorite handpiece.

Two Striper EURO follows the industry standard in shapes and nomenclature.



Euro-shape convenience

Popular European-numbered shapes simplify the selection process when desired

Grit indication

Grit-indicating grooves located by the neck of the shank for easy identification

Two Striper® Performance

Cuts faster, lasts longer for multiple preps

Sharper diamonds that stay put


Proprietary P.B.S.® brazing process


Distributed by premier®


5-packs
Available in over 50 grits/shapes.
*856-018ZX available in 3-pack




Order through your authorized dealer.

FG Shank Length



	 132-008	 133-010	 134-014	 135-014	 247-009	 368-023	 379-018	 379-023	 392-016	 801-023
Super Coarse						2015737		2015741		
Coarse						2015735	2015738	2015740		2015745
Fine	2015725	2015727	2015729	2015731	2015733	2015736	2015739	2015742	2015744	
Very Fine	2015726	2015728	2015730	2015732	2015734			2015743		
Minor/Major/ Length (mm)	0.83 / 1.05 / 3.01	1.03 / 1.25 / 4.22	1.42 / 1.65 / 6.01	1.42 / 1.65 / 9.01	0.92 / 1.15 / 3.22	– / 2.30 / 5.00	– / 1.80 / 3.40	– / 2.30 / 4.20	0.33 / 1.60 / 5.00	– / 2.30 / 1.92

	 811-033	 811-037	 830-012	 835-010	 835-012	 837-012	 847-016	 850-014	 850-016	 856-012
Super Coarse	2015747	2015749								
Coarse	2015746	2015748	2015750	2015752	2015753	2015754	2015755	2015756	2015757	2015758
Medium			2015751							
Minor/Major/ Length (mm)	– / 3.30 / 4.30	– / 3.70 / 7.00	– / 1.20 / 2.70	– / 1.00 / 4.00	– / 1.20 / 4.00	– / 1.20 / 8.00	1.05 / 1.60 / 8.00	0.73 / 1.40 / 10.00	0.93 / 1.60 / 10.00	0.80 / 1.20 / 8.00

	 856-014	 856-016	 856-018	 856-021	 856L-016	 856L-018	 859-010	 859-014	 862-012	 909-040
Super Coarse		2015762	2015764		2015767	2015769				
Coarse	2015759	2015760	2015763	2015765	2015766	2015768			2015772	2015774
Medium										
Fine		2015761					2015770	2015771	2015773	
ZX			2013687*							
Minor/Major/ Length (mm)	0.87 / 1.40 / 8.00	1.08 / 1.60 / 8.00	1.29 / 1.80 / 8.00	1.60 / 2.10 / 8.00	1.01 / 1.60 / 9.00	1.21 / 1.80 / 9.00	0.41 / 1.10 / 10.00	0.41 / 1.40 / 10.00	0.45 / 1.20 / 8.00	– / 4.00 / 1.52

Two Striper®
AN abrasive technology BRAND

TSZTECH™

NEW!

Removing & Adjusting Zirconia Restorations

TSZtech™ ZX Diamond Burs with LTech™ Technology

Zirconia removal has long been a challenge, but we're changing that! LTech™ Technology offers a breakthrough in zirconia removal, featuring an advanced contact layer design that prevents diamond dislodgement, ensuring sustained cutting performance.


Paired with tough, engineered crystals that remain sharp and intact even when cutting ultra-hard materials, the TSZtech™ ZX with LTech™ Technology streamlines zirconia removal procedures, saving valuable time and enhancing both patient and dentist experiences.



- 30% faster cutting[3], longer lasting than other zirconia burs
- Consistent diamond layer height for maximum cutting efficiency
- Reduced vibration for a smoother, more controlled operation


---

Proprietary LTech™ Technology



Precision plating creates a more consistent diamond layer height and a more secure attachment layer, resulting in faster cutting and improved durability.

Traditional Electroplated Technology



Inconsistent diamond layer results in poor prep surface quality as the edges are gouging and crystals are plucked away.


Distributed by premier®

The full line of TSZtech™ Diamond Burs optimize zirconia procedures through advanced grit selection, allowing professionals to achieve desired results with ease.

ZX Featuring LTech™ Proprietary Technology, this bur is designed for fast, efficient zirconia removal, reducing chair time and improving the experience for both patient and dentist.

ZM Ideal for efficient removal and moderate amounts of adjustments, preparing surfaces for the ZF bur.[4]

ZF Perfect for small, precise final adjustments of zirconia surfaces, before final zirconia polishing.[4]



	260.8	285.5	703.8KR	125	770.8
ZX		2013680		2013681	2013682
ZM	2015670	2015671	2015672	2015673	2015674
ZF	2015660	2015661	2015662	2015663	2015664



	550.8	125N	L260.8	L125	L770.8
ZX	2013683	2013684		2013685	2013686
ZM	2015675	2015676	2015677	2015678	2015679
ZF	2015665	2015666	2015667	2015668	2015669


TSZtech ZF and ZM diamond burs are available in 5-packs.
TSZtech ZX diamond burs are available in 3-packs.

TSZTECH™

Laboratory Diamonds

A New Diamond Technology for Today's Advanced Materials - Lithium Disilicate and Zirconia

TSZtech Laboratory diamond burs are available in 1-packs.



	X590.12	X769.10
ZF	2012848	2012849
ZM	2012850	2012851

Intra-oral Adjustments



285.5ZF **Football**

Atraumatic Crown and Bridge Removal



770.8ZX **Round-End Taper**

Endodontic Access




125ZM **Round**


Two Striper®
AN abrasive technology BRAND

"Combined lab and clinical data showed highest rated and best performing diamonds were the TS2000™ and 770.8C Two Striper® Ultra-Premium Diamonds...."[5]

TS2000™




	2000.7	2000.8	2000.9	2000.10	2001.7	2001.8	2001.9	2001.10	2002.7	2002.8	2002.9	2002.10	2001.8KR	2001.8 S	2004.9 S	2005.8 S
	Flat-End Taper												KR Taper	Short Shank		
Coarse	2000007	2000008	2000009	2000010	2000017	2000018	2000019	2000020	2000027	2000028	2000029	2000030	2000021	2000143	2000141	2000142
Minor/Major Diameter	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.3 / 1.8*	*1.3 / 1.8*	*1.3 / 1.8*	*1.3 / 1.8*	*1.0/1.6*			




	2003.7	2003.8	2003.9	2003.10	2004.7	2004.8	2004.9	2004.10	2005.7
	Round-End Taper								
Coarse	2000037	2000038	2000039	2000040	2000047	2000048	2000049	2000050	2000057
Minor/Major Diameter	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.3 / 1.8*



	2005.8	2005.9	2005.10	2009.7	2009.8	2009.9	2010.8	2011.8	2014.5
	Round-End Taper			Flame					Football
Coarse	2000058	2000059	2000060	2000097	2000098	2000099	2000108	2000118	2000145
Minor/Major Diameter	*1.3 / 1.8*	*1.3 / 1.8*	*1.3 / 1.8*	*1.4*	*1.4*	*1.4*	*1.6*	*1.8*	*2.3*




TS2000 Flat-End Cylinder	2013.7	2013.8	2013.10	2013.8KR	2016.8	2016.10	2006.7GC	2006.8GC	2006.9GC
	Flat-End Taper			KR Cylinder	Bevel-End Cylinder		Gold-Plated Curettage (GCP™)		
Coarse	2000136	2000138	2000140	2000139	2000168	2000170	2000067	2000068	2000069
Minor/Major Diameter	*1.5*	*1.5*	*1.5*	*1.6*	*1.4*	*1.8*	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*



	2006.10GC	2007.7GC	2007.8GC	2007.9GC	2007.10GC	2008.7GC	2008.8GC	2008.9GC	2008.10GC
	Gold-Plated Curettage (GCP™)								
Coarse	2000070	2000077	2000078	2000079	2000080	2000087	2000088	2000089	2000090
Minor/Major Diameter	*0.9 / 1.4*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.3 / 1.8*	*1.3 / 1.8*	*1.3 / 1.8*	*1.3 / 1.8*



Two Diamonds In One




"Fastest Brand"[6]

Less chair-time for you and your patient.


TS2000™ shapes are listed on pages 7, 8, 9, 10, 12 &13.

TS2000™ provide rapid gross reduction and a finished margin without changing burs.

A finished margin is created as shown. You can easily create a chamfer, bevel, or flat shoulder by selecting one of the more than forty shapes available.



Actual Preparation




Close-up occlusal view of finished margin

Distributed by premier®

ShortCut™ *Short Shank Diamonds*

Available in 5-packs only.

Assorted 5-Pack 2015502



	260.8 S	285.5 S	702.8 S	767.8 S	770.8 S
Coarse	2015901	2015900	2015909	2015907	2015905

microprep™




	MP89 S
Medium	2015968

TS2000™




	2001.8 S	2004.9 S	2005.8 S
Coarse	2000143	2000141	2000142

Gold-Plated Curettage (GCP™)



	252.SB S	253.SB S	254.SB S
Medium	2015920	2015921	2015922

Round-End Taper




	767.8 S	770.8 S	780.9 S	782.8 S
Coarse	2015907	2015905	2015906	2015928

Football




	285.5 S
Coarse	2015900
Fine	2015923

Flame



	260.8 S	261.8 S
Coarse	2015901	2015925
Fine	2015924	

Flat-End Taper




	702.8 S	703.8 S
Coarse	2015909	2015908

KR Taper



	703.8 KR S
Coarse	2015926

KR Cylinder



	521.8 S
Coarse	2015990

Amalgam Removers



	1250 S
Medium	2015912

Bevel-End Cylinder



	250.8 S
Coarse	2015902

Flat-End Cylinder



	515.7 S
Coarse	2015910

KS Cylinder (Round Edge)




	575.7 S
Coarse	2015911

Interproximal Trimmer



	201.3 S
Fine	2015927

Operative




	2 S	6 S	169L S	245 S	330 S	556 S	557 S
Medium	2015917	2015918	2015919	2015913	2015914	2015916	2015915

Skooter™



	1 S
Medium	T1 S - 2015903
Very Fine	F1 S - 2015904



Shorter Shank For Easier Access




S Available as a short shank diamond.
See page 10 for complete list and ordering information.


B Available in Bulk Packs.
See page 19 for complete list and ordering information.

Two Striper®
AN abrasive technology BRAND


Chamfer Prep

- ✓ Round-end tapered or cylindrical diamonds create a chamfered finish line
- ✓ Recommended prep design for full coverage crowns and veneers
- ✓ Recommended prep design for PFM crowns and some ceramic crowns







Round-End Taper	747.6	767.5	767.7	L767.7	767.8	767.9	770.5	770.7	770.8	L770.8	770.9	770.10	776.4	777.8	780.4	780.7
Super Coarse*			2015706			2015705			2015701		2015716	2015703				
Coarse	2015058	2015452	2015453	2015652	2015454	2015455	2015456	2015458	2015457	2015650	2015460	2015459	2015461	2015972	2015462	2015463
Medium														2015930		
Fine			2015483	2015653	2015484	2015485			2015487	2015651		2015489	2015491			
Minor/Major Diameter	*0.8 / 1.2*	*1.1 / 1.6*	*1.1 / 1.8*	*1.1 / 1.8*	*1.0 / 1.8*	*0.9 / 1.8*	*1.3 / 1.8*	*1.2 / 1.8*	*1.2 / 1.8*	*1.3 / 1.8*	*1.1 / 1.8*	*1.1 / 1.8*	*2.1 / 2.6*	*0.5 / 1.0*	*1.7 / 2.1*	*1.4 / 2.1*



Round-End Taper	780.8	780.9	781.6	781.8	781.10	782.6	782.8	782.10	783.6	783.8	783.10	784.6	784.8	784.10	785.4	785.7
Super Coarse*		2015702					2015709	2015713			2015721			2015715		
Coarse	2015464	2015467	2015465	2015478	2015466	2015423	2015424	2015425	2015426	2015427	2015468	2015429	2015430	2015431	2015470	2015469
Fine	2015494	2015497		2015498	2015496	2015433	2015434	2015435			2015438		2015440	2015441		2015499
Very Fine																2015529
Minor/Major Diameter	*1.4 / 2.1*	*1.4 / 2.1*	*0.8 / 1.1*	*0.8 / 1.2*	*0.8 / 1.3*	*1.1 / 1.4*	*1.1 / 1.5*	*1.1 / 1.6*	*1.3 / 1.6*	*1.3 / 1.7*	*1.3 / 1.8*	*1.6 / 1.9*	*1.6 / 2.0*	*1.6 / 2.1*	*2.2 / 2.6*	*1.5 / 2.4*




Round-End Taper	790.8	799.6.5	799.11
Super Coarse*	2015712		
Coarse	2015471	2015473	2015475
Medium		2015480	
Fine	2015501	2015503	2015505
Very Fine			2015535
Minor/Major Diameter	*2.0 / 2.6*	*0.8 / 1.5*	*0.8 / 1.8*



Round-Edge Cylinder	585.5	587.8
Coarse	2015305	2015304
Minor/Major Diameter	*1.3*	*1.8*



Parallel Cylinder Round-End	550.8	573.6	588.10
Coarse	2015311	2015310	2015312
Minor/Major Diameter	*1.4*	*1.2*	*1.3*




KS System	574.7	575.7	585.8	586.8	287.4	276.6.2	780.9	866
Coarse	2015953	2015302	2015303	2015307	2015981	2015898	2015467	2015546
Fine	2015954	2015322	2015306		2015982		2015497	
Minor/Major Diameter	*1.0*	*1.2*	*1.4*	*1.6*	*2.3*	*6.2*	*1.4 / 2.1*	*5.7*


TS2000™




TS2000 Round-End Taper	2003.7	2003.8	2003.9	2003.10	2004.7
Coarse	2000037	2000038	2000039	2000040	2000047
Minor/Major Diameter	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*1.1 / 1.6*



TS2000 Round-End Taper	2004.8	2004.9	2004.10	2005.7	2005.8
Coarse	2000048	2000049	2000050	2000057	2000058
Minor/Major Diameter	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.3 / 1.8*	*1.3 / 1.8*



TS2000 Round-End Taper	2005.9	2005.10
Coarse	2000059	2000060
Minor/Major Diameter	*1.3 / 1.8*	*1.3 / 1.8*





TS2000 Flame	2009.7	2009.8	2009.9	2010.8	2011.8
Coarse	2000097	2000098	2000099	2000108	2000118
Minor/Major Diameter	*1.4*	*1.4*	*1.4*	*1.6*	*1.8*

Distributed by premier®


Deep Chamfer / Curettage / Sub-gingival Preparation

- ✓ Provides clearly defined margins for CAD/CAM preps
- ✓ Popular for gingival curettage during crown preparation when margins are sub-gingival






Curettage "K"	763.8	763.10	764.8	764.10
Coarse	2015977	2015979	2015978	2015980
Minor/Major Diameter	*0.7 / 1.2*	*0.6 / 1.2*	*0.8 / 1.4*	*0.7 / 1.4*



Modified Beveled Cylinder	510.6	510.8	511.8	511.10
Coarse	2015973	2015974	2015975	2015976
Minor/Major Diameter	*1.0*	*1.0*	*1.2*	*1.2*


GCP™ Cutting Length

253.S**A**	= 7mm
253.S**B**	= 8mm
253.S**C**	= 10mm



Gold-Plated Curettage (GCP™)	252.SA	252.SB	252.SC	253.SA	253.SB	253.SC	254.SA	254.SB	254.SC	257.SB	258.SB
Medium	2015603	2015605	2015607	2015609	2015611	2015613	2015615	2015617	2015619	2015623	2015627
Minor/Major Diameter	*0.4 / 1.4*	*0.4 / 1.4*	*0.4 / 1.4*	*0.4 / 1.6*	*0.4 / 1.6*	*0.4 / 1.6*	*0.4 / 1.8*	*0.4 / 1.8*	*0.4 / 2.2*	*0.4 / 2.0*	*0.4 / 2.2*


TS2000™




Gold-Plated TS2000 Curettage (GCP™)	2006.7GC	2006.8GC	2006.9GC	2006.10GC	2007.7GC	2007.8GC	2007.9GC	2007.10GC	2008.7GC	2008.8GC	2008.9GC	2008.10GC
Coarse	2000067	2000068	2000069	2000070	2000077	2000078	2000079	2000080	2000087	2000088	2000089	2000090
Minor/Major Diameter	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.3 / 1.8*	*1.3 / 1.8*	*1.3 / 1.8*	*1.3 / 1.8*

All-Ceramic Guide Pin

- ✓ Guide pin 0.5mm chamfer protects biological width
- ✓ Controlled tooth preparation with uniform cutting depth
- ✓ Eliminate lip or "j" margin defects






Guide Pin	*Chamfer* SE740.8	*Modified Chamfer* SE259.8
Medium	2015966	2015964
Fine	2015967	2015965
Minor/Major Diameter	1.5 / 2.1	1.6 / 2.1

Assorted Crown Pack

(Chamfer/Shoulder)

Assorted 5-pack with our most popular Crown and Bridge Diamonds

2015500



770.8C **703.8C** **250.8C** **260.8C** **285.5C**

(S) Available as a short shank diamond.
See page 10 for complete list and ordering information.


*Two Striper® BLACK diamond.

(B) Available in Bulk Packs.
See page 19 for complete list and ordering information.

Two Striper®
AN abrasive technology BRAND


Shoulder Prep/Butt-Joint Margin

- ✓ Choose flat-end tapered or flat-end cylindrical diamonds to create shoulder preparations with a 90˚ margin






Flat-End Taper	700.5	700.6	700.8	700.9	700.11	701.5	701.7	L701.7	701.9	L701.9	702.8	703.8	703.9	703.10	704.9	708.3
Super Coarse*											2015708	2015704				
Coarse	2015332	2015333	2015334	2015335	2015337	2015338	2015339	2015656	2015341	2015657	2015343	2015345	2015346	2015347	2015348	2015350
Medium																2015445
Fine				2015365	2015367		2015369		2015371		2015373	2015375		2015377		
Very Fine				2015395	2015397						2015403					
Minor/Major Diameter	*0.6 / 1.1*	*0.7 / 1.2*	*0.6 / 1.8*	*0.6 / 1.4*	*0.7 / 1.8*	*0.8 / 1.2*	*1.1 / 1.8*	*1.1 / 1.8*	*0.9 / 1.8*	*0.9 / 1.8*	*0.9 / 1.4*	*1.2 / 1.8*	*1.1 / 1.8*	*1.1 / 1.8*	*1.3 / 1.8*	*1.3 / 2.1*



Flat-End Taper	708.4	721.6	721.8	721.10	722.6	722.8	722.10	723.6	723.8	723.10	724.8	724.10	725.10	726.10	727.10
Coarse	2015351	2015354	2015355	2015356	2015357	2015358	2015359	2015255	2015256	2015257	2015259	2015260	2015446	2015448	2015450
Fine				2015386			2015389		2015262				2015447		
Minor/Major Diameter	*1.3 / 2.0*	*0.8 / 1.1*	*0.8 / 1.2*	*0.8 / 1.3*	*1.0 / 1.4*	*1.0 / 1.5*	*1.0 / 1.6*	*1.3 / 1.6*	*1.3 / 1.7*	*1.3 / 1.8*	*1.5 / 2.0*	*1.5 / 2.1*	*0.9 / 1.8*	*1.1 / 2.0*	*1.4 / 2.3*




Flat-End Taper	797.11	798.10
Coarse		2015442
Medium	2015039	
Fine	2015049	2015443
Minor/Major Diameter	*0.6 / 1.2*	*0.5 / 1.5*




Flat-End Cylinder	513.4	513.5	514.3	514.4	514.5	514.7	515.5	515.7	515.8	516.7	520.4
Coarse	2015271	2015272	2015273	2015274	2015276	2015278	2015275	2015277	2015279	2015985	2015281
Medium	2015288										
Fine			2015293	2015294			2015295	2015297			
Minor/Major Diameter	*1.0*	*1.0*	*1.1*	*1.1*	*1.1*	*1.1*	*1.2*	*1.2*	*1.3*	*1.6*	*1.4*

TS2000™



TS2000 Flat-End Taper	2000.7	2000.8	2000.9	2000.10	2001.7	2001.8	2001.9	2001.10	2002.7	2002.8	2002.9	2002.10
Coarse	2000007	2000008	2000009	2000010	2000017	2000018	2000019	2000020	2000027	2000028	2000029	2000030
Minor/Major Diameter	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*0.9 / 1.4*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.1 / 1.6*	*1.3 / 1.8*	*1.3 / 1.8*	*1.3 / 1.8*	*1.3 / 1.8*




TS2000 Flat-End Cylinder	2013.7	2013.8	2013.10
Coarse	2000136	2000138	2000140
Minor/Major Diameter	*1.5*	*1.5*	*1.5*

Distributed by premier®


KR Modified Shoulder Prep

- ✓ Creates a margin with rounded internal angles
- ✓ Flat-end tapered or parallel diamond instruments featuring a rounded (beveled) edge





KR Cylinder	515.6	521.4	521.6
Coarse	2015986	2015987	2015989
Fine			2015988
Minor/Major Diameter	*1.2*	*1.4*	*1.4*



KR Taper	703.8 KR	708.4 KR	712.3 KR	722.8 KR	723.6 KR
Super Coarse*	2015718				
Coarse	2015352	2015951	2015992	2015360	2015996
Fine	2015353	2015949	2015995	2015658	2015997
Minor/Major Diameter	*1.2 / 1.8*	*1.3 / 2.0*	*2.0 / 2.2*	*1.0 / 1.5*	*1.0 / 1.5*




TS2000 KR Cylinder	2013.8KR
Coarse	2000139
Minor/Major Diameter	*1.6*
	837KRT-016




TS2000 KR Taper	2001.8KR
Coarse	2000021
Minor/Major Diameter	*1.0/1.6*
	847KRT-016

Bevel Shoulder Prep

- ✓ Recommended for porcelain-fused-to-metal or cast gold crowns
- ✓ Establish a 45˚ - 60˚ angle
- ✓ Place bevel finish line





Bevel-End Cylinder	244.10	248.8	250.8	250.9	250.11	251.8	255.8
Coarse	2015064	2015065	2015067	2015066	2015068	2015069	2015073
Medium							
Fine		2015085	2015087	2015086		2015089	2015093
Minor/Major Diameter	*0.5 / 1.3*	*0.4 / 1.1*	*0.4 / 1.2*	*0.4 / 1.2*	*0.4 / 1.2*	*0.4 / 1.2*	*0.4 / 1.8*



TS2000 Bevel-End Cylinder	2016.8	2016.10
Coarse	2000168	2000170
Minor/Major Diameter	*1.4*	*1.8*

(S) Available as a short shank diamond.
See page 10 for complete list and ordering information.


(B) Available in Bulk Packs.
See page 19 for complete list and ordering information.


*Two Striper® BLACK diamond.

Two Striper®
AN abrasive technology BRAND

Flame/Feather Edge Prep/Break Contact

- ✓ Can be used for full-contour Zirconia
- ✓ Popular design for cast-gold crowns
- ✓ Margin design for tight anterior crowns and tipped or perio-involved teeth
- ✓ Thin-shaped diamond for breaking interproximal contact






Flame	242.6	243.6	253.8	253.10	256.8	256.9	260.3	260.6.5	260.8	L260.8	260.10	261.8	262.6.5
Super Coarse*				2015719									
Coarse	2015059	2015061	2015070	2015071	2015062	2015075	2015121	2015122	2015123	2015133	2015125	2015124	2015983
Medium									2015956				
Fine	2015060		2015090	2015091		2015095		2015152	2015153	2015163	2015155	2015154	
Very Fine								2015182	2015183		2015185		
Minor/Major Diameter	*0.9/1.4*	*1.2/1.8*	*1.0/1.6*	*1.0/1.6*	*1.2/1.8*	*1.2/1.8*	*0.9/1.0*	*0.4/1.3*	*0.4/1.4*	*0.4/1.4*	*0.5/1.4*	*0.4/1.1*	*1.0/1.3*




Flame	262.8	263.8	265.6.5	265.8	L265.8	270.5	270.6.5	270.9	275.9	290.2
Coarse	2015126	2015127	2015129	2015131	2015654	2015135	2015137	2015139	2015955	
Medium										2015148
Fine		2015157	2015159	2015161	2015655	2015165	2015167	2015169		
Minor/Major Diameter	*1.1/1.7*	*0.4/1.6*	*0.5/1.4*	*0.4/1.4*	*0.4/1.4*	*0.5/1.6*	*0.4/1.9*	*0.5/1.9*	*0.4/2.2*	*1.6*




Interproximal Trimmer	200.3	201.3
Fine	2015042	2015043
Very Fine	2015052	2015053
Minor/Major Diameter	*0.3/0.8*	*0.3/0.8*

Break Contact/ Interproximal Reduction



Cone	207.7	207.10	209.6	209.10	608.9	797.11
Coarse	2015036	2015149	2015035	2015037	2015038	
Medium						2015039
Fine	2015046		2015045	2015047	2015048	2015049
Very Fine				2015057		
Minor/Major Diameter	*0.4/1.4*	*0.5/1.9*	*0.4/1.6*	*0.5/2.3*	*0.4/1.4*	*0.6/1.2*




Break Contact	203.5	700.9	700.11	777.8	798.10	799.11
Coarse		2015335	2015337	2015972	2015442	2015475
Medium				2015930		
Fine	2015054	2015365	2015367		2015443	2015505
Very Fine		2015395	2015397			2015535
Minor/Major Diameter	*0.3/1.5*	*0.6/1.4*	*0.7/1.8*	*0.5/1.0*	*0.5/1.5*	*0.8/1.8*

Distributed by premier®


Occlusal & Lingual Reduction (Bulk)

- ✓ An anatomically prepared occlusal surface provides adequate clearance and allows for uniform thickness of the restoration
- ✓ Football and wheel-shape diamonds reproduce lingual concavity in crown preparation





Football	283.4	285.5	290.4	292.3
Super Coarse*	2015710	2015700		
Coarse	2015142	2015143	2015144	2015145
Fine	2015172	2015173	2015174	2015175
Very Fine	2015202	2015203	2015204	2015205
Minor/Major Diameter	*2.5*	*0.7 / 2.3*	*0.7 / 1.9*	*0.4 / 1.6*



Depth - mm	1.0	1.5	1.5 "L" 1.5mm longer shank	1.7	2.0	2.0 "L" 1.5mm longer shank	2.4	2.8
Round	**115**	**120**	**L120**	**125**	**130**	**L130**	**135**	**138**
Super Coarse*							2015717	
Coarse	2015003	2015005	2015055	2015006	2015008	2015078	2015007	2015009
Fine	2015013	2015015		2015016	2015018		2015017	
Very Fine		2015025						
Minor/Major Diameter	*1.2*	*1.6*	*1.6*	*1.8*	*2.3*	*2.3*	*2.6*	*3.0*



Barrel	230	234	240
Super Coarse*	2015707	2015720	
Coarse	2015063	2015076	2015077
Fine	2015083		
Minor/Major Diameter	*1.3 / 3.4*	*1.5 / 4.1*	*1.2 / 4.8*




KS5	276.6.2
Coarse	2015898
Minor/Major Diameter	*3.1*



Egg	287.4
Coarse	2015981
Fine	2015982
Very Fine	2015998
Minor/Major Diameter	*2.3*



Occlusal Contour	204.3.5
Fine	2015952
Minor/Major Diameter	*3.10*




Wheel	860	862	863	866
Super Coarse*			2015714	
Coarse	2015543	2015544	2015545	2015546
Fine	2015553		2015555	
Minor/Major Diameter	*3.5*	*4.1*	*4.7*	*5.7*





TS2000 Football	2014.5
Coarse	2000145
Minor/Major Diameter	*2.3*

Veneer Prep



D.C. Depth-Cut	DC.5	DC.75	DC1.0	DC1.5
	Depth of cut (in mm): 0.5	0.75	1.0	1.5
Coarse	2015776	2015785	2015802	2015803
Minor/Major Diameter	*1.2*	*1.2*	*1.2*	*1.2*





Depth Marker	DC0.3	DC0.5
Medium	2015970	2015971
Minor/Major Diameter	*1.6*	*2.1*


(S) Available as a short shank diamond.
See page 10 for complete list and ordering information.

*Two Striper® BLACK diamond.

(B) Available in Bulk Packs.
See page 19 for complete list and ordering information.

Two Striper®
AN abrasive technology BRAND

Inlay/Onlay



- ✓ Flat-end diamonds with a rounded corner create rounded internal angles for ceramic/composite inlay and onlay preparation
- ✓ Tapered instruments establish proper divergence of axial walls

Flat-End Cylinder	513.4	513.5	514.3	514.4	514.5	514.7	515.5	L515.5 "L" 1.5mm longer shank	515.7 S	515.8	516.7	520.4
Super Coarse*									2015711			
Coarse	2015271	2015272	2015273	2015274	2015276	2015278	2015275	2015797	2015277	2015279	2015985	2015281
Medium	2015288											
Fine			2015293	2015294			2015295		2015297			
Minor/Major Diameter	*1.0*	*1.0*	*1.1*	*1.1*	*1.1*	*1.1*	*1.2*	*1.2*	*1.2*	*1.3*	*1.6*	*1.4*

KS Cylinder (Round Edge)	KS0 574.7	KS1 S 575.7	585.5	KS2 585.8	KS3 586.8	587.8
Coarse	2015953	2015302	2015305	2015303	2015307	2015304
Fine	2015954	2015322		2015306		
Minor/Major Diameter	*1.0*	*1.2*	*1.3*	*1.4*	*1.6*	*1.8*

Modified Shoulder

KR Cylinder	515.6	521.4	521.6	521.8 S
Coarse	2015986	2015987	2015989	2015990
Fine			2015988	
Minor/Major Diameter	*1.2*	*1.4*	*1.4*	*1.4*

Modified Shoulder

KR Taper	703.8 KR	708.4 KR	712.3 KR	722.8 KR	723.6 KR
Coarse	2015352	2015951	2015992	2015360	2015996
Fine	2015353		2015995	2015658	2015997
Minor/Major Diameter	*1.2 / 1.8*	*1.3 / 2.0*	*2.0 / 4.0*	*1.0 / 1.5*	*1.0 / 1.5*

TS2000™

TS2000 Flat-End Cylinder	2013.7	2013.8	2013.10
Coarse	2000136	2000138	2000140
Minor/Major Diameter	*1.5*	*1.5*	*1.5*

TS2000 KR Cylinder	2013.8KR
Coarse	2000139
Minor/Major Diameter	*1.6*

Modified Shoulder

Specialty

TGE™

- Refine shoulder
- Safe-sided

TGE™ (Tissue Guard End-Cutting)	1.0	1.2	1.4	1.6
Fine	2015836	2015838	2015840	2015842
Very Fine		2015839	2015841	
Minor/Major Diameter	*1.0*	*1.2*	*1.4*	*1.6*

Endo

Endo*	1C RA	291.4 RA
Coarse		2013197
Medium	2013194	

**Available in single-pack only.*

RA Latch Surgical 2.35mm x 26mm

Safe-End

Safe-End	S79
Fine	2015835
Minor/Major Diameter	*1.0/1.6*


Gross Reduction

Crown Cut™	SC5	SC8	SC10	ST6	ST8	ST11
Coarse	2015573	2015575	2015577	2015581	2015583	2015585
Minor/Major Diameter	*1.7*	*1.7*	*1.8*	*1.4/1.8*	*1.4/1.8*	*1.4/1.8*

Distributed by premier®


Operative

- ✓ Featuring diamond versions of popular carbides
- ✓ Cut cooler and faster with less patient trauma
- ✓ Inverted cones for removing defective restorations






Inverted Cone	310.1	315.1.75	317.4	318.5	320.2	324.1	390.3
Coarse	2015213	2015790	2015984	2015991	2015217	2015219	2015931
Fine	2015233						
Minor/Major Diameter	*0.9 / 1.3*	*1.1 / 1.6*	*1.6*	*1.8*	*1.3 / 1.7*	*1.1 / 2.1*	*1.2*



Operative	1/4	1	2	6	35	169L	170L	245	330	556	557	558	701
Medium	2013203	2013205	2013201	2013211	2013221	2013231	2013241	2013251	2013261	2013271	2013281	2013291	2013301
Minor/Major Diameter	*0.7*	*0.8*	*1.1*	*1.9*	*1.0*	*0.6 / 1.1*	*0.7 / 1.2*	*0.9 / 1.1*	*0.8 / 0.9*	*0.8 / 0.9*	*0.9 / 1.0*	*1.2 / 1.3*	*1.0 / 1.4*
	H1-005	H1-008	H1-010	H1-018	H2-010	H23L-009	H23L-010	H245-008	H24-008	H31-009	H31-010	H31-012	H33-012
Head Length mm	0.7	0.6	0.8	1.3	0.9	5.2	5.8	2.7	1.6	3.3	3.9	4.1	4.1



Pear	351.3.75	352.4	360.2	361.2.5	362.3	L362.3	364.5	365.4
Coarse	2015229	2015230	2015794	2015796	2015792	2015799	2015227	2015225
Fine					2015793			2015245
Minor/Major Diameter	*1.1*	*1.3*	*0.8 / 1.0*	*1.0 / 1.2*	*1.1 / 1.3*	*1.1 / 1.3*	*1.0 / 1.4*	*1.2 / 1.5*




Amalgam Removers	1210	1220	1240	1250
Medium	2015591	2015593	2015598	2015599
Minor/Major Diameter	*0.9 / 1.3*	*0.9 / 1.6*	*0.9 / 1.1*	*0.9 / 1.3*
Head Length mm	1.3	1.4	2.0	3.0

Two Striper® MFS® Finishing Diamonds




REALITY Four Star Award

MF1 *(45µ)*
MF2 *(20µ)*
MF3 *(10µ)*



MFS® Finishing Diamonds	135	201	257	265	285	292	799
MF1 *(45µ)*	2015850	2015853	2015856	2015859	2015862	2015865	2015868
MF2 *(20µ)*	2015851	2015854	2015857	2015860	2015863	2015866	2015869
MF3 *(10µ)*	2015852	2015855	2015858	2015861	2015864	2015867	2015870
Minor/Major Diameter	*2.4*	*0.2/0.6*	*1.2/1.8*	*0.3/1.2*	*2.0*	*1.4*	*0.5/1.6*
Head Length mm	2.2	3	8	6	5	3	11

microprep™






MicroPrep™	MP89	MP38	MP30R	MP30	MP53	MP53A
Medium	2015957	2015959	2015960	2015961	2015962	2015963
Fine	2015958					
Minor/Major Diameter	*0.7*	*0.7*	*0.9*	*1.2*	*1.4*	*1.4*



Skooter™	0	1	2
Medium	**TO** - 2015188	**T1** - 2015190	**T2** - 2015192
Very Fine	**FO** - 2015189	**F1** - 2015191	**F2** - 2015193
Minor/Major Diameter	*2.3*	*2.8*	*3.0*
Reference #	905-023	905-028	905-031



S Available as a short shank diamond.
See page 10 for complete list and ordering information.

*Two Striper® BLACK diamond.


B Available in Bulk Packs.
See page 19 for complete list and ordering information.

Two Striper®
AN abrasive technology BRAND


Two Striper® Bulk Packs

Your favorite Two Striper® diamonds are available in economic packs of 50!


- ✓ Cuts faster
  – *less chair time, happier patients*
- ✓ Sharper diamonds that stay put
  – *Proprietary P.B.S.® brazing process*
- ✓ Lasts longer for multiple preps
  – *Use fewer burs than 'bargain' diamonds*



P.B.S.® Brazing Process



	245	330	557	2005.8C	230C	250.8C	251.8C	260.10VF	260.8C	260.8F	260.8VF	261.8C	283.4C	285.5C	285.5F
	Medium	Medium	Medium	Coarse	Coarse	Coarse	Coarse	Very Fine	Coarse	Fine	Very Fine	Coarse	Coarse	Coarse	Fine
5-Pack	2013251	2013261	2013281	2000058	2015063	2015067	2015069	2015185	2015123	2015153	2015183	2015124	2015142	2015143	2015173
Bulk 50-Pack	2050245	2050330	2050557	2005058	2050063	2050067	2050069	2050185	2050123	2050153	2050183	2050124	2050142	2050143	2050173
Minor/Major Diameter	*0.9 / 1.1*	*0.8 / 0.9*	*0.9 /1.0*	*1.3 /1.8*	*1.3 /3.4*	*0.4 /1.2*	*0.4 /1.2*	*0.5 /1.4*	*0.4 /1.4*	*0.4 /1.4*	*0.4 /1.4*	*0.4 /1.1*	*n-a /2.5*	*0.7 /2.3*	*0.7 /2.3*



	285.5SC	285.5VF	515.7C	702.8C	703.8C	767.7C	767.9C	770.10C	770.8C	770.8F	770.8SC	780.9C	782.8C	790.8C	863C
	Super Coarse*	Very Fine	Coarse	Coarse	Coarse	Coarse	Coarse	Coarse	Coarse	Fine	Super Coarse*	Coarse	Coarse	Coarse	Coarse
5-Pack	2015700	2015203	2015277	2015343	2015345	2015453	2015455	2015459	2015457	2015487	2015701	2015467	2015424	2015471	2015545
Bulk 50-Pack	2050700	2050203	2050277	2050343	2050345	2050453	2050455	2050459	2050457	2050487	2050701	2050467	2050424	2050471	2050545
Minor/Major Diameter	*0.7 /2.3*	*0.7 /2.3*	*n-a /1.2*	*0.9 /1.4*	*1.2 /1.8*	*1.1 /1.8*	*0.9 /1.8*	*1.1 /1.8*	*1.3 /1.8*	*1.3 /1.8*	*1.3 /1.8*	*1.4 /2.1*	*1.1 /1.5*	*2.0 /2.6*	*n-a /4.7*


Distributed by premier®

*Two Striper® BLACK diamond.

Milling Tools


Premier® Dental for all your Planmeca PlanMill diamond needs!

P.B.S.® bonding maximizes the exposure of diamond cutting surfaces - especially at the tips and upper circumference of the diamond instrument. Crystals, placed evenly and precisely in a uniform matrix are fused to a surgical grade stainless steel shank.




Premier Two Striper Tapered Milling Tool Surface - 500X

PlanMill and PlanMill 40 are trademarks of Planmeca.



Tapered Ellipsoidal Conical



Tapered 13.2 Conical 12.6

PlanMill 40

Can also be used in 30 S, 40 S, & PlanMill 35

2016000 **Tapered** (6 per box)
2016001 **Ellipsoidal** (6 per box)
2016002 **Conical** (6 per box)
2016004 **Assorted** (2 of each)

PlanMill S Mills


2017000 **Tapered 13.2** (6 per box)
2017002 **Conical 12.6** (6 per box)
2017006 **Assorted**
(1) Tapered 13.2
(1) Conical 12.6

Two Striper®
AN abrasive technology BRAND

SilverStreak™

Diamond Model Trimming Wheel

The SilverStreak™ trims models faster, runs quieter and is substantially longer lasting than traditional wheels. It has an extra coarse diamond cutting surface which is patterned for aggressive, yet cool cutting. It is light-weight, durable and perfectly balanced for concentric operation.



Thickness - 6.4mm (1/4″)
Center Hole Diameter - 25.4mm (1″)

2013401 12″ Extra Coarse
2013405 10″ Extra Coarse

Thin-Flex®

Flexible Diamond Discs with True Edge Cutting

Designed for carving ceramic and composite materials. Thin-Flex® discs are excellent for contouring and shaping all surfaces including embrasures. Thin-Flex discs have diamond crystals wrapped around the edge of the disc to avoid "black marks" during carving procedures.

Thin-Flex discs may be ordered individually or in twin-packs. Twin-packs include two discs and a mandrel. The heavy duty stainless steel HP lab mandrels feature an enlarged neck area for safer operation. The screw head and washer system guarantees true running and reduces metal fatigue.



Single-sided cross-section



Double-sided cross-section





Single-Sided - 1/pkg.
2012633 926-7 *(0.10mm, 45 µ)*
2012637 928-7 *(0.24mm, 60 µ)*

Single-Sided - 2/pkg. with Mandrel
2012634 926-7 *(0.10mm, 45 µ)*
2012638 928-7 *(0.24mm, 60 µ)*

Double-Sided - 1/pkg.
2012653 927-7 *(0.15mm, 45 µ)*
2012657 929-7 *(0.28mm, 60 µ)*

Double-Sided - 2/pkg. with Mandrel
2012654 927-7 *(0.15mm, 45 µ)*
2012658 929-7 *(0.28mm, 60 µ)*

Recommended Speed - *20,000 rpm*
Outside Diameter - *22.2mm (7/8″)*
Center Hole Diameter - *1.6mm (1/16″)*



Mandrel Refills

Heavy Duty Stainless Steel

Package of 12
2014011 HP Shank


.25mm thread length

Distributed by premier®


Laboratory HP Diamonds

P.B.S.® Bonded - ISO Type 2 Shank 2.35mm

One per pack.



	X118	X120	X135	X150	X261.5	X265.10	X285.8.5	X324.1	X350.4	X535.7	X575.7
Coarse	2012700	2012701	2012703	2012705	2012707	2012708	2012709	2012711	2012713	2012719	
Fine		2012731	2012733			2012787	2012739	2012741	2012743	2012749	2012758
Reference #	801-018	801-016	801-026	801-036	860-014	863-014	862-026	805-022	807-030	836-027	830-026
Minor/Major Diameter	*1.8*	*1.6*	*2.6*	*3.6*	*0.4 / 1.4*	*0.5 / 1.4*	*0.3 / 2.6*	*2.2*	*3.0*	*2.7*	*2.6*




	X580.6	X590.12	X701.9	X702.10	X703.12	X767.7	X769.10	X772.10	X776.9	X884	X889
Coarse	2012720	2012682	2012723	2012759	2012724	2012725	2012762	2012765	2012726	2012770	2012727
Medium		2012681									
Fine	2012750	2012680	2012753			2012755	2012730	2012736	2012756		2012757
Reference #	836-052	842R-012	847-016	848-020	848L-025	855-018	850-018	850-025	856-040	825-046	825-067
Minor/Major Diameter	*5.2*	*1.2*	*0.8 / 1.8*	*1.2 / 2.0*	*1.4 / 2.5*	*1.0 / 1.8*	*0.7 / 1.8*	*1.4 / 2.5*	*2.9 / 4.0*	*4.6*	*6.7*

TSZTECH™

Laboratory Diamonds

A New Diamond Technology for Today's Advanced Materials
- Lithium Disilicate and Zirconia

One per pack.



	X590.12	X769.10
ZF	2012848	2012849
ZM	2012850	2012851


Two Striper®
AN abrasive technology BRAND

Two Striper® Perforated Diamond Finishing Strips

Two Striper® Perforated Diamond Finishing Strips are designed for interproximal reduction, shaping, and contouring. The perforations assist in debris removal. Available in two options:

Perforated, single-sided strips, with safe center

Available in 10-packs



Perforated Strip

✓ Perforated Strip – *Features an uncoated safe center for easy access through the contact point*



2013132	300F	
2013142	300T	
2013152	300UF	
2013162	350F	
2013172	350T	
2013182	350UF	
2013192	ASST	Assorted Package of 10 – Includes (2) 300 Medium, (2) 300 Fine, (1) 300 Ultra-Fine, (2) 350 Medium, (2) 350 Fine, (1) 350 Ultra Fine

Perforated, single-sided strips, with serrated center (combo-strips)

Available in 10-packs



Combo-Strip

✓ Perforated Combo-Strip – *Includes a serrated center to facilitate access in tight spaces between teeth*



2013161	350F-S	
2013171	350T-S	
2013181	350F-S	
2013191	ASST-S	Assorted Package of 10 – Includes (3) 350 Medium, (4) 350 Fine, (3) 350 Ultra-Fine

Length: 15.25mm

Width
300 = 2.5mm
350 = 3.75mm

Grit	Thickness
T = 60µ	0.15mm
F = 45µ	0.127mm
UF = 20µ	0.1mm

Compo-Disc®

Hand-Held Diamond Finishing Instruments

For extremely tight contacts, use a Compo-Disc® to gain access. Safe-sided, edge and face-cutting.

Thickness: 0.1mm • Diameter: 22.2mm • Grit: 35µ - 45µ




2013106 Compo-Discs (Pkg. of 2 discs)

Distributed by premier®


Two Striper® Diamond Finishing Strips

Hand-Held Stainless Steel Diamond Finishing Strips

Two Striper® Diamond Finishing Strips are thin, flexible, one-sided diamond abrasive strips available in two widths and three grits. They have an uncoated safe center for easy access through the contact point and are ideal for finishing Class III, IV and veneer restorations.







Standard strips, single-sided, with safe center

Available in 6-packs

	**2013107**	100F
	**2013111**	100T
	**2013115**	100UF
	**2013109**	150T
	**2013113**	150F
	**2013117**	150UF
Assorted Package of 6	**2013127**	ASST

Length: 15.25mm

Width
100 = 2.5mm
150 = 3.75mm

Grit	Thickness
T = 60µ	0.15mm
F = 45µ	0.127mm
UF = 20µ	0.1mm

Finishing Notes

- Use Two Striper® Diamond Finishing Strips in sequence; Blue, Red, Yellow.
- Trimming (Blue) may be used separately to roughen enamel.
- Finishing in a wet field is preferred. If light staining appears, flush area and use new section of strip.

Two Striper® 200T2

For Orthodontic Reproximation

Simultaneously strip away enamel from adjacent proximal surfaces.
Flexible, Double-sided with Safe Center

Length: 15.25mm • Width: 4mm • Grit: 60µ

Standard strips, double-sided, with safe center

Available in 3-packs



2013180 Two Striper® 200T2 (Ortho) (Pkg. of 3 strips)

Faster! Super-Charged Diamond Polishing

Single-paste System for Ceramics and Composites


Diamond Twist SCL™

Extra-Oral Polishing Kit

Glaze-like, multi-surface polishing - featuring advance open-weave fibra pre-polishers and 100% wool-paste applicators - for a super-charged mirror polish.

- ✓ Pre-polish with Fibra Points (barrel, wheel or taper); open-weave design generates significantly less heat
- ✓ Final polish with Diamond Twist Super-Charged Polishing Paste applied with felt barrel, wheel or point for a glaze-like polish

REALITY Five Star Award ★★★★★ REALITY'S CHOICES



2019001
6g Diamond Twist SCL Polishing Paste (unflavored)



2019032
Fibra Wheel, Thin (6/pack)



2019033
Fibra Wheel, Thick (4/pack)



2019031
Fibra Point, Taper



2019030
Fibra Point, Barrel



2019051
SurGrip Mandrel



2019021
Felt Wheel, Thin (10/pack)



2019020
Felt Wheel, Thick (5/pack)



2019022
Felt Point



Diamond Twist SCO™

Intra-Oral Polishing Kit

After permanent seating and adjustments, polish crowns, inlays, onlays and nano-hybrid composites intra-orally in one easy step using Diamond Twist Super-Charged Polishing Paste.

- ✓ Use the flat or taper polishing brushes (autoclavable) for posterior occlusal surfaces
- ✓ Use the Micro-Cloth disks for polishing and final high-gloss buffing
- ✓ Twist syringe avoids cross-contamination and dispenses just the right amount

REALITY Four Star Award ★★★★



2019003
3g Diamond Twist SCO Polishing Paste (mint flavored)



2019005
Tapered Polishing Brushes (25/pack)



2019004
Flat Polishing Brushes (25/pack)



2019091
FG Mandrels (4/pack)
not included in kit



2019090
RA Mandrels (4/pack)



2019060
1/2" Micro-Cloth Disks Dark Blue (25/pack)



2019061
5/16" Micro-Cloth Disks Light Blue (25/pack)



Recommended for Zirconia and Lithium Disilicate

Distributed by premier®

Luminescence®

Single-Gel Diamond Polishing

Polish all surfaces in half the time!

Luminescence is a patented, single-gel diamond polishing system. In less than two minutes, you can polish the surface of all restorative materials: composite, glass ionomer, compomer, amalgam, precious metal and enamel. Even porcelain achieves a lustrous surface in only minutes due to the optimal concentration of micron-sized diamond particles.

2014200 Luminescence Intro Kit *(3g syringe of Luminescence gel, 2 RA Mandrels, 50 applicators)*

2014202 Luminescence Gel (3g)

2014204 RA Mandrels (2)


2014206 Applicators (50)




"Optimum Polishing results on the CAD/CAM Vita Mark II ceramic."[7]

REALITY Four Star Award ★★★★

"Adjustment finishing with a Two Striper MF2 or MF3 alone, or as a final finishing stage after using a more coarse finishing diamond, followed by a 60-second polishing sequence with Luminescence diamond gel, was an effective procedure for optimum polishing results on the CAD-CAM Vita Mark II ceramic."



Vita Mark II Ceramic
Pre-finish surface reference



Vita Mark II
Final surface appearance after MF2 Diamond finish and Luminescence application

Luminescence *PLUS*™

Diamond Polishing Paste for Patient Comfort + All-Surface Esthetics





- ✓ All-Surface Polishing – *composite + porcelain + enamel + glass ionomer + precious metals + amalgam*
- ✓ Single-Grit/Single-Step Paste – *saves chair-time*
- ✓ Increase Patient Comfort – *contains potassium nitrate topical desensitizer*
- ✓ Convenient – *screw-tube dispenser + improved viscosity eliminates waste*


2014208 Luminescence Plus™! Gel (3g) *(Not available in Canada)*

Poli-Pro Disks™

Composite Polishing System

Faster Polishing! / No Gouging!
15% more surface area than metal-centered disks!

- ✓ Safe-Center - protects composite surface
- ✓ Faster - complete abrasive coverage of disk surface
- ✓ Thin and flexible - conforms to tooth anatomy
- ✓ Systematic finishing and polishing - single-use color-coded disks
- ✓ Two sizes: 1/2" for greater efficiency, 5/16" for control in gingival areas



All refill packs contain 40 disks.

1/2" Coarse (12mm) **Refill #2019071**	5/16" Coarse (8mm) **Refill #2019072**	1/2" Medium (12mm) **Refill #2019073**	1/2" Medium (12mm) **Refill #2019074**	5/16" Medium (8mm) **Refill #2019075**
1/2" Fine (12mm) **Refill #2019076**	5/16" Fine (8mm) **Refill #2019077**	1/2" Extra-Fine (12mm) **Refill #2019078**	5/16" Extra-Fine (8mm) **Refill #2019079**	RA Mandrels (4/pack) **Refill #2019090**
				FG Mandrels (4/pack) **Refill #2019091**

Coarse - 80 microns (grit up)

Medium - 40 microns (grit up)

Medium - 40 microns (grit down)

Fine - 30 microns (grit both sides)

Extra-Fine - 9 microns (grit both sides)



2019070 Intro Kit
(180 disks, 4 mandrels)

2019099 Starter Kit
(35 disks, 1 mandrel)

Priva Rotary Instruments®

Diamond Burs

- ✓ SmartEdge® Technology - *Intentionally positioned, uniformly shaped engineered diamond crystals maximize diamond-to-tooth contact*
- ✓ Fast Cutting - *Outstanding procedural efficiency*
- ✓ Reduced Vibration - *Greater control*
- ✓ Helps Prevent Diamond Dislodgement – *Provides consistent performance*
- ✓ Efficient Tooth Surface Preparation - *Creates a clean, consistent finish*

NEW!






SmartEdge® TECHNOLOGY

Priva Rotary Instruments® Diamond Burs, featuring SmartEdge® Technology, are high-performance electroplated diamond burs designed to enhance your rotary portfolio.


SmartEdge® optimizes rotational accuracy by creating an innovative cutting surface of intentionally positioned, uniformly shaped engineered diamond crystals – maximizing diamond-to-tooth contact. This technology helps to enable fast cutting for outstanding procedural efficiency and to reduce vibration (chatter) providing greater control. Additionally, SmartEdge® helps provide efficient tooth surface preparation, and prevent diamond dislodgement.



Maximize diamond-to-tooth contact - helps deliver efficient tooth surface preparation!




Recommended Crown Prep Setups




Chamfer Prep

- ✓ Round-end tapered or cylindrical diamonds create a chamfered finish line
- ✓ Recommended prep design for full coverage crowns and veneers
- ✓ Recommended prep design for PFM crowns and some ceramic crowns




Axial Reduction	P856016C*					P850014C**
Margin Refinement		P856016F*				P850014F**
Contact Break/Refinement					P859014C**	P862012F*
Occlusal Reduction			P379023C	P379023F		




Shoulder Prep

- ✓ Choose flat-end tapered or flat-end cylindrical diamonds to create shoulder preparations with a 90˚ margin




Axial Reduction	P847016C*				
Margin Refinement		P847016M*			
Contact Break/Refinement				P859014C**	P862012F*
Occlusal Reduction			P811037C		



Modified Shoulder Prep

- ✓ Creates a margin with rounded internal angles
- ✓ Flat-end tapered or parallel diamond instruments featuring a rounded (beveled) edge



Axial Reduction	P847KR016C*					
Margin Refinement		P856012F*				
Contact Break/Refinement					P859014C**	P862012F*
Occlusal Reduction			P811033C	P379023F		

* 8.0mm Abrasive Length / ** 10.0mm Abrasive Length



	245-009	330-008	368-016	368-018	368-023	369-025	379-016	379-018	379-023	392-016	557-010
Extra Coarse					P368023XC			P379018XC	P379023XC		
Coarse			P368016C	P368018C	P368023C	P369025C	P379016C	P379018C	P379023C		
Medium	P245009M	P330008M		P368018M	P368023M		P379016M	P379018M	P379023M		P557010M
Fine			P368016F	P368018F	P368023F		P379016F	P379018F	P379023F	P392016F	
Extra Fine				P368018XF	P368023XF			P379018XF	P379023XF		
Min.D/Maj.D/Length (mm)	– / 0.9 / 2.7	– / 0.8 / 1.6	– / 1.6 / 3.0	– / 1.8 / 4.5	– / 2.3 / 5.0	– / 2.5 / 5.5	– / 1.6 / 3.4	– / 1.8 / 3.4	– / 2.3 / 4.2	0.3 / 1.6 / 5.0	– / 1.0 / 4.2



	801-010	801-012	801-014	801-016	801-018	801-021	801-023	811-033	811-037	830-009	830-010
Extra Coarse			P801014XC	P801016XC	P801018XC			P811033XC	P811037XC		
Coarse		P801012C	P801014C	P801016C	P801018C	P801021C	P801023C	P811033C	P811037C		P830010C
Medium	P801010M	P801012M	P801014M	P801016M	P801018M		P801023M			P830009M	P830010M
Fine				P801016F	P801018F		P801023F				
Min.D/Maj.D/Length (mm)	– / 1.0 / –	– / 1.2 / –	– / 1.4 / –	– / 1.6 / –	– / 1.8 / –	– / 2.1 / –	– / 2.3 / –	– / 3.3 / 4.3	– / 3.7 / 7.0	– / 0.9 / 2.7	– / 1.0 / 2.7



	830-012	830L-012	830L-014	835-008	835-010	835-012	835KR-010	835KR-012	836-012	836-014	836KR-012
Extra Coarse					P835010XC	P835012XC					
Coarse	P830012C	P830L012C	P830L014C		P835010C	P835012C	P835KR010C	P835KR012C	P836012C	P836014C	P836KR012C
Medium	P830012M			P835008M	P835010M	P835012M					
Min.D/Maj.D/Length (mm)	– / 1.2 / 2.7	– / 1.2 / 4.0	– / 1.4 / 5.0	– / 0.8 / 3.0	– / 1.0 / 4.0	– / 1.2 / 4.0	– / 1.0 / 4.0	– / 1.2 / 4.0	– / 1.2 / 6.0	– / 1.4 / 6.0	– / 1.2 / 6.0




	837-012	837-014	837-018	837KR-012	846KR-016	847-014	847-016	847-018	847KR-016	848-016	848-018
Extra Coarse	P837012XC	P837014XC							P847KR016XC	P848016XC	P848018XC
Coarse	P837012C	P837014C	P837018C	P837KR012C	P846KR016C	P847014C	P847016C	P847018C	P847KR016C	P848016C	P848018C
Medium		P837014M					P847016M				
Fine	P837012F										
Min.D/Maj.D/Length (mm)	– / 1.2 / 8.0	– / 1.4 / 8.0	– / 1.8 / 8.0	– / 1.2 / 8.0	1.1 / 1.6 / 6.0	0.9 / 1.4 / 8.0	1.0 / 1.6 / 8.0	1.3 / 1.8 / 8.0	1.1 / 1.6 / 8.0	0.9 / 1.6 / 10.0	1.1 / 1.8 / 10.0




	850-012	850-014	850-016	850-018	852-010	852-012	855-012	855-014	855-016	856-012	856-014
Extra Coarse		P850014XC	P850016XC	P850018XC			P855012XC	P855014XC		P856012XC	P856014XC
Coarse	P850012C	P850014C	P850016C	P850018C		P852012C	P855012C	P855014C	P855016C	P856012C	P856014C
Medium	P850012M	P850014M	P850016M	P850018M	P852010M					P856012M	P856014M
Fine	P850012F	P850014F	P850016F	P850018F	P852010F	P852012F				P856012F	P856014F
Min.D/Maj.D/Length (mm)	0.6 / 1.2 / 10.0	0.7 / 1.4 / 10.0	0.9 / 1.6 / 10.0	1.2/ 1.8 / 10.0	0.5 / 1.0 / 4.4	0.5 / 1.2 / 6.0	0.8 / 1.2 / 6.0	1.0 / 1.4 / 6.0	1.3 / 1.6 / 6.0	0.8 / 1.2 / 8.0	0.9 / 1.4 / 8.0



	856-016	856-018	856-021	856-025	856L-014	856L-016	856L-018	856L-020	858-010	858-012	858-014
Extra Coarse	P856016XC	P856018XC	P856021XC	P856025XC	P856L014XC	P856L016XC	P856L018XC	P856L020XC			P858014XC
Coarse	P856016C	P856018C	P856021C	P856025C	P856L014C	P856L016C	P856L018C	P856L020C	P858010C	P858012C	P858014C
Medium	P856016M					P856L016M					
Fine	P856016F	P856018F				P856L016F			P858010F	P858012F	P858014F
Min.D/Maj.D/Length (mm)	1.1 / 1.6 / 8.0	1.3 / 1.8 / 8.0	1.6 / 2.1 / 8.0	2.0 / 2.5 / 8.0	0.8 / 1.4 / 9.0	1.0 / 1.6 / 9.0	1.2 / 1.8 / 9.0	1.4 / 2.0 / 9.0	0.5 / 1.0 / 8.0	0.5 / 1.2 / 8.0	0.5 / 1.4 / 8.0



	858-018	859-010	859-012	859-014	862-010	862-012	862-014	862-016
Extra Coarse						P862012XC	P862014XC	P862016XC
Coarse	P858018C	P859010C	P859012C	P859014C	P862010C	P862012C		P862016C
Medium		P859010M	P859012M	P859014M		P862012M	P862014M	
Fine		P859010F	P859012F	P859014F	P862010F	P862012F		
Extra Fine		P859010XF	P859012XF	P859014XF				
Min.D/Maj.D/Length (mm)	0.5 / 1.8 / 8.0	0.5 / 1.0/ 10.0	0.5 / 1.2 / 10.0	0.5 / 1.4 / 10.0	0.5 / 1.0 / 8.0	0.5/ 1.2 / 8.0	0.5 / 1.4 / 8.0	0.5 / 1.6 / 8.0



	863-010	863-012	878-012	878-014	878K-012	878K-014	879-014	879-016
Extra Coarse						P878K014XC	P879014XC	P879016XC
Coarse		P863012C	P878012C	P878014C	P878K012C	P878K014C	P879014C	P879016C
Fine	P863010F	P863012F						
Extra Fine		P863012XF						
Min.D/Maj.D/Length (mm)	0.5 / 1.0 / 10.0	0.5/ 1.2 / 10.0	– / 1.2 / 8.0	– / 1.4 / 8.0	0.7 / 1.2 / 8.0	0.9 / 1.4 / 8.0	– / 1.4 / 10.0	– / 1.6 / 10.0



	880-012	880-014	881-012	888-012	889-009	909-037	909-040
Extra Coarse						P909037XC	P909040XC
Coarse	P880012C	P880014C	P881012C	P888012C		P909037C	P909040C
Fine				P888012F	P889009F		
Min.D/Maj.D/Length (mm)	– / 1.2 / 6.0	– / 1.4 / 6.0	– / 1.2 / 8.0	0.5 / 1.2 / 8.0	0.4 / 0.9 / 3.5	– / 3.7 / 1.3	– / 4.0 / 1.7



Available in 5-packs

Solo Diamond®

Reliable Results. Pre-Sterilized.

Now that's Efficiency!

Solo Diamonds are engineered to be sharp and provide reliable results. What's more, each bur is pre-sterilized for efficiency & reduced cross-contamination risk.

- ✓ Sharp, high-quality diamonds
- ✓ Pre-sterilized = convenience + peace of mind
- ✓ Simplifies & eliminates unnecessary bur blocks
- ✓ Wide variety of clinical shapes & diamond grits
- ✓ 5-year pouch seal integrity validation






Instrument Code		
Coarse	151µ	862012 C
Medium	126µ	862012 M
Fine	64µ	862012 F
Extra Fine	46µ	862012 EF
Diameter/Length/Tip Diameter		*in mm*
Shape Number		1512.8
S - short-shank		862012C S

Price Levels
Standard
Large Specialty Finishing
Spiral


Pricing available from authorized dealer

Modified KR Taper




Coarse	846KR016C	845KR018C	847KR012C	847KR018C
Medium				
Fine	846KR016F		847KR012F	847KR018F
D/L/TD	*1.6/8.0/1.1*	*1.8/4.0/1.4*	*1.2/8.0/0.09*	*1.8/8.0/1.35*
Shape #	0816.8	0818.4	0812.8	0818.8

Flat-End Taper




Coarse	847014C	847L014C	847016C	848016C	847018C	848018C
Medium					847018M	
Fine	847014F		847016F	848016F		
D/L/TD	*1.4/8/1.0*	*1.4/9/0.7*	*1.6/8/1.2*	*1.6/10/1.1*	*1.8/8/1.4*	*1.8/10/1.4*
Shape #	0914.8	0914.9	0916.8	0916.10	0918.8	0918.10

Pointed Cone




Coarse	858014C	859014C
Medium		859014M
Fine	858014F	859014F
D/L	*1.4/8.0*	*1.4/10.0*
Shape #	1314.8	1314.10

Flame




Coarse	862010C	860012CS	862012C	863012C	862016C	862L018C
Medium	862010M		862012M			
Fine	862010F		862012F	863012F	862016F	862L018F
D/L	*1.0/8*	*1.2/6.0*	*1.2/8.0*	*1.2/10.0*	*1.6/8.0*	*1.8/9.0*
Shape #	1510.8	1512.6	1512.8	1512.10	1516.8	1518.9

Pear



Coarse			
Medium	830008MS	830008M	830012M
Fine			
D/L	*0.8/1.6*	*0.8/1.6*	*1.2/2.7*
Shape #	0508	0508	0512

Egg




Coarse	379018C	379023C
Medium		
Fine	379018F	379023F
D/L	*1.8/3.4*	*2.3/4.2*
Shape #	1908	1900

Football




Coarse	368016C	368020C	368023C
Medium		368020M	368023M
Fine	368016F	368020F	368023F
D/L	*1.6/3.2*	*2.0/5.0*	*2.3/5.2*
Shape #	1916	1920	1923

Spiral




Coarse	837T016C	856T014C	856T016C	850T016C	856T018C	856LT018C	368T023C
D/L/TD	*1.6/8*	*1.4/8/0.9*	*1.6/8/1.1*	*1.6/10/1.1*	*1.8/7/1.4*	*1.8/9/1.4*	*2.3/5.2*
Shape #	8716.8	8114.8	8116.8	8116.10	8118.7	8118.9	8923

Round




Coarse			801016C	801018C	801023C
Medium	801010M	801012M	801016M	801018M	801023M
Fine					801023F
D/L	*1.0/0.8*	*1.2/1.0*	*1.6/1.4*	*1.8/1.6*	*2.3/2.2*
Shape #	0110	0112	0116	0118	0123

Inverted Cone



Coarse		805016C
Medium	805012M	
Fine		
D/L	*1.2/0.9*	*1.6/1.4*
Shape #	0312	0316

Flat-End Cylinder



Coarse		835010C	837012C	837016C
Medium	835010MS	835010M		
Fine			837012F	
D/L	*1.0/4.0*	*1.0/4.0*	*1.2/7.0*	*1.6/8.0*
Shape #	0710MS	0710	0712.7	0716.8

Interproximal




Coarse	
Medium	132008M
Fine	132008F
D/L/TD	*0.8/3.0/0.3*
Shape #	1310.3

Curettage K




Coarse	878K012C	879K012C	878K014C	879K014C	878K016C	878K018C	879K018C
Medium				879K014M			879K018M
Fine	878K012F		878K014F				
D/L/TD	*1.2/8/0.8*	*1.2/10.0/0.8*	*1.4/8.0/1.0*	*1.4/10.0/1.1*	*1.6/8.0/1.2*	*1.8/8.0/1.4*	*1.8/10/1.4*
Shape #	1712.8	1712.10	1714.8	1714.10	1716.8	1718.8	1718.10

Wheel




Coarse	909037C	909040C
D/L	*3.5/1.3*	*4.2/1.5*
Shape #	2035	2042

Barrel



Coarse	811033C	811037C
D/L	*3.3/4.4*	*3.7/7.0*
Shape #	2133	2137

Gross Reduction



Coarse	6051017C	6055018C
D/L/TD	*1.7/8.0*	*1.8/8.0/1.3*
Shape #	2217	2218

Round-End Taper



Coarse	856025C
D/L/TD	*2.4/7/1.8*
Shape #	2424

Large Inverted Cone




Coarse	807018C
Medium	
Fine	
D/L	*1.8/5.0*
Shape #	0318.5

Double Inverted Cone




Coarse	
Medium	813014M
Fine	
D/L	*1.4/2.2*
Shape #	0416

Operative




Coarse			
Medium	H245008M	H330008M	H557010M
Fine			
D/L/TD	*1.1/2.8/0.9*	*0.9/1.6/0.8*	*1.0/4.2/0.9*
Shape #	FG245	FG330	FG557

Pit & Fissure




Coarse	
Medium	
Fine	8889MFS
D/L/TD	*0.8/3.0/0.4*
Shape #	1300FS

Round-End Taper




Coarse	856012C	850012C	856014C	850014C	855016CS	856015C	856016C	850016C	856018C	856L018C	856021C		
Medium			856014M	850014M			856016M	850016M		856L018M		850018M	
Fine			856014F	850014F			856016F	850016F					
D/L/TD	*1.2/7/0.7*	*1.2/10/0.7*	*1.4/8/0.9*	*1.4/10/0.9*	*1.6/6/1.1*	*1.6/8/1.1*	*1.8/8/1.1*	*1.6/10/1.1*	*1.8/7/1.2*	*1.8/9/1.2*	*2.1/8/1.4*	*1.8/10/1.2*	*2.2/8.5/1.8*
Shape #	1112.7	1112.10	1114.8	1114.10	1116.6	1115.8	1116.8	1116.10	1118.7	1118.9	1121.8	1118.10	1122.8

KS Cylinder




Coarse	KS0010C	KS1012C	KS2014C
Medium			
Fine		KS1012F	
D/L	*1.0/7.0*	*1.2/7.0*	*1.4/8.0*
Shape #	1210.7	1212.7	1214.8

Modified KR Cylinder




Coarse		837KR018C
Medium		
Fine		
D/L	*1.4/8.0*	*1.8/8.0*
Shape #	0614.8	0618.8

Modified Beveled Cylinder



Coarse	878012C	878014C
Medium		
Fine		
D/L	*1.2/8.0*	*1.4/8.0*
Shape #	1800.8	1804.8

Beveled Cylinder



Coarse	885012C	886012C
Medium		
Fine		
D/L	*1.2/8.0*	*1.2/10.0*
Shape #	1812.8	1812.10

Modified KR Taper Large




Coarse	845KR025C
D/L/TD	*2.5/4.0/2.0*
Shape #	2525

Acorn



Fine	905028F
D/L	*2.8/4.0*
Shape #	3028

Finishing



Extra Fine	132008EF	134014EF	135014EF	862012EF	379023EF	368023EF
D/L	*0.8/3.0*	*1.4/8*	*1.4/10.0*	*1.2/8.0*	*2.3/4.2*	*2.3/5.2*
Shape #	3310.3	3314.8	3314.10	3512.8	3900	3923



The Diamond Difference

Two Striper® | Priva Rotary Instruments® | Solo Diamond®

Find my Bur!

Easily convert many diamond brands to
Two Striper® or Solo Diamond®.
To use this versatile web-based tool, go to

findmybur.com

or scan the QR code.



premier® Inspired solutions for daily dentistry®

Premier® Dental Products Company • 1710 Romano Drive, Plymouth Meeting, PA 19462, U.S.A.
Phone 888-670-6100 • 610-239-6000 • Fax 610-239-6171 • premierdentalco.com • info@premierdentalco.com

® Premier is a registered trademark of Premier Dental Products Company

©2026 Premier Dental Co. | All Rights Reserved.

8621257 0826024 Rev0 NKP5M